


JeioTech

2025 - 2026 Catalog

Build a New Lab!

Lab Total Solution



Jeio Tech – Your trusted Lab Companion, providing solutions for scientific success. Dedicated to supporting researchers in achieving meaningful results.

Pursuing innovation in research and manufacturing technology,
we aim to provide a safe and reliable laboratory environment.

To Whom

Anyone who makes
the world a better place.



What

Ensure a safe and reliable
laboratory environment.



How

Provide high quality
products and services.



Your Lab Companion, Jeio Tech





Advancing Laboratories with High-Quality Equipment

- Manufacture of laboratory equipment in accordance with the IEC (International Electrotechnical Commission, such as CE, UL, etc.)
- Laboratory furniture compliant with European Standards (Fume Hood, Biosafety Cabinet, Fire Safety Cabinet, Acid/Corrosive Cabinet)
- Utilization of CFC-free refrigerants to reduce global warming (RoHS, Eco-Friendly Refrigerants)
- Jeio Tech supports researchers in the scientific field by donating laboratory equipment (Seeding Labs, Supporting Lab-equipment in developing countries)



2024-2025 NEW PRODUCTS

Biosafety Cabinet, Class II A2 type

For further details
Page 196

EN 12469 certified! Safe from biological hazards

- Protect user, samples, and environment
- Equipped with quiet and durable ECM motor
- Smart technology enhances user convenience



Watch introduction video

JB-15A




EN 12469 certified, a secure lab environment




Microbial Challenge Test Completed

Low noise level of 55 dB like in office




Ultra-quiet, measuring 12 dB quieter over 5 tests
(67 dB: noise of a car, 55 dB: office noise level)

Durable and energy-efficient ECM motor



Comparison of the motor control depending on voltage change (AC vs ECM)

Smart technology enhances user convenience



Using LC Bluetooth helps to prevent crosscontamination (Option)

Natural Convection Oven

For further details
Page 74

Enhanced temperature uniformity by Multi-path airflow



- > Heat transfers in multiple directions, enhancing temperature distribution
- > Natural convection with a wide temperature range: Amb.+15 to 250°C (Amb.+27 to 482°F)
- > Two models available: Advanced (ON4-V) and Standard (ON-S)
- > Available per model: 4 volume sizes and with/without window

30L model for use inside a fume hood



Effectively removes fumes through upper vent hole

FDA 21 CFR Part 11 Compliance



Integrated with 21 CFR Part 11 compliant software (LC DataKeeper), suitable for GMP audits (option)

Space efficiency by double stacking



Double stackable with optional stacking set (-05 / -15 models only)

Smart technology enhances user convenience



Real-time monitoring and control via mobile app (option)

2024-2025 NEW PRODUCTS

Laboratory Incubator

For further details
Page 90

4 models variations

IB4



> Convection Options

- Forced (Advanced: IB4-V, Standard: IB4-S)
- Air-jacket (Advanced: IB4-AV, Standard: IB4-AS)


> Capacity available in 4 sizes: Approx. 25L, 60L, 110L, and 150L

> Temperature range: Amb. +5 to 80°C (Amb. +9 to 176°F)

Microbial incubation is easy and convenient

Difficult to control incubation time due to dry media?

- Opti-flow™ adjusts the fan speed, thereby minimizing evaporation and increasing moisture retention time of the culture media
- Optimizing wind speed at level 2 minimizes media drying within two days



Intuitive Touchscreen Controller



Easy to operate with a high-visibility touchscreen controller

Data check in real-time from a distance



Real-time monitoring and control via mobile app (option)

FDA 21 CFR Part 11 Compliance



Integrated with 21 CFR Part 11 compliant software (LC DataKeeper), suitable for GMP audits (option)

Shaking Heating Bath

For further details
Page 60


More precise and stable shaking and temperature control

BS2-10




- > Precise rpm (Max. 200rpm) and temperature (Max. 100°C) control with superior accuracy
- > Excellent temperature variation due to the shaking motion ($50^{\circ}\text{C} \pm 0.1^{\circ}\text{C}$)
- > Minimized vibration under high loads by stable shaking ability
- > Extended height for larger sample volumes (Max. 2.8L flask)

Transparent gable-roof shape cover



Transparent cover allows easy observation and the gable shape prevents condensed water from falling onto the specimen

3-point calibration



Calibration maintains constant performance

Spring wire rack



Excellent corrosion resistance as made from stainless steel (1ea: standard)

Real-time mobile monitoring




Monitor and control the temperature and rpm via mobile (option)

LC DataKeeper



Data integrity software compliant with the FDA 21 CFR Part 11



A solution to prepare for GMP audits

LC DataKeeper is a Jeio Tech's proprietary software that complies with FDA 21 CFR Part 11 regulations following the ALCOA Plus guidelines. It supports **digital records, digital signs, and user access management**, including audit trails for the equipment.

Jeio Tech's equipment can be connected to LC DataKeeper to digitally record all experimental data. Only data approved by administrators through digital signatures can be recorded and exported, and permissions can be set differently by user. In addition, data can be viewed in real-time and automatic backup prevents data loss.



Screenshot of LC DataKeeper connected to equipment



Digitally signed validation document

Key factors for data integrity



User-specific Access level restrictions

Able to grant different levels of permission to access software functions by registered user.



User-specific data log recording

Support for audit trails and recording equipment usage history on a per-user basis.



Application of encryption modules for data security

Application of AES-256 two-way encryption method to prevent tampering with the original data.



Support for digital signatures by the chief administrator

The final approver can add notes and provide digital signatures for actions such as recording, deleting, and exporting research data.



Real-time data viewer (Graphs supported)

Real-time graphs, audit trails, event logs, etc. available for user-generated research data.



Support for data server or HDD backups

Automatic/manual backups to the PC to prevent loss of recorded data.

Audit trails through digital signatures and issuing reports

- Permission for equipment connection is granted only to administrators.
- Thorough audit trails supported by always requiring digital signatures for creation, recording, and conclusion of any experiment on connected equipment.
- Reports are issued in the form of electronic documents containing digital signatures for a concluded experiment to ensure reliability.



Trustworthy data with digital signature

Thorough security to prevent forgery and tampering

- Preventing third-party access by implementing AES-256 two-way encryption, the advanced encryption standard.
- Setting and assigning permissions for program function access, and maintaining transparent process quality by verifying user-specific logs.

User Permission	ADMIN	EDITOR	AUTHOR
Add/Delete equipment	<input type="radio"/>	-	-
System configuration	<input type="radio"/>	-	-
System backup and recovery	<input type="radio"/>	-	-
User management	<input type="radio"/>	-	-
Start/Stop equipment recording	<input type="radio"/>	<input type="radio"/>	-
Create/Delete experiment	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Start/Stop experiment	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Data export	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Digital signature	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Settings	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Version info, License management	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>


Different permission granted by user type

Data backup and recovery

- Setting an automatic backup interval to ensure the recorded data is preserved without loss.
- Basic version: manual backup, Pro version: automatic backup.

Convenient equipment connection and experimentation

- Capable of managing up to 100 connected equipment for the Pro version.
- Intuitively grasp the experiment progress through real-time graphs.
- More precise data analysis through zoom-in/out function.
- Graph extraction as a PNG image file when completing data through experiment initiation and conclusion.
- Data Storage Interval Setting (seconds)
 - You can set the interval at which data is recorded within the range of 1 to 999 seconds.
 - This setting is only available for Admin accounts.



Intuitive monitoring of experiment progress with real-time graphs

LC DataKeeper supported Equipment

We will continue to develop and expand the types of equipment that can be connected.




Specification

License	Basic	Pro
Duration of use	permanent	permanent
Max. number of equipment to connect	4 ea	100 ea
Manual backup and recovery	<input type="radio"/>	<input type="radio"/>
Data storage	<input type="radio"/>	<input type="radio"/>
Event/log management	<input type="radio"/>	<input type="radio"/>
Warning record function	<input type="radio"/>	<input type="radio"/>
CSV file extraction	<input type="radio"/>	<input type="radio"/>
User management	<input type="radio"/>	<input type="radio"/>
PDF Extraction for digital signatures	-	<input type="radio"/>
Automatic backup and recovery	-	<input type="radio"/>
Validation document provision	-	<input type="radio"/>
Cat. No.	SSA11001	SSA11002

※ The trial version can be used for 30 days and offers the same as the Basic version.

LC Connected (Mobile monitoring system)



LC Connected

Wherever! Whenever! Real-time monitoring and remote control of lab equipment anytime, anywhere

A combination of research/lab and network! Efficient management of lab research equipment in the same network

LC Connected is a smart mobile monitoring system that is introduced into modern laboratories for the IoT era when everything is connected, and new value is created.



[View the introduction video](#)

Whether you are on a business trip, dining out, or in a meeting, you will receive a notification to your mobile device in real-time if any error of the equipment is detected.

In addition, temperature / humidity / rpm control and Run / Stop functions can be controlled at any time, and from any location. When you also utilize for the data backup service, you can manage the stored data precisely on the central server through the LC Connected.




Reliable data management through backup of equipment history.

LC Connected System

LC GreenBox is a device that connects the lab's research equipment to a single network, allowing the experimenter to operate the equipment more efficiently. Experimenters can check the condition of their devices anytime, from any location via the mobile app. Each piece of equipment connected to the LC GreenBox creates an LC Connected system by sending data to the Cloud Server and transmitting information to the experimenter's mobile app.




LC GreenBox



LC Connected Flow diagram

- Real-time monitoring of equipment connected to LC Connected possible through the mobile application.
(Supported from Android 6.0 Marshmallow version and above. iOS 12.0 and above)
- Each lab member can access connected devices for which they have access permission through their personal mobile app.
- One LC GreenBox can communicate with up to four devices via an RS-232 port and can be expanded via multiple LC GreenBoxes.
- Access to the Cloud Server via Wi-Fi or wired LAN allows for a high degree of freedom in the lab.
- Real-time alarms for events concerning connected equipment.
(Temperature deviation, long-term door opening, power failure, and etc.)



Real-time monitoring through mobile devices




Connection up to 4 devices

※ LC Connected applied product list

Oven (OF4, OV4, OF3, OF, OV), Bath (RW3, CW3, VB), Incubated Shaker (ISS, IST, ISF), Shaker (OS), Incubator (IL3, IB, IL, ILP), Refrigerator (CLG3, PSR3), Power Supply (EPS), Temp. & Humid. Chamber (TH3), Heating & Cooling Chamber (TC3)

Description	Dimension (W x D x H, mm/inch)	Cat. No.
LC GreenBox	156 x 94 x 34 / 6.1 x 3.7 x 1.3	AAHQ1011K



Power failure notification

Highlight Products

Forced Convection Oven (OF4)

For further details
Page 66

Suitable for various applications such as drying laboratory instruments, dry heat sterilization, heat resistance testing, and burn-in tests across fields like biotechnology, chemistry, and biology. Operating from amb. +10°C to a max. 250°C and forced convection with adjustable fan speed to accommodate different samples and purposes. It enables optimum temperature uniformity throughout the chamber.



Vacuum Oven (OV4)

For further details
Page 78

Reliable and excellent temperature uniformity with variation less than $\pm 1.5^\circ\text{C}$ at 100°C and 5-inch TFT touch controller that enables intuitive control of the device.



Refrigerated Heating Bath Circulator (RW3)

For further details
Page 28

It offers precise and reliable temperature control within an extensive range of -35°C to 150°C. Ideal for experiments requiring minimal temperature deviation. Available in programmable and standard models, we offer eight variations for each type, tailored to different temperature ranges and capacities.



Incubator (IB4)

For further details
Page 90-96

Forced convection
with excellent temperature uniformity maintained through forced convection using an internal fan.

Natural convection
that minimizes scattering and contamination of samples using air-jacket type heat transfer function.



Incubated Shaker, Stackable

For further details
Page 114

The height of platform of the upper portion of triple stackable units is only 150 cm, and thus is convenient to remove and attach samples by pulling out the platform easily using a rail mechanism.



Temperature & Humidity Chamber (TH3)

For further details
Page 286

It is equipped with a dedicated controller that enhances user convenience, and the improved water supply line makes it possible to replenish water at the optimum time, making it excellent environmental testing equipment.



Safety Storage & Cabinet

For further details
Page 220

Storage of Flammable / Combustible Materials

In case of an unexpected fire that may occur at any time, flammable and combustible materials must be stored in a suitable location. Jeio Tech's flammable hazardous material storage cabinets to be certified the strictest EN 14470-1 standard, is the best choice.



Storage of corrosive substances such as strong acids and strong bases

The structure is optimized for long-term safe storage of corrosive substances such as strong acids and strong bases used in experiments and processes. It is composed of a solid steel structure and corrosion-resistant interior material, and provides differentiated safety performance due to its space specific exhaust structure.



Reagent cabinet with built-in circulation filter

It is a reagent cabinet in which air flow is uniformly formed and the harmful gas is adsorbed through the internal filter.

It is safe from fine particles and gas as it is equipped with Jeio Tech's exclusive carbon filter and HEPA filter.



Physically / Chemically sturdy cabinet

It is a reagent cabinet consisting of a polypropylene drawer and a stable structure with a rigid steel frame coated with epoxy. It is a window door type cabinet available to connect the extraction module.



Safe storage of high-pressure cylinders

Cylindrical high-pressure containers must be prevented from falling and must be properly protected from unexpected external impact. High-pressure cylinders can be safely stored in Jeio Tech's exclusive solid steel storage cabinets for peace of mind.



Contents

JEIOtech Research / Laboratory Equipment

Laboratory Equipment

— Bath Circulator

- 28 Refrigerated & Heating bath Circulator (Program, General)
- 32 Heating Bath Circulator (Program, General)

— Chiller (Recirculating Cooler)

- 40 Chiller (Recirculating Cooler), General
- 42 Chiller (Recirculating Cooler), Low temp. general
- 44 Chiller (Recirculating Cooler), Compact
- 46 Chiller (Recirculating Cooler), High temp. precision

— Heating Bath

- 54 Heating Bath (General, Agitator)
- 57 Cold Trap Bath
- 58 Heating Bath, Visible
- 60 Heating Bath, Shaking

— Oven & Furnace

- 66 Forced Convection Oven (Program, Advanced)
- 70 Forced Convection Oven, Large
- 72 Forced Convection Oven, Standard
- 73 Forced Convection Oven, Multi-Chamber
- 74 Natural Convection Oven (Advanced, Standard)
- 78 Vacuum Oven
- 81 Vacuum Oven System Station
- 82 Muffle Furnace

— Incubator

- 90 Incubator (Air-Jacket), Advanced
- 92 Incubator (Air-Jacket), Standard
- 94 Incubator (Forced Convection), Advanced
- 96 Incubator (Forced Convection), Standard
- 98 Incubator (Forced Convection), Mini
- 99 Incubator (Forced Convection), Multi-Chamber
- 100 Low Temp. Incubator, General
- 102 Low Temp. Incubator (Forced Convection), Personal
- 103 Low Temp. Incubator (Forced Convection), Multi Chamber
- 104 Plant Growth Chamber

— Incubated Shaker

- 112 Incubated Shaker, Stackable
- 114 Incubated Shaker, Stackable (Multi)
- 116 Incubated Shaker, Tabletop
- 118 Incubated Shaker, Floor

— Shaker

- 126 Shaker (Dual Action, General, Program)
- 130 Shaker, Compact
- 131 Microplate Shaker
- 132 Rocking / Waving Shaker
- 133 Funnel Shaker
- 134 Rotator

— Stirrer & Mixer

- 142 Hotplate & Magnetic Stirrer (Digital, Analog, Multi Type)
- 145 Hotplate, Digital
- 146 Round Flask & Vial Heating Block
- 148 Magnetic Stirrer (Digital, Analog, Multi Type)
- 158 Overhead Stirrer (Prestige Touch, High Performance, Easy Control)
- 168 Vortex Mixer

— Life Science Equipment

- 172 Heating Block
- 174 Heating & Cooling Block
- 175 Heating Shaker
- 178 Horizontal Electrophoresis System
- 180 Vertical Electrophoresis System
- 182 Power Supply for Electrophoresis System

— Autoclave

- 186 Autoclave

— Pump

- 188 Pump

— Ultrasonic Cleaner

- 190 Ultrasonic Cleaner

Biosafety Cabinet & Clean Bench, Fume Hood

— Biosafety Cabinet & Clean Bench

- 196 Biosafety Cabinet
- 200 Clean Bench (Advanced, Basic)
- 203 PCR Workstation
- 204 UV Sterilization Cabinet

— Fume Hood

- 208 Ductless Fume Hood
- 210 Fume Hood, Small
- 211 Molded Fume Hood, Small type (for observation)
- 212 Mobile Fume Extractor, Local exhaust type
- 213 Mobile Fume Extractor, Air purifying type
- 216 Extraction Arm Hood



Safety Storage & Cabinet

— Safety Storage & Cabinet

- 226 Fire Safety Storage Cabinet (Type 30, 90)
- 234 Safety Waste Disposal System
- 235 Acid/Corrosive Cabinet
- 240 Filtering Storage Cabinet
- 243 Gas Cylinder Cabinet
- 244 Multipurpose Cabinet

— Dry Cabinet & Desiccator

- 248 Dry Cabinet (Auto desiccator)
- 250 Vacuum Desiccator (Cubic, Round, Cylinder)
- 252 Desiccator

Refrigerator & Freezer for Laboratory/ Medical

— Laboratory Refrigerator & Freezer

- 260 Laboratory Refrigerator
- 262 Laboratory Freezer

— Medical Refrigerator

- 264 Pharmaceutical Refrigerator

Lab Essentials

- 270 Plastic Ware (Liquid Handling, Storage, Connection)

JEIOtech Environmental Reliability Test Equipment

Environmental Reliability Test Equipment

— Temp. & Humid. Chamber

- 286 Temp. & Humid Chamber, Tabletop type (-40°C, -20°C, 0°C)
- 290 Temp. & Humid. Chamber, Ultra and Low temp. type (-70°C, -40°C)
- 295 Temp. & Humid. Chamber, General type (-5°C)
- 299 Temp. & Humid. Chamber, Dual type
- 302 Temp. & Humid. Chamber, Basic type
- 305 Pharmaceutical Stability Test Chamber
- 310 Temp. & Humid. Chamber, Walk-in type
- 311 Temp. & Humid. Chamber, Customized

— Heating & Cooling Chamber

- 316 Heating & Cooling Chamber, Table top (-40°C, -20°C) Horizontal
- 320 Heating & Cooling Chamber, Ultra Low Temp. (-70°C) Horizontal
- 324 Heating & Cooling Chamber, Low Temp. (-35°C) Horizontal
- 328 Heating & Cooling Chamber, Low Temp. (-35°C) Vertical
- 332 Heating & Cooling Chamber, Low Temp. (-25°C) Horizontal
- 336 Heating & Cooling Chamber, Ambient temp. (-5°C) Horizontal
- 340 Heating & Cooling Chamber, General (-20°C, 0°C) Horizontal

— Heating Chamber

- 346 Heating Chamber



Products at a Glance

JEIOtech Research / Laboratory Equipment

Laboratory Equipment



Bath Circulator 24



Chiller (Recirculating Cooler) 36



Heating Bath 54



Cold Trap Bath 57



Heating Bath, Visible 58



Heating Bath, Shaking 60



Forced (Natural) Convection Oven 66 (74)



Forced Convection Oven, Multi-Chamber 73



Vacuum Oven 78



Furnace 82



Incubator 86



Plant Growth Chamber 104



Incubated Shaker, Stackable 114



Incubated Shaker, Tabletop 116



Incubated Shaker Floor 118



Shaker 126

JEIOtech Research / Laboratory Equipment



Microplate Shaker **131**



Rocking / Waving Shaker **132**



Funnel Shaker **133**



Rotator **134**



Hotplate & Magnetic Stirrer **142**



Hotplate **145**



Magnetic Stirrer **148**



Overhead Stirrer **156**



Vortex Mixer **166**



Heating & Cooling Block **174**



Heating Shaker **175**



Electrophoresis System (horizontal/vertical) **178**



Power Supply for Electrophoresis System **182**



Autoclaves **184**



Electrical Aspirator & Vacuum pump **188**



Ultrasonic Cleaner **190**

Products at a Glance

JEIOtech Research / Laboratory Equipment

Biosafety Cabinet & Clean Bench, Fume Hood



Biosafety Cabinet 196



Clean Bench 200



PCR Workstation 203



UV Sterilization Cabinet 204



Ductless Fume Hood 208



Fume Hood, Small 210



Mobile Fume Extractor, Local Exhaust Type 212



Extraction Arm Hood 216

Safety Storage & Cabinet



Fire Safety Storage Cabinet 226



Acid/Corrosive Cabinet 236



Filtering Storage Cabinet 240



Gas Cylinder Cabinet 243



Multipurpose Cabinet 244



Dry Cabinet (Auto desiccator) 248



Vacuum Desiccator 250



Desiccator 252

JEIOtech Research / Laboratory Equipment

Refrigerator & Freezer for Laboratory/ Medical



Laboratory Refrigerator **260**



Laboratory Freezer **262**



Pharmaceutical Refrigerator **264**

Lab Essentials



Plastic Ware **270**



Safety Product **275**

JEIOtech Environmental Reliability Test Equipment



Temp. & Humid. Chamber Tabletop type **286**



Temp. & Humid. Chamber General type **295**



Temp. & Humid. Chamber Dual type **299**



Pharmaceutical Stability Chamber **305**



Heating & Cooling Chamber Tabletop type **316**



Heating & Cooling Chamber Ultra low/Low/Ambient Temp. **320**



Heating & Cooling Chamber General type **340**



Heating Chamber **346**

Features at a glance

Certificate / Patent / Common Features



Registered safety patent based on Jeio Tech's proprietary technology.



1 year warranty free A/S.



Digital Input / Output for external device interlock control.



Registered patent based on Jeio Tech's proprietary technology.



2 year warranty free A/S.



It memorizes the state of power failure and operates automatically when power is restored.



Satisfaction with the requirements of European regulations and guidelines.



3 year warranty free A/S.



More precise temperature control through 3-point temperature calibration.



Certification of conformity of electricity, gas, etc. in accordance with the safety standards of Canada and USA.



Intuitive operation with color touch display.



More precise temperature control through temperature calibration.



Recognized as an excellent design product by the Ministry of Commerce, Industry, and Energy.



Real-time equipment monitoring and control system using mobile app.



Warning alarm in case of deviation of temperature control.



Compliance with Guidelines for Restriction of Hazardous Substances use.



Possible to set end time or start time of device operation.



Over temperature protection function.



Eco-friendly product that has passed environmental impact assessment tests.



Control of temperature / humidity / rpm / time, etc. through dedicated program.



Electronic over temperature protection system.



5 year warranty for direct drive motor.



PC communication via RS-232 / RS-485 / USB port.



Conventional over temperature protection, a kind of backup device against electronic over temperature protection failure.




Ensuring the reliability of electronic records and signatures to support GMP audit.



Real-time monitoring to check status of device.



Data storage using USB ports.



Stirrer / Hotplate



Bluetooth connection available for mobile control.



Possible to check and control temperature of sample by connecting to device with an external sensor.



Possible to connect an external controller.



Maintains possible maximum orbital motion without stopping, even if overload occurs.



If the surface temperature of the hot plate exceeds 50 degrees, the warning lamp turns on.

Storage / Cabinet



Door locking system with key.



Excellence in acid/alkaline storage with a corrosion-resistant structure and exhaust system.



30 minute fire resistance confirmed.



60 minute fire resistance confirmed.



90 minute fire resistance confirmed.



User safety by blocking potentially harmful UV light.



Leakage minimized using an integrated mold.

Shaker / Incubated Shaker



Orbital motion either clockwise or counterclockwise direction can be set.



Orbital motion.



Rocking motion.



Waving motion.



Vertical reciprocal motion.



Reciprocal motion in inclined up to 20°.



Selectable shaking motion either orbital or reciprocal motion.



Mixing through 360 degree rotation.



More accurate speed control with rpm calibration.

Bath / Chiller / Incubator / Temperature & Humidity Chamber



Indicator to ascertain water level.



Alarm if the water is low level.

Hood / Clean Bench



Differentiated diffusing muffler produces good laminar flow.



Equipped with UV Light.



Exhaust via duct connection.



Purifies and discharges harmful gas through internal filter.



If sash door is opened, automatically UV light is off and fluorescent light and fan are operated.



Equipped with the same controller inside and outside for convenience.

Refrigerator / Freezer



Certificate of medical devices management system for the manufacturers of medical devices.



Automatic defrosting does not need the work of unnecessary frost elimination.



In case of abnormal operation of the device, a voice message alarm is transmitted via telephone.



More convenient security system with electronic lock function.



Real-time temperature monitoring by thermal line recorder.



Laboratory Equipment



Laboratory Equipment



Bath Circulator



General Application

Internal/external temperature control: Spectrophotometer, refractometer, polarimeter, viscometer, fermenter, electrophoresis equipment, column chromatography, rotary evaporator etc.

	Description	Temp. Range (°C/°F)	Max. Temp. Stability (°C/°F)	Volume (L / cu ft)	Model	Page
Refrigerated & Heating	Program <ul style="list-style-type: none"> • 10 Programs & 10 Steps • Intuitive operation with color touch display 	-25 to 150 /-13 to 302 (RW3-25P) -35 to 150 /-31 to 302 (RW3-35P)	±0.05 / 0.09 at -10 °C	5, 10, 20, 30 / 0.2, 0.4, 0.7, 1.1	RW3-P	28
	General <ul style="list-style-type: none"> • Clear VFD & simple operation • Applies smart refrigeration system 	-25 to 150 /-13 to 302 (RW3-25) -35 to 150 /-31 to 302 (RW3-35)	±0.05 / 0.09 at -10°C	5, 10, 20, 30 / 0.2, 0.4, 0.7, 1.1	RW3	30
Heating	Program <ul style="list-style-type: none"> • 10 Programs & 10 Steps • Intuitive operation with color touch display 	Amb. +10 to 150 / Amb. +18 to 302	±0.05 / 0.09 at 50°C	5, 10, 20, 30 / 0.2, 0.4, 0.7, 1.1	CW3-P	32
	General <ul style="list-style-type: none"> • Clear VFD & simple operation • Detection and notification of temperature deviation 	Amb. +10 to 150 / Amb. +18 to 302	±0.05 / 0.09 at 50°C	5, 10, 20, 30 / 0.2, 0.4, 0.7, 1.1	CW3	32

※ The contents of the above and the contents of this catalog may differ depending on the specific model and conditions of use. For the information about the features and specifications that applying to each models, please check the information on the corresponding page of each models.

	Satisfaction with the requirements of European regulations and guidelines.		Real-time equipment monitoring and control system using mobile app.		Registered safety patent based on Jeio Tech's proprietary technology.
	Eco-friendly product that has passed environmental impact assessment tests.		Digital Input / Output for external device interlock control.		Intuitive operation with color touch display.
	Warning alarm in case of deviation of temperature control.		Possible to set end time or start time of device operation.		Control of temperature / humidity / rpm / time, etc. through dedicated program.
	PC communication via RS-232 / RS-485 / USB port.		More precise temperature control through 3-point temperature calibration.		Over temperature protection function.
	Alarm if the water is low level.		It memorizes the state of power failure and operates automatically when power is restored.		2 year warranty free A/S.



Innovative control system & design Application of specialized exclusive control system

> Selection of 2 control systems

Program Controller Type

- Program operation controller, 10 programs. (10 steps maximum)
- Intuitive operation with color touch display.

General Controller Type

- Controller for constant temperature value control and temperature deviation notification. (deviation limit/alarm delay setting)
- Precise and stable control of set temperature.

> Sharper screen display

Program type is applied to color touch display, general type is applied to clear VFD.

> Provision of remote equipment monitoring service (option)

Real-time monitoring of equipment operation status using smartphone. (when purchased LC GreenBox)

> More stable temperature control (RW3/CW3)

- Depending on the type of fluid used in the bath, three different auto-tuning values ($\Delta T1-3$) are saved.
- Auto-tuning is activated after 3 operating temperatures are inputted for each auto-tuning.

> Optimized temperature control

Auto-tuning function enables optimal control even when mainly used temperature or ambient conditions are frequently changed.

> 3-point temperature calibration function (RW3/CW3)

Temperature calibration of up to 3 points provides more reliable temperature control.

> Computer Interface (USB, RS-232)

Convenient remote control and data management on PC via built-in dedicated software.



Color touch display controller.



Clear VFD controller.



Provision of remote equipment monitoring service.



Provision of convenient and useful functions for actual use Applies Jeio Tech's patented safety system

- > **Convenient timer function provided**
Wait On/Wait Off is run after the set time elapses.
(Max. 999 hours and 59 minutes)
- > **Smart cooling control system**
According to the temperature, the flow of refrigerant is controlled.
The cooling speed from high temperature to room temperature is faster than previous product.
- > **Convenient solution discharge**
Drain valve is included as standard, making it both easy and convenient to discharge the solution in the bath.
- > **Stable performance of circulation pump**
Stable circulation pump with max. 0.4 bar pressure makes it easy to adjust pump speed in 5 steps.
- > **Top-grade protection system (Jeio Tech's patent)**
The safety device is configured to operate preferentially even if the malfunction is occurred in the device control unit.
- > **Shutdown function in case of over temperature**
Safe structure with shutdown function using over temperature limiter in case of overheating.
- > **Solution level detection**
Detects the level of solution in the bath and stops the device and notifying when the level is low.
- > **Prepared for unexpected power failure**
Auto-run function capable of selecting reoperation when power is restored after sudden power failure.

※ Some of the above contents are limited to specific models.



Sophisticated refrigerant volume control system.



Removable freezing condenser grill structure. (RW3)



Convenient solution discharge (CW3) with drain valve located on the front.


Refrigerated & Heating Bath Circulator

Program type

Convenient program control function

Program Type

- Possible to set and save 10 programs.
- Up to 10 steps can be set in the one program.
- Set time up to 99 hours and 59 minutes per step.
- Includes useful scheduler management functions such as program operation range and repeat interval setting.
- Program repetition possible up to 99 times.
- Provides temperature ramp control in program operation.
- Digital Input / Output function makes it easy to connect external equipment.
- Intuitive operation with color touch display.



RW3-0525P

with PPS (poly phenylene sulfide)
Bath Cover (standard)

Structural Functional Features

- Internal circulation and internal/external circulation of solution.
- Built-in bath cover with handle.
- Internal bath with seamless integral structure. (except RW3-3025P, 3035P models)
- Stainless steel material inside and at the top of the bath is excellent in terms of corrosion resistance and clean maintenance.
- Drain valve for easy draining of bath solution.
- Easy-to-clean refrigerator condenser filter.
- Equipped with easy-to-move/install caster.
- 8 models provided according to temperature range, capacity, etc.
- 5-steps pump speed adjustment with simple operation.
- Possible to check and control temperature of sample by connecting to device with an external sensor. (only available option when ordering)

Outstanding Safety

- Top-rated overheating protection system. (Registration KR 10-0397583)
- Over temperature limiting function. (Over Temperature Limiter)
- Prevention of mal-operation with controller lock function.
- Check and alert of temperature sensor abnormality.
- If an error occurs during operation, the details of the error are displayed on the screen and an alarm is sounded.
- Safe structure in which the heater, sensor, and pump inside the bath are separated by covering plates.


Use Convenience Features

- Temperature calibration of up to 3 points provides more reliable temperature control.
- Optimized control with auto-tuning.
- Wait On/Off timer. (up to 99 hours 59 minutes)
- Microprocessor PID method for precise temperature control.
- USB/RS-232 connection and software provide convenient computer operation and data management.
- Notification of low solution. (low level limiter)
- Automatic restart of operation when power is restored after power failure.
- Monitor via mobile app anytime, anywhere with LC Connected. (Requires LC GreenBox, option)

Pump Capacity



Pump Step	1	2	3	4	5
Max. pressure (psi)	0.5	1.1	2.4	4.9	5.9
Max. flow rate (L/min)	5.6	8.0	12.8	18.8	28.0

※ Bath fluid : Water





Heat-up Time, Cool-down Time



Specification

Model		RW3-0525P	RW3-1025P	RW3-2025P	RW3-3025P
Bath volume (L / cu ft)		5 / 0.2	10 / 0.4	20 / 0.7	30 / 1.1
Temperature	Working temp. range (°C / °F)	-25 to 150 / -13 to 302			
	Temp. stability ¹⁾ (±°C / °F)	0.05 / 0.09	0.05 / 0.09	0.05 / 0.09	0.05 / 0.09
Cooling capacity	at +20°C (+68°F), W	320	380	550	740
	at 0°C (+32°F), W	230	270	420	630
	at -20°C (-4°F), W	80	150	300	470
Dimension	Bath opening / depth (W x D x H, mm / inch)	150 x 99 x 160 / 5.9 x 3.9 x 6.3	211 x 154 x 180 / 8.3 x 6.1 x 7.1	290 x 214 x 200 / 11.4 x 8.4 x 7.9	300 x 264 x 230 / 11.8 x 10.4 x 9.1
	Exterior (W x D x H, mm / inch)	308 x 480 x 733 / 12.1 x 18.9 x 28.9	366 x 532 x 780 / 14.4 x 20.9 x 30.7	386 x 602 x 825 / 15.2 x 23.7 x 32.5	446 x 662 x 895 / 17.6 x 26.1 x 35.2
	Net weight (kg / lbs)	37 / 81.6	44 / 97	60 / 132.2	67 / 147.7
Electrical requirements (230V, A)		50/60Hz, 11	50/60Hz, 11	50/60Hz, 12.5	50Hz, 13.5 60Hz, 13.5
Cat. No.		AAH57642K	AAH57652K	AAH57662K	AAH57672K AAH57671K
Cat. No. (EU only)		AAH57645K	AAH57655K	AAH57665K	AAH57675K AAH57671K
Electrical requirements (120V, A)		60Hz, 12.6	60Hz, 12.6		
Cat. No.		AAH57643U	AAH57653U	-	-

Model		RW3-0535P	RW3-1035P	RW3-2035P	RW3-3035P
Bath volume (L / cu ft)		5 / 0.2	10 / 0.4	20 / 0.7	30 / 1.1
Temperature	Working temp. range (°C / °F)	-35 to 150 / -31 to 302			
	Temp. stability ¹⁾ (±°C / °F)	0.05 / 0.09	0.05 / 0.09	0.05 / 0.09	0.05 / 0.09
Cooling capacity	at +20°C (+68°F), W	280	450	800	880
	at 0°C (+32°F), W	150	330	660	700
	at -20°C (-4°F), W	80	150	400	440
Dimension	Bath opening / depth (W x D x H, mm / inch)	150 x 99 x 160 / 5.9 x 3.9 x 6.3	211 x 154 x 180 / 8.3 x 6.1 x 7.1	290 x 214 x 200 / 11.4 x 8.4 x 7.9	300 x 264 x 230 / 11.8 x 10.4 x 9.1
	Exterior (W x D x H, mm / inch)	308 x 480 x 733 / 12.1 x 18.9 x 28.9	366 x 532 x 780 / 14.4 x 20.9 x 30.7	386 x 602 x 825 / 15.2 x 23.7 x 32.5	446 x 662 x 895 / 17.6 x 26.1 x 35.2
	Net weight (kg / lbs)	38 / 83.8	44 / 97	60 / 132.2	67 / 147.7
Electrical requirements (230V, A)		50/60Hz, 12.5	50/60Hz, 12.5	50Hz, 13.5 60Hz, 13.5	50Hz, 14.5 60Hz, 14.5
Cat. No.		AAH57742K	AAH57752K	AAH57762K AAH57761K	AAH57772K AAH57771K
Cat. No. (EU only)		AAH57745K	AAH57755K	AAH57765K AAH57761K	AAH57775K AAH57771K

※ Technical data (according to DIN 12876)

※ In case using silicon oil, you need to exchange other material tube (inlet and out late part), such as copper tube, a teflon tube instead of a silicon tube in fitting.

1) Bath fluid : Ethanol at -10°C / 14°F

Refrigerated & Heating Bath Circulator

General type

Stable internal/external temperature control

General Type

- Alert function when deviation from control temperature upper limit and lower limit occurs.
- Easy and convenient to use with simple operation.
- Clear VFD panel for clear information display.
- Applies smart refrigeration system using Jeio Tech's technology.
- Optimized control reduces overall power consumption.
- Easily removable condenser grill structure.
- Equipped with advanced filter with flat plate structure in refrigeration condenser.
- Compact structure design.

RW3-2025

with PPS (poly phenylene sulfide)
Bath Cover (standard)



Structural Functional Features

- Internal circulation and internal/external circulation of solution.
- Refrigerator stopped automatically when it rises too high of temperature, saving power, extending freezer life, and reducing noise generation.
- High-efficiency refrigeration through application of electronic expansion valve.
- Built-in bath cover with handle.
- Internal bath with seamless integral structure. (except RW3-3025, 3035 models)
- Stainless steel material inside and at the top of the bath is excellent in terms of corrosion resistance and clean maintenance.
- Left and right handles for convenient handling.
- Drain valve for easy draining of bath solution.
- Easy-to-clean refrigerator condenser filter.
- Equipped with easy-to-move/install caster.
- 8 models provided according to temperature range, capacity, etc.
- 5-steps pump speed adjustment with simple operation.
- Possible to check and control temperature of sample by connecting to device with an external sensor. (only available option when ordering)

Use Convenience Features

- Temperature calibration of up to 3 points provides more reliable temperature control.
- Optimized control with auto-tuning.
- Save and use 3 frequently used temperatures.
- Wait On/Off timer. (up to 999 hours 59 minutes)
- Microprocessor PID method for precise temperature control.
- USB/RS-232 connection and software provide convenient computer operation and data management.
- Notification of low solution. (low level limiter)
- Automatic restart of operation when power is restored after power failure.
- Monitor via mobile app anytime, anywhere with LC Connected. (Requires LC GreenBox, option)


Outstanding Safety

- Top-rated overheating protection system. (Registration KR 10-0397583)
- Over temperature limiting function.
- Prevention of mal-operation with controller lock function.
- Check and alert of temperature sensor abnormality.
- If an error occurs during operation, the details of the error are displayed on the screen and an alarm is sounded.
- Safe structure in which the heater, sensor, and pump inside the bath are separated by covering plates.

Pump Capacity

Pump Step	1	2	3	4	5
Max. pressure (psi)	0.5	1.1	2.4	4.9	5.9
Max. flow rate (L/min)	5.6	8.0	12.8	18.8	28.0

※ Bath fluid : Water





Single-Action refrigeration condenser grill structure allows for detaching and attaching without additional tools.



The drain valve located on the side is simple and convenient to use as it is ball valve type.

Specification

Model		RW3-0525	RW3-1025	RW3-2025	RW3-3025
Bath volume (L / cu ft)		5 / 0.2	10 / 0.4	20 / 0.7	30 / 1.1
Temperature	Working temp. range (°C / °F)	-25 to 150 / -13 to 302			
	Temp. stability ¹⁾ (±°C / °F)	0.05 / 0.09	0.05 / 0.09	0.05 / 0.09	0.05 / 0.09
Cooling capacity	at 20°C (+68°F), W	400	560	850	980
	at 0°C (+32°F), W	270	360	610	700
	at -20°C (-4°F), W	140	180	310	465
Dimension	Bath opening / depth (W x D x H, mm / inch)	150 x 99 x 160 / 5.9 x 3.9 x 6.3	211 x 154 x 180 / 8.3 x 6.1 x 7.1	290 x 214 x 200 / 11.4 x 8.4 x 7.9	300 x 264 x 230 / 11.8 x 10.4 x 9.1
	Exterior (W x D x H, mm / inch)	308 x 480 x 733 / 12.1 x 18.9 x 28.9	366 x 532 x 780 / 14.4 x 20.9 x 30.7	386 x 602 x 825 / 15.2 x 23.7 x 32.5	446 x 662 x 895 / 17.6 x 26.1 x 35.2
	Net weight (kg / lbs)	37 / 81.6	44 / 97	60 / 132.3	67 / 147.7
Electrical requirements (230V, A)		50Hz, 11	60Hz, 11	50/60Hz, 12.5	50Hz, 13.5
Cat. No.		AAH57042K	AAH57041K	AAH57052K	AAH57051K
Cat. No. (EU only)		AAH57045K	AAH57041K	AAH57055K	AAH57051K
Electrical requirements (120V, A)		60Hz, 12.6	60Hz, 12.6		
Cat. No.		AAH57043U	AAH57053U	-	-

Model		RW3-0535	RW3-1035	RW3-2035	RW3-3035
Bath volume (L / cu ft)		5 / 0.2	10 / 0.4	20 / 0.7	30 / 1.1
Temperature	Working temp. range (°C / °F)	-35 to 150 / -31 to 302			
	Temp. stability ¹⁾ (±°C / °F)	0.05 / 0.09	0.05 / 0.09	0.05 / 0.09	0.05 / 0.09
Cooling capacity	at 20°C (+68°F), W	400	560	850	980
	at 0°C (+32°F), W	270	360	610	700
	at -20°C (-4°F), W	140	180	310	465
Dimension	Bath opening / depth (W x D x H, mm / inch)	150 x 99 x 160 / 5.9 x 3.9 x 6.3	211 x 154 x 180 / 8.3 x 6.1 x 7.1	290 x 214 x 200 / 11.4 x 8.4 x 7.9	300 x 264 x 230 / 11.8 x 10.4 x 9.1
	Exterior (W x D x H, mm / inch)	308 x 480 x 733 / 12.1 x 18.9 x 28.9	366 x 532 x 780 / 14.4 x 20.9 x 30.7	386 x 602 x 825 / 15.2 x 23.7 x 32.5	446 x 662 x 895 / 17.6 x 26.1 x 35.2
	Net weight (kg / lbs)	38 / 83.8	45 / 99.2	60 / 132.3	68 / 149.9
Electrical requirements (230V, 50Hz, A)		12.5	12.5	13.5	13.5
Cat. No.		AAH57142K	AAH57152K	AAH57162K	AAH57172K
Cat. No. (EU only)		AAH57145K	AAH57155K	AAH57165K	AAH57175K
Electrical requirements (230V, 60Hz, A)		12.5	12.5	13.5	13.5
Cat. No.		AAH57141K	AAH57151K	AAH57161K	AAH57171K

※ Technical data (according to DIN 12876)

※ In case using silicon oil, you need to exchange other material tube (inlet and out late part), such as copper tube, a teflon tube instead of a silicon tube in fitting.

1) Bath fluid : Ethanol at -10°C / 14°F

Heating Bath Circulator

Suitable for precise temperature testing above ambient temperature

CW3-05P

with PPS (poly phenylene sulfide)
Bath Cover (standard)



CW3-05

with PPS (poly phenylene sulfide)
Bath Cover (standard)



Program Type Features

- 10 Programs & 10 steps program control
- Set time up to 99 hours and 59 minutes per step.
- Intuitive operation with color touch display.
- Digital Input / Output function makes it easy to connect external equipment.

Structural Functional Features

- Minimizes evaporation of internal solution at high temperatures and low-heat conductivity bath cover is provided as standard.
- Minimized usage space with compact design.
- Excellent durability even at high temperatures and repeated use.
- Internal circulation and internal/external circulation of solution.
- Easy to move and manage with the body handle.
- Stainless steel material inside and at the top of the bath is excellent in terms of corrosion resistance and clean maintenance.
- Drain valve for easy draining of bath solution.
- Includes rubber feet for stable product operation.
- 5-steps pump speed adjustment with simple operation.

Outstanding Safety

- Top-rated overheating protection system. (Registration KR 10-0397583)
- Over temperature limiting function.
- Prevention of mal-operation with controller lock function.
- Check and alert of temperature sensor abnormality.
- If an error occurs during operation, the details of the error are displayed on the screen and an alarm is sounded.
- Safe structure in which the heater, sensor, and pump inside the bath are separated by covering plates.

General Type Features

- Alert function when deviation from control temperature upper limit and lower limit occurs.
- Easy and convenient to use with simple operation.
- Clear VFD panel for clear information display.
- Convenient drain valve structure on the front panel.
- Built-in USB connection and dedicated software.


Use Convenience Features

- Temperature calibration of up to 3 points provides more reliable temperature control.
- Optimized control with auto-tuning.
- Save and use 3 frequently used temperatures.
- Wait On/Off timer. (up to 999 hours 59 minutes)
- Microprocessor PID method for precise temperature control.
- Convenient remote control and data management through built-in dedicated software.
- Notification of low solution. (low level limiter)
- Automatic restart of operation when power is restored after power failure.
- Monitor via mobile app anytime, anywhere with LC Connected. (Requires LC GreenBox, option)

Pump Capacity

Pump Step	1	2	3	4	5
Max. pressure (psi)	0.5	1.1	2.4	4.9	5.9
Max. flow rate (L/min)	5.6	8.0	12.8	18.8	28.0

※ Bath fluid : Water





Specification

Model		CW3-05P	CW3-10P	CW3-20P	CW3-30P
Bath volume (L / cu ft)		5 / 0.2	10 / 0.4	20 / 0.7	30 / 1.1
Temperature	Working temp. range (°C / °F)	Amb. +10 to 150 / Amb. +18 to 302	Amb. +10 to 150 / Amb. +18 to 302	Amb. +10 to 150 / Amb. +18 to 302	Amb. +10 to 150 / Amb. +18 to 302
	Temp. stability ¹⁾ (±°C / °F)	0.05 / 0.09	0.05 / 0.09	0.05 / 0.09	0.08 / 0.162
	Heat up time to 70°C / 158 °F (min)	24	31	43	55
Dimension	Bath opening / depth (W x D x H, mm / inch)	150 x 99 x 160 / 5.9 x 3.9 x 6.3	211 x 154 x 180 / 8.3 x 6.1 x 7.1	290 x 214 x 200 / 11.4 x 8.4 x 7.9	300 x 264 x 230 / 11.8 x 10.4 x 9.1
	Exterior (W x D x H, mm / inch)	228 x 395 x 424 / 9.0 x 15.6 x 16.7	289 x 450 x 444 / 11.4 x 17.7 x 17.5	368 x 510 x 464 / 14.5 x 20.1 x 18.3	378 x 560 x 494 / 14.9 x 22 x 19.4
	Net weight (kg / lbs)	12.5 / 27.6	15.0 / 33.1	17.5 / 38.6	19.5 / 43
Electrical requirements (230V, 50/60Hz, A)		8.5	8.5	8.5	8.5
Cat. No.		AAH52741K	AAH52751K	AAH52761K	AAH52771K
Electrical requirements (120V, 60Hz, A)		8.6	8.6	8.6	8.6
Cat. No.		AAH52742U	AAH52752U	AAH52762U	AAH52772U

Model		CW3-05	CW3-10	CW3-20	CW3-30
Bath volume (L / cu ft)		5 / 0.2	10 / 0.4	20 / 0.7	30 / 1.1
Temperature	Working temp. range (°C / °F)	Amb. +10 to 150 / Amb. +18 to 302	Amb. +10 to 150 / Amb. +18 to 302	Amb. +10 to 150 / Amb. +18 to 302	Amb. +10 to 150 / Amb. +18 to 302
	Temp. stability ¹⁾ (±°C / °F)	0.05 / 0.09	0.05 / 0.09	0.05 / 0.09	0.08 / 0.162
	Heat up time to 70°C / 158 °F (min)	15	23	43	55
Dimension	Bath opening / depth (W x D x H, mm / inch)	150 x 99 x 160 / 5.9 x 3.9 x 6.3	211 x 154 x 180 / 8.3 x 6.1 x 7.1	290 x 214 x 200 / 11.4 x 8.4 x 7.9	300 x 264 x 230 / 11.8 x 10.4 x 9.1
	Exterior (W x D x H, mm / inch)	228 x 395 x 424 / 9.0 x 15.6 x 16.7	289 x 450 x 444 / 11.4 x 17.7 x 17.5	368 x 510 x 464 / 14.5 x 20.1 x 18.3	378 x 560 x 494 / 14.9 x 22 x 19.4
	Net weight (kg / lbs)	12.5 / 27.6	15.0 / 33.1	17.5 / 38.6	19.5 / 43
Electrical requirements (230V, 50/60Hz, A)		8.5	8.5	8.5	8.5
Cat. No.		AAH52341K	AAH52351K	AAH52361K	AAH52371K
Electrical requirements (120V, 60Hz, A)		8.6	8.6	8.6	8.6
Cat. No.		AAH52342U	AAH52352U	AAH52362U	AAH52372U

※ The lowest temperature that can be controlled depends on changes in ambient temperature and sample temperature. Therefore, please contact the distributor of our products before purchasing for technical consultation.

※ Technical data (according to DIN 12876)

※ In case using silicon oil, you need to exchange other material tube (inlet and out late part), such as copper tube, a teflon tube instead of a silicon tube in fitting.

1) Bath fluid : Water at 50°C / 122°F

Accessories Page 34 Spring Wire Rack, Test Tube Rack, Silicone Tube, Tubing Insulation, Clamp, Cover, LC GreenBox

Accessories

Stainless Steel Cover

Sturdy, corrosion-resistant, stainless steel material.

Cat. No.	Model
AAA51531	RW3-0525(P), RW3-0535(P), CW3-05(P)
AAA51532	RW3-1025(P), RW3-1035(P), CW3-10(P)
AAA51533	RW3-2025(P), RW3-2035(P), CW3-20(P)
AAA51534	RW3-3025(P), RW3-3035(P), CW3-30(P)



Spring Wire Rack

It has a lattice-shaped double spring structure, which makes it easy to fix various kinds of experiment containers and the height can also be easily adjusted.

Cat. No.	Model
AAA51522	RW3-2025(P), 2035(P)
AAA51523	RW3-3025(P), 3035(P)
AAA51522	CW3-20(P)
AAA51523	CW3-30(P)



Test Tube Rack

Test tubes can be easily fixed in the bath and various models are available according to test tube size. Structure allows for easy adjustment of the degree of tilt and lock.

Cat. No.	Description (mm / inch)	Max. Quantity of Racks
AAA44581	ø8 / 0.3 x 86 tubes	CW3-20(P): 2ea CW3-30(P): 2ea RW3-2025(P): 2ea RW3-2035(P): 2ea RW3-3025(P): 2ea RW3-3035(P): 2ea
AAA44582	ø10 / 0.4 x 86 tubes	
AAA44583	ø12 / 0.5 x 58 tubes	
AAA44585	ø16 / 0.6 x 32 tubes	
AAA44586	ø25 / 1.0 x 19 tubes	
AAA41531	ø8 / 0.3 x 50 tubes	CW3-10(P): 1ea RW3-1025(P): 1ea RW3-1035(P): 1ea
AAA41532	ø10 / 0.4 x 50 tubes	
AAA41533	ø12 / 0.5 x 33 tubes	
AAA41535	ø16 / 0.6 x 16 tubes	
AAA41539	ø25 / 1.0 x 10 tubes	



External Sensor (for RW3)

Possible to check and control temperature of sample by connecting to device with an external sensor. (only available option when ordering)



Jeio Tech



How to Control Temperature using External Sensor

Cat. No.	Description
AAA51571	External sensor board + External Sensor (3 m)
AAA51572	External sensor board + External Sensor (5 m)

※ If you purchased the RW3 initially without an external sensor, it cannot be added to the product later.



LC GreenBox

Real-time monitoring and control of device operation via mobile app, LC Connected. Up to 4 devices can be connected at once.

Cat. No.	Dimension (W x D x H, mm / inch)
AAHQ1011K	156 x 94 x 34 / 6.1 x 3.7 x 1.3



Silicone Tube

The operating temperature range is -50 to 250°C / -58 to 482°F.
Inner diameter of 10 mm,
Outside diameter of 16 mm.

Cat. No.	Length (m / inch)
00FAA0001110	10 / 393.7



Tubing Insulation

EPDM material insulation prevents heat loss when using external circulation.

Cat. No.	Description
00FEA0001453	I.D. 25 mm (13T)



Tube Clamp

Attachment that stably holds the tube connection.

Cat. No.	Description
00HXE0001099	3/4"
00HXE0001100	5/8"



Bath Circulator

VS.

Chiller



	Bath Circulator (ex. RW3-3035)	Chiller (Recirculating Cooler) (ex. HX-25H)
Application	Preserves precise temperature of external devices where heat removal capacity is not required largely.	Maintains temperature by stably removing large heat generated from external devices.
Cooling Capacity	Heat removal capacity is small. 880W at 20°C / 68°F.	Heat removal capacity is large. 2,800W at 20°C / 68°F.
Temperature Range	Temperature range is wide. -35 to 150°C / -31 to 302°F.	Temperature range is narrow. 3 to 40°C / 37 to 104°F.
Pump Pressure	Circulation pressure is weak. Max. 0.4 bar.	Circulation pressure is strong. Max. 6 bar.

Jeio Tech's Chiller Information [Page 36](#)

Chiller for Rotary Evaporator

- Chiller optimized for rotary evaporator.
- Easy to install and use with compact design.
- Pump In/Out port is located at the upper part, so it is convenient to connect with evaporator.
- Convenient bath fluid injection on the top.
- LCD screen for easy checking of operation status.
- Includes low-level solution alert notification.

Model	RC-05
Temp. Range (°C / °F)	-10 to 30 / 14 to 86
Temp. Stability (±°C / °F)	1.0 / 1.8 at 15°C
Cooling Capacity	580W at 20°C / 68°F
230V, 50/60Hz, A	4
Cat. No.	ACH651012K
Cat. No. (EU only)	ACH651015K



Chiller [Page 44](#)

Chiller (Recirculating Cooler)



General Application

- Sample thawing, and control of temperature of diverse experimental equipment.
- Analytical equipment requiring cooling. (rotary evaporator, vacuum oven, spectrophotometer, etc.)
- Plant, medical equipment, metalwork laser, injection machine.

	Description	Max. Temp. Stability (±°C/°F)	Temp. Range (°C/°F)	Max. Cooling Capacity (kW)	Model	Page
Ambient temp. type	General <ul style="list-style-type: none"> • Superior durability and performance • Verified according to international standards 	±1 / 1.8	3 to 40 / 37.4 to 104	8.1 at 20°C	HX	40
Low temp. type	General <ul style="list-style-type: none"> • Applies clear VFD display • Includes temperature calibration function 	±1 / 1.8	-20 to 40 / -4 to 104	7.1 at 20°C	HL	42
	Small <ul style="list-style-type: none"> • Optimization of rotary evaporator use • Compact design 	±1 / 1.8	-10 to 30 / 14 to 86	0.5 at 20°C	RC	44
High / Low temp. type	Precision <ul style="list-style-type: none"> • Wide temperature zone control • 0.1°C precision control 	±0.1 / 0.2	-20 to 80 / -4 to 176	6.5 at 20°C	HH	46

※ The contents of the above and the contents of this catalog may differ depending on the specific model and conditions of use. For the information about the features and specifications that applying to each models, please check the information on the corresponding page of each models.



Satisfaction with the requirements of European regulations and guidelines.



Real-time equipment monitoring and control system using mobile app.



Eco-friendly product that has passed environmental impact assessment tests.



Alarm if the water is low level.



Over temperature protection function.



More precise temperature control through temperature calibration.



PC communication via RS-232 / RS-485 / USB port.



2 year warranty free A/S.

Chiller (Recirculating Cooler)



Innovative refrigeration system and dedicated controller for chiller Improved durability and stable control

> Innovative refrigeration system

Increased durability of refrigerant adjustment-type refrigeration system by minimizing compressor on/off operation.

> Verified by objective criteria

Designed and tested according to stringent international standards.
Proven performance specifications and high reliability.

> Heat absorption ability of excellent reproducibility

Convenient to select and use the optimal chiller with the heat absorption ability of each model verified by objective criteria.

> Variety of products optimized for their specified purpose of use

Available in 4 types and 24 models as a standard type according to temperature range, precision, and endothermic/pump ability.

> Chiller dedicated controller

Allows more stable control by using dedicated controller developed with Jeio Tech's proprietary technology to optimize for chiller characteristics.

> Temperature calibration function


Controller with temperature calibration function allows for stable implementation of more accurate temperature control.

> On/Off of freezer and pump

The refrigeration system and pump can each be turned On/Off separately.
This is useful for maintenance and experiment setting, etc.

> Computer interface

Convenient device operation and data processing function using computer connection and dedicated software.



Bath inlet / outlet port.



Pressure Gauge.



Convenient bath fluid level checking.



Useful features and structural features in actual use

Easy maintenance

> Provision of remote equipment monitoring service (option)

Real-time monitoring of equipment operation status using smartphone.
(when purchased LC GreenBox) (HH model)

> Convenient bath fluid level check

Convenient to check the level of bath fluid at the front of the instrument.
Level Indicator with LED backlight.

> Check pump pressure on the front of the instrument

Pump pressure can be confirmed from the front of the machine, so it is always convenient to check the operation status.

> Casters for easy installation

Convenient casters installed for easy moving and installation during initial installation and during use.

> Over temperature warning function

Over temperature warning function that generates an alarm when over temperature occurs.

> Low bath fluid warning

If the bath fluid is low, an alarm is generated so that fluid can be supplied when appropriate.

> Application support priority

Chiller operation continues to protect user application even when abnormal temperature or bath fluid is detected.

> Refrigeration system protection device

Refrigerator overload and over temperature protection.

Automatic stop when refrigerant pressure (high/low pressure) abnormality occurs.

※ Some of above contents are limited to specific models.



Bath filling inlet.



Includes drain valve.



Easy-to-install casters.

Chiller (Recirculating Cooler)

HX general type

Verified performance and durability according to international standards.

Structural Functional Features

- Designed and tested according to strict international standards.
- Provides proven performance specifications and high reliability.
- Optimized slim and compact design minimizes installation space.
- The conventional simple refrigerant compressor control method can cause the compressor life to be shortened by turning the refrigeration system on/off from time to time to control the temperature. However, the Jeio Tech's Chiller is a unique refrigerant control type cooling system completed using our proprietary technology development. Compressor On/Off operation is minimized, and durability of the cooling system is greatly increased.
- Double service valve prevents refrigerant leakage.
- Easy installation with air-cooled type integrated refrigeration system.
- Easy to clean as it is easy to remove refrigerator condenser grill, making it convenient to maintain efficiency of refrigeration.
- Installation of refrigeration system using environmentally friendly refrigerant.
- Pump pressure can be confirmed from the front of the machine, so it is always convenient to check the operation status.
- Pressure-sensitive user application protection (H Type) due to by-pass function to reduce pump pressure.
- Bath fluid inlet is wide, making installation and maintenance convenient.
- Equipped with removable casters for easy movement and installation.



HX-65H

HX-45H



Use Convenience Features

- Bright LCD display provides excellent visual perception.
- Identify all operations and operating conditions from the front panel.
- Simultaneously displays the set value and the current value, and settings can be changed during equipment use, making it convenient.
- Highly-reliable control through calibration.
- The refrigeration compressor and pump can be each turned On/Off.
- Level indicator with LED backlight for easy checking of bath fluid level.
- Easy to maintain as a drain valve is included.


Outstanding Safety

- Over temperature alarm.
- Provides a warning when abnormal temperature or bath fluid shortage is detected, and keeps the chiller operating continuously to protect the user application.
- Over-current and short circuit protection of device.
- Auto stop when operating current of refrigerator is overloaded.




Pumping capacities

Magnetic pump



※ Application model : HX-20, 25

Centrifugal pump



※ Application model : HX-20H, 25H, 35H, 45H, 55H

Specification

Model	HX-20	HX-25	HX-20H	HX-25H	HX-35H	HX-45H	HX-55H
Interior dimensions							
Max. filling volume (L / cu ft)	35 / 1.24	35 / 1.24	35 / 1.24	35 / 1.24	45 / 1.59	45 / 1.59	45 / 1.59
Refrigerator capacity (HP)	3/4	1	3/4	1	1.5	2	3
For tubing dia. (mm / inch)	20 / (3/4")	20 / (3/4")	20 / (3/4")	20 / (3/4")	20 / (3/4")	20 / (3/4")	20 / (3/4")
Filling inlet (Ø, mm / inch)	50 / 1.9	50 / 1.9	50 / 1.9	50 / 1.9	50 / 1.9	50 / 1.9	50 / 1.9
Exterior dimensions							
Width (mm / inch)	620 / 24.4	620 / 24.4	620 / 24.4	620 / 24.4	745 / 29.3	745 / 29.3	745 / 29.3
Depth (with drain valve), (mm / inch)	715 (785) / 28.1(30.9)	715 (785) / 28.1(30.9)	715 (785) / 28.1(30.9)	715 (785) / 28.1(30.9)	735 (800) / 28.9 (31.5)	735 (800) / 28.9 (31.5)	735 (800) / 28.9 (31.5)
Height (with lid), (mm / inch)	955 (990) / 37.6 (39.0)	955 (990) / 37.6 (39.0)	955 (990) / 37.6 (39.0)	955 (990) / 37.6 (39.0)	1060 (1095) / 41.7 (43.1)	1060 (1095) / 41.7 (43.1)	1060 (1095) / 41.7 (43.1)
Weight (kg / lbs)	93 / 205	101 / 223	97 / 214	105 / 231	120 / 265	130 / 287	140 / 309
Temperature data							
Working temperature range (°C / °F)	3 to 40 / 37.4 to 104	3 to 40 / 37.4 to 104	3 to 40 / 37.4 to 104				
Temperature stability at 15°C (±°C / °F)	1 / 1.8	1 / 1.8	1 / 1.8	1 / 1.8	1 / 1.8	1 / 1.8	1 / 1.8
Cooling capacity, Max							
at 20°C (kW)	2.2	2.8	2.2	2.8	4.2	5.5	8.1
at 10°C (kW)	1.6	2	1.6	2	3.6	4.2	5.9
at 5°C (kW)	1.2	1.6	1.2	1.6	2.8	3.4	4.8
Pump data							
Max. flow rate (L / min, gal / min)	40 / 10.6	40 / 10.6	32 / 8.5	32 / 8.5	32 / 8.5	32 / 8.5	32 / 8.5
Max. pressure (bar / psi)	2 / 29.0	2 / 29.0	6 / 87.0	6 / 87.0	6 / 87.0	6 / 87.0	6 / 87.0
Electrical data & Ordering information							
230V, 50Hz, A	7.0	7.6	10.5	11.0	17.0		
Cat. No.	AAH64012K	AAH64022K	AAH64112K	AAH64122K	AAH64132K	-	-
Cat. No. (EU only)	AAH64015K	AAH64025K	AAH64115K	AAH64125K	AAH64135K	-	-
230V, 60Hz, A	8.5	9.0	9.0	9.5	15.0	19.0	
Cat. No.	AAH64011K	AAH64021K	AAH64111K	AAH64121K	AAH64132K	AAH64141	-
120V, 60Hz, A	17.0	22.0					
Cat. No.	AAH64013U	AAH64023U	-	-	-	-	-
380V, 50Hz, A						7.0	8.5
Cat. No.	-	-	-	-	-	AAH64148K	AAH64158K
Cat. No. (EU only)	-	-	-	-	-	AAH64140K	AAH64150K
380V, 60Hz, A						10.0	11.5
Cat. No.	-	-	-	-	-	AAH64149K	AAH64159K

- ※ According to DIN 12876
- ※ Temperature stability / Cooling capacity / Pump data : Water
- ※ Above specification value is recorded by 230V 60Hz. (HX-45H, 55H is recorded by 380V 60Hz)
- ※ Product performance may be affected by ambient temperature.

Chiller (Recirculating Cooler)

HL low temperature general type

Low temperature-type general model controlled down to -20°C.

Structural Functional Features

- Low-temperature chiller controlled from -20°C to 40°C has proven performance specifications and high reliability.
- Optimized slim and compact design minimizes installation space.
- The conventional simple refrigerant compressor control method can cause the compressor life to be shortened by turning the refrigeration system on/off from time to time to control the temperature. However, the Jeio Tech's Chiller is a unique refrigerant control type cooling system completed using our proprietary technology development. Compressor On/Off operation is minimized, and durability of the cooling system is greatly increased.
- Double service valve prevents refrigerant leakage.
- Easy installation with air-cooled type integrated refrigeration system.
- Easy to clean as it is easy to remove refrigerator condenser grill, making it convenient to maintain efficiency of refrigeration.
- Installation of refrigeration system using environmentally friendly refrigerant.
- Pump pressure can be confirmed from the front of the machine, so it is always convenient to check the operation status.
- Bath fluid inlet is wide, making installation and maintenance convenient.
- Equipped with removable casters for easy movement and installation.



HL-25H

HL-20H



Use Convenience Features

- Clear VFD display provides excellent visual perception.
- Identify all operations and operating conditions from the front panel.
- Simultaneously displays the set value and the current value, and settings can be changed during equipment use, making it convenient.
- Highly-reliable control through calibration.
- The refrigeration compressor and pump can be each turned On/Off.
- Level indicator with LED backlight for easy checking of bath fluid level.
- Easy to maintain as a drain valve is included.


Outstanding Safety

- Over temperature alarm.
- Provides a warning when abnormal temperature or bath fluid shortage is detected, and keeps the chiller operating continuously to protect the user application.
- Over-current and short circuit protection of device.
- Automatic stop in case of over temperature of the compressor.
- Malfunction prevented by controller lock function.




Pumping capacities

SS housing ABS pump




※ Application model : HL - 05, 10, 15, 20

Submersible pump



※ Application model : HL - 15H, 20H, 25H, 35H



※ Application model : HL - 45H, 55H

Specification

Model	HL-05	HL-10	HL-15	HL-20	HL-15H	HL-20H	HL-25H	HL-35H	HL-45H	HL-55H
Interior dimensions										
Max. filling volume (L / cu ft)	8 / 0.28	8 / 0.28	14 / 0.49	14 / 0.49	14 / 0.49	14 / 0.49	25 / 0.88	25 / 0.88	39 / 1.38	39 / 1.38
Refrigerator capacity (HP)	1/5	1/3	1/2	3/4	1/2	3/4	1	1.5	2	3
For tubing dia. (mm / inch)	20 / (3/4")	20 / (3/4")	20 / (3/4")	20 / (3/4")	20 / (3/4")	20 / (3/4")	20 / (3/4")	20 / (3/4")	20 / (3/4")	20 / (3/4")
Filling inlet (Ø, mm / inch)	120 / 4.7	120 / 4.7	120 / 4.7	120 / 4.7	120 / 4.7	120 / 4.7	120 / 4.7	120 / 4.7	120 / 4.7	120 / 4.7
Exterior dimensions										
Width (mm / inch)	405 / 15.9	405 / 15.9	515 / 20.3	515 / 20.3	515 / 20.3	515 / 20.3	550 / 21.7	550 / 21.7	605 / 23.8	605 / 23.8
Depth (with drain valve), (mm / inch)	550 (620) / 21.7 (24.4)	550 (620) / 21.7 (24.4)	645 (715) / 25.4 (28.1)	645 (715) / 25.4 (28.1)	645 (715) / 25.4 (28.1)	645 (715) / 25.4 (28.1)	830 (900) / 32.7 (35.4)	830 (900) / 32.7 (35.4)	975 (1045) / 38.4 (41.1)	975 (1045) / 38.4 (41.1)
Height (with lid), (mm / inch)	675 (710) / 26.6 (10.6)	675 (710) / 26.6 (10.6)	800 (835) / 31.5 (32.9)	800 (835) / 31.5 (32.9)	800 (835) / 31.5 (32.9)	800 (835) / 31.5 (32.9)	1105 (1140) / 43.5 (44.9)	1105 (1140) / 43.5 (44.9)	1265 (1300) / 49.8 (51.2)	1265 (1300) / 49.8 (51.2)
Weight (kg / lbs)	62.9 / 138.7	64.9 / 143.1	86.5 / 190.7	87.4 / 192.7	91.8 / 202.4	92.7 / 204.4	141.3 / 311.5	146.3 / 322.5	171±10 / 377±20.1	176±10 / 388±20.1
Temperature data										
Working temperature range (°C / °F)	-20 to 40 / -4 to 104	-20 to 40 / -4 to 104	-20 to 40 / -4 to 104	-20 to 40 / -4 to 104						
Temperature stability at 15°C (±°C / °F)	1 / 1.8	1 / 1.8	1 / 1.8	1 / 1.8	1 / 1.8	1 / 1.8	1 / 1.8	1 / 1.8	1 / 1.8	1 / 1.8
Cooling capacity, Max										
at 20°C (kW)	0.7	0.76	1.7	1.85	1.7	1.85	2.8	3.5	6.5	7.1
at 10°C (kW)	0.6	0.65	1.35	1.6	1.35	1.6	2.2	2.5	4.5	6
at 0°C (kW)	0.4	0.5	0.87	1.2	0.87	1.2	1.3	1.8	3	4.1
at -10°C (kW)	0.2	0.36	0.67	0.87	0.67	0.87	0.9	1.1	2.1	2.5
at -20°C (kW)	0.06	0.16	0.32	0.45	0.32	0.45	0.45	0.6	1.2	1.5
Pump data										
Max. flow rate (L / min, gal / min)	60 / 15.9	60 / 15.9	60 / 15.9	60 / 15.9	40 / 10.6	40 / 10.6	40 / 10.6	40 / 10.6	70 / 18.5	70 / 18.5
Max. pressure (bar / psi)	1.4 / 20.3	1.4 / 20.3	1.4 / 20.3	1.4 / 20.3	3.3 / 47.9	3.3 / 47.9	3.3 / 47.9	3.3 / 47.9	6 / 87.0	6 / 87.0
Electrical data & Ordering information										
230V, 50Hz, A	4.5	5.0	6.0	7.0	6.0	7.0	7.5	12.5	-	-
Cat. No.	AAH65002K	AAH65012K	AAH65022K	AAH65032K	AAH65122K	AAH65132K	AAH65142K	AAH65152K	-	-
Cat. No. (EU only)	AAH65005K	AAH65015K	AAH65025K	AAH65035K	AAH65125K	AAH65135K	AAH65145K	AAH65155K	-	-
230V, 60Hz, A	5.0	5.5	6.5	7.5	6.5	7.5	9.5	13.5	-	-
Cat. No.	AAH65001K	AAH65011K	AAH65021K	AAH65031K	AAH65121K	AAH65131K	AAH65141K	AAH65151K	-	-
120V, 60Hz, A	10.0	11.0	13.0	15.0	-	-	-	-	-	-
Cat. No.	AAH65003U	AAH65013U	AAH65023U	AAH65033U	-	-	-	-	-	-
380V, 50Hz, A	-	-	-	-	-	-	-	-	5.0	6.0
Cat. No.	-	-	-	-	-	-	-	-	AAH65168K	AAH65178K
Cat. No. (EU only)	-	-	-	-	-	-	-	-	AAH65160K	AAH65170K
380V, 60Hz, A	-	-	-	-	-	-	-	-	5.5	7.0
Cat. No.	-	-	-	-	-	-	-	-	AAH65169K	AAH65179K

※ According to DIN 12876

※ Temperature stability / Pump data : Water, Cooling capacity : Ethanol

※ Above specification value is recorded by 230V 60Hz. (HL-45H, 55H is recorded by 380V 60Hz)

※ Product performance may be affected by ambient temperature.

Chiller (Recirculating Cooler)

Compact type

Optimized for use with Rotary Evaporators

▣ Structural Functional Features

- Optimized for use with rotary evaporator.
- Pump In/Out port is located at the top, making it easy to connect and disassemble with the evaporator set.
- Optimized slim and compact design minimizes installation space.
- The conventional simple refrigerant compressor control method can cause the compressor life to be shortened by turning the refrigeration system on/off from time to time to control the temperature. However, the Jeio Tech's Chiller is a unique refrigerant control type cooling system completed using our proprietary technology development. Compressor On/Off operation is minimized, and durability of the cooling system is greatly increased.
- Double service valve prevents refrigerant leakage.
- Easy installation with air-cooled type integrated refrigeration system.
- Easy to clean as it is easy to remove refrigerator condenser grill, making it convenient to maintain efficiency of refrigeration.
- Installation of refrigeration system using environmentally friendly refrigerant.
- Pump pressure can be confirmed from the front of the machine, so it is always convenient to check the operation status.
- Bath fluid inlet is wide, making installation and maintenance convenient.



RC-05

▣ Use Convenience Features

- Bright LCD display provides excellent visual perception.
- Identify all operations and operating conditions from the front panel.
- Simultaneously displays the set value and the current value, and settings can be changed during equipment use, making it convenient.
- Highly-reliable control through calibration.
- The refrigeration compressor and pump can each be turned On/Off.
- Level indicator with LED backlight for easy checking of bath fluid level.
- Easy to maintain as a drain valve is included.




▣ Outstanding Safety

- Over temperature alarm.
- Provides a warning when abnormal temperature or bath fluid shortage is detected, and keeps the chiller operating continuously to protect the user application.
- Over-current and short circuit protection of device.
- Malfunction prevented by controller lock function.
- Automatic stop in case of over temperature of the compressor

Pumping capacities

Magnetic pump



※ Application model : RC-05



RC-05

Rotary Evaporator

Specification

Model	RC-05
Interior dimensions	
Max. filling volume (L / cu ft)	5 / 0.18
Refrigerator capacity (HP)	1/3
For tubing dia. (mm / inch)	9.5 / (3/8")
Filling inlet (Ø, mm / inch)	37 / 1.46
Exterior dimensions	
Width (mm / inch)	300 / 11.8
Depth (with drain valve), (mm / inch)	550 / 21.7
Height (with lid), (mm / inch)	530 / 20.9
Weight (kg / lbs)	45 / 99.2
Temperature data	
Working temperature range (°C / °F)	-10 to 30 / 14 to 86
Temperature stability at 15°C (±°C / °F)	1 / 1.8
Cooling capacity, Max	
at 20°C (W)	580
at 10°C (W)	450
at 0°C (W)	350
at -10°C (W)	270
at -20°C (W)	110
Pump data	
Max. flow rate (L / min, gal / min)	26 / 6.87
Max. pressure (bar / psi)	0.7 / 10.15
Electrical requirements (230V, 50/60Hz, A)	
Cat. No.	ACH651012K
Cat. No. (EU only)	ACH651015K

※ According to DIN 12876
 ※ Temperature stability / Pump data: Water, Cooling capacity: Ethanol
 ※ Above specification value is recorded by 230V 60Hz.

Custom-made chillers produced according to customer needs

Customized design, production and installation services with optimized performance for your application.

◆ Expert technical consultation

We provide professional consultation based on our accumulated skills and experience centered on actual customer requirements.

◆ Optimized proposal

We provide optimized proposals that comprehensively consider the intended use and required performance / conditions.



Chiller (Recirculating Cooler)

HH high-temperature precision type

Precise control with $\pm 0.1^{\circ}\text{C}$ stability up to 80°C

Structural Functional Features

- It is high temperature type able to control from -20°C to 80°C , and with precise temperature control-type chiller with $\pm 0.1^{\circ}\text{C}$ stability, has proven performance specifications and high reliability.
- Optimized slim and compact design minimizes installation space.
- The conventional simple refrigerant compressor control method can cause the compressor life to be shortened by turning the refrigeration system on/off from time to time to control the temperature. However, the Jeio Tech's Chiller is a unique refrigerant control type cooling system completed using our proprietary technology development. Compressor On/Off operation is minimized, and durability of the cooling system is greatly increased.
- Double service valve prevents refrigerant leakage.
- Easy installation with air-cooled type integrated refrigeration system.
- Easy to clean as it is easy to remove refrigerator condenser grill, making it convenient to maintain efficiency of refrigeration.
- Installation of refrigeration system using environmentally friendly refrigerant.
- Pump pressure can be confirmed from the front of the machine, so it is always convenient to check the operation status.
- Bath fluid inlet is wide, making installation and maintenance convenient.
- Equipped with removable casters for easy movement and installation.



Use Convenience Features

- Can be controlled by connecting external temperature sensor. (option)
- Convenient control and easy data processing with computer connection and dedicated software.
- Alarm output port provided as standard.
- Identify all operations and operating conditions from the front panel.
- Simultaneously displays the set value and the current value, and settings can be changed during equipment use, making it convenient.
- Highly-reliable control through calibration.
- The refrigeration compressor and pump can be each turned On/Off.
- Level indicator with LED backlight for easy checking of bath fluid level.
- Easy to maintain as a drain valve is included.
- Over temperature alarm.


Outstanding Safety

- Over temperature alarm.
- Provides a warning when abnormal temperature or bath fluid shortage is detected, and keeps the chiller operating continuously to protect the user application.
- Over-current and short circuit protection of device.
- Malfunction prevented by controller lock function.
- Automatic stop in case of over temperature of the compressor.




Pumping capacities

Submersible pump



※ Application model : HH - 15, 20, 25, 35



※ Application model : HH - 45H, 55H

Specification

Model	HH-15	HH-20	HH-25	HH-35	HH-45H	HH-55H
Interior dimensions						
Max. filling volume (L / cu ft)	14 / 0.49	14 / 0.49	25 / 0.88	25 / 0.88	39 / 1.38	39 / 1.38
Refrigerator capacity (HP)	1/2	3/4	1	1.5	2	3
For tubing dia. (mm / inch)	20 / (3/4")	20 / (3/4")	20 / (3/4")	20 / (3/4")	20 / (3/4")	20 / (3/4")
Filling inlet (Ø, mm / inch)	120 / 4.7	120 / 4.7	120 / 4.7	120 / 4.7	120 / 4.7	120 / 4.7
Exterior dimensions						
Width (mm / inch)	515 / 20.3	515 / 20.3	550 / 21.7	550 / 21.7	605 / 23.8	605 / 23.8
Depth (with drain valve), (mm / inch)	645 (715) / 25.4 (28.1)	645 (715) / 25.4 (28.1)	830 (900) / 32.7 (35.4)	830 (900) / 32.7 (35.4)	975 (1045) / 38.4 (41.4)	975 (1045) / 38.4 (41.4)
Height (with lid), (mm / inch)	800 (835) / 31.5 (32.9)	800 (835) / 31.5 (32.9)	975 (1045) / 38.4 (41.1)	975 (1045) / 38.4 (41.1)	1265 (1300) / 49.8 (51.2)	1265 (1300) / 49.8 (51.2)
Weight (kg / lbs)	92 / 202.8	93 / 205	142 / 313	147 / 324	171±10 / 377±20.1	176±10 / 388±20.1
Temperature data						
Working temperature range (°C / °F)	-20 to 80 / -4 to 176	-20 to 80 / -4 to 176	-20 to 80 / -4 to 176	-20 to 80 / -4 to 176	-20 to 80 / -4 to 176	-20 to 80 / -4 to 176
Temperature stability at 15°C (±°C / °F)	0.1 / 0.2	0.1 / 0.2	0.1 / 0.2	0.1 / 0.2	0.1 / 0.2	0.1 / 0.2
Heating capacity (kW)	2	2	4	4	8.4	8.4
Cooling capacity, Max						
at 20°C (kW)	1.3	1.6	2.7	3	5	6.5
at 0°C (kW)	0.9	1.2	1.3	1.8	3	4.1
at -20°C (kW)	0.2	0.3	0.5	0.6	1	1.3
Pump data						
Max. flow rate (L / min, gal / min)	40 / 10.6	40 / 10.6	40 / 10.6	40 / 10.6	70 / 18.5	70 / 18.5
Max. pressure (bar / psi)	3.3 / 47.9	3.3 / 47.9	3.3 / 47.9	3.3 / 47.9	6 / 87.02	6 / 87.02
Electrical data & Ordering information						
230V, 50Hz, A	15.5	15.5	25.0	30.0		
Cat. No.	AAH67012K	AAH67022K	AAH67032K	AAH67042K	-	-
Cat. No. (EU only)	AAH67015K	AAH67025K	AAH67035K	AAH67045K	-	-
230V, 60Hz, A	15.5	16	27	31		
Cat. No.	AAH67012K	AAH67021K	AAH67031K	AAH67041K	-	-
380V, 50Hz, A					17.8	19.3
Cat. No.	-	-	-	-	AAH67158K	AAH67168K
Cat. No. (EU only)	-	-	-	-	AAH67150K	AAH67160K
380V, 60Hz, A					18.3	19.8
Cat. No.	-	-	-	-	AAH67159K	AAH67169K

※ According to DIN 12876

※ Temperature stability / Pump data : Water, Cooling capacity : Ethanol

※ Above specification value is recorded by 230V 60Hz. (HH-45H, 55H is recorded by 380V 60Hz)

※ Product performance may be affected by ambient temperature.

※ Under the condition of no-load operation, cooling capacity should be lower than heating capacity.

Accessories Page 48 Fitting, Connectors, Adapters, Tubing, Tube Clamp, External Sensor

Accessories

Proven reliability & reproducibility

Made and tested according to international standards

- ◆ Meets safety requirements, including electrical safety standards for chillers.
- ◆ Conduction of device performance tests required as per the standards to provide reliable product performance.
- ◆ Verified test results provided as specification to allow for selection of the suitable model for the user.



Barbed Fittings

Cat. No.	Description
00HXE0001066	3/4" Male to barbed fitting for tubing 1" inner dia
00HXE0001067	3/4" Male to barbed fitting for tubing 3/4" inner dia
00HXE0001068	3/4" Male to barbed fitting for tubing 5/8" inner dia
00HXE0001069	3/4" Male to barbed fitting for tubing 1/2" inner dia
00HXE0001070	3/4" Male to barbed fitting for tubing 3/8" inner dia
AAA64501	1" Barbed fittings set with Ball valve
AAA64502	3/4" Barbed fittings set with Ball valve
AAA64503	5/8" Barbed fittings set with Ball valve
AAA64504	1/2" Barbed fittings set with Ball valve
AAA64505	3/8" Barbed fittings set with Ball valve
AAA64506	1" Barbed fittings set with Gate valve
AAA64507	3/4" Barbed fittings set with Gate valve
AAA64508	5/8" Barbed fittings set with Gate valve
AAA64509	1/2" Barbed fittings set with Gate valve
AAA64510	3/8" Barbed fittings set with Gate valve



Barbed Fittings

- Suitable for connecting tubes of various sizes.
- Pipe fittings, tube fittings, hose fittings, etc. can be fitted without tools.

Connectors / Adapters

Cat. No.	Description
00HXE0001105	3/4" Male to 3/4" Female
00HXE0001106	3/4" Male to 5/8" Female
00HXE0001107	3/4" Male to 1/2" Female
00HXE0001108	3/4" Male to 3/8" Female
00HXE0001109	3/4" Male to 3/4" Male
00HXE0001110	3/4" Male to 5/8" Male
00HXE0001111	3/4" Male to 1/2" Male
00HXE0001112	3/4" Male to 3/8" Male



Connectors / Adapters

- Used when connecting tubes and devices.

One Touch Adapters Set

Cat. No.	Description
00HXE0001075	12 mm One touch adapter ID 9.0 mm, OD 12.0 mm
00HXE0001076	10 mm One touch adapter ID 6.5 mm, OD 10.0 mm
00HXE0001077	8 mm One touch adapter ID 5.5 mm, OD 8.0 mm
00HXE0001078	6 mm One touch adapter ID 4.0 mm, OD 6.0 mm
AAA64531	12 mm Adapters set with Ball valve
AAA64532	10 mm Adapters set with Ball valve
AAA64533	8 mm Adapters set with Ball valve
AAA64534	6 mm Adapters set with Ball valve



One Touch Adapters Set

- Easy to connect or replace tubing.

※ Temperature Range : -15°C to 60°C

Flexible Fittings

Cat. No.	Description
00GBA0008589	3/4" Male to 3/4" Flexible fitting
00HXE0001072	3/4" Male to 1/2" Female with 1/2" Male to 1/2" Flexible fitting
AAA64521	3/4" Flexible fittings set with Ball valve
AAA64522	1/2" Flexible fittings set with Ball valve
AAA64523	3/4" Flexible fittings set with Gate valve
AAA64524	1/2" Flexible fittings set with Gate valve



Flexible Fittings

- Tubing is made of stainless steel material.

Tubings

Cat. No.	Description	
00HXE0001079	1 m 12 mm Urethane tubing	ID 9.0 mm, OD 12.0 mm
00HXE0001080	1 m 10 mm Urethane tubing	ID 6.5 mm, OD 10.0 mm
00HXE0001081	1 m 8 mm Urethane tubing	ID 5.5 mm, OD 8.0 mm
00HXE0001082	1 m 6 mm Urethane tubing	ID 4.0 mm, OD 6.0 mm
00HXE0001083	1 m 1" PVC tubing with thread	ID 25.0 mm, OD 31.0 mm
00HXE0001084	1 m 3/4" PVC tubing with thread	ID 19.0 mm, OD 24.0 mm
00HXE0001085	1 m 5/8" PVC tubing with thread	ID 16.0 mm, OD 20.5 mm
00HXE0001086	1 m 1/2" PVC tubing with thread	ID 12.0 mm, OD 16.0 mm
00HXE0001087	1 m 3/8" PVC tubing with thread	ID 10.0 mm, OD 14.0 mm
00HXE0001088	1 m 1" PVC tubing with wire	ID 25.0 mm, OD 33.0 mm
00HXE0001089	1 m 3/4" PVC tubing with wire	ID 19.0 mm, OD 26.0 mm
00HXE0001090	1 m 5/8" PVC tubing with wire	ID 15.0 mm, OD 22.0 mm
00HXE0001091	1 m 1/2" PVC tubing with wire	ID 12.0 mm, OD 18.0 mm
00HXE0001092	1 m 3/8" PVC tubing with wire	ID 9.0 mm, OD 15.0 mm

※ Temperature Range : -15°C to 60°C



Tubings

- Excellent flexibility compared to PVC tubing.
- Excellent heat and corrosion resistance.



Tubing Insulations

- Made from EPDM material, which has excellent flame resistance / oxidation resistance.
- Excellent flexibility at low temperatures.

Tubing Insulations

Cat. No.	Description
00FEA0001453	EPDM Insulation, 25 mm inner dia. (13T)

Flexible Tubing

Cat. No.	Description
00HXE0001128	3/4" Flexible tubing
00HXE0001074	1/2" Flexible tubing



Tube Clamps

- Clamp is made from stainless steel material less susceptible to rust and corrosion.

Tube Clamps

Cat. No.	Description
00HXE0001098	1" Clamp
00HXE0001099	3/4" Clamp
00HXE0001100	5/8" Clamp
00HXE0001101	1/2" Clamp
00HXE0001102	3/8" Clamp



Distributing Fittings

- Multiple tubes can be connected simultaneously.

Distributing Fittings

Cat. No.	Description
AAA64541	3-Way distributing barbed fittings set (1/4")
AAA64542	3-Way distributing barbed fittings set (3/8")

External Sensor (for HH)

Cat. No.	Description
00CFA0001946	Pt 100 sensor with 3m cable



External Sensor

- Possible to check and control temperature of sample by connecting to device with an external sensor



Jeio Tech



How to Control Temperature using External Sensor

Heating Bath



General Application

General experiments: Determination of coliform, dissolution of samples, microbial experiments, bacteria experiments, cell cultures, etc.

	Description	Temp. Range (°C/°F)	Max. Temp. Stability (°C/°F)	Volume (L / cu ft)	Model	Page
Heating Bath	General <ul style="list-style-type: none"> Temperature calibration and timer function. Protection against low fluid level and dry running. 	Amb. +7 to 100 / Amb. +12.6 to 212	0.3 / 0.54 at 50°C	3.5, 11.5, 20 / 0.1, 0.4, 0.7	BW3-G	54
	Agitator <ul style="list-style-type: none"> Improved distribution by agitation. Stirring speed control with the scale knob. 	Amb. +5 to 100 / Amb. +9 to 212	0.1 / 0.18 at 50°C	3.5, 11.5, 20 / 0.1, 0.4, 0.7	BW-H	55
	Visible <ul style="list-style-type: none"> Intuitive viewing window for convenient internal observation Suitable for viscosity measurement as well. 	Amb. +5 to 150 / Amb. +9 to 302	0.1 / 0.18 at 50°C	25, 40 / 0.9, 1.4	VB-G	58
	Shaking <ul style="list-style-type: none"> Provides reliable reciprocating motion Convenient structure for sample handling 	Amb. +5 to 100 / Amb. +9 to 212	0.1 / 0.18 at 50°C	17, 26, 37, 55 / 0.6, 0.9, 1.3, 1.9	BS2	60
Others	Cold Trap <ul style="list-style-type: none"> Trap for water, organic solvents, etc. Possible to use in cooling bath. 	Down to -40 / -40	N/A	4.8 / 0.17	CTB	57

※ The contents of the above and the contents of this catalog may differ depending on the specific model and conditions of use. For the information about the features and specifications that applying to each models, please check the information on the corresponding page of each models.

	Satisfaction with the requirements of European regulations and guidelines.		Real-time equipment monitoring and control system using mobile app.		Registered safety patent based on Jeio Tech's proprietary technology.
	Compliance with Guidelines for Restriction of Hazardous Substances use.		Eco-friendly product that has passed environmental impact assessment tests.		Possible to set end time or start time of device operation.
	PC communication via RS-232 / RS-485 / USB port.		More precise temperature control through 3-point temperature calibration.		More precise temperature control through temperature calibration.
	Over temperature protection function.		Alarm if the water is low level.		It memorizes the state of power failure and operates automatically when power is restored.
	2 year warranty free A/S.				


Heating Bath



Constant temperature bath with proven performance and useful functions applying Jeio Tech's patented safety system

- > **More stable temperature control**
Optimized Microprocessor PID control for more accurate temperature control.
- > **Optimized temperature controls**
Auto-tuning function enabling optimal control even when room temperature or set value is changed frequently.
- > **Temperature calibration function**
Optimized Microprocessor PID control for more accurate temperature control.
- > **Convenient timer function**
Dual wait on/off timer modes.
(1 min. to 999 hr. 59 min.)
- > **Top-grade overheating prevention system (Jeio Tech's patent)**
Even if a failure occurs in the device control unit, the safety function is activated so that it is preferentially operated.
- > **Over temperature shutdown function**
Safety structure with shutdown function using over temperature limiter in case of overheating.
- > **Detection of solution level (BW3, VB, BS2)**
Detecting the level of solution in the bath and stops the device and notifies the user when the fluid level is low.
- > **Convenient solution drainage (BW3, VB, BS2)**
Drain Valve is included as standard, making it easy and convenient to discharge solution in the bath.


※ Some of the above contents are limited to specific models.



Simultaneously displays set temperature and actual temperature.



Safe use with over temperature protection function.



Warning alarm when water level is low.



Provides optimal choice for diverse experiment needs
Model series according to function classification and usage



General



Agitator



Visible




Shaking

Main Feature	Representative water bath with advanced features and specifications	Uniform heating through magnetic stirring	Convenient for internal observation and easy to measure viscosity	Precision water bath with 200rpm shaking function
Technical Feature	Auto-tuning	Auto-tuning	Auto-tuning	Auto-tuning
	Calibration	Calibration	Calibration	Calibration
	Timer	-	Timer	Timer
	-	-	Communication	Communication
	Drain Valve	-	Drain Valve	Drain Valve
Safety	Over Temp.	Over Temp.	Over Temp.	Over Temp.
	Water Level	-	Water Level	Water Level
	Controller Lock	Controller Lock	Controller Lock	Controller Lock

Heating Bath

General type



Representative water bath with advanced features and specifications.

Structural Functional Features

- Drain valve for easy draining of bath solution. (except for BW3-05G model)
- Convenient structure for sample handling with device control part located on the inside.
- Convenient for carrying and moving with handles on the left and right sides.
- The internal bath is integral structure. It is made of stainless steel and has excellent durability.
- Compact design for excellent space utilization.

Use Convenience Features

- Clear VFD panel for clear information display.
- Microprocessor PID method for precise temperature control.
- Optimized temperature control with temperature auto-tuning.
- 3-point calibration function minimizes temperature deviation.
- Wait On/Off timer. (up to 999 hours 59 minutes)
- Save and use 3 frequently used temperatures.
- Automatic restart of operation when power is restored after sudden power failure.

Outstanding Safety

- Top-rated overheating prevention system. (registration KR 10-0397583)
- Over temperature limit function.
- Safe structure in which the heater, sensor, and pump inside the bath are separated by covering plates.
- Malfunction prevented by controller lock function.
- Self-diagnosis function to check abnormal condition of temperature sensor.



BW3-05G



Larger and brighter clear VFD panels intuitively provide device operation information.

Specification

Model		BW3-05G	BW3-10G	BW3-20G
Bath volume (L / cu ft)		3.5 / 0.1	11.5 / 0.4	20 / 0.7
Temperature	Working temp. range (°C / °F)	Amb. +7 to 100 / Amb. +12.6 to 212	Amb. +7 to 100 / Amb. +12.6 to 212	Amb. +7 to 100 / Amb. +12.6 to 212
	Temp. stability at 50°C (±°C / °F)	0.4 / 0.72	0.3 / 0.54	0.3 / 0.54
	Heat up time to 70°C (min)	34	36	44
Dimension	Bath opening / depth (W x L, D) (mm / inch)	240 x 136, 150 / 9.4 x 5.4, 5.9	300 x 240, 200 / 11.8 x 9.5, 7.9	498 x 300, 200 / 19.6 x 11.8, 7.9
	Exterior (W x L x H) (mm / inch)	307 x 200 x 264 / 12.1 x 7.9 x 10.4	360 x 300 x 294 / 14.2 x 11.8 x 11.6	564 x 356 x 294 / 22.2 x 14 x 11.6
	Net weight (kg / lbs)	6.5 / 14.3	10 / 22	18.5 / 40.8
Electrical requirements (230V, 50/60Hz, A)		3	4.3	8.7
Cat. No.		AAH41305K	AAH41315K	AAH41325K
Electrical requirements (120V, 60Hz, A)		5.8	8.3	8.3
Cat. No.		AAH41306U	AAH41316U	AAH41326U

※ The lowest controllable temperature depends on variations in room and sample temperatures. Please contact technical support for consultation before purchasing.

※ Technical data (according to DIN 12876)

※ Bath fluid : Water

Accessories Page 56 Cover, Spring Wire Rack, Half Shelf Adjuster, Test Tube Rack

Heating Bath Agitator type



Uniform heating through magnetic stirring

Structural Functional Features

- Incorporates agitation function by magnetic stirring to mix solution uniformly and improve temperature uniformity.
- Compact design for excellent space utilization.
- The internal bath is integral structure. It is made of stainless steel and has excellent durability.
- The control part at the bottom front is tilted, making it easier to set up the product and check information.
- Convenient for carrying and moving with handles on the left and right sides.

Use Convenience Features

- Microprocessor PID method for precise temperature control.
- Optimized temperature control with temperature auto-tuning.
- Calibration function minimizes temperature deviation.
- Save and use 3 frequently used temperatures.
- Possible to select whether to use the agitation function or not, and scale-type agitation speed control makes it easy to operate.
- Automatic restart of operation when power is restored after sudden power failure.


Outstanding Safety

- Top-rated overheating prevention system. (registration 10-0397583)
- Over temperature limit function.
- Malfunction prevented by controller lock function.
- Safe structure in which the heater, sensor, and pump inside the bath are separated by covering plates.
- If an error occurs during operation, the details of the error are shown on the display and an alarm is sounded.



BW-10H

with Gable Cover & Test Tube Rack (Option)



Circulation of solution by an agitator

Specification

Model		BW-05H	BW-10H	BW-20H
Bath volume (L / cu ft)		3.5 / 0.1	11.5 / 0.4	20 / 0.7
Temperature	Working temp. range (°C / °F)	Amb. +5 to 100 / Amb. +9 to 212	Amb. +5 to 100 / Amb. +9 to 212	Amb. +5 to 100 / Amb. +9 to 212
	Temp. stability at 50°C (±°C / °F)	0.1 / 0.18	0.1 / 0.18	0.1 / 0.18
	Heat up time to 70°C (min)	25	45	50
Dimension	Bath opening / depth (W x L, D) (mm / inch)	240 x 136, 150 / 9.4 x 5.4, 5.9	300 x 240, 200 / 11.8 x 9.5, 7.9	498 x 300, 200 / 19.6 x 11.8, 7.9
	Exterior (W x L x H) (mm / inch)	307 x 216 x 266 / 12.1 x 8.5 x 10.5	364 x 316 x 318 / 14.3 x 14.4 x 12.5	564 x 372 x 318 / 22.2 x 14.6 x 12.5
	Net weight (kg / lbs)	6.5 / 14.3	9.5 / 21	15.5 / 34.2
Electrical requirements (230V, 50/60Hz, A)		3.0	4.3	8.7
Cat. No.		AAH45115K	AAH45125K	AAH45135K
Electrical requirements (120V, 60Hz, A)		5.8	8.3	8.3
Cat. No.		AAH45113U	AAH45123U	AAH45133U

- ※ The lowest controllable temperature depends on variations in room and sample temperatures. Please contact technical support for consultation before purchasing.
- ※ Technical data (according to DIN 12876)
- ※ Bath fluid : Water, Using agitator

Accessories Page 56 Cover, Spring Wire Rack, Half Shelf Adjuster, Test Tube Rack

Accessories



Gable Type

- Prevents condensed water from falling onto specimen with gable roof shape.
- Transparency making it convenient internal observation.




Flat Type

- Made from stainless steel with high corrosion resistance even with long-term use.
- Convenient handle.



Open Ring Type

- Select the ring size according to flask size.
- Multiple rings in one hole allow for adjustment.



Spring Wire Rack

- It has a lattice-shaped double spring structure, which makes it easy to fix various kinds of experiment containers and the height can also be easily adjusted.
- Excellent corrosion resistance as made from stainless steel.



Half Shelf Adjuster

- Easy-to-use pedestal structure when using low-height experiment container.
- Convenient height adjustment.
- Convenient for handling with handles on either side.



Test Tube Rack

- Fix test tubes of various sizes.
- Excellent durability as made from stainless steel.
- Structure allowing tilt adjustment.

Model		BW3-05G	BW3-10G	BW3-20G	BW-05H	BW-10H	BW-20H	
Cover	Gable Type	AAA45531	AAA45532	AAA45533	AAA45531	AAA45532	AAA45533	
	Flat Type	AAA45541	AAA45542	AAA45543	AAA45541	AAA45542	AAA45543	
	Open Ring (mm / inch)	Ø80 / Ø3.1 x 2ea	AAA45521	-	-	AAA45521	-	-
		Ø80 / Ø3.1 x 4ea	-	AAA45522	-	-	AAA45522	-
		Ø80 / Ø3.1 x 8ea	-	-	AAA45523	-	-	AAA45523
		Ø110 / Ø4.3 x 1ea	AAA45524	-	-	AAA45524	-	-
		Ø110 / Ø4.3 x 2ea	-	AAA45525	-	-	AAA45525	-
		Ø110 / Ø4.3 x 6ea	-	-	AAA45527	-	-	AAA45527
		Ø185 / Ø7.3 x 1ea	-	AAA45526	-	-	AAA45526	-
Ø185 / Ø7.3 x 2ea	-	-	AAA45528	-	-	AAA45528		
Rack ¹⁾	Spring Wire Rack	-	AAA45502	AAA45503	-	AAA45502	AAA45503	
	Half Shelf Adjuster	AAA45511	AAA45512	AAA45513	AAA45511	AAA45512	AAA45513	
	Test Tube Rack (mm / inch)	Ø8 / Ø0.3 x 50 tubes	AAA41531	-	-	AAA41531	-	-
		Ø10 / Ø0.4 x 50 tubes	AAA41532	-	-	AAA41532	-	-
		Ø12 / Ø0.5 x 33 tubes	AAA41533	-	-	AAA41533	-	-
		Ø16 / Ø0.6 x 16 tubes	AAA41535	-	-	AAA41535	-	-
		Ø25 / Ø1.0 x 10 tubes	AAA41539	-	-	AAA41539	-	-
		Ø8 / Ø0.3 x 176 tubes	-	AAA45551	AAA45551	-	AAA45551	AAA45551
		Ø10 / Ø0.4 x 176 tubes	-	AAA45552	AAA45552	-	AAA45552	AAA45552
		Ø12 / Ø0.5 x 84 tubes	-	AAA45554	AAA45554	-	AAA45554	AAA45554
		Ø16 / Ø0.6 x 68 tubes	-	AAA45556	AAA45556	-	AAA45556	AAA45556
		Ø25 / Ø1.0 x 33 tubes	-	AAA45561	AAA45561	-	AAA45561	AAA45561
		Max.mounting quantity	1ea	1ea	2ea	1ea	1ea	2ea

1) These accessories cannot be use with flat cover. Please use gable cover.

Cold Trap Bath



Collect organic solvent at as low as -40°C / -40°F

Structural Functional Features

- Capture moisture, organic solvents, etc. at as low as -40°C / -40°F.
- 2 glass traps provided as standard.
- Glass trap is simple to install and fix, and it's easy to handle collected liquid/solvents.
- Can be used in baths that requires rapid low temperature.
- Compact design makes installation convenient.
- Easy to clean condenser filter to maintain refrigeration efficiency with detachable condenser cover.

Use Convenience Features

- Digitally displays the internal temperature at the front of the device.
- Provides the ability to dissolve and discharge ice or solid solvents inside.



CTB-10

Accessories

Glass Trap

- Ability to capture trapped solvents or moisture.

Cat. No.	Model
00FCA0001110	Glass Trap (1ea)



Specification

Model		CTB-10
Temperature range (°C / °F)		Down to -40 / -40
Material	Bath external	Stainless steel, 1.2t, Double painted & baked
	Bath up desk	Stainless steel, 1.5t
	Insulation	Polyurethane foam (30mm)
	Refrigerator	1/2 HP, R-404A (LBP)
Dimension	Volume (L / cu ft)	4.8 / 0.17
	Bath (Ø x H, mm / inch)	210 x 204 / 8.3 x 8
	Glass trap (Ø x H, mm / inch)	60 x 262 / 2.4 x 103
	Exterior (W x D x H) (mm / inch)	360 x 465 x 555 / 14.1 x 18.3 x 21.9
	Net weight (kg / lbs)	46 / 101.4
Electrical requirements (230V, 50/60Hz, A)		4.5
Cat. No.		AAH62012K
Cat. No. (EU only)		AAH62015K

Heating Bath

Visible type

Convenient for internal observation and easy to measure viscosity

Structural Functional Features

- Excellent temperature variation with internal circulation pump.
- Includes 5-step pump circulation speed adjustment function.
- Includes cooling line connecting with external cooling equipment.
- Bath cover provided as standard for easy viscometer installation.
- Convenient to use thermometer hole in bath cover.
- Drain valve for easy draining of bath solution.
- Convenient structure for sample handling due to the dropped location of the internal bath and device control part.

Use Convenience Features

- Temperature adjustment control with 0.01°C resolution.
- Microprocessor PID method for precise temperature control.
- Optimized temperature control with temperature auto-tuning.
- Calibration function minimizes temperature deviation.
- Wait On/Off timer. (up to 99 hours 59 minutes)
- Automatic restart of operation when power is restored after sudden power failure.
- Save and use 3 frequently used temperatures.
- Low level limiter.
- Operation and data management with computer connection.
- Monitor via mobile app anytime, anywhere with LC Connected. (Requires LC GreenBox, option)



VB-40G

with Bath Cover (standard)

Outstanding Safety

- Top-rated overheating prevention system. (registration KR 10-0397583)
- Over temperature limit function.
- Malfunction prevented by controller lock function.
- Self-diagnosis function to check abnormal condition of temperature sensor.

Specification

Model		VB-25G	VB-40G
Bath volume (L / cu ft)		25 / 0.9	40 / 1.4
Temperature	Working temp. range (°C / °F)	Amb. +5 to 150 / Amb. +9 to 302	Amb. +5 to 150 / Amb. +9 to 302
	Temp. stability at 50°C (±°C / °F)	0.1 / 0.18	0.1 / 0.18
	Heat up time to 70°C (min)	35	55
Dimension	Bath opening / depth (W x L, D) (mm / inch)	145 x 265, 330 / 5.7 x 10.4, 13	295 x 265, 330 / 11.6 x 10.4, 13
	Cover lids (mm / inch)	Φ51 x 5ea / Φ2 x 5ea	Φ51 x 8ea / Φ2 x 8ea
	Hole for thermometer (mm / inch)	Φ12.6 x 1ea	Φ0.5 x 1ea
	Exterior (W x L x H) (mm / inch)	492 x 374 x 536 / 19.4 x 14.7 x 21.1	642 x 374 x 536 / 25.3 x 14.7 x 21.1
	Net weight (kg / lbs)	28 / 61.7	37 / 81.6
Electrical requirements (230V, 50/60Hz, A)		13	13
Cat. No.		AAH43302K	AAH43402K

※ The lowest controllable temperature depends on variations in room and sample temperatures. Please contact technical support for consultation before purchasing.



Superior Internal Observation

- The front and back parts of the bath are made of transparent glass and can be easily used for experiments that require internal observation, viscosity measurement, etc.
- The front/rear view windows are made of reinforced glass with double-sided structure on each side and have excellent insulation and tightness.



Convenient Viscometer Fixation

- Hole and hole cover for viscometer installation can be mounted on stainless steel bath cover.
- Viscometer holder (option) to fix the viscometer to the bath to perform viscosity measurement.



Accessories

Viscometer Holder

- Holder that can fix to the hole in the bath cover by fixing the viscometer.
- Convenient structure for detachment and mount of viscometer.

Cat. No.	Max. Mountable Holders	
	VB-25G	VB-40G
AAA43501	5	8



Bath Circulator

- External circulation required for maintaining low temperature when using transparent constant temperature water bath at temperatures under high temperature.
- Connects to the cooling line of the transparent constant temperature water bath.
- For more information on the Bath Circulator, see [Page 24](#)



LC GreenBox

- Real-time monitoring and control of device operation via mobile app, LC Connected.
- Up to 4 devices can be connected.



Description	Dimension (W x D x H, mm / inch)	Cat. No.
LC GreenBox	156 x 94 x 34 / 6.1 x 3.7 x 1.3	AAHQ1011K

Shaking Heating Bath

More precise and stable shaking and temperature control.

Structural Functional Features

- Precise rpm(Max 200rpm) and temperature(Max 100°C) control with superior accuracy.
- Excellent temperature distribution through stirring motion. (50°C±0.1°C)
- Better performance with 3-points temperature calibration and 1-point shaking calibration.
- Stable shaking ability minimizes vibration under high loads.
- PID control optimal for the bath.
- Bath structure with a smooth and circulatory flow of water.
- Extended height for larger sample volumes (maximum 2.8L of flasks).



BS2-10

Use Convenience Features

- Real-time temperature and rpm control through mobile monitoring system(option)
- Scheduling optimization through an automatic on/off timer (up to a maximum of 999 hours and 59 minutes).
- Easy operation and monitoring with a touch color display.
- Drain port of high accessibility located on the side.
- Real-time monitoring and control of RPM even while in operation.

	Uncalibrated	Calibrated
Calibration 1	30.0 °C	31.2 °C
Calibration 2	40.0 °C	39.6 °C
Calibration 3	50.0 °C	49.0 °C

3-point calibration



Visual monitoring and control of RPM during operation

Outstanding Safety

- Audio and visual alarms ring when there is not enough water to prevent overheating.
- Alarm when temperature and speed deviation occur or abnormal condition of temperature sensor.
- Over-current and short circuit protection of device.



Real-time monitoring and control via mobile (option)

Inner Gable Cover

- Minimizes evaporation by using the mounted gable cover.
- Gable roof-shape interior prevents condensed water from falling onto specimen
- Open and close smoothly by applying torque hinge technology.



Spring Wire Rack

- Excellent corrosion resistance as made from stainless steel (1ea: standard)
- Easy to maintain as structure makes for convenient attachment to body.
- It has a lattice-shaped double spring structure, which makes it easy to fix various kinds of experiment containers and the height can also be easily adjusted.

Cat. No.	Model
AAA44571	BS2-10
AAA44572	BS2-15
AAA44573	BS2-20
AAA44574	BS2-30





Specification

Model		BS2-10	BS2-15	BS2-20	BS2-30
Bath volume (L / cu ft)		17 / 0.6	26 / 0.9	37 / 1.3	55 / 1.9
Temperature	Working temp. range (°C / °F)	Amb. +5 to 100 / Amb. +9 to 212	Amb. +5 to 100 / Amb. +9 to 212	Amb. +5 to 100 / Amb. +9 to 212	Amb. +5 to 100 / Amb. +9 to 212
	Temp. stability at 50°C (± °C / °F)	0.1 / 0.18	0.1 / 0.18	0.1 / 0.18	0.1 / 0.18
Shaking system	Motion type	Reciprocal	Reciprocal	Reciprocal	Reciprocal
	Frequency (rpm)	30 to 200	30 to 200	30 to 200	30 to 200
	Stroke length (mm / inch)	25 / 0.98	25 / 0.98	30 / 1.18	30 / 1.18
Dimension	Bath opening, depth (W x L, D) (mm / inch)	143 x 269, 192 / 5.6 x 10.6, 7.6	263 x 269, 192 / 10.4 x 10.6, 7.6	415 x 269, 192 / 16.3 x 10.6, 7.6	565 x 319, 192 / 22.2 x 12.6, 7.6
	Exterior (W x L x H) (mm / inch)	474 x 417 x 426 / 18.7 x 16.4 x 16.8	594 x 417 x 426 / 23.5 x 16.4 x 16.8	746 x 417 x 426 / 29.4 x 16.4 x 16.8	882 x 467 x 426 / 34.7 x 18.4 x 16.8
	Net weight (kg / lbs)	23 / 51	27 / 60	31 / 68	38 / 84
Electrical requirements (230V/60Hz, A)		3.7	4.6	8.9	8.9
Cat. No.		AAH44105K	AAH44155K	AAH44205K	AAH44305K
Electrical requirements (120V, 60Hz, A)		7.1	8.8	-	-
Cat. No.		AAH44103U	AAH44153U	-	-

- ※ The lowest controllable temperature depends on variations in room and sample temperatures. Please contact technical support for consultation before purchasing.
- ※ Technical data (according to DIN 12876)
- ※ Bath fluid : Water

Accessories

Flask Clamp

- Clamp being able to hold an erlenmeyer flask.
- Excellent durability as made of stainless steel.
- Equipped with flask clamps of the desired size/ quantity on the universal platform and used inside the bath.

Cat. No.	Flask Clamp	Max. Mountable Flask Clamps			
		BS2-10	BS2-15	BS2-20	BS2-30
AAA23550	50 ml	9	12	20	32
AAA23551	100 ml	3	9	15	24
AAA23552	250 ml	2	4	8	12
AAA23556	300 ml	2	4	8	12
AAA23553	500 ml	2	4	6	12
AAA23554	1 L	1	2	4	6
AAA23555	2 L ¹⁾	-	1	2	3
AAA23557	2.8 L ¹⁾	-	1	1	2

1) 2 L and 2.8 L flasks can be used with the bath cover opened.



Universal Platform

- Accessories for mounting flask clamp, test tube rack, etc. inside bath.
- Can be used instead of the spring part in the spring wire rack.

Cat. No.	Model
AAA44521	BS2-10
AAA44522	BS2-15
AAA44523	BS2-20
AAA44524	BS2-30



Dedicated Platform

- Set of 250 ml/500 ml flask clamps mounted on Universal Platform.
- Can be used instead of the spring part in the spring wire rack.

Cat. No.	Flask Clamps	Model
AAA44531	250 ml x 2ea	BS2-10
AAA44532	500 ml x 2ea	
AAA44541	250 ml x 4ea	BS2-15
AAA44542	500 ml x 4ea	
AAA44551	250 ml x 8ea	BS2-20
AAA44552	500 ml x 6ea	
AAA44561	250 ml x 12ea	BS2-30
AAA44562	500 ml x 12ea	



Test Tube Rack

- Fix test tubes of various sizes.
- Excellent durability as made of stainless steel.
- Structure allows tilt adjustment.
- Can be used instead of the spring part in the spring wire rack.

Cat. No.	Description (mm / inch)	Max. Quantity of Racks
AAA44581	ø8 / ø0.3 x 86 Tubes	BS2-10 : 1ea BS2-15 : 2ea BS2-20 : 3ea BS2-30 : 4ea
AAA44582	ø10 / ø0.4 x 86 Tubes	
AAA44583	ø12 / ø0.5 x 58 Tubes	
AAA44585	ø16 / ø0.6 x 32 Tubes	
AAA44586	ø25 / ø1.0 x 19 Tubes	



Stainless Steel Cover

- Sturdy, corrosion-resistant, stainless steel material.

Cat. No.	Model
AAA44600	BS2-10
AAA44601	BS2-15
AAA44602	BS2-20
AAA44603	BS2-30

LC GreenBox

- Real-time monitoring the state of a device via mobile app (LC Connected) under an internet connection.
- Available to connect 4 devices in max for 1 unit of LC GreenBox.

Cat. No.	Dimension (W x D x H, mm / inch)
AAHQ1011K	156 x 94 x 34 / 6.14 x 3.7 x 1.34

Oven & Furnace



- | | | |
|--|--|---|
|  <p>Satisfaction with the requirements of European regulations and guidelines.</p> |  <p>Registered safety patent based on Jeio Tech's proprietary technology.</p> |  <p>Compliance with Guidelines for Restriction of Hazardous Substances use.</p> |
|  <p>Ensuring the reliability of electronic records and signatures to support GMP audit.</p> |  <p>Intuitive operation with color touch display.</p> |  <p>Digital Input / Output for external device interlock control.</p> |
|  <p>Real-time equipment monitoring and control system using mobile app.</p> |  <p>Real-time monitoring to check status of device.</p> |  <p>More precise temperature control through 3-point temperature calibration.</p> |
|  <p>More precise temperature control through temperature calibration.</p> |  <p>Data storage using USB ports.</p> |  <p>Electronic over temperature protection system.</p> |
|  <p>Over temperature protection function.</p> |  <p>Warning alarm in case of deviation of temperature control.</p> |  <p>Possible to set end time or start time of device operation.</p> |
|  <p>Control of temperature / humidity / rpm / time, etc. through dedicated program.</p> |  <p>PC communication via RS-232 / RS-485 / USB port.</p> |  <p>It memorizes the state of power failure and operates automatically when power is restored.</p> |
|  <p>2 year warranty free A/S.</p> | | |

General Application

Forced convection oven: Aging and curing, agricultural genetics, degradation/digestion of proteins and starches, universal/general drying/heating, QC/batch testing, stability testing, sterilization.

Natural convection oven : Annealing, plant tissue drying, drug metabolism, purified chemical precipitation/fine chemical precipitation, universal/general drying/heating, serum protein analysis, sterilization, etc.

Vacuum oven: Agricultural genetics, annealing, desiccating, outgassing, vacuum embedding, volatile resin, polymer test, etc.

Furnace: New material development, superconducting material heating, ceramic sintering, heat treatment reaction, melting of glass and metal, casting moulding, organic matter weight analysis and coal ash measurement, heat resistance test, etc.

	Description	Temp. Range (°C/°F)	Max. Temp. Fluctuation (°C/°F)	Volume (L / cu ft)	Model	Page
Forced Convection Ovens	Program	<ul style="list-style-type: none"> • 10 programs & 10 steps • Easy operation with color touch 	Amb. +10 to 250 / Amb. +18 to 482	±0.3 / 0.54 at 100°C	26, 61, 112, 147 / 0.9, 2.1, 3.9, 5.1	OF4-P 66
	Advanced	<ul style="list-style-type: none"> • Most representative forced convection type • Includes safety device as standard 	Amb. +10 to 250 / Amb. +18 to 482	±0.3 / 0.54 at 100°C	26, 61, 112, 147 / 0.9, 2.1, 3.9, 5.1	OF4-V 68
	Large	<ul style="list-style-type: none"> • Large capacity standard type oven • Uniform and precise temperature control 	Amb. +15 to 200 / Amb. +27 to 392 Amb. +15 to 300 / Amb. +27 to 572	±0.2 / 0.36 at 100°C	314, 450, 760 / 11.1, 15.9, 26.8	OF3 OF3-H 70
	Standard	<ul style="list-style-type: none"> • Stable product that is faithful to the basics • Includes safety device as standard 	Amb. +10 to 250 / Amb. +18 to 482	±0.3 / 0.54 at 100°C	26, 61, 112, 147 / 0.9, 2.1, 3.9, 5.1	OF4-S 72
	Multi	<ul style="list-style-type: none"> • Independent temperature control by chamber • 4-in-1 	Amb. +10 to 200 / Amb. +18 to 392	±0.5 / 0.9 at 100°C	240 / 8.5 (60/2.1 x 4 Chambers)	OF-02G -4C 73
Natural Convection Ovens	Advanced	<ul style="list-style-type: none"> • Most representative natural convection type • Suitable for drying powder samples 	Amb. +15 to 250 / Amb. +27 to 482	±0.5 / 0.9 at 100°C	26, 58, 110, 147 / 0.92, 2.05, 3.86, 5.19	ON4-V 74
	Standard	<ul style="list-style-type: none"> • Stable product that is faithful to the basics • Includes safety device as standard 	Amb. +15 to 250 / Amb. +27 to 482	±0.5 / 0.9 at 100°C	26, 58, 110, 147 / 0.92, 2.05, 3.86, 5.19	ON4-S 76
Vacuum Oven	Vacuum	<ul style="list-style-type: none"> • Product optimized for vacuum drying • Structure separating vacuum/release line 	Amb. +15 to 250 / Amb. +27 to 482	±0.2 / 0.36 at 100°C	28, 65 / 1.0, 2.3	OV4 78
Furnace	Furnace	<ul style="list-style-type: none"> • High temperature test up to 1100°C • Convenient program control 	Max. 1100 / 2012	±1.0 / 1.8 at 1100°C	3.0, 4.5, 12.0 / 0.1, 0.2, 0.4	MF2-G 82

※ The contents of the above and the contents of this catalog may differ depending on the specific model and conditions of use. For the information about the features and specifications that applying to each models, please check the information on the corresponding page of each models.



Oven Series with proven performance and durability Minimization of energy loss and provision of precise control

> **Variety of products optimized for specific purpose of use**

Optimal model provided based on heat transfer method, control temperature range, chamber capacity and program control.

> **Verified according to international standards**

Fluctuation, variation, heating time, and recovery time are tested according to DIN 12880.

> **Minimization of energy loss**

Superior insulation and enclosed design minimize energy loss and reduce operating costs.

> **Provision of remote equipment monitoring service (option)**

Real-time monitoring of equipment operation status using smartphone. (when purchased LC GreenBox)

> **Stable and accurate temperature control**

Optimized microprocessor PID control for accurate temperature control.

> **3-point temperature calibration function (OF4 / ON4 / OV4)**

Temperature calibration of up to 3 points allows for more reliable temperature control.

> **Optimized temperature control**

Auto-tuning function enables optimal control even when temperature or ambient conditions are changed frequently.

> **Computer interface control**

Convenient remote control and data management through built-in dedicated software. (applicable models)



Inside of forced convection oven



Hands-free door system



Provision of remote equipment monitoring service



Has functions that are convenient for actual use and applies Jeio Tech's patented safety system

- > **Convenient timer function provided**
Device is run (Wait On/Wait Off) after the set time elapses. (Max. 999 hours and 59 minutes)
- > **Controller lock function**
Lock function for preventing change of controller settings prevents erroneous operation due to arbitrary changing.
- > **Convenient temperature save and operation**
Conveniently save and use frequently used temperatures. (up to 3 presets allowed)
- > **Safe for use at maximum temperature**
According to international standards, even when operated at maximum temperature, the surface temperature does not exceed 51°C.
- > **Top-grade overheating prevention system (Jeio Tech's patent)**
Even if a failure occurs in the device control unit, the safety device is structured so that it is preferentially operated.
- > **Overheat temperature cut-off function**
Safety structure with shutdown function using Over Temperature Limiter in case of overheating.
- > **Door open detection**
Stops the device and provides notification/alarms when the door is open for an extended period of time.
- > **Prepared for unexpected power failure**
Auto-run function capable of selecting re-operation when power is restored after sudden power failure.

※ Some of the above contents are limited to specific models.



5-inch touch display (OF4, ON4, OV4, MF2)



General controller with clear VFD panel (OF3)



Double-stacked for space efficiency (Option) (OF4, ON4)

Forced Convection Oven

Program type

Convenient program control function

Program Type

- User interface that provides a quick view of temperature changes (heating, stable, and cooling) and fan speed also.
- Set various test condition using 10 programs and 10 steps, which can be repeated up to 99 times.
- Setting time up to 99 hours 59 minutes per step.
- Start temperature, temperature and time for each step, and Opti-flow™ all adjustable as desired.
- Scheduler function to operate at scheduled times.



Use Convenience Features

- Opti-flow™, a function to flexibly adjust the fan speed according to the specimen type and test conditions.
- Scheduler function to operate at scheduled times.
- Test data is automatically saved and easily downloaded via USB easy to edit and use. (in CSV format)
- Real-time graphs that allow users to view events and alarms at a glance so that the test progress can be easily and intuitively monitored.
- Ramping rates (°C/min) flexibly adjustable for each test performed
- RS-232/USB port and compatible software provided to enable easy device control and data management using PC.
- 36 latest events, such as errors and door opened or closed, automatically recorded, helping user obtain feedback on experimental results.
- Monitor and control via mobile app anytime, anywhere with LC Connected. (Requires LC GreenBox, option)
- Integrated with 21 CFR Part 11 compliant LC DataKeeper, suitable for GMP audits (option)

OF4-10PW


with Wire Shelves 2ea (standard)

Outstanding Safety


- Top-rated overheating protection system. (registration KR 10-0397583)
- BMS communication port that helps monitor the device temperature on the building management system, enabling quick response in case of abnormality. (option)
- Device surface not overheating even in operation at the highest temperature. (70°C or below, at 250°C)
- Over temperature limiter function
- Triple-toughened viewing window for safe and convenient observation even at high temperatures.
- Self-diagnostics function to check abnormal condition of temperature sens



Easy and intuitive program setting interface




Fan speed can be adjusted manually according




Structural Functional Features

- Wide temperature range of up to 250°C.
- Excellent visibility and easy control with a 5-inch touch screen
- Airflow optimized for uniform heat transfer.
- Superior insulation and enclosed design minimize energy loss and reduce operating costs.
- Hands-free™ door system to allow users to easily open the door while carrying samples with both hands.
- Double stackable with stacking set (option) maximizes space utilization (-05 / -15 models only)
- Vent holes for gas discharge at the top.
- Compact models designed to be used in a fume hood available. (OF4-03)



Graph generation to visualize both stored and real-time data



Scheduler function to operate at scheduled times



Real-time monitoring and control via mobile (option)

Specification

Model	Solid door (without window)	OF4-03P	OF4-05P	OF4-10P	OF4-15P
	Window door	OF4-03PW	OF4-05PW	OF4-10PW	OF4-15PW
Chamber volume (L / cu ft)		26 / 0.92	61 / 2.15	112 / 3.96	147 / 5.19
Temperature¹⁾	Range (°C / °F)	Amb. + 10 to 250 / Amb. + 18 to 482			
	Fluctuation at 100°C (±°C / ±°F)	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54
	Variation at 100°C (±°C / ±°F)	1.4 / 2.52	1.9 / 3.42	1.7 / 3.06	1.5 / 2.70
	Recovery time at 100°C (min)	4	4	4	4
Control panel	5 inch TFT Color LCD				
Communication interface	USB-B, RS-232				
Function	Program setting	10 program, 10 step, start temp.			
	Scheduler	Set day, time, and program No.			
	Recently event record	36 ea			
	Fan speed adjustable	Three-level wind speed settings			
	Recommended temperature	Three frequently used temperature settings suggested			
	Graph	Available checking the stored and real-time data as graphs			
	Data storage	Save the test data in CSV format			
Dimensions	Temp. ramping rate (°C/min. / °F/min.)	0 to 13 / 0 to 23.4	0 to 9 / 0 to 16.2	0 to 6 / 0 to 10.8	0 to 5 / 0 to 9
	Interior (W x D x H, mm / inch)	300 x 210 x 400 / 11.8 x 8.3 x 15.8	380 x 330 x 480 / 15.0 x 13.0 x 18.9	485 x 350 x 655 / 19.1 x 13.8 x 25.8	550 x 376 x 710 / 21.7 x 14.8 x 28.0
	Exterior (W x D x H, mm / inch)	478 x 527 x 604 / 18.8 x 20.8 x 23.8	558 x 648 x 684 / 22.0 x 25.5 x 26.9	663 x 677 x 859 / 26.1 x 26.7 x 33.8	728 x 704 x 914 / 28.7 x 27.7 x 36.0
	Net weight (kg / lbs)	33 / 72.8	53 / 116.8	63 / 138.9	73 / 160.9
Shelves	Quantity of shelves (Standard / max.)	2 / 4	2 / 5	2 / 8	2 / 9
	Max. Load per shelf (kg / lbs)	20 / 44.09			
Number of air change per hour at 100°C²⁾	123	22	18	11	
Electrical requirements (230V, 50/60Hz, A)	6.4				
Cat. No.	Solid door (without window)	AAH12915K	AAH12925K	AAH12935K	AAH12945K
	Window door	AAH12955K	AAH12965K	AAH12975K	AAH12985K
Electrical requirements (120V, 60Hz, A)	12.3				
Cat. No.	Solid door (without window)	AAH12916U	AAH12926U	AAH12936U	AAH12946U
	Window door	AAH12956U	AAH12966U	AAH12976U	AAH12986U

1) Data measured at fan speed 3, with possible increase of Fluctuation and Variation values at fan speeds 1 and 2 (According to DIN 12880, pre-2013)
 2) Data measured at fan speed 3, with possible decreases in Air change at fan speeds 1 and 2. (According to ASTM D 5374, reapproved 1999)

Forced Convection Oven

Advanced type

Improved functionality and reliable temperature control performance

Excellent Functionality

- Wide temperature range of up to 250°C.
- Excellent visibility and easy control with a 5-inch touch screen
- Airflow optimized for uniform heat transfer.
- Superior insulation and enclosed design minimize energy loss and reduce operating costs.
- Hands-free™ door system to allow users to easily open the door while carrying samples with both hands.
- Double stackable with stacking set (option) maximizes space utilization. (-05 / -15 models only)
- Vent holes for gas discharge at the top.
- Compact models designed to be used in a fume hood available. (OF4-03)



Use Convenience Features

- Opti-flow™, a function to flexibly adjust the fan speed according to the specimen type and test conditions.
- Three-point temperature calibration that ensures high-temperature accuracy over a wide temperature range.
- Ramping rates (°C/min) flexibly adjustable for each test performed
- 36 latest events, such as errors and door opened or closed, automatically recorded, helping user obtain feedback on experimental results.
- Monitor and control via mobile app anytime, anywhere with LC Connected. (Requires LC GreenBox, option)

OF4-10VW


with Wire Shelves 2ea (standard)

Convenient checking for the test progress and data


- RS-232/USB port and compatible software provided to enable easy device control and data management using PC.
- Test data is automatically saved and easily downloaded via USB; easy to edit and use. (in CSV format)
- Real-time graphs that allow users to view events and alarms at a glance.
- Integrated with 21 CFR Part 11 compliant LC DataKeeper, suitable for GMP audits (option)



5-Inch touch display



Fan speed can be adjusted manually according



Outstanding Safety


- Top-rated overheating protection system. (registration KR 10-0397583)
- BMS communication port that helps monitor the device temperature on the building management system, enabling quick response in case of abnormality. (option)
- Device surface not overheating even in operation at the highest temperature. (70°C or below, at 250°C)
- Over temperature limiter function
- Triple-toughened viewing window for safe and convenient observation even at high temperature.(OF4-VW)
- Self-diagnostics function to check abnormal condition of temperature sensor.



Fumes inside the chamber effectively removed



Data storage using USB ports



Safe surface of the device even in operation at the highest temperature

Specification

Model	Solid door (without window)	OF4-03V	OF4-05V	OF4-10V	OF4-15V
	Window door	OF4-03VW	OF4-05VW	OF4-10VW	OF4-15VW
Chamber volume (L / cu ft)		26 / 0.92	61 / 2.15	112 / 3.96	147 / 5.19
Temperature¹⁾	Range (°C / °F)	Amb. + 10 to 250 / Amb. + 18 to 482			
	Fluctuation at 100°C (±°C / ±°F)	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54
	Variation at 100°C (±°C / ±°F)	1.4 / 2.52	1.9 / 3.42	1.7 / 3.06	1.5 / 2.70
	Heating time to 100°C (min)	8	10	14	15
	Recovery time at 100°C (min)	4	4	4	4
Control panel	5 inch TFT Color LCD				
Communication interface	USB-B, RS-232				
Function	Recently event record	36 ea			
	Fan speed adjustable	Three-level wind speed settings			
	Recommended temperature	Three frequently used temperature settings suggested			
	Graph	Available checking the stored and real-time data as graphs			
	Data storage	Save the test data in CSV format			
	Temp. ramping rate (°C/min. / °F/min.)	0 to 13 / 0 to 23.4	0 to 9 / 0 to 16.2	0 to 6 / 0 to 10.8	0 to 5 / 0 to 9
Dimensions	Interior (W x D x H, mm / inch)	300 x 210 x 400 / 11.8 x 8.3 x 15.8	380 x 330 x 480 / 15.0 x 13.0 x 18.9	485 x 350 x 655 / 19.1 x 13.8 x 25.8	550 x 376 x 710 / 21.7 x 14.8 x 28.0
	Exterior (W x D x H, mm / inch)	478 x 527 x 604 / 18.8 x 20.8 x 23.8	558 x 648 x 684 / 22.0 x 25.5 x 26.9	663 x 677 x 859 / 26.1 x 26.7 x 33.8	728 x 704 x 914 / 28.7 x 27.8 x 36.0
	Net weight (kg / lbs)	33 / 72.8	53 / 116.8	63 / 138.9	73 / 160.9
Shelves	Quantity of shelves (standard / max.)	2 / 4	2 / 5	2 / 8	2 / 9
	Max. Load per shelf (kg / lbs)	20 / 44.09			
Number of air change per hour at 100°C²⁾	123	22	18	11	
Electrical requirements (230V, 50/60Hz, A)	6.4				
Cat. No.	Solid door (without window)	AAH12715K	AAH12725K	AAH12735K	AAH12745K
	Window door	AAH12815K	AAH12825K	AAH12835K	AAH12845K
Electrical requirements (120V, 60Hz, A)	12.3				
Cat. No.	Solid door (without window)	AAH12716U	AAH12726U	AAH12736U	AAH12746U
	Window door	AAH12816U	AAH12826U	AAH12836U	AAH12846U

1) Data measured at fan speed 3, with possible increase of Fluctuation and Variation values at fan speeds 1 and 2 (According to DIN 12880, pre-2013)
 2) Data measured at fan speed 3, with possible decreases in Air change at fan speeds 1 and 2. (According to ASTM D 5374, reapproved 1999)

Forced Convection Oven

Large type

Up to 300°C and program/general control selection

Structural Functional Features

- Optimized model provision and selection based on temperature range, capacity, controller, etc.
- Air flow optimized for uniform heat transfer ensures stable overall temperature uniformity.
- Includes vent hole base with cover for gas exhaust, cable connection, etc. on the side.
- Includes castors making it easy to move and install.

Use Convenience Features

- Monitor via mobile app anytime, anywhere with LC Connected. (Requires LC GreenBox, option)
- Temperature calibration of up to 3 points allows for more reliable temperature control.
- Optimized control with temperature auto-tuning.
- Microprocessor PID method for precise temperature control.

Outstanding Safety

- Top-rated overheating protection system. (registration KR 10-0397583)
- Over temperature limit function.
- Safe product surface, even when operated at maximum temperature.
- Triple-toughened viewing window for safe and convenient observation even at high temperatures. (OF3-W)
- Malfunction prevented by controller lock function.
- Self-diagnostics function to check abnormal condition of temperature sensor.

Program Type

- Possible to set and save 10 programs.
- Up to 10 steps can be set per program.
- Set time up to 99 hours and 59 minutes per step.
- Includes useful scheduler management functions such as program operation range and repeat interval setting.
- Program repetition possible up to 99 times.
- Digital Input / Output function makes it easy to connect external equipment.



General Type

- 4.7-inch large VFD panel for clear information display.
- Easy and convenient to use with simple operation.
- Notifies when deviation from control temperature upper limit and lower limit occurs.
- USB/RS-232 connection and software provide convenient computer operation and data managing.
- Save and use 3 frequently used temperatures.



OF3-45PW

with Wire Shelves 2ea (standard)



General Controller

Simultaneous display of set temperature and current temperature.



Program Controller

- Color touch display
- Intuitive operation environment.



Specification

		Max. 200°C		
Model	General (Without Window)	OF3-30	OF3-45 (CE Completed)	OF3-75 (CE Completed)
	General (Window Door)	OF3-30W	OF3-45W (CE Completed)	OF3-75W (CE Completed)
	Program (Without Window)	OF3-30P	OF3-45P (CE Completed)	OF3-75P (CE Completed)
	Program (Window Door)	OF3-30PW	OF3-45PW (CE Completed)	OF3-75PW (CE Completed)
Chamber volume (L / cu ft)		314 / 11.1	450 / 15.9	760 / 26.8
Temperature	Range (°C / °F)	Amb. +15 to 200 / Amb.+27 to 392	Amb. +15 to 200 / Amb.+27 to 392	Amb. +15 to 200 / Amb.+27 to 392
	Fluctuation at 100°C (±°C / °F)	0.2 / 0.36	0.2 / 0.36	0.2 / 0.36
	Variation at 100°C (±°C / °F)	3.5 / 6.3	3.5 / 6.3	4.5 / 8.1
	Heating time to 100°C (min.)	20	25	25
	Recovery time at 100°C (min.)	7	7	7
Dimension	Interior (W x D x H, mm / inch)	645 x 750 x 650 / 25.4 x 29.5 x 25.6	645 x 750 x 930 / 25.4 x 29.5 x 36.6	745 x 850 x 1200 / 29.3 x 33.5 x 47.2
	Exterior (W x D x H, mm / inch)	850 x 1100 x 1235 / 33.5 x 43.3 x 48.6	850 x 1100 x 1525 / 33.5 x 43.3 x 60	950 x 1200 x 1785 / 37.4 x 47.2 x 70.3
	Net weight (Kg / lbs)	200 / 440.9	240 / 529.1	300 / 661.4
Shelves	No. of wire shelf (standard / max.)	2 / 5	3 / 8	3 / 11
	Max. load per shelf (kg / lbs)	25 / 55.1	25 / 55.1	25 / 55.1
Electrical requirements (230V, 50/60Hz, A)		1P, 13.2		
Cat. No.	General (Without Window)	AAH18215K	-	-
	General (Window Door)	AAH18315K	-	-
	Program (Without Window)	AAH18615K	-	-
	Program (Window Door)	AAH18715K	-	-
Electrical requirements (380V, 50/60Hz, A)		-	3P, 6.7	3P, 9.0
Cat. No.	General (Without Window)	-	AAH182214K	AAH182314K
	General (Window Door)	-	AAH183214K	AAH183314K
	Program (Without Window)	-	AAH186214K	AAH186314K
	Program (Window Door)	-	AAH187214K	AAH187314K

※ Technical data (according to DIN 12880 before 2013)

		Max. 300°C		
Model	General (Without Window)	OF3-30H (CE Completed)	OF3-45H (CE Completed)	OF3-75H (CE Completed)
	General (Window Door)	OF3-30HW (CE Completed)	OF3-45HW (CE Completed)	OF3-75HW (CE Completed)
	Program (Without Window)	OF3-30HP (CE Completed)	OF3-45HP (CE Completed)	OF3-75HP (CE Completed)
	Program (Window Door)	OF3-30HPW (CE Completed)	OF3-45HPW (CE Completed)	OF3-75HPW (CE Completed)
Chamber volume (L / cu ft)		314 / 11.1	450 / 15.9	760 / 26.8
Temperature	Range (°C / °F)	Amb. +15 to 300 / Amb.+27 to 572	Amb. +15 to 300 / Amb.+27 to 572	Amb. +15 to 300 / Amb.+27 to 572
	Fluctuation at 100°C (±°C / °F)	0.2 / 0.36	0.2 / 0.36	0.2 / 0.36
	Variation at 100°C (±°C / °F)	3.5 / 6.3	3.5 / 6.3	4.5 / 8.1
	Heating time to 100°C (min.)	15	20	20
	Recovery time at 100°C (min.)	7	7	7
Dimension	Interior (W x D x H, mm / inch)	645 x 750 x 650 / 25.4 x 29.5 x 25.6	645 x 750 x 930 / 25.4 x 29.5 x 36.6	745 x 850 x 1200 / 29.3 x 33.5 x 47.2
	Exterior (W x D x H, mm / inch)	850 x 1100 x 1235 / 33.5 x 43.3 x 48.6	850 x 1100 x 1525 / 33.5 x 43.3 x 60	950 x 1200 x 1785 / 37.4 x 47.2 x 70.3
	Net weight (Kg / lbs)	200 / 440.9	240 / 529.1	300 / 661.4
Shelves	No. of wire shelf (standard / max.)	2 / 5	3 / 8	3 / 11
	Max. load per shelf (kg / lbs)	25 / 55.1	25 / 55.1	25 / 55.1
Electrical requirements (380V, 50/60Hz, A)		3P, 6.7	3P, 9.0	3P, 10.8
Cat. No.	General (Without Window)	AAH192114K	AAH192214K	AAH192314K
	General (Window Door)	AAH193114K	AAH193214K	AAH193314K
	Program (Without Window)	AAH196114K	AAH196214K	AAH196314K
	Program (Window Door)	AAH197114K	AAH197214K	AAH197314K

※ The lowest controllable temperature depends on variations in room and sample temperatures. Please contact technical support for consultation before purchasing.


※ Technical data (according to DIN 12880 before 2013)

Oven with uniform horizontal airflow (up to 350°C), check [Page 346](#) for industrial heating chamber



Accessories [Page 84](#) Wire Shelves, Perforated Shelves, Recorder, LC GreenBox

Forced Convection Oven Standard type



Provide the same temperature uniformity performance as advanced/programmable types

Stable Temperature Control

- Wide temperature range of up to 250°C.
- Airflow optimized for uniform heat transfer.
- Superior insulation and enclosed design minimize energy loss and reduce operating costs.
- Three-point temperature calibration that ensures high-temperature accuracy over a wide temperature range.
- Available to set the upper and lower temperature limit, alarm activated to quick response when occurring the deviation.

Use Convenience Features

- Recording the 12 most recent events, such as start/stop of operation, the operating time, and error occurrence.
- Monitor and control via mobile app anytime, anywhere with LC Connected. (Requires LC GreenBox, option)
- BMS communication port that helps monitor the device temperature on the building management system. (option)

Outstanding Safety

- Top-rated overheating protection system. (registration KR 10-0397583)
- Device surface not overheating even in operation at the highest temperature. (70°C or below, at 250°C)
- Over temperature limiter function
- Triple-toughened viewing window for safe and convenient observation even at high temperatures.



OF4-10SW

with Wire Shelves 2ea (standard)



3.5 Inch touch display

Specification

Model	Solid door (without window)	OF4-03S	OF4-05S	OF4-10S	OF4-15S
	Window door	OF4-03SW	OF4-05SW	OF4-10SW	OF4-15SW
Chamber volume (L / cu ft)		26 / 0.92	61 / 2.15	112 / 3.96	147 / 5.19
Temperature¹⁾	Range (°C / °F)	Amb. +10 to 250 / Amb. +18 to 482			
	Fluctuation at 100°C (±°C / ±°F)	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54
	Variation at 100°C (±°C / ±°F)	1.4 / 2.52	1.9 / 3.42	1.7 / 3.06	1.5 / 2.70
	Heating time to 100°C (min)	8	10	14	15
	Recovery time at 100°C (min)	4	4	4	4
Control panel		3.5 inch TFT Color LCD			
Communication interface		USB-B, RS-232 (option)			
Function	Recently event record	12 ea			
Dimensions	Interior (W x D x H, mm / inch)	300 x 210 x 400 / 11.8 x 8.3 x 15.8	380 x 330 x 480 / 15.0 x 13.0 x 18.9	485 x 350 x 655 / 19.1 x 13.8 x 25.8	550 x 376 x 710 / 21.7 x 14.8 x 28.0
	Exterior (W x D x H, mm / inch)	478 x 527 x 604 / 18.8 x 20.8 x 23.8	558 x 648 x 684 / 22.0 x 25.5 x 26.9	663 x 677 x 859 / 26.1 x 26.7 x 33.8	728 x 704 x 914 / 28.7 x 27.7 x 36.0
	Net weight (kg / lbs)	33/72.8	53/116.8	63/138.9	73/160.9
Shelves	Quantity of shelves (Standard / max.)	2/4	2/5	2/8	2/9
	Max. Load per shelf (kg / lbs)	20 / 44.09			
Number of air change per hour at 100°C²⁾		123	22	18	11
Electrical requirements (230V, 50/60Hz, A)		6.4			
Cat. No.	Solid door (without window)	AAH14115K	AAH14125K	AAH14135K	AAH14145K
	Window door	AAH14215K	AAH14225K	AAH14235K	AAH14245K
Electrical requirements (120V, 60Hz, A)		12.3			
Cat. No.	Solid door (without window)	AAH14116U	AAH14126U	AAH14136U	AAH14146U
	Window door	AAH14216U	AAH14226U	AAH14236U	AAH14246U

Accessories Page 84 Wire Shelves, Perforated Shelves, Stacking set, BMS Port, LC GreenBox, Communication port, LC DataKeeper

Forced Convection Oven Multi type



Independent temperature control and optimal space efficiency per chamber

Structural Functional Features

- Airflow optimized for uniform heat transfer.
- Superior insulation and enclosed design minimize energy loss and reduce operating costs.
- Easy to open and close the door, soft handle minimizes shock to the sample.
- Stainless steel interior and shelves are excellent in terms of corrosion resistance and clean maintenance.
- Easy to clean as inner edges are curved.
- Includes vent hole base with cover for gas exhaust, cable connection, etc. on the side.

Use Convenience Features

- Calibration function minimizes temperature deviation.
- Optimized control with temperature auto-tuning.
- Save and use 3 frequently used temperatures.
- Wait On/Off timer. (up to 99 hours 59 minutes)
- Microprocessor PID method for precise temperature control.
- Automatic restart of operation when power is restored after sudden power failure.


Outstanding Safety

- Top-rated overheating protection system. (registration KR 10-0397583)
- Over temperature limit function.
- Safe product surface, even when operated at maximum temperature.
- Malfunction prevented by controller lock function.
- Self-diagnostics function to check abnormal condition of temperature sensor.



OF-02G-4C

with Wire Shelves 2ea (standard)



Specification

Model		OF-02G-4C
Chamber volume (L / cu ft)		240 / 8.5 (60 / 2.1 x 4 chambers)
Temperature	Range (°C / °F)	Amb. +10 to 200 / Amb. +18 to 392
	Fluctuation at 100°C (±°C / °F)	0.5 / 0.9
	Variation at 100°C (±°C / °F)	1 / 1.8
	Heating time to 100°C (min.)	15
	Recovery time at 100°C (min.)	2
Dimension	Interior (W x D x H, mm / inch)	400 x 360 x 420 / 15.7 x 14.2 x 16.5 (each chamber)
	Exterior (W x D x H, mm / inch)	1170 x 640 x 1360 / 46.1 x 25.2 x 53.5
	Net weight (Kg / lbs)	170 / 374.8
Shelves	Quantity of shelves (standard / max.)	2 / 4 per chamber
	Max. Load per shelf (Kg / lbs)	26 / 57.3
Electrical requirements (230V, 50/60Hz, A)		16.8
Cat. No.		AAH1A015K
Electrical requirements (120V, 60Hz, A)		33.4
Cat. No.		AAH1A016U

※ The lowest controllable temperature depends on variations in room and sample temperatures. Please contact technical support for consultation before purchasing.

※ Technical data (according to DIN 12880, before 2013)

Accessories Page 84 Wire Shelves, Perforated Shelves


Natural Convection Oven, Advanced type



ON4-10V



ON4-10VW



Airflow

Intuitive Touch Controller

- Intuitive interface with 5-inch display
- Events and alerts in a single view with real-time graph
- Three frequent temperatures automatically suggested
- 'Wait on / Wait off' timer function



Large intuitive display




Graphs for real-time process monitoring

Convenient Structural Functions

- Designed with a depth of less than 750 mm, making it easy and stable on a lab table
- Uniform and stable temperature distribution with multi-path airflow improves reliability of experiment data



Depth of less than 750 mm allows stable positioning on lab table



Uniform temperature control with multi-path airflow

Supervise Equipment with Real-time Monitoring

- Monitor the temperature on BMS (Building Management System) for quick responses to abnormalities. (option)
- Monitor via mobile app anytime, anywhere with LC Connected. (Requires LC GreenBox, option)



Timely responses via BMS-connected monitoring (BMS port: option)



Real-time monitoring via mobile devices (LC GreenBox: option)



Convenient Experiment and Data Checking

- Support RS-232/USB port and software making easy operation and data management from PC
- Test data are automatically stored and easily exported via USB; easy to edit data (in CSV format)
- Logs 36 recent events, including errors, opened door, etc
- Integrated with 21 CFR Part 11 compliant LC DataKeeper, suitable for GMP audits (option)



Data storage using USB ports




Store recent 36 events

Reliable Temperature Control Performance

- Wide temperature range of up to 250°C
- Temperature accuracy for wide range by 3-point calibration
- 「Ramp Control™」, function that controls the rate of temperature rise by setting ramping rate(°C/min) per experiment
- Uniform temperature distribution through airflow optimized for uniform heat transfer (within ±3.6 °C at 100°C)



Accurate temperature control with 3-point calibration



Ramp Control™ per experiment condition

Specification

Model	Solid door (without window)	ON4-03V	ON4-05V	ON4-10V	ON4-15V
	Window door	ON4-03VW	ON4-05VW	ON4-10VW	ON4-15VW
Chamber volume (L / cu ft)		26 / 0.92	58 / 2.05	110 / 3.88	147 / 5.19
Temperature	Range (°C / °F)	Amb.+15 ~ 250°C / Amb.+27 ~ 482°F			
	Fluctuation at 100°C (°C / °F)	±0.4 / ±0.72	±0.5 / ±0.9	±0.5 / ±0.9	±0.5 / ±0.9
	Variation at 100°C(°C / °F)	±2.8 / ±5.04	±2.6 / ±4.68	±3.6 / ±6.48	±3.5 / ±6.3
	Heating time to 100°C (min.) ¹⁾	25	39	22	21
	Recovery time at 100°C (min.)	17	19	10	9
Control panel	5 inch TFT Color LCD				
Communication interface	USB-B, RS-232				
Function	Recently event record	36 ea			
	Recommended temperature	Suggest three frequently used temperature			
	Graph	Available viewing stored and real-time graphs			
	Data storage	Save your experiment data (in CSV format)			
	Temp. ramping rate (°C/min. °F/min)	0 to 2.2 0 to 3.96	0 to 2 0 to 3.6	0 to 4.4 0 to 7.92	0 to 2.6 0 to 4.68
Dimensions	Interior (W x D x H, mm / inch)	273 x 270 x 350 / 10.75 x 10.63 x 13.78	340 x 390 x 422 / 13.4 x 15.35 x 16.61	441 x 420 x 597 / 17.36 x 16.54 x 23.5	500 x 446 x 659 / 19.69 x 17.56 x 25.95
	Exterior (W x D x H, mm / inch)	478 x 527 x 604 / 18.8 x 20.75 x 23.78	558 x 648 x 684 / 21.97 x 25.5 x 26.93	663 x 677 x 859 / 26.1 x 26.65 x 33.82	728 x 704 x 914 / 28.66 x 27.7 x 36.0
	Net weight (kg / lbs)	33 / 72.8	53 / 116.8	63 / 138.9	73 / 160.9
Shelves	Quantity of shelves (standard/max.)	2 / 4	2 / 5	2 / 8	2 / 9
	Max. Load per shelf (kg / lbs)	20 / 44.09			
Number of air change per hour at 100°C	16				
Electrical requirements	230V, 50/60Hz, A	3.2	3.5	5.7	7.4
Cat. No.	Solid door (without window)	AAH111105K	AAH111115K	AAH111125K	AAH111135K
	Window door	AAH112105K	AAH112115K	AAH112125K	AAH112135K
Electrical requirements	120V, 60Hz, A	6	6.7	10.9	14.2
Cat. No.	Solid door (without window)	AAH111106U	AAH111116U	AAH111126U	AAH111136U
	Window door	AAH112106U	AAH112116U	AAH112126U	AAH112136U

1) Heating time is defined as the time taken to reach 98% of the set temperature.
 ※ Technical data according to DIN 12880


Natural Convection Oven, Standard type



ON4-10S



ON4-10SW




Intuitive Touch Controller

- Intuitive interface with 3.5-inch display
- 'Wait on / Wait off' timer function



3.5-inch touch display




Wait on & Wait off timer

Reliable Temperature Control Performance

- Wide temperature range of up to 250°C
- Temperature accuracy for wide range by 3-point calibration
- The upper and lower temperature limits can be set, enabling alerts for deviations to trigger quick actions
- Uniform temperature distribution through airflow optimized for uniform heat transfer (within $\pm 3.6^{\circ}\text{C}$ at 100°C)



Accurate temperature control with 3-point calibration



Settings for temperature deviation

Convenient Experiment and Data Checking

- Support RS-232/USB port and software making easy operation and data management from PC (option)
- Monitor via mobile app anytime, anywhere with LC Connected. (Requires LC GreenBox, option)
- Monitor the temperature on BMS (Building Management System) for quick responses to abnormalities. (option)
- Integrated with 21 CFR Part 11 compliant LC DataKeeper, suitable for GMP audits (option)



Manipulation via PC (option)



Real-time monitoring via mobile devices (LC GreenBox: option)



Ensuring User Safety

- Effectively removes fumes from the inside of the chamber via the top vent hole (ON4-03)
- Safe surface even when operating at max temperatures (less than 70°C at 250°C)
- Implemented with the top-rated overheating prevention system (Patented in Korea)
- Over Temperature Limiter
- Complaint to IEC Class II protection standards and temperature protection in line with DIN 12880 Class 3.1



Removes fume from the inside of the chamber



Safe surface even at maximum temperature

Convenient Structure and Functions

- Logs the last 12 events, including run times, stops, and errors
- Double stackable with stacking set (option) maximizes space utilization (-05 / -15 models only)
- 「Hands-free™」 door allows both hands free while holding samples



Logs the recent 12 events



Increased space utilization by double stacking (stacking set : option)
(※ When double stacking, the minimum control temperature of the upper chamber may be limited.)

Specification

Model	Solid door (without window)	ON4-03S	ON4-05S	ON4-10S	ON4-15S
	Window door	ON4-03SW	ON4-05SW	ON4-10SW	ON4-15SW
Chamber volume (L / cu ft)		26 / 0.92	58 / 2.05	110 / 3.88	147 / 5.19
Temperature	Range (°C / °F)	Amb.+15 ~ 250°C / Amb.+27 ~ 482°F			
	Fluctuation at 100°C (°C / °F)	±0.4 / ±0.72	±0.5 / ±0.9	±0.5 / ±0.9	±0.5 / ±0.9
	Variation at 100°C (°C / °F)	±2.8 / ±5.04	±2.6 / ±4.68	±3.6 / ±6.48	±3.5 / ±6.3
	Heating time to 100°C (min.) ¹⁾	25	39	22	21
	Recovery time at 100°C (min.)	17	19	10	9
Control panel		3.5 inch TFT Color LCD			
Communication interface		USB-B, RS-232 (option)			
Function	Recently event record	12 ea			
Dimensions	Interior (W x D x H, mm / inch)	273 x 270 x 350 / 10.75 x 10.63 x 13.78	340 x 390 x 422 / 13.4 x 15.35 x 16.61	441 x 420 x 597 / 17.36 x 16.54 x 23.5	500 x 446 x 659 / 19.69 x 17.56 x 25.95
	Exterior (W x D x H, mm / inch)	478 x 527 x 604 / 18.8 x 20.75 x 23.78	558 x 648 x 684 / 21.97 x 25.5 x 26.93	663 x 677 x 859 / 26.1 x 26.65 x 33.82	728 x 704 x 914 / 28.66 x 27.7 x 36.0
	Net weight (kg / lbs)	33 / 72.8	53 / 116.8	63 / 138.9	73 / 160.9
Shelves	Quantity of shelves (standard/max.)	2 / 4	2 / 5	2 / 8	2 / 9
	Max. Load per shelf (kg / lbs)	20 / 44.09			
Number of air change per hour at 100°C		16			
Electrical requirements	230V, 50/60Hz, A	3.2	3.5	5.7	7.4
Cat. No.	Solid door (without window)	AAH15045K	AAH15065K	AAH15075K	AAH15085K
	Window door	AAH15145K	AAH15165K	AAH15175K	AAH15185K
Electrical requirements	120V, 60Hz, A	6	6.7	10.9	14.2
Cat. No.	Solid door (without window)	AAH15046U	AAH15066U	AAH15076U	AAH15086U
	Window door	AAH15146U	AAH15166U	AAH15176U	AAH15186U

1) Heating time is defined as the time taken to reach 98% of the set temperature.
※ Technical data according to DIN 12880

Vacuum Oven

Excellent temperature uniformity, allowing testing at any position

Precise temperature control

- Wide temperature range of up to 250 °C
- Temperature variation up to $\pm 1.5^{\circ}\text{C}$ at 100°C according to DIN 12880:2013
- Actual experiment area is expanded by excellent temperature uniformity
- 3-point calibration improves temperature accuracy in the entire area

Structural Functions

- Excellent corrosion resistance by stainless steel and aluminum shelves
- Safety circuits protect from overcurrent and over temperature
- Vacuum line and release line are separated and convenient to use
- Viton gasket when drying highly reactive chemicals (option)
- In vacuum, tempered glass adheres tightly to the gasket, maintaining an airtight seal


Outstanding User Safety

- Door buffer system that maintains stable vacuum condition. In case of explosion, internal pressure is discharged via the upper/lower exhaust ports
- BMS port is provided as standard to monitor within a building management system (option)
- Polycarbonate cover is provided as standard in case of tempered glass breakage
- Top-rated overheating protection system (Patented in Korea)
- The surface remains safe (< 60°C) at high temperature




OV4-30

2 shelves included (standard)



Buffer system door structure



Variation less than $\pm 1.5^{\circ}\text{C}$ at 100°C

System Settings > Calibration		
Actual Temperature		17.0 °C
	Uncalibrated	Calibrated
Temp 1	60.2 °C	60.0 °C
Temp 2	100.7 °C	100.0 °C
Temp 3	151.0 °C	150.0 °C

3-point calibration for accuracy



Safe from burns with surface below 60°C even at max temperature



User Convenience

- Wide LCD Color touch display controller
- Monitor via mobile app anytime, anywhere with LC Connected. (Requires LC GreenBox, option)
- Test data are automatically stored and easily exported with USB in CSV format (except for OV-S)
- Events/alerts in a single view with real-time graph (except for OV-S)
- Able to limit speed of temperature rise with [Ramp Control™] (except for OV-S)
- Wait On/Off timer (up to 99hrs 59mins)
- Auto recommendation of 3 temperatures based on user history
- Integrated with 21 CFR Part 11 compliant LC DataKeeper, suitable for GMP audits (option)




Intuitive touch controller
(OV4 : 5 inch / OV4-S : 3.5 inch)



Real-time monitoring via mobile devices (LC GreenBox: option)



Timely responses via BMS-connected Monitoring (BMS port: option)



Application of replacing internal air of chamber with inert gas

Event	Date	Time	Details
Event	2020-03-15	07:00:00	Basic control started
Error	2020-03-15	07:01:00	Warning: Over Temp.
Event	2020-03-15	07:02:00	Basic control started
Error	2020-03-15	07:03:00	Warning: Over Temp.

Store recent 36 events (OV-S: 12 events)

No.	Date_Name	Check
1	20200423_123400.csv	<input checked="" type="checkbox"/>
2	20200423_123410.csv	<input checked="" type="checkbox"/>
3	20200423_123420.csv	<input checked="" type="checkbox"/>
4	20200423_123430.csv	<input checked="" type="checkbox"/>
5	20200423_123440.csv	<input checked="" type="checkbox"/>
6	20200423_123450.csv	<input checked="" type="checkbox"/>

Automatic storage of experimental data in the internal memory




Data storage using USB ports



Graphs for real-time process monitoring

Vacuum Oven Use Configuration

- It is recommended to use the cold trap bath on the vacuum line to protect the weak vacuum pump against water/solvent or to supply dry air to the sample.



Specification

Model		OV4-30	OV4-30S*	OV4-65	OV4-65S*
Control panel		5 inch	3.5 inch	5 inch	3.5 inch
Chamber volume (L / cu ft)		28 / 1		65 / 2.3	
Temperature	Range (°C / °F)	Amb. +15 to 250 / Amb. +27 to 482		Amb. +15 to 250 / Amb. +27 to 482	
	Fluctuation at 100°C (±°C / °F)	0.2 / 0.36		0.2 / 0.36	
	Variation at 100°C (±°C / °F)	1.5 / 2.7		1.5 / 2.7	
	Heating time to 100°C (min.)	70		90	
Vacuum gauge range, analog (Mpa)		0 to 0.1		0 to 0.1	
Dimension	Interior (W x D x H, mm / inch)	302 x 305 x 302 / 11.9 x 12 x 11.9		402 x 405 x 402 / 15.8 x 15.9 x 15.8	
	Exterior (W x D x H, mm / inch)	537 x 522 x 655 / 21.1 x 20.6 x 25.8		636 x 622 x 755 / 25.0 x 24.5 x 29.7	
	Net weight (Kg / lbs)	68 / 150		107 / 235.9	
Shelves	Quantity of shelves (standard / max.)	2 / 3		3 / 4	
	Max. Load per shelf (Kg / lbs)	20 / 44.1		20 / 44.1	
Material		Body: Stainless Steel, Shelf: Aluminum		Body: Stainless Steel, Shelf: Aluminum	
Electrical requirements (230V, 50/60Hz, A)		5.7		6.7	
Cat. No.		AAH13315K	AAH13335K	AAH13325K	AAH13345K
Electrical requirements (120V, 60Hz, A)		10.2		11.7	
Cat. No.		AAH13316U	AAH13336U	AAH13326U	AAH13346U

※ The lowest controllable temperature depends on variations in room and sample temperatures. Please contact technical support for consultation before purchasing.

※ Technical data (according to DIN 12880, after 2013)

* The OV4-30S and OV4-65S exclude some convenient features. (real-time graphs, data storage, etc.)



Jeio Tech



Jeio Tech Releases
New Vacuum Oven OV4 series



Key Measurement of the
Chamber-type Performance



How to choose the right shelf
for the OV4 series

Accessories Page 84 BMS Port, Silicone Gasket, Viton Gasket, Shelf, Shelf Spacer, LC GreenBox, LC DataKeeper

Vacuum Oven System Station

Dedicated Table for Neat Organization of Vacuum Drying

Compact System Structure

- Compact system placing vacuum oven on top, cold trap bath and vacuum pump below.
- Neat arrangement of power cords and vacuum hoses using ring bolts and cable clamps.

Convenient Use and Maintenance

- Two outlets for cold trap bath and vacuum pump and convenient on/off control of using front-panel switches.
- Easy removal for cleaning and replacing glass traps of cold trap bath.
- Easy installation and movement with casters.



Vacuum Oven System Station (OVT)



Neat hose arrangement



Convenient caster for moving/fixing



Easy replacement and cleaning




Convenient power On/Off

Composition Example

Table	OVT
Vacuum Oven	OV4-30
Cold Trap Bath	CTB-20 with glass trap
Vacuum Pump	G-50DA, 1.3 ⁻¹ Pa

※ The above configuration is an example using Jeio Tech products. Please purchase the model that suits your needs.



Specification

Cat. No.	Description
AAA13551	Vacuum Oven Table (OVT), 910 x 710 x 890 mm

blog


Jeio Tech



Read more about OVT

Muffle Furnace

Stable high-temperature experiments at 1100 °C Convenient and intuitive touch controller



Up to 1100°C
(4-side molding type heater)

MF2-12G



Enhanced user convenience of touch controller

- Intuitive use with a large 5-inch TFT touch controller.
- Auto 3 temperature recommendation based on operating history without any recording or storage.
- Timer(wait on, wait off) and Ramp(temperature increase rate) function
- Monitor via mobile app anytime, anywhere with LC Connected. (Requires LC GreenBox, option)
- Recording the 36 most recent events, such as start/stop of operation and error occurrence.





Large 5-inch touch display



Real-time monitoring and control via mobile

Easy program settings intuitively without a manual (MF2-G : program type only)

- Available to set up to 10 programs and 10 steps.
- Available to set up to 99 repetitions of program from 0.
- Function of setting the start temperature to begin an experiment at an accurate temperature the user wants.
- Possible to skip to the next step or hold the current step during program operation.
- Available to set a deviation of temperature or time to keep going a program operation even if the temperature is not reached to the target accurately.
- A useful scheduler function available to select and schedule the desired day, time, and program for repetitive experimental settings.




Setting a program for an ashing process
(Waste Processing Standard ES 06305.2)



Highly reliable temperature control performance

- Possible to set the rate of temperature change (°C/min) (Ramp function) to control it.
- Minimizes heat loss through the use of durable ceramic insulation materials.
- Excellent interior sealability due to shock absorbers on the doors and a sturdy double-door structure.
- Stable temperature control through digital PID control.
- High-precision temperature control over a wide range of temperatures through a 3-point calibration process.



Temperature ramp control (Fastest, 5°C/min, 10°C/min) at the set temperature of 800 °C.



Precision temperature control with a 3-point calibration process

Various safety functions ensure user safety

- A safety door lock function that opens only below the door lock temperature set by user.
- A built-in electronic overheat protection system that prevents the equipment from being damaged due to overheating.
- Front surface of the product is safe even when operating at the maximum temperature.
- The equipment stops working when the door is opened for safe use




Safety door lock function




Electronic overheat protection system

Interior Dimensions (Unit: mm / inch)

Model	W	H	H'
MF2-12G(F)	130 / 5.1	90 / 3.5	-
MF2-22G(F)	150 / 5.9	100 / 3.9	-
MF2-32G(F)	250 / 9.8	150 / 5.9	175 / 6.8



MF2-12/22G(F)
MF2-12/22/32G(F)



MF2-32G(F)

Specification

Model	General Type	MF2-12GF	MF2-22GF	MF2-32GF
	Program Type	MF2-12G	MF2-22G	MF2-32G
Chamber volume (L / cu ft)		3 / 0.11	4.5 / 0.16	12 / 0.42
Temperature				
Range (°C / °F)		Amb.+15 to 1100 / Amb.+27 to 2012		
Continuous operating temperature (°C/°F)*		Amb.+15 to 960 / Amb.+27 to 1760		
Fluctuation (±°C / °F) (at max. temp.)		1.0 / 1.8		
Dimensions				
Interior (W x D x H, mm / inch)		130 x 250 x 90 / 5.1 x 9.8 x 3.5	150 x 300 x 100 / 5.9 x 11.8 x 3.9	250 x 300 x 150 (175) ¹⁾ / 9.8 x 11.8 x 5.9 (6.8)
Exterior (W x D x H, mm / inch)		395 x 500 x 590 / 15.6 x 19.7 x 23.2	415 x 550 x 600 / 16.3 x 21.7 x 23.6	525 x 550 x 690 / 20.7 x 21.7 x 27.2
Weight (kg / lbs)		42 / 92.6	50 / 110.2	70 / 154.3
Electrical requirements (230V, 60Hz, A)		9	13.6	20.4
Insulating material		Alumina ceramic		
Heater		4-side molding wire		
Cat. No.	General Type	AAHP2112K	AAHP2122K	AAHP2132K
	Program Type	AAHP2111K	AAHP2121K	AAHP2131K

* We recommend operating the equipment within the range of continuous operating temperatures for ensuring the stability and integrity of the equipment during its lifetime.

1) The shape of the top inside is an arch, and the height in parentheses is a measured value to the highest point.

Accessories for Oven & Furnace

Shelves



Wire Shelf (Oven)

- Excellent corrosion resistance as made of stainless steel.
- Excellent air permeability and uniform temperature uniformity.



Perforated Shelf (Oven)

- Porous structure made of stainless steel.
- Solid structure makes it suitable for high load samples.



Ceramic Shelf (Furnace)

- Plate for placing a crucible.
- Made of Mullite for thermal resistance.

Shelves for Forced/Natural Convection Oven

Model	Wire shelves		Perforated shelves	
	Cat. No.	Quantity (standard/max.)	Cat. No.	Quantity (max.)
OF4-03P/V/S, OF4-03PW/VW/SW	00FDA0009739	2 / 4	AAA125341	4
OF4-05P/V/S, OF4-05PW/VW/SW	00FDA0009738	2 / 5	AAA125342	5
OF4-10P/V/S, OF4-10PW/VW/SW	00FDA0009737	2 / 8	AAA125343	8
OF4-15P/V/S, OF4-15PW/VW/SW	00FDA0009736	2 / 9	AAA125344	9
OF3-30/30W/30P/30PW	00RTD0001196	2 / 5	AAA18501	5
OF3-30H/30HW/30HP/30HPW	00RTD0001196	2 / 5	AAA18501	5
OF3-45/45W/45P/45PW	00RTD0001196	3 / 11	AAA18501	8
OF3-45H/45HW/45HP/45HPW	00RTD0001196	3 / 11	AAA18501	8
OF3-75/75W/75P/75PW	00RTD0001197	3 / 11	AAA18502	11
OF3-75H/75HW/75HP/75HPW	00RTD0001197	3 / 11	AAA18502	11
OF-02G-4C	00EDA0008136	8 / 16	AAA12531	16
ON4-03V/S, ON4-03VW/SW	00FDA0009739	2 / 4	AAA125341	4
ON4-05V/S, ON4-05VW/SW	00FDA0009949	2 / 5	AAA125362	5
ON4-10V/S, ON4-10VW/SW	00FDA0009950	2 / 8	AAA125363	8
ON4-15V/S, ON4-15VW/SW	00FDA0009951	2 / 9	AAA125364	9

※ Brackets for each shelf are included.

Shelves for Vacuum Oven

Description	Position	Material	Cat. No.	
			OV4-30/30S	OV4-65/65S
Bottom shelf, Basic ²⁾		Aluminum	AAA13519	AAA13520
Extra shelf, Optional ^{1), 4)}		Stainless Steel	AAA13529	AAA13530
		Aluminum	AAA13511	AAA13513
Lower shelf, Basic ²⁾		Stainless Steel	AAA13521	AAA13523
		Aluminum	AAA13517	AAA13518
Upper shelf, Basic ^{2), 3)}		Stainless Steel	AAA13527	AAA13528
		Aluminum	AAA13512	AAA13514
Multiple shelf, Optional ^{1), 4)}		Stainless Steel	AAA13522	AAA13524
		Aluminum	AAA13531	AAA13533
Shelf spacer, Basic ²⁾		Aluminum	AAA13532	AAA13534
		Stainless Steel	AAA13515	AAA13516
			Stainless Steel	AAA13525

1) You can add an extra shelf or multiple shelf that can be placed at the bottom to experiment with more samples same time, and the value of Variation and Fluctuation may increase slightly.

2) The shelf and shelf spacer are included as bundled accessories, and you can purchase when need to replace them.

3) For OV4-65 models, 2 upper shelves are basic provided.

4) Maximum number of shelf for OV4-30 model : 4 extra shelves / 5 multiple shelves and for OV4-65 model : 4 extra shelves / 8 multiple shelves.

Shelves for Furnace

Model	Dimension (W x D x H, mm)	Ceramic Shelf
MF2-12 series	120 x 230 x 12t	00FDA0009706
MF2-22 series	140 x 280 x 12t	00FDA0009707
MF2-32 series	260 x 260 x 12t	00FDA0009708

Recorders & Brackets



Recorder (Dot)

- 6-channel graph-type temperature recorder.
- Adjustable recording range, speed, etc.



Recorder (Thermal Line)

- Graph/Digital mode.
- Thermal-type temperature recorder.



Stacking Set (-05 / -15 models only)

- Double stackable
- Firm stacking allows efficient utilization of space (※ When double stacking, the minimum control temperature of the upper chamber may be limited.)

Model	Description	Cat. No.
OF3-30/45/75/30H/45H/75H Series	Recorder (Dot)	AAA8T500
	Recorder (Thermal line)	AAAE1503
MF2 series	Recorder (Thermal line)	AAAE1503
OF4/ON4-05	Stacking set	AAA125547
OF4/ON4-15	Stacking set	AAA125549

Door Gaskets for Vacuum Oven

Description	Cat. No.		
	OV4-30/30S	OV4-65/65S	
Silicone Gasket	Included as standard. Silicon gasket for sealing the door.	00FDA0009912	00FDA0009913
Viton Gasket	High chemical resistance	00FDA0009914	00FDA0009915

Communication



BMS Port

- Real-time monitoring at a central monitoring facility.
- Easier to ensure safety, especially in labs where multiple units are operated.



Communication Port

- Control PC connection and data by dedicated software.
- Included as default in OF4-P/V and optional for OF4-S.



LC GreenBox

(Mobile Monitoring System)

- Remote monitoring and control via mobile application.
- One LC GreenBox can connect up to 4 devices.
- Connect using RS-232 port on the communication port.



LC DataKeeper

- High-integrity software compliant with FDA 21 CFR Part 11
- Integrated with Jeio Tech equipment to undergo GMP audits

Ports	OF4-P/V, ON4-P/V	OF4-S, ON4-S	OV4
BMS Port	AAA125534		AAA125558
Communication Port¹⁾	Standard	AAA125533	-
LC GreenBox	AAHQ1011K		
LC DataKeeper Basic²⁾	SSA11001		
LC DataKeeper Pro²⁾	SSA11002		

1) Communication port is provided as a standard in advanced (-V) model and is an optional add-on for the basic (-S) model for connecting LC DataKeeper, LC GreenBox, BMS, and PC.

2) Trial version is available for 30 days, offering the same features as the Basic version, and can be requested through the website.

Incubator



Satisfaction with the requirements of European regulations and guidelines.



Eco-friendly product that has passed environmental impact assessment tests.



More precise temperature control through temperature calibration.



Possible to set end time or start time of device operation.



Control of temperature / humidity / rpm / time, etc. through dedicated program.



2 year warranty free A/S.



Real-time equipment monitoring and control system using mobile app.



Ensuring the reliability of electronic records and signatures to support GMP audit.



Warning alarm in case of deviation of temperature control.



Indicator to ascertain water level.



PC communication via RS-232 / RS-485 / USB port.



Registered safety patent based on Jeio Tech's proprietary technology.



More precise temperature control through 3-point temperature calibration.



Over temperature protection function.



Alarm if the water is low level.



It memorizes the state of power failure and operates automatically when power is restored.

General Application

Organic bacteria, organic drugs, tissue cultures in the areas of antibody testing, clinical studies, hematology, microbiology, and others.

	Description	Temp. Range (°C/°F)	Max. Temp. Fluctuation (±°C/°F)	Volume (L / cu ft)	Model	Page
Incubator	Air-jacket, Advanced (Natural Convection) <ul style="list-style-type: none"> • Ultimate temperature control • Advanced air-jacket technology 	Amb.+5 ~ 80°C / Amb.+9 ~ 176°F	±0.15 / ±0.27 at 37°C	25, 60, 114, 150 (L) / 0.88, 2.12, 4.03, 5.31 (cu ft)	IB4-AV	90
	Air-jacket, Standard (Natural Convection) <ul style="list-style-type: none"> • Simple yet effective • Film Heater incubator 	Amb.+5 ~ 80°C / Amb.+9 ~ 176°F	±0.15 / ±0.27 at 37°C	25, 60, 114, 150 (L) / 0.88, 2.12, 4.03, 5.31 (cu ft)	IB4-AS	92
	Advanced (Forced Convection) <ul style="list-style-type: none"> • Consistent temperature with forced air circulation • Prevent media drying with adjustable airflow 	Amb.+5 ~ 80°C / Amb.+9 ~ 176°F ¹⁾	±0.15 / ±0.27 at 37°C	25, 60, 111, 146 (L) / 0.88, 2.12, 3.92, 5.16 (cu ft)	IB4-V	94
	Standard (Forced Convection) <ul style="list-style-type: none"> • User-friendly interface with built-in forced circulation • Uniform temperature distribution 	Amb.+5 ~ 80°C / Amb.+9 ~ 176°F ¹⁾	±0.15 / ±0.27 at 37°C	25, 60, 111, 146 (L) / 0.88, 2.12, 3.92, 5.16 (cu ft)	IB4-S	96
	Mini (Forced Convection) <ul style="list-style-type: none"> • Space-saving design • Easy access with bottom-hinged door 	Amb. +5 to 65 / Amb. +9 to 149	0.5 at 37°C / 0.9 at 98.6°F	10 / 0.4	IM	98
	Multi (Forced Convection) <ul style="list-style-type: none"> • Independent temperature control for each chamber • 2-in-1, 4-in-1 models 	Amb. +5 to 70 / Amb. +9 to 158	0.1 at 37°C / 0.18 at 98.6°F	120 / 4.2 (60 L / 2.1 cu ft x 2 chamber) 240 / 8.5 (60 L / 2.1 cu ft x 4 chamber)	IB-02G-2C IB-02G-4C	99
Low Temperature Incubator	Air-Jacket (Natural Convection) <ul style="list-style-type: none"> • Minimization of culture medium dry-out • Air-jacket type incubator 	4 to 70 / 39.2 to 158	0.2 at 37°C / 0.36 at 98.6°F	156, 254 / 5.5, 9.0	IL3-A	100
	General (Forced Convection) <ul style="list-style-type: none"> • Optimized for BOD experiments • Forced convection method 	0 to 80 / 32 to 176	0.1 at 37°C / 0.18 at 98.6°F	150, 242 / 5.3, 8.5	IL3	100
	Personal (Forced Convection) <ul style="list-style-type: none"> • Application of peltier cooling device system • Structure minimizes noise/vibration 	5 to 40 / +41 to 104 (at 20°C / 68°F) 10 to 40 / +50 to 104 (at 25°C / 77°F) 15 to 40 / +59 to 104 (at 30°C / 86°F)	0.1 at 25°C / 0.18 at 77°F	14, 48 / 0.5, 1.7	ILP	102
	Multi (Forced Convection) <ul style="list-style-type: none"> • Independent temperature control for each chamber • 2-in-1, 4-in-1 models 	0 to 60 / 32 to 140	0.1 at 25°C / 0.18 at 77°F	300 / 10.6 (150 L / 5.3 cu ft x 2 chamber) 600 / 21.2 (150 L / 5.3 cu ft x 4 chamber)	IL-11-2C IL-11-4C	103
Plant Growth Chamber	Growth Chamber <ul style="list-style-type: none"> • 10-step program control • Temperature/Humidity auto-tuning 	5 to 50 / 41 to 122 (Lamp off) 10 to 50 / 50 to 122 (Lamp on) 20 to 50 / 68 to 122 (with humidity)	0.1 at 25°C / 0.18 at 77°F	300, 1000 / 10.6, 35.3	GC	104

※ The contents of the above and the contents of this catalog may differ depending on the specific model and conditions of use. For the information about the features and specifications that applying to each model, please check the information on the corresponding page of each model.

1) IB4-V and IB4-S (forced convection type) can control down to amb +5°C (amb +9°F) when the water tray is filled with water. (Water tray: option)



An incubator series that provides optimal culture environment Application of dedicated controller for precise control

> Optimized selection model

Optimal model provided based on heat transfer method and airflow, capacity, performance and function, controller type, etc.

> Verified according to international standards

Fluctuation, variation, heating time, and recovery time are tested according to DIN 12880.

> Glass door for contamination prevention

Internal glass door allows for convenient observation of sample without temperature change and prevents inflow of external contaminants.

> Clean use and maintenance

Internal stainless steel structure allows for clean maintenance and excellent corrosion resistance.

> More stable temperature control

Optimized Microprocessor PID control for more accurate temperature control.

> 3-point temperature calibration function (IB4, IL3)

Temperature calibration of up to 3 points provides more reliable temperature control.

> Optimized temperature control

Auto-tuning function enables optimal control even when ambient temperature or conditions are changed frequently.

> Computer interface control

Convenient remote control and data management through built-in dedicated software. (IB4, IB4-A)



Sealed glass door for contamination prevention.



Large intuitive touch display controller



Easy shelf height adjustment.



Functions convenient for actual use and Jeio Tech's patented safety system

- > **Provide convenient timer function**
Wait On/Wait Off is run after the set time elapses. (Max. 999 hours and 59 minutes)
- > **Controller lock function**
Lock function for preventing change of controller settings prevents erroneous operation due to arbitrary changing.
- > **Convenient temperature value save and operation**
Conveniently save and used frequently used temperatures value. (up to 3 presets allowed)
- > **Cable port provided as standard**
38mm diameter cable port and cover for convenient utilization of culture conditions. (IB4, IB-multi models)
- > **Top-grade overheating prevention system (Jeio Tech's patent)**
Even if a failure occurs in the device control unit, the safety device is structured so that it is preferentially operated.
- > **Over heating limit function in case of accidental over heating.**
Safety structure with shutdown function using over temperature limiter in case of overheating.
- > **Support for remote equipment monitoring services (option)**
Real-time monitoring of equipment operating status using smart phones.
- > **Prepared for unexpected power failure**
Auto-run function capable of selecting reoperation when power is restored after power failure.

※ Some of the above contents are limited to specific models.




Various models for selecting the desired structure.



Cable port provided as standard on the side.

Incubator (Air-jacket) Advanced type



IB4-15AV

Intuitive Large Touch Controller

- Intuitive interface with 5-inch display
- Events and alerts in a single view with real-time graph
- Three frequent temperatures automatically suggested
- 'Wait on / Wait off' timer function



Large intuitive display



Graphs for real-time process monitoring

Unique Structural Functions

- Using Water tray delays media drying, making incubation time longer (option)
- Double stackable with stacking set (option) maximizes space utilization



Water tray that minimizes media drying (option)



Increased space utilization by double stacking (stacking set : option)
*Double stacking example: (top) forced convection oven (bottom) forced convection incubator

Supervise Equipment with Real-time Monitoring

- Monitor the temperature on BMS (Building Management System) for quick responses to abnormalities. (option)
- Monitor via mobile app anytime, anywhere with LC Connected. (Requires LC GreenBox, option)



Timely responses via BMS-connected monitoring (BMS port: option)



Real-time monitoring via mobile devices (LC GreenBox: option)



Convenient Experiment and Data Checking

- Support RS-232/USB port and software making easy operation and data management from PC
- Test data are automatically stored and easily exported via USB; easy to edit data (in CSV format)
- Logs 36 recent events, including errors, opened door, etc
- Integrated with 21 CFR Part 11 compliant LC DataKeeper, suitable for GMP audits (option)



Data storage using USB ports




Store recent 36 events

Reliable Temperature Control Performance

- Temperature accuracy for wide range by 3-point calibration
- Ramp Control™, function that controls the rate of temperature rise by setting ramping rate(°C/min) per experiment
- Uniform temperature distribution through airflow optimized for uniform heat transfer (within ±0.55°C at 37°C)



Accurate temperature control with 3-point calibration



Ramp Control™ per experiment condition


Specification

Model		IB4-03AV	IB4-05AV	IB4-10AV	IB4-15AV
Chamber volume (L / cu ft)		25 / 0.88	60 / 2.19	114 / 4.03	150 / 5.3
Temperature	Range (°C / °F)	Amb.+5 ~ 80°C / Amb.+9 ~ 176°F			
	Fluctuation at 37°C (°C / °F)	±0.1 / ±0.18	±0.1 / ±0.18	±0.15 / ±0.27	±0.1 / ±0.18
	Variation at 37°C (°C / °F)	±0.31 / ±0.56	±0.48 / ±0.86	±0.46 / ±0.83	±0.55 / ±0.99
	Heating time to 37°C (min.)	44	40	46	49
	Recovery time at 37°C (min.)	5	9	9	9
Control panel		5 inch TFT Color LCD			
Communication interface		USB-B, RS-232			
Function	Recently event record	36 ea			
	Recommended temperature	Suggest three frequently used temperature			
	Graph	Available viewing stored and real-time graphs			
	Data storage	Save your experiment data (in CSV format)			
	Temp. ramping rate (°C/min. °F/min)	0 ~ 2 0 ~ 3.6			
Dimensions	Interior (W x D x H, mm / inch)	300 x 210 x 400 / 11.81 x 8.27 x 15.75	380 x 330 x 480 / 14.96 x 12.99 x 18.9	485 x 350 x 655 / 19.1 x 13.78 x 25.79	550 x 376 x 710 / 21.62 x 14.8 x 27.95
	Exterior (W x D x H, mm / inch)	478 x 527 x 604 / 18.8 x 20.75 x 23.78	558 x 648 x 684 / 21.97 x 25.5 x 26.93	663 x 677 x 859 / 26.1 x 26.65 x 33.82	728 x 704 x 914 / 28.66 x 27.7 x 36.0
	Net weight (kg / lbs)	31 / 68.34	40 / 88.18	55 / 121.25	60 / 132.28
Shelves	Quantity of shelves (standard/max.)	2 / 4	2 / 5	2 / 8	2 / 9
	Max. Load per shelf (kg / lbs)	20 / 44.09			
Electrical requirements (230V, 50/60Hz, A)		2.5	3.5	3.5	3.5
Cat. No.		AAH21755K	AAH21765K	AAH21775K	AAH21785K
Electrical requirements (120V, 60Hz, A)		4.8	6.7	6.7	6.7
Cat. No.		AAH21753U	AAH21763U	AAH21773U	AAH21783U

※ Technical data according to DIN 12880

Accessories Page 107 Wire Shelf, Perforated Shelf, Stacking set, BMS Port, LC GreenBox, LC DataKeeper

Incubator (Air-jacket) Standard type




IB4-15AS

Intuitive Touch Controller

- Intuitive interface with 3.5-inch display
- 'Wait on / Wait off' timer function




3.5-inch touch display



Wait on & Wait off timer

Reliable Temperature Control Performance

- Temperature accuracy for wide range by 3-point calibration
- The upper and lower temperature limits can be set, enabling alerts for deviations to trigger quick actions
- Uniform temperature distribution through airflow optimized for uniform heat transfer (within $\pm 0.55^{\circ}\text{C}$ at 37°C)



Accurate temperature control with 3-point calibration



Settings for temperature deviation

Convenient Experiment and Data Checking

- Support RS-232/USB port and software making easy operation and data management from PC (option)
- Monitor via mobile app anytime, anywhere with LC Connected. (Requires LC GreenBox, option)
- Monitor the temperature on BMS (Building Management System) for quick responses to abnormalities. (option)
- Integrated with 21 CFR Part 11 compliant LC DataKeeper, suitable for GMP audits (option)



Manipulation via PC (option)



Real-time monitoring via mobile devices (LC GreenBox: option)



Enhanced Safety and Convenience

- Implemented with the highest-level overheating prevention system (Patented in Korea)
- Over Temperature Limiter
- Complaint to IEC Class II protection standards and temperature protection in line with DIN 12880 Class 3.1
- Logs the last 12 events, including run times, stops, and errors
- Double stackable with stacking set (option) maximizes space utilization



Logs the recent 12 events



Increased space utilization by double stacking (stacking set : option)
*Double stacking example: (top) forced convection oven (bottom) forced convection incubator

Multiple Features to Minimize Incubation Failure

- Using Water tray delays media drying, making incubation time longer (option)
- Automatically resume the control of previous operations when restored within 30 seconds from blackouts
- Controller can be locked with the button on the top, preventing any actions from unintended pressing
- Abnormalities are notified with audiovisual alerts to trigger quick actions



Water tray that minimizes media drying (option)




Audiovisual alerts to notify abnormalities

Specification

Model		IB4-03AS	IB4-05AS	IB4-10AS	IB4-15AS
Chamber volume (L / cu ft)		25 / 0.88	60 / 2.19	114 / 4.03	150 / 5.3
Temperature	Range (°C / °F)	Amb.+5 ~ 80°C / Amb.+9 ~ 176°F			
	Fluctuation at 37°C (°C / °F)	±0.1 / ±0.18	±0.1 / ±0.18	±0.15 / ±0.27	±0.1 / ±0.18
	Variation at 37°C (°C / °F)	±0.31 / ±0.56	±0.48 / ±0.86	±0.46 / ±0.83	±0.55 / ±0.99
	Heating time to 37°C (min.)	44	40	46	49
	Recovery time at 37°C (min.)	5	9	9	9
Control panel		3.5 inch TFT Color LCD			
Communication interface		USB-B, RS-232 (option)			
Function	Recently event record	12 ea			
Dimensions	Interior (W x D x H, mm / inch)	300 x 210 x 400 / 11.81 x 8.27 x 15.75	380 x 330 x 480 / 14.96 x 12.99 x 18.9	485 x 350 x 655 / 19.1 x 13.78 x 25.79	550 x 376 x 710 / 21.62 x 14.8 x 27.95
	Exterior (W x D x H, mm / inch)	478 x 527 x 604 / 18.8 x 20.75 x 23.78	558 x 648 x 684 / 21.97 x 25.5 x 26.93	663 x 677 x 859 / 26.1 x 26.65 x 33.82	728 x 704 x 914 / 28.66 x 27.7 x 36.0
	Net weight (kg / lbs)	31 / 68.34	40 / 88.18	55 / 121.25	60 / 132.28
Shelves	Quantity of shelves (standard/max.)	2 / 4	2 / 5	2 / 8	2 / 9
	Max. Load per shelf (kg / lbs)	20 / 44.09			
Electrical requirements (230V, 50/60Hz, A)		2.5	3.5	3.5	3.5
Cat. No.		AAH21555K	AAH21565K	AAH21575K	AAH21585K
Electrical requirements (120V, 60Hz, A)		4.8	6.7	6.7	6.7
Cat. No.		AAH21553U	AAH21563U	AAH21573U	AAH21583U

※ Technical data according to DIN 12880

Incubator (Forced Convection), Advanced type



Airflow

IB4-15V

Intuitive Large Touch Controller

- Intuitive interface with 5-inch display
- Events and alerts in a single view with real-time graph
- Three frequent temperatures automatically suggested
- 'Wait on / Wait off' timer function




Large intuitive display



Graphs for real-time process monitoring

Unique Structural Functions

- Opti-flow™, the fan speed control function based on sample and testing conditions
- Double stackable with stacking set (option) maximizes space utilization



Minimize media drying with Opti-flow™



Increased space utilization by double stacking (stacking set : option)
*Double stacking example: (top) forced convection oven (bottom) forced convection incubator

Supervise Equipment with Real-time Monitoring

- Monitor the temperature on BMS (Building Management System) for quick responses to abnormalities (option)
- Monitor via mobile app anytime, anywhere with LC Connected. (Requires LC GreenBox, option)



Timely responses via BMS-connected monitoring (BMS port: option)



Real-time monitoring via mobile devices (LC GreenBox: option)



Convenient Experiment and Data Checking

- Support RS-232/USB port and software making easy operation and data management from PC
- Test data are automatically stored and easily exported via USB; easy to edit data (in CSV format)
- Logs 36 recent events, including errors, opened door, etc
- Integrated with 21 CFR Part 11 compliant LC DataKeeper, suitable for GMP audits (option)



Data storage using USB ports




Store recent 36 events

Reliable Temperature Control Performance

- Temperature accuracy for wide range by 3-point calibration
- Ramp Control™, function that controls the rate of temperature rise by setting ramping rate(°C/min) per experiment
- Uniform temperature distribution through airflow optimized for uniform heat transfer (within ±0.41°C at 37°C)



Accurate temperature control with 3-point calibration



Ramp Control™ per experiment condition


Specification

Model		IB4-03V	IB4-05V	IB4-10V	IB4-15V
Chamber volume (L / cu ft)		25 / 0.88	60 / 2.12	111 / 3.92	146 / 5.16
Temperature	Range (°C / °F)	Amb.+5 ¹⁾ ~ 80°C / Amb.+9 ~ 176°F			
	Fluctuation at 37°C (°C / °F)	±0.1 / ±0.18	±0.15 / ±0.27	±0.1 / ±0.18	±0.1 / ±0.18
	Variation at 37°C (°C / °F)	±0.23 / ±0.41	±0.41 / ±0.74	±0.39 / ±0.70	±0.35 / ±0.63
	Heating time to 37°C (min.)	5	4	8	7
	Recovery time at 37°C (min.)	5	4	4	3
Control panel		5 inch TFT Color LCD			
Communication interface		USB-B, RS-232			
Function	Recently event record	36 ea			
	Recommended temperature	Suggest three frequently used temperature			
	Graph	Available viewing stored and real-time graphs			
	Data storage	Save your experiment data (in CSV format)			
	Temp. ramping rate (°C/min. °F/min)	0 ~ 8 0 ~ 14.4	0 ~ 6 0 ~ 10.8	0 ~ 5 0 ~ 9	0 ~ 4 0 ~ 7.2
Dimensions	Interior (W x D x H, mm / inch)	300 x 210 x 400 / 11.81 x 8.27 x 15.75	380 x 330 x 480 / 14.96 x 12.99 x 18.9	485 x 350 x 655 / 19.1 x 13.78 x 25.79	550 x 376 x 710 / 21.62 x 14.8 x 27.95
	Exterior (W x D x H, mm / inch)	478 x 527 x 604 / 18.8 x 20.75 x 23.78	558 x 648 x 684 / 21.97 x 25.5 x 26.93	663 x 677 x 859 / 26.1 x 26.65 x 33.82	728 x 704 x 914 / 28.66 x 27.7 x 36.0
	Net weight (kg / lbs)	32 / 70.55	41 / 90.39	56 / 123.46	62 / 136.69
Shelves	Quantity of shelves (standard / max.)	2 / 4	2 / 5	2 / 8	2 / 9
	Max. Load per shelf (kg / lbs)	20 / 44.09			
Electrical requirements (230V, 50/60Hz, A)		3			
Cat. No.		AAH21655K	AAH21665K	AAH21675K	AAH21685K
Electrical requirements (120V, 60Hz, A)		5.8			
Cat. No.		AAH21653U	AAH21663U	AAH21673U	AAH21683U

1) Able to control down to amb +5°C(amb +9°F) when the water tray is filled with water. (Water tray: option)

※ Technical data according to DIN 12880

Incubator (Forced Convection), Standard type




IB4-15S

Intuitive Touch Controller

- Intuitive interface with 3.5-inch display
- 'Wait on / Wait off' timer function




3.5-inch touch display



Wait on & Wait off timer

Reliable Temperature Control Performance

- Temperature accuracy for wide range by 3-point calibration
- The upper and lower temperature limits can be set, enabling alerts for deviations to trigger quick actions
- Uniform temperature distribution through airflow optimized for uniform heat transfer (within $\pm 0.41^{\circ}\text{C}$ at 37°C)



Accurate temperature control with 3-point calibration



Settings for temperature deviation

Convenient Experiment and Data Checking

- Support RS-232/USB port and software making easy operation and data management from PC (option)
- Monitor via mobile app anytime, anywhere with LC Connected. (Requires LC GreenBox, option)
- Monitor the temperature on BMS (Building Management System) for quick responses to abnormalities. (option)
- Integrated with 21 CFR Part 11 compliant LC DataKeeper, suitable for GMP audits (option)



Manipulation via PC (option)



Real-time monitoring via mobile devices (LC GreenBox: option)



Enhanced Safety and Convenience

- Implemented with the highest-level overheating prevention system (Patented in Korea)
- Over Temperature Limiter
- Complaint to IEC Class II protection standards and temperature protection in line with DIN 12880 Class 3.1
- Logs the last 12 events, including run times, stops, and errors
- Double stackable with stacking set (option) maximizes space utilization



Logs the recent 12 events



Increased space utilization by double stacking (stacking set : option)
*Double stacking example: (top) forced convection oven (bottom) forced convection incubator

Multiple Features to Minimize Incubation Failure

- Using Water tray delays media drying, making incubation time longer (option)
- Alarm is activated when door remains open longer than the set time to maintain internal temperature
- Automatically resume the control of previous operations when restored within 30 seconds from blackouts
- Controller can be locked with the button on the top, preventing any actions from unintended pressing



Water tray that minimizes media drying (option)



Audiovisual alerts to notify abnormalities

Specification

Model		IB4-03S	IB4-05S	IB4-10S	IB4-15S
Chamber volume (L / cu ft)		25 / 0.88	60 / 2.19	111 / 3.92	146 / 5.16
Temperature	Range (°C)	Amb.+5 ¹⁾ ~ 80°C / Amb.+9 ~ 176°F			
	Fluctuation at 37°C (°C / °F)	±0.1 / ±0.18	±0.15 / ±0.27	±0.1 / ±0.18	±0.1 / ±0.18
	Variation at 37°C (°C / °F)	±0.23 / ±0.41	±0.41 / ±0.74	±0.39 / ±0.70	±0.35 / ±0.63
	Heating time to 37°C (min.)	5	4	8	7
	Recovery time at 37°C (min.)	5	4	4	3
Control panel		3.5 inch TFT Color LCD			
Communication interface		USB-B, RS-232 (option)			
Function	Recently event record	12 ea			
Dimensions	Interior (W x D x H, mm / inch)	300 x 210 x 400 / 11.81 x 8.27 x 15.75	380 x 330 x 480 / 14.96 x 12.99 x 18.9	485 x 350 x 655 / 19.1 x 13.78 x 25.79	550 x 376 x 710 / 21.62 x 14.8 x 27.95
	Exterior (W x D x H, mm / inch)	478 x 527 x 604 / 18.8 x 20.75 x 23.78	558 x 648 x 684 / 21.97 x 25.5 x 26.93	663 x 677 x 859 / 26.1 x 26.65 x 33.82	728 x 704 x 914 / 28.66 x 27.7 x 36.0
	Net weight (kg / lbs)	32 / 70.55	41 / 90.39	56 / 123.46	62 / 136.69
Shelves	Quantity of shelves (standard/max.)	2 / 4	2 / 5	2 / 8	2 / 9
	Max. Load per shelf (kg / lbs)	20 / 44.09			
Electrical requirements (230V, 50/60Hz, A)		3			
Cat. No.		AAH21455K	AAH21465K	AAH21475K	AAH21485K
Electrical requirements (120V, 60Hz, A)		5.8			
Cat. No.		AAH21453U	AAH21463U	AAH21473U	AAH21483U

1) Able to control down to amb +5°C(amb +9°F) when the water tray is filled with water. (Water tray: option)

※ Technical data according to DIN 12880

Incubator (Forced Convection) Mini type



Small-sized incubator with compact design for ambient temperature use

Structural Functional Features

- Small incubator with compact design.
- Structure that circulates internal air for temperature control.
- Design with door that opens downward, minimizing the required space on the side of the device.
- Convenient internal observation door with a sight window.
- Stainless steel interior and shelves are excellent in terms of corrosion resistance and clean maintenance.




Use Convenience Features

- Accurate temperature control through temperature calibration and microprocessor PID.
- Clear VFD panel for convenient operation status checking.
- Easy operation with touch button and jog dial.
- Over-temperature and over-current limiter function.
- Provides notification/alarms when the door is open for an extended period of time.
- Alarm function when temperature sensor is disconnected.

IM-10

with Wire Shelves 2ea (standard)



※ When the chamber is operating with forced convection type, the amount of moisture evaporation of the medium is relatively high.

Specification

Model		IM-10
Chamber volume (L / cu ft)		10 / 0.4
Temperature	Range (°C / °F)	Amb. +5 to 65 / Amb. +9 to 149
	Fluctuation at 37°C (±°C / °F)	0.5 / 0.9
	Variation at 37°C (±°C / °F)	0.5 / 0.9
	Heating time to 37°C (min.)	5
	Recovery time at 37°C (min.)	2
Dimensions	Interior (W x D x H, mm / inch)	278 x 262 x 160 / 11 x 10.4 x 6.3
	Exterior (W x D x H, mm / inch)	400 x 410 x 264 / 15.7 x 16.1 x 10.4
	Net weight (Kg / lbs)	13.0 / 28.7
Shelves	Quantity of shelves (standard/max.)	2 / 4
Electrical requirements (230V, 50/60Hz, A)		1
Cat. No.		AAH24315K
Electrical requirements (120V, 60Hz, A)		2.1
Cat. No.		AAH24316U

※ The lowest controllable temperature depends on variations in room and sample temperatures. Please contact technical support for consultation before purchasing.

※ Technical data (according to DIN 12880 before 2013)

Accessories Page107 Wire Shelves, Perforated Shelves

Incubator (Forced Convection) Multi type



Independent control and effective space utilization per chamber

Structural Functional Features

- 2 / 4 chambers can be controlled independently.
- Easy to use with maximally space-efficient structure.
- Structure that circulates internal air for temperature control.
- Easy to open and close the door. Soft handle minimizes shock to the sample.
- Stainless steel interior and shelves are excellent in terms of corrosion resistance and clean maintenance.
- Easy to clean as inner edges are curved.
- Useful cable port and cover on the side.
- Equipped with easy-to-move/install caster.

Use Convenience Features

- Accurate temperature control through temperature calibration and microprocessor PID.
- Save and use 3 frequently used temperatures.
- Wait On/Off timer. (up to 99 hours 59 minutes)
- Mal-operation prevented by controller lock function.
- Over temperature limit function.
- Alarm function when temperature sensor is disconnected.
- Provides notification/alarms when the door is open for an extended period of time.




IB-02G-2C

with Wire Shelves
2ea per chamber (standard)

IB-02G-4C

with Wire Shelves
2ea per chamber (standard)



※ When the chamber is operating with forced convection type, the amount of moisture evaporation of the medium is relatively high.

Specification

Model		IB-02G-2C	IB-02G-4C
Chamber volume (L / cu ft)		120 / 4.2 (60 / 2.1 x 2 chambers)	240 / 8.5 (60 / 2.1 x 4 chambers)
Temperature	Range (°C / °F)	Amb. +5 to 70 / Amb. +9 to 158	Amb. +5 to 70 / Amb. +9 to 158
	Fluctuation at 37°C (±°C / °F)	0.1 / 0.18	0.1 / 0.18
	Variation at 37°C (±°C / °F)	0.6 / 1.08	0.6 / 1.08
	Heating time to 37°C (min.)	16	16
	Recovery time at 37°C (min.)	3	3
Dimension	Interior (W x D x H, mm / inch)	400 x 360 x 420 / 15.7 x 14.2 x 16.5 (each chamber)	400 x 360 x 420 / 15.7 x 14.2 x 16.5 (each chamber)
	Exterior (W x D x H, mm / inch)	570 x 640 x 1360 / 22.4 x 25.2 x 53.5	1170 x 640 x 1360 / 46.1 x 25.2 x 53.5
	Net weight (Kg / lbs)	110 / 242.5	168 / 370.4
Shelves	Quantity of shelves (standard/max.)	2 / 4 (each chamber)	2 / 4 (each chamber)
	Max. Load per shelf (Kg / lbs)	30 / 66	30 / 66
Electrical requirements (230V, 50/60Hz, A)		4.1	4.1 x 2ea
Cat. No		AAH2A115K	AAH2A015K
Electrical requirements (120V, 60Hz, A)		7.5	7.5 x 2ea
Cat. No		AAH2A116U	AAH2A016U

※ The lowest controllable temperature depends on variations in room and sample temperatures. Please contact technical support for consultation before purchasing.

※ Technical data (according to DIN 12880 before 2013)

Accessories Page 107 Wire Shelves, Perforated Shelves

Low Temperature Incubator (Forced Convection/Air-jacket) General type

Provides stable, low-temperature culture environment



IL3-15A

with Wire Shelves 2ea (standard)



IL3-25

with Wire Shelves 3ea (standard)

Structural Functional Features

- Compact structure design.
- Clear VFD panel for clear information display.
- The inner toughened-glass door is completely enclosed with silicone packing, making it easy to observe the sample without changing the temperature.
- External door can be opened and closed smoothly with magnetic packing to minimize sample impact.
- The sample space is independently kept clean and the structure protects against sources of contamination.
- Internal power outlet included as standard for easy operation of small products.
- Easily removable condenser grill structure.
- Stainless steel interior and shelves are excellent in terms of corrosion resistance and clean maintenance.
- Easy to clean as inner edges are curved.

Use Convenience Features

- Accurately controlling temperature based on temperature calibration /auto-tuning / microprocessor PID.
- Save and use 3 frequently used temperatures.
- Wait On/Off timer. (up to 999 hours 59 minutes)
- USB/RS-232 connection and software provide convenient computer operation and data managing.
- Operation and data management with computer connection.
- Alarm function when temperature sensor is disconnected.
- Monitor via mobile app anytime, anywhere with LC Connected. (Requires LC GreenBox, option)


Outstanding Safety

- Top-rated overheating protection system. (registration KR 10-0397583)
- Mechanical/ electronic overheat protection.
- Set the temperature deviation notification.
- Mal-operation prevented by controller lock function.




Temperature Variation

- Superior temperature uniformity with excellent temperature control capability.




Two Series Of Models According To Heat Transfer Method

- Forced Convection**
The structure allows for the inside air to be circulated and the temperature to be controlled, producing stable temperature control and good distribution. However, the amount of moisture evaporation of the culture is relatively high.
- Natural Convection (Air-jacket)**
The air-jacket structure heats the outside air in the space where the sample is placed and transfers the heat to the inside, and the moisture evaporation amount of the culture is relatively low.



Natural Convection (Air-jacket)



Forced Convection

Specification

Model		Forced Convection		Natural Convection (Air-jacket)	
		IL3-15	IL3-25	IL3-15A	IL3-25A
Chamber volume (L / cu ft)		150 / 5.3	242 / 8.5	156 / 5.5	254 / 8.97
Refrigerator (Hp)		1 / 6	1 / 6	1 / 6	1 / 6
Temperature	Range (°C / °F)	0 to 80 / 32 to 176	0 to 80 / 32 to 176	4 to 70 / 39.2 to 158	4 to 70 / 39.2 to 158
	Fluctuation at 37°C (±°C / °F)	0.1 / 0.18	0.1 / 0.18	0.2 / 0.36	0.2 / 0.36
	Variation at 37°C (±°C / °F)	0.8 / 1.4	1.0 / 1.8	0.7 / 1.3	1.0 / 1.8
	Heating time 20 to 37°C (min.)	43	45	45	50
	Cooling time 20 to 5°C (min.)	45	50	53	58
Recovery time at 37°C (min.)		4	4	5	5
Dimension	Interior (W x D x H, mm / inch)	600 x 500 x 500 / 24 x 20 x 20	500 x 520 x 930 / 20 x 20.5 x 36.6	600 x 520 x 500 / 24 x 20.5 x 20	500 x 540 x 940 / 20 x 21.3 x 37
	Exterior (W x D x H, mm / inch)	720 x 760 x 1070 / 28.3 x 29.9 x 42.1	620 x 770 x 1618 / 24.4 x 30.3 x 63.7	770 x 790 x 1115 / 30.3 x 31.1 x 43.9	670 x 805 x 1668 / 26.4 x 31.7 x 65.7
	Net weight (Kg / lbs)	100 / 220	135 / 297	108 / 238.1	145 / 319.7
Shelves	Quantity of shelves (standard / max.)	2 / 6	3 / 13	2 / 6	3 / 13
	Max. Load per shelf (Kg / lbs)	26 / 57.3	30 / 66	26 / 57.3	30 / 66
Electrical requirements (230V, 50/60Hz, A)		4.5	6.4	6.6	8
Cat. No.	AAH220412K	AAH220512K	AAH220612K	AAH220712K	
Electrical requirements (120V, 60Hz, A)		8.3	10.6	10.2	12.3
Cat. No.	AAH220462U	AAH220562U	AAH220632U	AAH220732U	

※ Technical data (according to DIN 12880 before 2013)

Accessories Page 107 Wire Shelves, Perforated Shelves, LC GreenBox

Low Temperature Incubator Personal type



Low temperature personal incubator

Structural Functional Features

- Small incubator with compact design.
- Minimizes vibration and power consumption by using electronic cooling element. (peltier module)
- Structure that circulates internal air for temperature control.
- The inner glass door is completely enclosed with silicone packing, making it easy to observe the sample without changing the temperature.
- External door can be opened and closed smoothly with magnetic packing to minimize sample impact.
- Stainless steel interior and shelves are excellent in terms of corrosion resistance and clean maintenance.



Use Convenience Features


- Accurately controlling temperature based on temperature calibration / auto-tuning / microprocessor PID.
- Save and use 3 frequently used temperatures.
- Wait On/Off timer. (up to 99 hours 59 minutes)
- Operation and data management with computer connection.
- Alarm function when temperature sensor is disconnected.
- Internal air circulation speed and defrost cycle control.
- Monitor via mobile app anytime, anywhere with LC Connected. (Requires LC GreenBox, option)

ILP-02 (stackable unit)

with Wire Shelf 1ea (standard)

ILP-12

with Wire Shelf 2ea (standard)



Outstanding Safety

- Top-rated overheating protection system. (registration KR 10-0397583)
- Over temperature limit function.
- Mal-operation prevented by controller lock function.

Convenient Program Control

- Various culture and incubating experiments with the 9 step program.
- Temperature/time settings for each step.
- Up to 99 hours and 59 minutes.
- Repeat operation up to 200 times.

Specification

Model		ILP-02	ILP-12
Chamber volume (L / cu ft)		14 / 0.5	48 / 1.7
Temperature	Range (°C / °F)	5 to 40 / +41 to 104 (at 20°C / 68°F) 10 to 40 / +50 to 104 (at 25°C / 77°F) 15 to 40 / +59 to 104 (at 30°C / 86°F)	5 to 40 / +41 to 104 (at 20°C / 68°F) 10 to 40 / +50 to 104 (at 25°C / 77°F) 15 to 40 / +59 to 104 (at 30°C / 86°F)
	Fluctuation at 25°C (±°C / °F)	0.1 / 0.18	0.1 / 0.18
	Variation at 25°C (±°C / °F)	0.5 / 0.9	0.4 / 0.72
	Heating time to 40°C (min.)	17	9
	Cooling time 40 to 5°C (min.)	60	60
Dimension	Recovery time at 25°C (min.)	1.5	1.5
	Interior (W x D x H, mm / inch)	315 x 200 x 230 / 12.4 x 7.9 x 9.1	334 x 334 x 430 / 13.1 x 13.1 x 17
	Exterior (W x D x H, mm / inch)	430 x 495 x 400 / 17 x 19.4 x 15.7	430 x 665 x 620 / 17 x 26.2 x 24.4
Shelves	Net weight (Kg / lbs)	27.5 / 60.6	49 / 108
	Quantity of shelves (standard / max.)	1 / 6	2 / 13
Electrical requirements (230V, 50/60Hz, A)	Max. Load per shelf (Kg / lbs)	15 / 33.0	20 / 44.1
Cat. No.	Electrical requirements (120V, 60Hz, A)	1.8	3.5
		AAH24032K	AAH24042K
		3.4	6.7
Cat. No.		AAH24033U	AAH24043U

※ Technical data (according to DIN 12880 before 2013)

Accessories Page 107 Wire Shelves, Perforated Shelves, LC GreenBox

Low Temperature Incubator Multi type (Forced Convection)



Independent control and efficient space utilization for each chamber

Structural Functional Features

- Each chamber can be controlled independently and easy to use with maximally space-efficient structure.
- Structure that circulates internal air for temperature control.
- The inner glass door is completely enclosed with silicone packing, making it easy to observe the sample without changing the temperature.
- External door can be opened and closed smoothly with magnetic packing to minimize sample impact.
- Stainless steel interior and shelves are excellent in terms of corrosion resistance and clean maintenance.
- Easy to clean as inner edges are curved.



IL-11-2C




IL-11-4C

Incubator

Use Convenience Features

- Accurately controlling temperature based on temperature calibration / auto-tuning / microprocessor PID.
- Save and use 3 frequently used temperatures.
- Wait On/Off timer. (up to 99 hours 59 minutes)
- Alarm function when temperature sensor is disconnected.
- Internal air circulation speed and defrost cycle control

※ When the chamber is operating with forced convection type, the amount of moisture evaporation of the medium is relatively high.



Outstanding Safety

- Top-rated overheating protection system. (registration KR 10-0397583)
- Over temperature limit function.
- Mal-operation prevented by controller lock function.

Convenient Program Control

- Various culture and incubating experiments with the 9 step program.
- Temperature/time settings for each step.
- Up to 99 hours and 59 minutes.
- Repeat operation up to 200 times.

Specification

Model		IL-11-2C		IL-11-4C	
Chamber volume (L / cu ft)		300 / 10.6 (150 / 5.3 x 2 chambers)		600 / 21.2 (150 / 5.3 x 4 chambers)	
Refrigerator (Hp)		1/8 x 2ea		1/8 x 4ea	
Temperature	Range (°C / °F)	0 to 60 / 32 to 140		0 to 60 / 32 to 140	
	Fluctuation at 25°C (±°C / °F)	0.1 / 0.18		0.1 / 0.18	
	Variation at 25°C (±°C / °F)	0.5 / 0.9		0.5 / 0.9	
	Heating time to 40°C (min.)	16		16	
	Cooling time 40 to 5°C (min.)	53		53	
Dimension	Recovery time at 40°C (min.)	2		2	
	Interior (W x D x H, mm / inch)	600 x 500 x 500 / 23.6 x 19.7 x 19.7 (each chamber)		600 x 500 x 500 / 23.6 x 19.7 x 19.7 (each chamber)	
	Exterior (W x D x H, mm / inch)	760 x 760 x 1690 / 29.9 x 29.9 x 66.1		1500 x 760 x 1690 / 59.1 x 29.9 x 66.1	
Shelves	Net weight (Kg / lbs)	180 / 397		350 / 772	
	Quantity of shelves (standard / max.)	2 / 12 (each chamber)		2 / 12 (each chamber)	
Electrical requirements (230V, A)	Max. Load per shelf (Kg / lbs)	26 / 57.3		26 / 57.3	
		50Hz, 9	60Hz, 9	50Hz, 9 x 2ea	60Hz, 9 x 2ea
Cat No		AAH28112K	AAH28111K	AAH28012K	AAH28011K
Electrical requirements (120V, 60Hz, A)		16.5		16.5 x 2ea	
Cat No		AAH28113U		AAH28013U	

※ Technical data (according to DIN 12880 before 2013)

※ This product is forced convection type. If the sample dries easily or is sensitive, we recommend using the natural convection type.

Accessories Page 107 Wire Shelves, Perforated Shelves

Plant Growth Chamber

Provision of optimum photosynthesis conditions for plant growth



GC-1000TLH



GC-300TLH

Optimal photosynthesis condition for plant growth.

- Uniform temperature and humidity. (TLH models)
- Reinforced Blue and Red spectrums lighting system for plant photosynthesis.
- Optional CO₂ sensor. (option)
- Programmable temperature, humidity, and lighting.
- Progressive temperature, humidity, illumination program control for optimization of plant growth environment for night and day.


Structural Functional Features

General control system


- 5°C to 50°C(lamp off) / 10°C to 50°C. (lamp on)
- Max. to 35,000Lux for GC-1000.
Max. to 20,000Lux for GC-300 Models.
- 40 to 80% RH for GC-300TLH. (at 20 to 35°C)
50 to 90% RH for GC-1000TLH. (at 20 to 35°C)
- Max. 5,000ppm CO₂ on/off system. (option)
- Microprocessor PID control / Temperature calibration / Automatic tuning.
- 10 step programmable temperature, humidity, and illumination profiles and repeatable steps of up to 999 cycles.
- Advanced performance of low and high temperature and humidity control, achieved by adapting the humidification tank exterior to the chamber, added benefit for long term test.

Illuminance control system

- Uniformed luminescence distribution.
- Broad distribution of side lamps for hastening the growth of plants.
- High intensity illumination of upper lamp for light efficient and low thermal load. (for GC-1000TLH/1000TL)
- Unique construction for minimization of heat increase from surrounding lamps.
 - Tempered glass door blocks heated air from lamps.
 - Designed to exhaust heated air through upper vent holes. (for GC-1000TLH)
- Stable and long lasting lighting through introduction of high frequency electronic ballast lamps.



GC-300




GC-1000

Outstanding Safety

- Automatic shut off after overheat alarm.
- Low and empty water level alarm.
- Power supply leakage breaker.
- Over current protection.
- Open door alarm.


Lamps setting

- Control illumination values with the below lamp setting arrangements.



[FRONT]

GC-300 lamp arrangement




[FRONT]

GC-1000 lamp arrangement

Use Convenience Features


- Ergonomically designed door construction.
 - Inner glass door with silicone and external door with magnetic sealing for dual airtight packing.
 - Well designed providing a smooth open for minimization of damage to plants.
- Wide inner tempered glass door for clear observation of plant growth without affecting inner chamber's environment.
- Tall plants can be grown inside of the chamber with adjustment of shelf level.
- Maintenance of water supply are easily performed with a front water tank. (GC-1000TLH)
 - Water level check indicated by water level bar in the tank.
- Water supply during operation is also available for long term test purposes.
- Adjustable water tank position. (for GC-300TLH)
- Detachable condenser air-filter for easy maintenance of refrigerating efficiency.
 - Maintenance of air filter no longer cumbersome with our detachable condenser air filter.
- Casters for easy mobility during installation or relocation.
- Eco-friendly CFC-free refrigerant use.



Specification


Model		GC-300TL	GC-300TLH	GC-1000TLH
Chamber volume (L / cu ft)		300 / 10.6	300 / 10.6	1000 / 35.3
Temperature	Range (°C / °F)	5 to 50 / 41 to 122 - Lamp off 10 to 50 / 50 to 122 - Lamp on 20 to 50 / 68 to 122 - with humidity		
	Fluctuation at 25°C (±°C / °F)	0.4 / 0.72 - without humidity		
	Variation at 25°C (±°C / °F)	1.0 / 1.8 - without humidity		
Humidity	Range (%RH)	-	40 to 80 at 20 to 35°C 70 to 90 at 36 to 50°C	50 to 90 at 20 to 35°C 60 to 90 at 36 to 50°C
	Fluctuation at 60%RH (±%RH)	-	3	3
Illumination	Range (Lux)	0 to 20000	0 to 20000	0 to 35000
	Control (steps)	10	10	11
	Fluorescent lamp (W)	32 x 18ea	32 x 18ea	32 x 16ea
	Metal lamp (W)	-	-	400 x 3ea
CO₂ (Option)	Range	750 to 5000	750 to 5000	750 to 5000
	Sensor	NDIR CO ₂ sensor	NDIR CO ₂ sensor	NDIR CO ₂ sensor
Dimensions	Interior (W x D x H, mm / inch)	510 x 540 x 1100 / 20 x 21.3 x 43.3	510 x 540 x 1100 / 20 x 21.3 x 43.3	1200 x 800 x 1080 / 47.2 x 31.5 x 42.5
	Exterior (W x D x H, mm / inch)	700 x 805 x 1900 / 27.6 x 31.7 x 74.8	700 x 805 x 1900 / 27.6 x 31.7 x 74.8	1410 x 1070 x 2150 / 55.5 x 42.1 x 84.6
	Net weight (Kg / lbs)	250 / 551.2	250 / 551.2	550 / 1212.5
Electrical requirements (230V, A)		50/60Hz, 12.5	50/60Hz, 16	50Hz, 23.2 60Hz, 30
Cat. No.		AAHA1012K	AAHA1022K	AAHA1032K AAHA1031K

Accessories for Incubator



Wire Shelf

- (2 shelves included as standard)
- Stainless Steel
- Uniform temperature distribution based on high air permeability



Perforated Shelf

- Stainless Steel
- Robust structure is ideal for heavy load samples



Recorder (Dot)

- 6-channel graph-type temperature recorder.
- Adjustable recording range, speed, etc.



CO₂ Control System

- CO₂ control system for optimal photosynthetic conditions (Sensor, PCB, Valve, etc)

Model	Wire Shelves		Perforated Shelves		Recorder	CO ₂ Control System
	Cat. No.	Quantity of shelves (standard/max.)	Cat. No.	Quantity of shelves (max.)	Cat. No.	Cat. No.
IB4-03	00FDA0009739	2 / 4	AAA125341	4	-	-
IB4-05	00FDA0009738	2 / 5	AAA125342	5	-	-
IB4-10	00FDA0009737	2 / 8	AAA125343	8	-	-
IB4-15	00FDA0009736	2 / 9	AAA125344	9	-	-
IM-10	00EDA0008244	2 / 4	AAA24511	4	-	-
IB-02G-2C/4C ¹⁾	00EDA0008136	2 / 4	AAA12531	4	-	-
ILP-02	00EDA0008223	1 / 6	AAA24501	6	-	-
ILP-12	00EDA0008221	2 / 13	AAA24502	13	-	-
IL3-15/15A	00EDA0008219	2 / 6	AAA22521	6	-	-
IL3-25/25A	00EDA0008220	3 / 13	AAA22522	13	-	-
IL-11-2C/4C ¹⁾	00EDA0008219	2 / 12	AAA12532	12	-	-
GC-300TL	00EDA0008220	3 / 14	AAA22522	14	AAAA1501	AAAA1521
GC-300TLH	00EDA0008220	3 / 14	AAA22522	14	AAAA1501	AAAA1521
GC-1000TLH	00EDA0008222	6 / 29	AAAA1512	29	AAAA1502	AAAA1521

1) Standard for each chamber



LC DataKeeper

- High-integrity software compliant with FDA 21 CFR Part 11
- Integrated with Jeio Tech equipment to undergo GMP audits



LC GreenBox

- Real-time monitoring and control via mobile app
- Connect up to 4 devices to one
- Communication port connectivity via a RS-232 port



BMS Port

- Real-time equipment monitoring from a central monitoring facility
- Useful for safety management at laboratories that operate multiple devices



Stacking Set (-05 / -15 models only)

- Double stackable
- Firm stacking allows efficient utilization of space
(※ When double stacking, the minimum control temperature of the upper chamber may be limited due to the heat generated by the lower chamber)

Model	IB4-03	IB4-05	IB4-10	IB4-15
LC DataKeeper-Basic ¹⁾			SSA11001	
LC DataKeeper-Pro ¹⁾			SSA11002	
LC GreenBox			AAHQ1011K	
BMS Port			AAA125534	
Stacking Set	-	AAA125547	-	AAA125549
Water Tray	AAA125443	AAA125444	AAA125445	AAA124554
Communication Port ²⁾			AAA125533	

1) Trial version is available for 30 days, offering the same features as the Basic version, and can be requested through the website.

2) Communication port is provided as a standard in advanced (-V) model and is an optional add-on for the basic (-S) model for connecting LC DataKeeper, LC GreenBox, BMS, and PC.

Incubated Shaker



General Application

Molecular biology, animal/plant cell culture, fermentation, genetics, new drug development, etc.

	Description	Temp. Range (°C/°F)	Motion Type	Max. Speed (rpm)	250 ml Flask Capacity	Model	Page	
Multi-type Stackable	Stackable Model <ul style="list-style-type: none"> • Double-stacked for space efficiency • Includes an electric outlet and wire shelves to use small equipment 	Amb. +5 to 80 / Amb. +9 to 176	Orbital	500	24	ISS	112	
		Amb. -20 (min. 4) to 80 / Amb. -36 (min. 7.2) to 176		500	24	ISS-R	112	
	Stackable Model (Multi) <ul style="list-style-type: none"> • Maximum triple-stacked for space efficiency • Powerful DD motor for large-capacity shaking 	Amb. +10 to 80 / Amb. +18 to 176	Orbital	400	39	ISS-7100	114	
		Amb. -20 (min. 4) to 80 / Amb. -36 (min. 7.2) to 176		400	39	ISS-7100R	114	
		Amb. +10 to 80 / Amb. +18 to 176		300	39	ISS-7200	114	
		Amb. -20 (min. 4) to 80 / Amb. -36 (min. 7.2) to 176		300	39	ISS-7200R	114	
	Single-type	Benchtop Model <ul style="list-style-type: none"> • Convenient internal observation • Optimized for use on test bench top 	Amb. +5 to 80 / Amb. +9 to 176	Orbital	500	20	IST	116
			Amb. -15 (min. 15) to 80 / Amb. -27 (min. 27) to 176		500	20	IST-R	116
Floor Standing Model <ul style="list-style-type: none"> • Optimized for large-capacity of culture • Large capacity up to 5ea x 6 L flasks at the same time 		Amb. +5 to 80 / Amb. +9 to 176	Orbital	500	39	ISF-7100	118	
		Amb. -20 (min. 4) to 80 / Amb. -36 (min. 7.2) to 176		500	39	ISF-7100R	118	
		Amb. +5 to 80 / Amb. +9 to 176		300	39	ISF-7200	118	
		Amb. -20 (min. 4) to 80 / Amb. -36 (min. 7.2) to 176		300	39	ISF-7200R	118	

※ The contents of the above and the contents of this catalog may differ depending on the specific model and conditions of use. For the information about the features and specifications that applying to each models, please check the information on the corresponding page of each models.


 Satisfaction with the requirements of European regulations and guidelines.	 Certification of conformity of electricity, gas, etc. in accordance with the safety standards of Canada and USA.	 Real-time equipment monitoring and control system using mobile app.
 Registered patent based on Jeio Tech's proprietary technology.	 Compliance with Guidelines for Restriction of Hazardous Substances Use.	 Eco-friendly product that has passed environmental impact assessment tests.
 Intuitive operation with color touch display.	 Orbital motion either clockwise or counterclockwise direction can be set.	 Possible to set end time or start time of device operation.
 Warning alarm in case of deviation of temperature control.	 Electronic over temperature protection system.	 Conventional over temperature protection, a kind of backup device against electronic over temperature protection failure.
 PC communication via RS-232 / RS-485 / USB port.	 It memorizes the state of power failure and operates automatically when power is restored.	 More precise temperature control through 3-point temperature calibration.
 5 year warranty for direct drive motor	 2 year warranty free A/S.	

Incubated Shaker



Reliable precision temperature control up to 80°C Innovatory temperature control system using Jeio Tech's technology

- > **Extended temperature range up to 80°C**
Broad range of applications, including culturing of thermophilic microorganisms at a temperature range up to 80°C, which is a wider range than previous models.
- > **Remote monitoring service (option)**
Monitor via mobile app anytime, anywhere with LC Connected (Requires LC GreenBox, option)
- > **More stable temperature control**
Optimized Microprocessor PID control for more accurate temperature control.
- > **Temperature calibration function**
Temperature calibration of up to 3 points allows for more reliable temperature control.
- > **Real-time monitoring of temperature deviation**
Prevents experiment errors through user notification when the actual temperature differs from the set temperature.
- > **Strong and stable shaking**
DD motor & triple cam system provide a powerful and stable shaking motion even under high loads. (ISF, ISS-7000) (registration KR 10-1678160)
Provide 5-year warranty for DD motor.
- > **Electronic over temperature protection system**
Independent precision temperature sensor and controller provides more secure and reliable over temperature protection. (ISF, ISS)
- > **Conventional over-temperature protection system**
Backup device against electronic system failure. Implementation of over-temperature protection even for electronic errors.



USB, RS-232 connection port provided.



Easy to attach and detach large-volume sample (ISF) thanks to the inclined structure. (registered design 30-2014-0060254)



Viewing window and LED lighting for internal observation.



Optimized model based on shape and performance Easy and intuitive operation with color touch screen display

- > **Achieve intuitive and convenient operation**
Color LCD touch screen with excellent touch feeling makes researcher easy operation and it shows a status of operating.
- > **One action control provides convenience**
Convenient operation environment with easy Start/Stop button for fast and simple operation and stop.
- > **Complete model lineup according to usage environment**
Three models (tabletop, stackable, floor) provide the ability to choose the optimized model for the user environment.
- > **Maximal utilization with three operation modes**
The three modes of temperature and shaking, temperature, or shaking, provide the ability to perform various experiments with a single device.
- > **Computer interface environment provided**
USB, RS-232 connection and dedicated software provide convenient device operation and data accumulation/management function.
- > **Convenient sample observation environment provided**
Low-power consumed LED illumination provided as standard provides a convenient environment for sample observation.
- > **Lock function of controller**
Lock function for preventing change of controller settings prevents erroneous operation due to arbitrary changing.
- > **Spill drain port for internal cleanliness and easy maintenance**
Even if the sample leaks inside the equipment, it can be cleaned through the spill drain port.

※ Some of the above contents are limited to specific models.



Color LCD touch screen.



Safe use with over temperature protection function.



Spill drain port for clean internal management.

Incubated Shaker Stackable type

Double stackable for space efficiency

Structural Functional Features

- Best space-efficient model with double stacking.
- Up to 80°C, 500rpm enables setting of various experimental conditions.
- Three operating modes: Temperature & Shaking, Temperature, or Shaking.
- Includes internal power outlet.
- Shaking structure with minimized vibration.

Use Convenience Features

- Accurate rpm and temperature control through temperature auto-tuning / microprocessor PID and calibration.
- Monitor via mobile app anytime, anywhere with LC Connected. (Requires LC GreenBox, option)
- Countdown timer (999 hours 59 minutes) of shaking operation.
- Intuitive operation with color touch display.
- High brightness / low-power LED illumination makes it easy to observe inside.
- Spill drain port for internal cleanliness and easy maintenance.


Outstanding Safety

- Electronic over temperature protection system. (advanced over temperature (A-OT) limiter)
Independent precision temperature sensor and controller provides more secure and reliable over temperature protection.
- Conventional over-temperature protection system. (B-OT, backup over temperature limiter)
A backup device against electronic system failure, mechanically preventing over temperature even when there are electronic errors.
- Soft start and stop shaking function.
- Malfunction prevented by controller lock function.
- Displays all information such as device operation and error occurrence.
- Overcurrent protection device restricts operation when objects are caught inside the system and shaking is not possible.



ISS-3075R

with Shelf 2ea (standard)
Universal Platform, Flask Clamp (option)
Magnetic Stirrer (MS-17G)



Use the stacking kit for double stacked (option)



Color touch display



Smart Self-diagnosis

Real-time self-diagnostics to check equipment for abnormal conditions.

- Alarm is generated by monitoring difference between actual rpm/temperature and set rpm/temperature.
- Automatic restart as it was operated after power recovery.



Specification

Model	ISS-3075	ISS-3075R	ISS-4075	ISS-4075R
Temperature				
Range (°C / °F) ¹⁾	Amb. +5 to 80 / Amb. +9 to 176	Amb. -20 (Min. 4) to 80 / Amb. -36 (Min. 7.2) to 176	Amb. +5 to 80 / Amb. +9 to 176	Amb. -20 (Min. 4) to 80 / Amb. -36 (Min. 7.2) to 176
Fluctuation at 37°C in flask (±°C / °F)	0.1 / 0.18	0.1 / 0.18	0.1 / 0.18	0.1 / 0.18
Variation at 37°C in flask (±°C / °F)	0.5 / 0.9	0.5 / 0.9	0.5 / 0.9	0.5 / 0.9
Refrigerator (Hp)	-	1/6	-	1/6
Shaking system				
Motion type	Orbital	Orbital	Orbital	Orbital
Speed range (rpm)	20 to 500 (stackable : 20 to 250)			
Max. speed per permissible load (rpm)	500 (max. 10 kg / 22.0 lbs) 400 (max. 15 kg / 33.1 lbs) 250 (max. 15 kg / 33.1 lbs) - Stacked	500 (max. 10 kg / 22.0 lbs) 400 (max. 15 kg / 33.1 lbs) 250 (max. 15 kg / 33.1 lbs) - Stacked	500 (max. 10 kg / 22.0 lbs) 400 (max. 21 kg / 46.3 lbs) 250 (max. 21 kg / 46.3 lbs) - Stacked	500 (max. 10 kg / 22.0 lbs) 400 (max. 21 kg / 46.3 lbs) 250 (max. 21 kg / 46.3 lbs) - Stacked
Amplitude size (mm / inch, dia.)	19.1 / 0.75	19.1 / 0.75	19.1 / 0.75	19.1 / 0.75
Speed stability ²⁾ (rpm, %)	±1	±1	±1	±1
Timer (count-down)	999hr 59min	999hr 59min	999hr 59min	999hr 59min
Dimension				
Chamber volume (L / cu ft)	80 / 2.8	80 / 2.8	150 / 5.3	150 / 5.3
Platform (W x D, mm / inch)	350 x 350 / 13.8 x 13.8	350 x 350 / 13.8 x 13.8	450 x 450 / 17.7 x 17.7	450 x 450 / 17.7 x 17.7
Interior (W x D x H, mm / inch)	440 x 440 x 418 / 17.3 x 17.3 x 16.5	440 x 440 x 418 / 17.3 x 17.3 x 16.5	540 x 540 x 518 / 21.3 x 21.3 x 20.4	540 x 540 x 518 / 21.3 x 21.3 x 20.4
Exterior (W x D x H, mm / inch)	578 x 815 x 854 / 22.8 x 32.1 x 33.6	578 x 815 x 854 / 22.8 x 32.1 x 33.6	678 x 895 x 934 / 26.7 x 35.2 x 36.8	678 x 895 x 934 / 26.7 x 35.2 x 36.8
Quantity of shelves (standard / max.)	2 / 8	2 / 8	2 / 11	2 / 11
Net weight (kg / lbs)	153 / 337.3	160 / 352.7	183 / 403.4	190 / 418.9
Electrical requirements (230V, 50/60Hz, A)	6.2	7.9	6.2	7.9
Cat. No.	AAH23432K	AAH23532K	AAH23632K	AAH23732K
Electrical requirements (120V, 60Hz, A)	9.5	13	9.5	13
Cat. No.	AAH23415U	AAH23515U	AAH23615U	AAH23715U

1) The temperature control range of the ISS-3075 and ISS-4075 models are from 5°C above the ambient temperature (15°C to 40°C) to 80°C. The temperature control range of the ISS-3075R and ISS-4075R models is from 20°C lower than the ambient temperature. (18°C to 35°C) of the device to 80°C. (However, the lowest control temperature is 4°C)

2) Actual rpm allowable error for set rpm is ±1rpm if setting is less than 100rpm, ±1% if setting is higher than 100rpm.

Accessories Page 120

Universal Platform, Flask Clamp, Funnel Clamp, Test Tube Rack, Microplate Holder, Lab Sticker, Spring Wire Rack, Dedicated Platform, Rubber Mat, Universal Attachment, Rectangular Foot Caster, Stacking Kit, Wire Shelves, LC GreenBox, Starter Kit

Incubated Shaker

Multi-stackable type

Maximal space-saving multi-stackable model

Structural Functional Features

- Large-capacity triple stackable incubated shaker.
- Up to 80°C, 400rpm enables experimental conditions for various applications.
- Three operating modes: Temperature & Shaking, Temperature, or Shaking.
- Shaking structure with minimized vibration.
- Convenient to change samples with glide-up door.
- Stable shaking of high-load samples. (registration KR 10-1678160)
- Applied a direct drive motor that has outstanding durability, a 5-year warranty for the motor.

Use Convenience Features

- Accurate rpm and temperature control through temperature auto-tuning / microprocessor PID and calibration.
- Monitor via mobile app anytime, anywhere with LC Connected. (Requires LC GreenBox, option)
- Countdown timer (999 hours 59 minutes) of shaking operation.
- Intuitive operation with color touch display.
- High brightness / low-power LED illumination makes it easy to observe inside.
- Spill drain port for internal cleanliness and easy maintenance.


Outstanding Safety

- Electronic over temperature protection system. (A-OT, advanced over temperature limiter)
Independent precision temperature sensor and controller provides more secure and reliable over temperature protection.
- Conventional over-temperature protection system. (B-OT, backup over temperature limiter)
A backup device against electronic system failure, mechanically preventing over temperature even when there are electronic errors.
- Soft start and stop shaking function.
- Malfunction prevented by controller lock function.
- Displays all information such as device operation and error occurrence.
- Overcurrent protection device restricts operation when objects are caught inside the system and shaking is not possible.



ISS-7100R

with Universal Platform, Flask Clamp (option)



View the introduction video



Door lift-up structure for safe and convenient sample handling.



Sliding platform (Push-pull Slide™, Patented in Korea)



Color touch display.



Smart Self-diagnosis

Real-time self-diagnostics to check equipment for abnormal conditions.


- Alarm is generated by monitoring difference between actual rpm/temperature and set rpm/temperature.
- Automatic restart as it was operated after power recovery.

DD Shaking System Patented

Powerful DD motor for large-capacity shaking
The triple cam system provides stable shaking with high load samples.

(Patented in Korea)

Multi-Stacking Options



Specification

Model	ISS-7100	ISS-7100R	ISS-7200	ISS-7200R
Temperature				
Range (°C / °F) ¹⁾	Amb. +10 to 80 / Amb. +18 to 176	Amb. -20 (Min. 4) to 80 / Amb. -36 (Min. 7.2) to 176	Amb. +10 to 80 / Amb. +18 to 176	Amb. -20 (Min. 4) to 80 / Amb. -36 (Min. 7.2) to 176
Fluctuation at 37°C in flask (±°C / °F)	0.1 / 0.18	0.1 / 0.18	0.1 / 0.18	0.1 / 0.18
Variation at 37°C in flask (±°C / °F)	0.6 / 1.08	0.5 / 0.90	0.6 / 1.08	0.5 / 0.90
Refrigerator (Hp)	-	1/6	-	1/6
Shaking system				
Motion type	Orbital	Orbital	Orbital	Orbital
Speed range (rpm)	30 to 400	30 to 400	30 to 300	30 to 300
Max. speed per permissible load (rpm)	1 unit only 400 (max. 20 kg / 44.1 lbs)	1 unit only 400 (max. 20 kg / 44.1 lbs)	1 unit only 300 (max. 20 kg / 44.1 lbs)	1 unit only 300 (max. 20 kg / 44.1 lbs)
	Stacked units 250 (max. 15 kg / 33.1 lbs) 200 (max. 20 kg / 44.1 lbs)	Stacked units 250 (max. 15 kg / 33.1 lbs) 200 (max. 20 kg / 44.1 lbs)	Stacked units 200 (max. 15 kg / 33.1 lbs) 150 (max. 20 kg / 44.1 lbs)	Stacked units 200 (max. 15 kg / 33.1 lbs) 150 (max. 20 kg / 44.1 lbs)
	Amplitude size (mm / inch, dia.)	25.4 / 1	25.4 / 1	50.8 / 2
Speed stability ²⁾ (rpm, %)	±1	±1	±1	±1
Timer (count-down)	999 hr 59 min.	999 hr 59 min.	999 hr 59 min.	999 hr 59 min.
Dimension (1 unit)				
Chamber volume (L / cu ft)	198 / 7.0	198 / 7.0	198 / 7.0	198 / 7.0
Platform (W x D, mm / inch)	755 x 520 / 29.7 x 20.5	755 x 520 / 29.7 x 20.5	755 x 520 / 29.7 x 20.5	755 x 520 / 29.7 x 20.5
Interior (W x D x H, mm / inch)	809 x 624 x 393 / 31.9 x 24.6 x 15.5	809 x 624 x 393 / 31.9 x 24.6 x 15.5	809 x 624 x 393 / 31.9 x 24.6 x 15.5	809 x 624 x 393 / 31.9 x 24.6 x 15.5
Exterior (W x D x H, mm / inch) ³⁾	1 unit ; 1307 x 996 x 755 / 51.5 x 39.2 x 29.7 1 unit+stand-I or 2 units stacked ; 1307 x 996 x 1355 / 51.5 x 39.2 x 53.3 2 units+stand-II stacked ; 1307 x 996 x 1655 / 51.5 x 39.2 x 65.2 3 units stacked ; 1307 x 996 x 1955 / 51.5 x 39.2 x 77.0			
Net weight (Kg / lbs)	300 / 661.4	320 / 705.5	300 / 661.4	320 / 705.5
Electrical requirements (230V, 50/60Hz, A)	3.7	5.6	3.7	5.6
Cat. No.	AAH238132K	AAH238232K	AAH239132K	AAH239232K
Electrical requirements (120V, 60Hz, A)	7.1	10	7.1	10
Cat. No.	AAH238115U	AAH238215U	AAH239115U	AAH239215U

1) The temperature control range of the ISS-7100 and ISS-7200 models are from 10°C above the ambient temperature of the equipment. (15°C to 40°C) to 80°C. The temperature control range of the ISS-7100R and ISS-7200R models are from 20°C below the ambient temperature (18°C to 35°C) of the device to 80°C. (However, the lowest control temperature is 4°C)

2) Actual rpm allowable error for set rpm is 1rpm if setting is less than 100rpm, ±1% if setting is higher than 100rpm.

3) Product height: 600mm, height of Stand I: 600mm, height of Stand II: 300mm, height of wheel attachment: 155mm.

Accessories Page 120

Universal Platform, Flask Clamp, Funnel Clamp, Test Tube Rack, Microplate Holder, Lab Sticker, Spring Wire Rack, Dedicated Platform, Rubber Mat, Universal Attachment, Stand, Stacking Kit, LC GreenBox, Starter Kit

Incubated Shaker Tabletop type

Provides convenient internal observation environment

Structural Functional Features

- Model optimized for use on experiment bench.
- Up to 80°C, 500rpm enables setting of various experimental conditions.
- Three operating modes: Temperature & Shaking, Temperature, or Shaking.
- Door structure allows for easy sample removal/attach at the front.
- Shaking structure with minimized vibration.
- Four models according to refrigeration system and the size of the platform.

Use Convenience Features

- Accurate rpm and temperature control through temperature auto-tuning / microprocessor PID and calibration.
- Monitor via mobile app anytime, anywhere with LC Connected. (Requires LC GreenBox, option)
- Countdown timer (999 hours 59 minutes) of shaking operation.
- Intuitive operation with color touch display.
- High brightness / low-power LED illumination makes it easy to observe inside.
- Spill drain port for internal cleanliness and easy maintenance.

Outstanding Safety

- Over temperature limit function.
- Soft start and stop shaking function.
- Malfunction prevented by controller lock function.
- Displays all information such as device operation and error occurrence.
- Overcurrent protection device restricts operation when objects are caught inside the system and shaking is not possible.



IST series

with Universal Platform, Flask Clamp (option)



Color touch display



Smart Self-diagnosis

Real-time self-diagnostics to check equipment for abnormal conditions.

- Alarm is generated by monitoring difference between actual rpm/temperature and set rpm/temperature.
- Automatic restart as it was operated after power recovery.



Specification

Model	IST-3075	IST-3075R	IST-4075	IST-4075R
Temperature				
Range (°C / °F) ¹⁾	Amb. +5 to 80 / Amb. +9 to 176	Amb. -15 (Min. 15) to 80 / Amb. -27 (Min. 27) to 176	Amb. +5 to 80 / Amb. +9 to 176	Amb. -15 (Min. 15) to 80 / Amb. -27 (Min. 27) to 176
Fluctuation at 37°C in flask (±°C / °F)	0.1 / 0.18	0.1 / 0.18	0.1 / 0.18	0.1 / 0.18
Variation at 37°C in flask (±°C / °F)	0.5 / 0.90	0.5 / 0.90	0.5 / 0.90	0.5 / 0.90
Refrigerator (Hp)	-	1/6	-	1/6
Shaking system				
Motion type	Orbital	Orbital	Orbital	Orbital
Speed range (rpm)	20 to 500	20 to 500	20 to 500	20 to 500
Max. speed per permissible load (rpm)	500 (max. 10 kg / 22.0 lbs) 400 (max. 15 kg / 33.1 lbs)	500 (max. 10 kg / 22.0 lbs) 400 (max. 15 kg / 33.1 lbs)	500 (max. 14 kg / 30.9 lbs) 400 (max. 21 kg / 46.3 lbs)	500 (max. 14 kg / 30.9 lbs) 400 (max. 21 kg / 46.3 lbs)
Amplitude size (mm / inch, dia.)	19.1 / 0.75	19.1 / 0.75	19.1 / 0.75	19.1 / 0.75
Speed stability ²⁾ (rpm, %)	±1	±1	±1	±1
Timer (count-down)	999 hr 59 min.			
Dimension				
Chamber volume (L / cu ft)	53 / 1.9	53 / 1.9	83 / 2.9	83 / 2.9
Platform (W x D, mm / inch)	350 x 350 / 13.8 x 13.8	350 x 350 / 13.8 x 13.8	450 x 450 / 17.7 x 17.7	450 x 450 / 17.7 x 17.7
Interior (W x D x H, mm / inch)	410 x 410 x 320 / 16.1 x 16.1 x 12.6	410 x 410 x 320 / 16.1 x 16.1 x 12.6	510 x 510 x 320 / 20.1 x 20.1 x 12.6	510 x 510 x 320 / 20.1 x 20.1 x 12.6
Exterior (W x D x H, mm / inch)	440 x 785 x 510 / 17.3 x 30.9 x 20.1	440 x 785 x 510 / 17.3 x 30.9 x 20.1	540 x 890 x 510 / 21.3 x 35 x 20.1	540 x 890 x 510 / 21.3 x 35 x 20.1
Net weight (Kg / lbs)	65 / 143.3	73 / 160.9	81 / 178.6	90 / 198.4
Electrical requirements (230V, 50/60Hz, A)				
Cat. No.	AAH23422K	AAH23522K	AAH23622K	AAH23722K
Electrical requirements (120V, 60Hz, A)				
Cat. No.	AAH23465U	AAH23565U	AAH23665U	AAH23765U

1) The temperature control range of the IST-3075 and IST-4075 models are from 5°C above the ambient temperature of the equipment. (15°C to 40°C) to 80°C The temperature control range of the IST-3075R and IST-4075R models are from 15°C below the ambient temperature (18°C to 35°C) of the device to 80°C. (However, the lowest control temperature is 15°C)

2) Actual rpm allowable error for set rpm is ±1rpm if setting is less than 100rpm, ±1% if setting is higher than 100rpm.

Accessories Page 120

Universal Platform, Flask Clamp, Funnel Clamp, Test Tube Rack, Microplate Holder, Lab Sticker, Spring Wire Rack, Dedicated Platform, Rubber Mat, Universal Attachment, LC GreenBox, Starter Kit

Incubated Shaker Floor type

Optimal solution for mass culturing

Structural Functional Features

- Large-capacity incubated shaker model.
- Up to 80°C, 500rpm enables setting of various experimental conditions.
- It is convenient for attaching and detaching a large amount of sample because structure of door is lowering from the rear to the front. (registration KR 30-0837007)
- Three operating modes: Temperature & Shaking, Temperature, or Shaking.
- Easy to move/fix the product.
- Stable shaking of high-load samples. (registration KR 10-1678160)
- Applied a direct drive motor that has outstanding durability, a 5-year warranty for the motor.

Use Convenience Features

- Accurate rpm and temperature control through temperature auto-tuning / microprocessor PID and calibration.
- Monitor via mobile app anytime, anywhere with LC Connected. (Requires LC GreenBox, option)
- Countdown timer (999 hours 59 minutes) of shaking operation.
- Spill drain port for internal cleanliness and easy maintenance.


Outstanding Safety

- Electronic over temperature protection system. (A-OT, advanced over temperature limiter)
Independent precision temperature sensor and controller provides more secure and reliable over temperature protection.
- Conventional over-temperature protection system. (B-OT, backup over temperature limiter)
A backup device against electronic system failure, mechanically preventing over temperature even when there are electronic errors.
- Soft start and stop shaking function.
- Malfunction prevented by controller lock function.
- Displays all information such as device operation and error occurrence.
- Overcurrent protection device restricts operation when objects are caught inside the system and shaking is not possible.



ISF series

with Optional Accessories



Easy to attach and detach large volume sample thanks to the inclined lid structure. (Registered design 30-2014-0060254)



Color touch display



Smart Self-diagnosis

Real-time self-diagnostics to check equipment for abnormal conditions.

- Alarm is generated by monitoring difference between actual rpm/temperature and set rpm/temperature.
- Automatic restart as it was operated after power recovery.



DD Shaking System Patented

Powerful DD motor for large-capacity shaking. The triple cam system provides stable shaking with high load samples.

(Patented in Korea)

Specification

Model	ISF-7100	ISF-7100R	ISF-7200	ISF-7200R
Temperature				
Range (°C / °F) ¹⁾	Amb. +5 to 80 / Amb. +9 to 176	Amb. -20 (Min. 4) to 80 / Amb. -36 (Min. 7.2) to 176	Amb. +5 to 80 / Amb. +9 to 176	Amb. -20 (Min. 4) to 80 / Amb. -36 (Min. 7.2) to 176
Fluctuation at 37°C in flask (±°C / °F)	0.1 / 0.18	0.1 / 0.18	0.1 / 0.18	0.1 / 0.18
Variation at 37°C in flask (±°C / °F)	0.5 / 0.90	0.5 / 0.90	0.5 / 0.90	0.5 / 0.90
Refrigerator (Hp)	-	1/6	-	1/6
Shaking system				
Motion type	Orbital	Orbital	Orbital	Orbital
Speed range (rpm)	30 to 500	30 to 500	30 to 300	30 to 300
Max. speed per permissible load (rpm)	500 (max. 15 kg / 33.1 lbs) 400 (max. 25 kg / 55.1 lbs) 350 (max. 35 kg / 77.2 lbs)	500 (max. 15 kg / 33.1 lbs) 400 (max. 25 kg / 55.1 lbs) 350 (max. 35 kg / 77.2 lbs)	300 (max. 20 kg / 44.1 lbs) 200 (max. 35 kg / 77.2 lbs)	300 (max. 20 kg / 44.1 lbs) 200 (max. 35 kg / 77.2 lbs)
Amplitude size (mm / inch, dia.)	25.4 / 1	25.4 / 1	50.8 / 2	50.8 / 2
Speed stability ²⁾ (rpm, %)	±1	±1	±1	±1
Timer (count-down)	999 hr 59 min.	999 hr 59 min.	999 hr 59 min.	999 hr 59 min.
Dimension				
Chamber volume (L / cu ft)	270 / 9.5	270 / 9.5	270 / 9.5	270 / 9.5
Platform (W x D, mm / inch)	755 x 520 / 29.7 x 20.5	755 x 520 / 29.7 x 20.5	755 x 520 / 29.7 x 20.5	755 x 520 / 29.7 x 20.5
Interior (W x D x H, mm / inch)	894 x 634 x 480 / 35.2 x 25 x 18.9	894 x 634 x 480 / 35.2 x 25 x 18.9	894 x 634 x 480 / 35.2 x 25 x 18.9	894 x 634 x 480 / 35.2 x 25 x 18.9
Exterior (W x D x H, mm / inch)	1128 x 854 x 1035 / 44.4 x 33.6 x 40.8	1128 x 854 x 1035 / 44.4 x 33.6 x 40.8	1128 x 854 x 1035 / 44.4 x 33.6 x 40.8	1128 x 854 x 1035 / 44.4 x 33.6 x 40.8
Net weight (Kg / lbs)	277 / 610.7	287 / 632.7	277 / 610.7	287 / 632.7
Electrical requirements (230V, 50/60Hz, A)				
Cat. No.	AAH23342K	AAH23541K	AAH23442K	AAH23641K
Electrical requirements (120V, 60Hz, A)				
Cat. No.	AAH23375U	AAH23575U	AAH23475U	AAH23675U

1) The temperature control range of the ISF-7100 and ISF-7200 models are from 5°C above the ambient temperature of the equipment. (15°C to 40°C) to 80°C The temperature control range of the ISF-7100R and ISF-7200R models are from 20°C below the ambient temperature (18°C to 35°C) of the device to 80°C. (However, the lowest control temperature is 4°C)

2) Actual rpm allowable error for set rpm is ±1rpm if setting is less than 100rpm, ±1% if setting is higher than 100rpm.

Accessories Page 120

Universal Platform, Flask Clamp, Funnel Clamp, Test Tube Rack, Microplate Holder, Lab Sticker, Spring Wire Rack, Dedicated Platform, Rubber Mat, Universal Attachment, LC GreenBox, Starter Kit

Accessories

Universal Platform



Universal Platform
Basic accessories required to use Erlenmeyer flasks, test tubes, separated funnels, etc.

Accessories for Universal Platform



Flask Clamp (SUS / Plastic)
Erlenmeyer flask mounting (50ml – 6L)
High-strength spring maintains excellent clamping.



Starter Kit
Consisting of plastic clamps that are commonly used on Universal Platforms. 50ml, 100ml, 250ml, 500ml, 1L, 2L.



Funnel Clamp
Separated funnel mounting (250 ml - 2L)
Easy detachment/attachment and stable clamping.



Test Tube Rack
Test tube mounting (Φ8 mm - Φ35 mm)
Easily adjustable angle.



Microplate Holder
Various options for Microplate Single, Tower, Flat A (Large), Flat B (Small)



Lab Sticker
Excellent adhesion as size is 200x200mm, allowing for application to petri dishes, culture bottles, etc.
Temp. Range (°C): 15 to 60 Max. Speed (rpm): 250

LC GreenBox



LC GreenBox
Monitor the operation status of the device in real time via the mobile app with a simple internet connection.

Set Accessories



Spring Wire Rack
Fixed with 2-stage grid springs to accommodate various sample container shapes.




Dedicated Platform
Accessories (250 ml, 500 ml) installed on platform with maximum quantity of flask clamps.



Rubber Mat
Convenient for use with culture bottles and petri dishes with rubber mat to prevent slipping on the bottom.



Universal Attachment
Rubber mat on the bottom and adjustable clamping rod make it a convenient structure for a wide range of applications.



Wire Shelf
Shelf and its bracket enabling sample loading.
(applied to ISS-3075(R)/4075(R) series)



Retractable Foot Caster
Device that makes it easy to move/fix.
(applied to ISS-3075(R)/4075(R) series)



Stand 1 & 2
When using 1 unit of the ISS-7000 model or when using more than 2 units stackable, different height of stands provide convenient access to a platform.



Stacking Kit
Fixing tool for stacked units.
Applied to ISS series models.

Description	IST-3075 Series	IST-4075 Series	ISS-3075 Series	ISS-4075 Series	ISF-7000 Series	ISS-7000 Series
Dedicated Platform (250 ml)	AAA23611-V1	AAA23623-V1	AAA23611-V1	AAA23621-V1	AAA23634	AAA23637
Dedicated Platform (500 ml)	AAA23612-V1	AAA23622-V1	AAA23612-V1	AAA23622-V1	AAA23635	AAA23636
Rubber Mat	AAA23531	AAA23532	AAA23531-V1	AAA23532-V1	AAA23534	AAA23537-V1
Universal Attachment	AAA23511-V1	AAA23512-V1	AAA23511-V1	AAA23512-V1	AAA23505	AAA23517-V1
Lab Sticker	AAA30551	AAA30551	AAA30551	AAA30551	AAA30551	AAA30551
Wire Shelf	-	-	00EDA0008235	00EDA0008233	-	-
Stand 1 with Stacking kit (600 mm height for stacking 1 unit)	-	-	-	-	-	AAA23571
Stand 2 with Stacking kit (300 mm height for stacking 2 unit)	-	-	-	-	-	AAA23572
Retractable Foot Caster	-	-	AAA23801	AAA23802	-	-
Stacking Kit	-	-	AAA23811	AAA23812	-	AAA23573

LC GreenBox

Description	Dimension (W x D x H, mm/inch)	IST-3075 Series	IST-4075 Series	ISS-3075 Series	ISS-4075 Series	ISF-7000 Series	ISS-7000 Series
LC GreenBox	156 x 94 x 34 / 6.1 x 3.7 x 1.3	AAHQ1011K					

Description		Model	IST-3075 Series	IST-4075 Series	ISS-3075 Series	ISS-4075 Series	ISF-7000 Series	ISS-7000 Series
		Platform Size (mm/inch)	350 x 350 / 13.8 x 13.8	450 x 450 / 17.7 x 17.7	350 x 350 / 13.8 x 13.8	450 x 450 / 17.7 x 17.7	755 x 520 / 29.7 x 20.5	755 x 520 / 29.7 x 20.5
Universal Platform		Cat. No.	AAA23501-V3	AAA23502-V3	AAA23501-V3	AAA23502-V3	AAA23503-V3	AAA23507-V1
Flask Clamp for Erlenmeyer Flask (Stainless Steel)	50 ml	AAA23550	36	49	36	49	88	88
	100 ml	AAA23551	24	39	24	39	68	68
	250 ml	AAA23552	14	20	14	24	39	39
	300 ml	AAA23556	14	20	14	24	39	39
	500 ml	AAA23553	9	16	9	16	28	28
	1 L	AAA23554	6	9	6	9	18	18
	2 L	AAA23555	4	5	4	6	11	11
	2.8 L	AAA23557	2	4	2	5	8	8
	4 L	AAA23558	-	-	-	4	6	6
	6 L	AAA23559	-	-	-	2	5	-
Flask Clamp for Erlenmeyer Flask (Plastic)¹⁾	50 ml	AAA30570	36	49	36	49	88	88
	100~125 ml	AAA30571	19	29	19	29	48	48
	200 ml	AAA30572	13	21	13	21	35	35
	250 ml	AAA30573	13	19	13	19	35	34
	300 ml	AAA30574	12	19	12	19	35	34
	500 ml ²⁾	AAA30575	9	12	9	12	24	21
	1 L ²⁾	AAA30576	4	9	5	9	12	11
	2 L ²⁾	AAA30577	2	4	4	5	8	6
Funnel Clamp	250 ml	AAA23562	4	6	4	6	15	15
	500 ml	AAA23563	2	4	3	4	10	10
	1 L	AAA23564	-	3	-	3	7	7
	2 L	AAA23565	-	-	-	-	5	5
Test Tube Rack³⁾	Ø8 mm (86ea)	AAA23581	2	4	2	4	7	7
	Ø10 mm (86ea)	AAA23582	2	4	2	4	7	7
	Ø12 mm (58ea)	AAA23583	2	4	2	4	7	7
	Ø14 mm (58ea)	AAA23584	2	4	2	4	7	7
	Ø17 mm (32ea)	AAA23585	2	4	2	4	7	7
	Ø25 mm (19ea)	AAA23586	2	4	2	4	7	7
	Ø30 mm (10ea)	AAA23587	2	4	2	4	7	7
	Ø35 mm (10ea)	AAA23594	2	4	2	4	7	7
Microplate Holder	Single	AAA23654	8	15	8	15	25	25
	Tower	AAA23651	5	8	6	10	21	21
	Flat A (Large)	AAA23652	-	2	-	2	4	4
	Flat B (Small)	AAA23653	2	2	2	3	6	6

1) In the case of plastic flask clamps mounting, deformation may occur at temperatures over 40°C. Stainless steel flask clamps are recommended when using over 40°C and high temperature environment.

2) Possible to use under 300rpm

3) Maximum quantity on each platform. (no angle condition for test tube rack), it is the maximum quantity that does not cause the rack to tilt.

Description	IST-3075 Series	IST-4075 Series	ISS-3075 Series	ISS-4075 Series	ISF-7000 Series	ISS-7000 Series
Spring Wire Rack (Cat. No.)	AAA23521-V1	AAA23522-V2	AAA31521-V1	AAA31522-V1	AAA23514	AAA23527-V1
Quantity of Erlenmeyer Flask	50 ml	16	25	16	25	45
	100 ml	9	16	9	16	28
	250 ml	4	9	4	9	18
	300 ml	4	9	4	9	18
	500 ml	4	5	4	5	15
	1 L	2	4	2	4	10
	2 L	1	2	1	2	6
	2.8 L	1	2	1	2	4

Starter Kit

Description	Composition of Clamp	Cat. No.
Starter Kit I	50 ml (2), 100 ml (2), 250 ml (4), 500 ml (4), 1L (2), 2L (2)	AAA30581
Starter Kit II	50 ml (3), 100 ml (3), 250 ml (6), 500 ml (6)	AAA30582

Shaker



General Application

General mixing, drug analysis, ELISA, enzyme immunoassay, dyeing, and bleaching of samples, protein synthesis.

	Description	Motion Type	Max. Speed (rpm)	Max. Platform Size (mm)	250 ml Flask Capacity	Model	Page
General type	Dual Action <ul style="list-style-type: none"> Orbital/Reciprocal motion selection Provides powerful shaking 	Orbital / Reciprocal	500	320 x 260	8	OS-2000	126
			500	350 x 350	14	OS-3000	126
			500	450 x 450	24	OS-4000	126
	General <ul style="list-style-type: none"> Stable shaking even with high load samples 	Orbital	500	755 x 520	39	OS-7100	126
			300	755 x 520	39	OS-7200	126
	Program <ul style="list-style-type: none"> Program control mode Mounts the large capacity of samples 	Orbital	500	885 x 520	102	SKC	128
Compact <ul style="list-style-type: none"> Easy to use as a small shaker Can be used inside of chamber 	Orbital	300	250 x 250	9	CMS	130	
Type by usages	Microplate <ul style="list-style-type: none"> Suitable for sample mixing for microplates, tubes, etc. 	Orbital	1200	269 x 185	N/A	CPS	131
	Rocking <ul style="list-style-type: none"> Easy to adjust the angle of motion Up to 99 hours timer function 	Rocking	100	276 x 276	N/A	CRS	132
	Waving <ul style="list-style-type: none"> Easy to adjust the angle of motion Can be used inside of chamber 	Waving	100	276 x 276	N/A	CWS	132
	Funnel <ul style="list-style-type: none"> Powerful up and down reciprocal motion Tilt adjustment. (0 to 20°) 	Reciprocal	300	Funnel 500 ml x 10 ea		RS-2	133
	Rotating <ul style="list-style-type: none"> Smooth and stable operation Easy attachment and detachment of accessories 	Rotating	70	1.5 ml Tube x 120ea (When double disks are stacked)		CRT	134

※ The contents of the above and the contents of this catalog may differ depending on the specific model and conditions of use. For the information about the features and specifications that applying to each models, please check the information on the corresponding page of each models.

 Satisfaction with the requirements of European regulations and guidelines.	 Certification of conformity of electricity, gas, etc. in accordance with the safety standards of Canada and USA.	 Real-time equipment monitoring and control system using mobile app.
 Registered patent based on Jeio Tech's proprietary technology.	 Recognized as an excellent design product by the Ministry of Commerce, Industry, and Energy.	 Compliance with Guidelines for Restriction of Hazardous Substances Use.
 Possible to set end time or start time of device operation.	 PC communication via RS-232 / RS-485 / USB port.	 Selectable shaking motion either orbital or reciprocal motion.
 Orbital motion.	 Rocking motion.	 Waving motion.
 Vertical reciprocal motion.	 Reciprocal motion in inclined up to 20°.	 Mixing through 360 degree rotation.
 It memorizes the state of power failure and operates automatically when power is restored.	 5 year warranty for direct drive motor.	 2 year warranty free A/S.



Shaking speed of 500rpm with the world highest levels Powerful Shaking & Flexible Amplitude

- > **Primarily 500rpm shaking speed in Korea**
Reliable shaking operation at high speed by allowing stable 500rpm control. (OS, SKC)
- > **Precise control technology allows for stable control**
Precise speed control in all shaking areas by adopting control system developed using Jeio Tech's own technology.
- > **Specialized shaking mechanism (Jeio Tech's patent)**
 - BLDC motor and belt system applied to small model enables compact and strong shaking ability.
 - The medium and large class model DD motor and triple cam system allow for stable shaking of high load samples.
- > **Increased efficiency with diverse shaking radius**
Four kinds of shaking radius (amplitude) are applied to provide a suitable shaking environment.
- > **Orbital/reciprocal shaking with dual action**
Rotational shaking and reciprocal shaking are selectively used in one shaker. (OS-2000, 3000, 4000)
- > **Intuitive and convenient operation**
 - Large clear VFD panel provides intuitive operation information for all devices.
 - Smooth and comfortable control knob for easy and convenient operation.

※ Some of the above contents are limited to specific models.



Intuitive checking of device operation information.



USB & RS-232 Support. (OS)



Various accessories provided according to the sample.



Provides functionality that is practically necessary for users Jeio Tech's Shaker with useful aggregated functions

- > **Convenient setting by changing rpm in real time**
Set the appropriate rpm while checking the actual speed of sample shaking.
- > **DD motor of outstanding durability**
Applied a direct drive motor that has outstanding durability, a 5-year warranty for the motor. (OS-7000, SKC)
- > **Remote monitoring service (only OS)**
Monitor and control equipment remotely through a mobile application anytime and anywhere. (option) (OS)
- > **Computer interface environment**
USB, RS-232 connection and dedicated software provide convenient device operation and data accumulation/management function. (OS)
- > **Count-down timer for accurate operation**
Count-down timer that can be set up to 999 hours and 59 minutes for convenient use.
- > **Prepared for unexpected power shutdown**
Auto-run function that automatically restarts to the normal state immediately in the event of sudden power failure. (OS)
- > **Real-time monitoring of rpm error**
Prevents experiment errors through user notification when the actual rpm differs from the set rpm. (OS)
- > **Prevention of malfunction with motor self-protection function**
Motor operation is controlled if overcurrent occurs in the circuit configuration with the motor self-protection function.





Heavy Duty Holder (Option) for large capacity shaking

- Shaking up to 50L carboy on OS-7100 / 7200 models.
- Select 20L / 50L Holder Sponge depending on the amount of shaking.

Controllable rpm according to water weight

Model	Holder sponge	Water weight (L)					
		0	10	20	30	40	50
OS-7100 (rpm)	20L x 2	300	280	180	160	160	-
	50L x 1	300	190	150	130	130	130
OS-7200 (rpm)	20L x 2	210	180	160	150	150	-
	50L x 1	210	140	120	110	100	100

Cat.No	Description
AAA3D501	Heavy duty holder for OS-7000 series
AAA3D502	Heavy duty holder sponge (Φ376mm x 1hole) for 50L Carboy
AAA3D503	Heavy duty holder sponge (Φ290mm x 2hole) for 20L Carboy
AAAJ3031	20L Carboy (260 x 317 x 539mm)



OS-7100

with Heavy duty holder (option) + 1 hole-Heavy duty holder sponge (option)



OS-7100

with Heavy duty holder (option) + 2-hole Heavy duty holder sponge (option)

Specification

Model	OS-2000	OS-3000	OS-4000	OS-7100	OS-7200
Shaking system					
Motion type	Dual Action (Orbital/Reciprocal)	Dual Action (Orbital/Reciprocal)	Dual Action (Orbital/Reciprocal)	Orbital	Orbital
Speed range (rpm)	20 to 500	20 to 500	20 to 500 (Orbital) / 20 to 450 (Reciprocal)	30 to 500	30 to 300
Max. speed per permissible load ¹⁾	500rpm (max. 10 kg) 300rpm (max. 16 kg)	500rpm (max. 10 kg) 300rpm (max. 18 kg)	Orbital: 500 (max. 14kg), 400 (max. 21 kg) Recipro.: 450 (max. 10kg), 300 (max. 21 kg)	500rpm (max. 10 kg) 400rpm (max. 30 kg)	300rpm (max. 23 kg) 250rpm (max. 35 kg)
Amplitude size (mm / inch, dia.)	19.1 / 0.8 (default) (12.7 / 0.5 available)	19.1 / 0.8 (default) (25.4 / 1 available)	19.1 / 0.75 (default) 25.4 / 1 (available)	25.4 / 1	50.8 / 2
Speed stability ²⁾ (rpm, %)	±1	±1	±1	±1	±1
Timer (count-down)	999hr 59min	999hr 59min	999hr 59min	999hr 59min	999hr 59min
Dimension					
Platform (W x D, mm / inch)	320 x 260 / 12.6 x 10.2	350 x 350 / 13.8 x 13.8	450 x 450 / 17.7 x 17.7	755 x 520 / 29.6 x 20.5	
Body (W x D x H, mm/inch)	354 x 388 x 146 / 13.9 x 15.3 x 5.7	409 x 508 x 141 / 16.1 x 20 x 5.6	510 x 632 x 156 / 20.1 x 24.9 x 6.1	755 x 627 x 151 / 29.7 x 24.7 x 5.9	
Net weight (Kg / lbs)	22 / 48.5	28 / 61.7	43 / 94.8	117 / 257.9	
Electrical requirements (230V, 50/60Hz, A)					
	0.3	0.3	0.3	0.4	
Cat. No.	AAH3D1115K	AAH3D1215K	AAH3D1315K	AAH3D1415K	AAH3D1515K
Electrical requirements (120V, 60Hz, A)					
	0.6	0.6	0.6	0.8	
Cat. No.	AAH3D1165U	AAH3D1265U	AAH3D1365U	AAH3D1465U	AAH3D1565U

1) In Orbital Shaking mode, Regarding the maximum controllable rpm and amplitude by weight, the measurement is performed based on the default value.
2) Actual rpm allowable error for set rpm (1rpm if setting is less than 100rpm, ±1% if setting is higher than 100rpm)

Accessories Page 136

Universal Platform, Spring Wire Rack, Dedicated Platforms, Rubber Mat, Universal Attachment, Lab Sticker, Microplate Tray, Stage Rack, Flask Clamp, Funnel Clamp, Test Tube Rack, Microplate Holder, LC Green Box, Starter Kit

Shaker

Program type

Program operation and maximum sample loading

▣ Structural Functional Features

- Standard 350rpm, Maximum 500rpm. (option)
- Application of DD motor & triple cam system stable shaking ability. (registration KR 10-1678160)
- Provides the widest platform size model.
- Up to 102 pieces of 250ml flasks can be mounted at once. (when using dual platform.)
- Dual platform mountable. (only flasks less than 250ml can be mounted on the bottom part)
- Provided Maximum 50.8mm shaking radius model.
- Convenient platform detachment system.
- Stable shaking with low center of gravity structure.
- Corrosion-resistant outer surface with powder coating and PP material.
- Structure prevents the sample from flowing into the inside.
- Model awarded Good Design certification.
- Applied a direct drive motor that has outstanding durability, a 5-year warranty for the motor.




SKC-7000 Series

with Universal Platform, Flask Clamp (option)

▣ Use Convenience Features

- Accurate shaking with PID feedback control.
- Provides Quick, Timer, and Program modes.
- Wait On/Off timer. (up to 99 hours 59 minutes)
- More reliable shaking simulation by clockwise / counterclockwise rotation.
- Stable control system even at low-speed agitation.
- Soft start and stop shaking function.
- Easy to use with fast acceleration/deceleration function.
- Possible to stop the equipment at the same position.
- Automatic rpm adjustment in case of instability.
- Equipment stoppage and alarm in case of abnormal shock.



- Displays sets rpm and actual rpm simultaneously.
- Provides all information on equipment operation.
- Convenient operation with shaking start/stop button.
- Convenient operation by changing rpm with control knob.
- Clear display. (clear VFD)




Outstanding Safety

- Alarm occurs when the current rpm differs from the set rpm.
- Over-current protection.

Differentiated Program Operation

- 20 patterns can be saved.
- 20 segments can be set per pattern.
- Time, rpm, and orbital direction can be set for each segment.
- Includes up to 99 repeat operations.



Specification

Model	SKC-7075	SKC-7100	SKC-7200
Shaking system			
Motion type	Orbital	Orbital	Orbital
Speed range (rpm) ¹⁾	30 to 500	30 to 500	30 to 500
Max. speed per permissible load ²⁾	350rpm (max. 13kg) 310rpm (max. 28kg)	340rpm (max. 13 kg) 260rpm (max. 28 kg)	250rpm (max. 13 kg) 200rpm (max. 28 kg)
Amplitude size (mm / inch, dia)	19.1 / 0.8	25.4 / 1	50.8 / 2
Speed stability ³⁾ (rpm, %)	±1	±1	±1
Timer (wait on/off)	99hr 59min	99hr 59min	99hr 59min
Dimension			
Platform (Small, W x D, mm / inch)	755 x 520 / 29.7 x 20.5		
Platform (Large, W x D, mm / inch)	885 x 520 / 34.8 x 20.5		
Exterior (W x D x H, mm / inch)	755 x 662 x 122 / 29.7 x 26.1 x 4.8		
Net weight (kg / lbs)	75 / 165.4	75 / 165.4	75 / 165.4
Electric requirements (230V, 50/60Hz, A)	0.4	0.4	0.4
Cat. No.	AAH30205K	AAH30215K	AAH30225K
Electric requirements (120V, 60Hz, A)	0.8	0.8	0.8
Cat. No.	AAH30206U	AAH30216U	AAH30226U

1) 350 to 500rpm applies when separately purchased options are mounted.

2) The maximum rpm by weight in Orbital Shaking mode.

3) Actual rpm allowable error for set rpm (1rpm if setting is less than 100rpm, ±1% if setting is higher than 100rpm)

Accessories Page 136

Universal Platform, Spring Wire Rack, Rubber Mat, Universal Attachment, Lab Sticker, Flask Clamp, Funnel Clamp, Test Tube Rack, Microplate Holder, Starter Kit

Shaker

Compact type



Compact Mini Shaker for easy use

Structural Functional Features

- Max. rpm of 300 is suitable for life science experiments sensitive to low rpm such as hybridization and immunocytochemistry.
- Orbital motion with shaking radius of 13 mm.
- Stable control up to 4 kg.
- It can be used mounted nine 250ml flasks.
- Count-down timer. (up to 99 hours and 59 minutes)
- Displays the elapsed time.
- Various accessories are available.

Use Convenience Features

- It can be used in a temperature chamber such as cold room or incubator because it can be used at ambient temperatures between 2°C and 60°C.
- Small size and light weight provide good mobility and usability in various laboratory environments and experiment requirements.
- Accurate shaking control according to sample weight with PID feedback control.
- Easy to operate with clear VFD display, touch keys and control knobs.
- Displays the set rpm and actual rpm simultaneously, and displays timer operation information.



CMS-350

with Universal Platform, Flask Clamp (option)

Outstanding Safety

- Self-checks for errors.
- Body made from PP material with excellent chemical resistance.
- Protected from overcurrent.

Specification

Model	CMS-350	
Shaking system		
Motion type	Orbital	
Speed range (rpm)	50 to 300	
Rpm control at 100rpm (\pm rpm)	5	
Orbit size (mm / inch, dia)	13	
Timer	Run time : 1min ~ 99hr 59min	
Dimension		
Platform (W x D, mm / inch)	250 x 250 / 9.8 x 9.8	
Body ¹⁾ (W x D x H, mm / inch)	249 x 325 x 109 / 9.8 x 12.8 x 4.3	
Net weight (kg / lbs)	8 / 17.6	
Electrical requirements		
AC 100 ~ 240V, 50/60Hz, 0.2A Max.		
Cat. No.	KR Plug	AAH3B1115K
	US Plug	AAH3B1115U

※ When operating at room temperature after use at low temperature, it may be affected by condensation. We recommend that you turn on the power after leaving at room temperature for a while.

1) Main body dimension without accessories

Accessories Page 136

Universal Platform, Spring Wire Rack, Rubber Mat, Universal Attachment, Lab Sticker, Flask Clamp, Funnel Clamp, Test Tube Rack, Microplate Holder, Starter Kit

Microplate Shaker



Exclusive Shaker for microplates and tubes

Structural Functional Features

- Suitable for sample mixing with microplate, micro-tube and conical-tube at a maximum rpm of 1200.
- Orbital motion with shaking radius of 3 mm.
- Smooth shaking start control.
- Count-down timer. (up to 99 hours and 59 minutes)
- Displays the elapsed time.
- Various tube rack accessories.

Use Convenience Features

- It can be used in a temperature chamber such as cold room or incubator because it can be used at ambient temperatures between 2°C and 60°C.
- Small size and light weight provide good mobility and usability in various laboratory environments and experiment requirements.
- Accurate shaking control according to sample weight with PID feedback control.
- Easy to operate with clear VFD display, touch keys and control knobs.
- Displays the set rpm and actual rpm simultaneously, and displays timer operation information.



CPS-350

Outstanding Safety

- Self-checks for errors.
- Body made from PP material with excellent chemical resistance.
- Protected from overcurrent.

Specification

Model	CPS-350	
Shaking system		
Motion type	Orbital	
Speed range (rpm) with 4 Microplates	150 to 1200	
Orbit size (mm / inch, dia)	3 / 0.12	
Speed control	Feedback control	
Timer	Run time : 1min ~ 99hr 59min	
Dimension		
Tray (W x D) (mm /inch)	269 x 185 / 10.6 x 7.3	
Body ¹⁾ (W x D x H) (mm /inch)	269 x 325 x 128 / 10.6 x 12.8 x 5.0	
Net weight (kg / lbs)	7.7 / 17	
Electrical requirements		
AC 100 ~ 240V, 50/60Hz, 0.2A Max.		
Cat. No.	KR Plug	AAH3C1115K
	US Plug	AAH3C1115U

※ When operating at room temperature after use at low temperature, it may be affected by condensation. We recommend that you turn on the power after leaving at room temperature for a while.

1) Including tray

Accessories



1.5 ml Micro-tube Rack



15 ml Conical Tube Rack



50 ml Conical Tube Rack

Cat. No.	Description	Max. mountable tube
AAAJ2311	1.5 ml Micro-tube rack	80 (5 x 16)
00CPS0000011	15 ml Conical tube rack	30 (3 x 10)
00CPS0000012	50 ml Conical tube rack	12 (2 x 6)

Rocking / Waving Shaker



Optimized rocking / waving shaker

Structural Functional Features

- Easily adjustable operating angle to optimize according to application.
- Provides speeds up to 100rpm.
- Provides stable control at low speeds.
- 6 micro-plates can be mounted.
- Convenient tray structure for easy mounting and removal.
- Count-down timer. (up to 99 hours and 59 minutes)
- Displays the elapsed time.

Use Convenience Features

- It can be used in a temperature chamber such as cold room or incubator because it can be used at ambient temperatures between 2°C and 60°C.
- Small size and light weight provide good mobility and usability in various laboratory environments and experiment requirements.
- Accurate shaking control according to sample weight with PID feedback control.
- Easy to operate with clear VFD display, touch keys and control knobs.
- Displays the set rpm and actual rpm simultaneously, and displays timer operation information.

Outstanding Safety

- Self-checks for errors.
- Body made from PP material with excellent chemical resistance.
- Protected from overcurrent.



CWS-350

with Dimpled Mat (standard)



Rocking
(CRS-350 Model)



Waving
(CWS-350 Model)


Specification

Model	CRS-350	CWS-350	
Shaking system			
Motion type	Rocking	Waving	
Speed range (rpm)	5 to 100	8 to 100	
Rpm control at 100rpm (± rpm)	5	5	
Adjustable tilt (°)	7 ~ 13	0 ~ 13	
Timer	Run time : 1min ~ 99hr	Run time : 1min ~ 99hr	
Dimension			
Tray (W x D, mm / inch)	276 x 276 / 10.87 x 10.87	276 x 276 / 10.87 x 10.87	
Body ¹⁾ (W x D x H, mm / inch)	276 x 334 x 164 / 10.8 x 13.1 x 6.4	276 x 334 x 164 / 10.8 x 13.1 x 6.4	
Net weight (kg / lbs)	3.7 / 8.16	3.7 / 8.16	
Electrical requirements			
	AC 100 ~ 240V, 50/60Hz, 0.1A	AC 100 ~ 240V, 50/60Hz, 0.1A	
Cat. No.	KR Plug	AAH391115K	AAH381115K
	US Plug	AAH391115U	AAH381115U

※ When operating at room temperature after use at low temperature, it may be affected by condensation. We recommend that you turn on the power after leaving at room temperature for a while.

1) Main body dimension without accessories

Accessories



Rubber Mat



Dimpled Mat



Universal Harness

Cat. No.	Description	Composition
AAA38501	Rubber mat	Tray, Sticky pad
AAA38502	Dimpled mat	Tray, Silicone pad
AAA38503	Universal harness	Tray, Sticky pad, Rubber Band (6ea)

Funnel Shaker



Powerful up-down reciprocal motion

Structural Functional Features

- Vertical reciprocation at a maximum of 300 rpm.
- 40 mm stroke length.
- Adjustable shake tilt (0~20°) to achieve high mixing efficiency at low rpm.
- Count-down timer. (up to 999 hours and 59 minutes)
- 5kg on one side, 10kg in total.
- Convenient holder height adjustment.
- Stable device structure at high speeds.

Use Convenience Features

- Digitally displays up-down reciprocal motion and it's easy to adjust speed with dial.
- Easy to detach separated funnel, maximizing convenience for repeat users.
- Adjustable height of funnel holder makes it possible to mount funnels of various sizes. (up to 2L)
- 4 holders are provided as basic, 10 holders can be mounted in maximum when it adds. (option, but the available mounted quantities are depends on the funnel size.)
- Monitor via mobile app anytime, anywhere with LC Connected. (Requires LC GreenBox, option)

Outstanding Safety

- Alarms occur when there is a deviation between set rpm and actual rpm.
- Over-current protection.
- Optimized balance design to use in stable even when samples are installed at one side or tilted reciprocal motion.
- Automatic restart in the event of sudden power loss.



RS-2

4 funnel holders are included in the product package.



Up to 5 holders may be installed on each side.



Tilt angle adjustment from 0 to 20° for maximum mixing efficiency.

Specification

Model		RS-2
Shaking system		
Motion type	Vertical & Inclined reciprocating (0 to 20°)	
Speed range (rpm)	Vertical 30 to 300 / Inclined 30 to 250	
Speed accuracy at 100rpm (± rpm)	1	
Amplitude, recip. stroke length (mm / inch)	40 / 1.6	
Timer	999hr 59 min	
Dimension		
Exterior (W x D x H, mm / inch)	600 x 460 x 470 / 23.6 x 18.1 x 18.5	
Minimum required space considering 20° tilt	900 x 460 x 670 / 35.4 x 18.1 x 26.3	
Max. load per each side (kg / lbs)	5 / 11	
Net weight (kg / lbs)	52 / 114.6	
Electrical data & Ordering No.		
100 ~ 240V, 50/60Hz, A	0.9 Max.	
Cat. No.	KR plug	AAH32025K
	US plug	AAH32023U

Maximum number of installed funnels (per side)

Funnel Capacity	Round type	Cylindrical type
Less than 500 mL*	3-5ea	3-5ea
500 mL	3ea	3ea
1L	3ea	3ea
2L	2ea	2ea

* The number of installation may vary according to the funnel sizes.

Accessories

Cat. No.	Description
AAA32502-V1	Funnel holder
AAHQ1011K	LC GreenBox

Rotator

Optimal model configuration with powerful performance

Structural Functional Features

- Easily adjustable angle up to 90°.
- Provides a wide range of rpm speeds. (2 to 70rpm)
- Provides smooth and stable control even at low speeds.
- Stable control up to 4 kg.
- Count-down timer. (up to 99 hours and 59 minutes)
- Displays the elapsed time.
- The structure of minimizing vibration and movement.

Use Convenience Features

- Possible to combine various accessories in one body.
- Easy disk removal.
- Up to 60 pieces mountable (total 120 pieces) based on 1.5ml tube per one disk when 2-stage disk loading.
- Easy to operate with clear VFD display, touch keys and control knobs.
- Displays set rpm and actual rpm simultaneously.

Outstanding Safety

- Self-checks for errors.
- Body made from PP material with excellent chemical resistance.
- Over-current protection.



CRT-350

with Disk, 1.5 ml Tube Holder (option)



Adjustable angle up to 90°



Easy angle adjustment

Specification

Model	CRT-350
Shaking system	
Speed range (rpm)	2 to 70
Angle (°)	0 to 90
Timer	Run time : 1min ~ 99hr 59min
Dimension	
Exterior with disk (option) (W x D x H) (mm / inch)	320 x 325 x 463 / 12.6 x 12.8 x 18.2
Net weight with disk (option) (kg / lbs)	7.5 / 16.5
Electrical requirements	AC 100 to 240V, 50/60Hz, 0.2A Max.
Cat. No.	
KR plug	AAH3E115K
US plug	AAH3E115U



Number of tubes used at the same time when double disks are stacked.

Tube	Quantity
1.5 ml	120
5 ml	36
15 ml	36
50 ml	24



Accessories



1.5 ml (Ø 10 mm)
Tube holder



5 ml (Ø 12 mm)
Tube holder




15 ml (Ø 16.5 mm)
Tube holder



50 ml (Ø 29 mm)
Tube holder



Disk for flask clamp



Flask clamp (plastic)



Flask clamp (stainless steel)



Post (use for double disk loading)

Tube Holder

Description	Max. Mountable Tube	Cat. No.
1.5 ml (Ø 10 mm)	60	AAA3E501
5 ml (Ø 12 mm)	18	AAA3E502
15 ml (Ø 16.5 mm)	18	AAA3E503
50 ml (Ø 29 mm)	12	AAA3E504

Erlenmeyer Flask Clamp Holder

Description	Max. Mountable Tube		Cat. No.
Disk for flask clamp	N/A		AAA3E511
Flask clamp (Stainless steel)	50 ml	13	AAA23550
	100 ml	8	AAA23551
	300 ml	6	AAA23556
Flask clamp (Plastic)	50 ml	13	AAA30570
	100 ml~125 ml	8	AAA30571
	300 ml	6	AAA30574

Post for double stacking

Description	Cat. No.
Post	AAA3E521



Accessories for Shaker

Universal Platform



Universal Platform

Basic accessories required to use separated funnels, erlenmeyer flask, test tube, etc.



Dual Universal Platform

Maximum sample loading possible with two-stage structure. (SKC-7000 series model exclusively)

Accessories for Universal Platform



Flask Clamp (STS / Plastic)

Erlenmeyer flask mounting (50ml - 6L)
High-strength spring maintains excellent clamping.



Starter Kit

Consisting of plastic clamps that are commonly used on Universal Platforms. 50ml, 100ml, 250ml, 500ml, 1L, 2L.



Funnel Clamp

Separated funnel mounting (250ml - 2L)
Easy detachment and stable clamping.



Test Tube Rack

Test tube mounting (Ø8 mm - Ø35 mm)
Easily adjustable angle.



Microplate Holder

Various configurations for microplate
Single, Tower, Flat A (Large), Flat B (Small)



Lab Sticker

Excellent adhesion at 200mm x 200mm size, applied to petri dishes, culture bottles, etc.
Temp. Range (°C): 15 to 60 Max. Speed (rpm): 250

LC Connected



LC GreenBox

Real-time monitoring of device operation via mobile app.
Up to 4 devices can be connected.

Set Accessories



Spring Wire Rack, Half Spring Wire Rack

Available to load various sample container shapes.
For SKC-7000, half type can be used that providing both a universal platform and spring wire rack.



Dedicated Platform

Accessories (250 ml, 500 ml) installed on platform with maximum quantity of flask clamps.



Rubber Mat

Convenient for use with culture bottles and petri dishes with rubber plate to prevent slipping on the floor.



Universal Attachment

Floor rubber plate and adjustable clamping rod make it a convenient structure for a wide range of applications.



Microplate Tray (Dual Stacking Tray)

Dual stacking structure for microplate allows for excellent space utilization.



Microplate Tray (Large Tray)

The tray for microplates with anti-slip rubber plate on the bottom. (OS-2000, CMS-350)



Stage Rack

The rack maximizes space efficiency by stacking up to four shakers.
(OS-3000/4000/7000)



Weight Increaser

Accessories allowing for use of SKC series shaker up to 500rpm.

Starter Kit

Description	Composition of Clamp	Cat. No.
Starter Kit I	50 ml (2), 100 ml (2), 250 ml (4), 500 ml (4), 1L (2), 2L (2)	AAA30581
Starter Kit II	50 ml (3), 100 ml (3), 250 ml (6), 500 ml (6)	AAA30582

Description	OS-2000	OS-3000	OS-4000	OS-7000 Series	SKC-7000 Series		CMS-350
Spring Wire Rack (Half)	-	-	-	-	AAA30525 ¹⁾	AAA30526 ²⁾	-
Dedicated Platforms (250 mL)	-	AAA23611-V1	AAA23621-V1	AAA23634	-	-	-
Dedicated Platforms (500 mL)	-	AAA23613-V1	AAA23622-V1	AAA23635	-	-	-
Rubber Mat	AAA3A531-V1	AAA23531-V1	AAA23532-V1	AAA23534	AAA30533 ¹⁾	AAA30534 ²⁾	AAA38501
Universal Attachment	AAA3A511-V1	AAA23511-V1	AAA23512-V1	AAA23505	AAA30513 ¹⁾	AAA30514 ²⁾	AAA3B521
Lab Sticker	AAA30551	AAA30551	AAA30551	AAA30551	AAA30551	AAA30551	AAA30551
Microplate Tray	Dual Stacking	AAA3A541	-	-	-	-	AAA3B531
	Large Tray	AAA3A542	-	-	-	-	AAA3B532
Stage Rack	4 Stages	-	AAA31541-V1	AAA31542-V1	AAA31702-V1	-	-
	3 Stages	-	-	-	AAA31703-V1	-	-
Weight Increaser	-	-	-	-	AAA30561	AAA30561	-

1) Small Size
2) Large Size

Description	Model	OS-2000	OS-3000	OS-4000	OS-7000 Series	SKC-7000 series		CMS-350	
	Platform Size (mm / inch)	320 x 260 / 13.8 x 10.2	350 x 350 / 13.8 x 13.8	450 x 450 / 17.7 x 17.7	755 x 520 / 29.7 x 20.5	755 x 520 / 29.7 x 20.5 (Small)	885 x 520 / 34.8 x 20.5 (Large)	276 x 276 / 10.9 x 10.9	
Universal Platform	Cat. No.	AAA3A501-V2	AAA23501-V3	AAA23502-V3	AAA23503-V3	AAA30502¹⁾ AAA30541²⁾	AAA30504¹⁾ AAA30542²⁾	AAA3B501	
Flask Clamp for Erlenmeyer Flask (Stainless steel)⁴⁾	50 mL	AAA23550	20	36	49	88	108 (216) ³⁾	126 (252) ³⁾	16
	100 mL	AAA23551	16	28	39	68	72 (144) ³⁾	84 (168) ³⁾	9
	250 mL	AAA23552	8	14	24	39	42 (84) ³⁾	51 (102) ³⁾	9
	300 mL	AAA23556	8	14	24	39	36	40	9
	500 mL	AAA23553	6	10	16	28	27	30	5
	1 L	AAA23554	3	6	9	18	18	20	4
	2 L	AAA23555	-	4	6	11	12	13	1
	2.8 L	AAA23557	-	3	5	8	9	11	-
	4 L	AAA23558	-	1	4	6	8	9	-
6 L	AAA23559	-	1	3	5	6	7	-	
Flask Clamp for Erlenmeyer Flask (Plastic)⁵⁾	50 mL	AAA30570	20	36	49	88	72	80	16
	100~125 mL	AAA30571	12	19	29	48	54	60	9
	200 mL	AAA30572	8	13	21	35	30	32	9
	250 mL	AAA30573	8	13	19	35	30	32	9
	300 mL	AAA30574	6	12	19	35	30	32	9
	500 mL ⁶⁾	AAA30575	6	9	12	24	24	28	5
	1 L ⁶⁾	AAA30576	3	5	9	12	12	15	4
	2 L ⁶⁾	AAA30577	-	4	5	8	6	9	1
Funnel Clamp	250 mL	AAA23562	4	4	6	15	18	20	3
	500 mL	AAA23563	3	3	4	10	11	14	2
	1 L	AAA23564	3	3	3	7	9	12	2
	2 L	AAA23565	-	2	2	5	7	8	-
Test Tube Rack⁴⁾	Ø8 mm (86ea)	AAA23581	2	2	4	7	9	12	2
	Ø10 mm (86ea)	AAA23582	2	2	4	7	9	12	2
	Ø12 mm (58ea)	AAA23583	2	2	4	7	9	12	2
	Ø14 mm (58ea)	AAA23584	2	2	4	7	9	12	2
	Ø17 mm (32ea)	AAA23585	2	2	4	7	9	12	2
	Ø25 mm (19ea)	AAA23586	2	2	4	7	9	12	2
	Ø30 mm (10ea)	AAA23587	2	2	4	7	9	12	2
	Ø35 mm (10ea)	AAA23594	2	2	4	7	9	12	2
Microplate Holder⁴⁾	Single	AAA23654	6	8	15	25	25	28	6
	Tower	AAA23651	4	6	10	21	22	25	6
	Flat A (Large)	AAA23652	-	-	2	4	4	5	-
	Flat B (Small)	AAA23653	1	2	3	6	6	7	-

1) Single platform installed in SKC-7000

2) Dual platform installed in SKC-7000, only flasks less than 250ml can be installed on the lower platform. (Height of lower platform is about 175mm)

3) The numbers in parentheses indicate the maximum number of 250ml flasks that can be installed when Dual Platform is mounted.

4) Maximum quantity that can be mounted on each platform (For test tube rack, it is the maximum quantity when the rack is not to be tilted)

5) In the case of plastic flask clamps mounting, deformation may occur at temperatures over 40°C. Stainless steel flask clamps are recommended when using in a high temperature environment.

6) Possible to use under 300rpm

Spring Wire Rack		OS-2000	OS-3000	OS-4000	OS-7000 Series	SKC-7000 Series		CMS-350
						Small	Large	
Spring Wire Rack (Cat. No.)		AAA3A521	AAA31521-V1	AAA31522-V1	AAA23514	AAA30523	AAA30524	AAA3B511
Quantity of Erlenmeyer Flask	50 mL	15	16	25	45	45	50	15
	100 mL	9	9	16	28	28	32	9
	250 mL	4	4	9	18	18	21	4
	300 mL	4	4	9	18	18	21	4
	500 mL	4	4	5	15	15	18	4
	1 L	2	2	4	10	10	12	2
	2 L	-	1	2	6	6	7	-
2.8 L	-	1	2	4	4	5	-	

LC GreenBox

Description	OS-2000	OS-3000	OS-4000	OS-7000 Series
LC GreenBox				AAHQ1011K

Stirrer & Mixer



Satisfaction with the requirements of European regulations and guidelines.



GOOD DESIGN

Recognized as an excellent design product by the Ministry of Commerce, Industry, and Energy.



Over Temp.

Over temperature protection function.



Program

Control of temperature / humidity / rpm / time, etc. through dedicated program.



Timer

Possible to set end time or start time of device operation.



CSA US

Certification of conformity of electricity, gas, etc. in accordance with the safety standards of Canada and USA.



RoHS Compliant

Compliance with Guidelines for Restriction of Hazardous Substances Use.



Non-stop Best Effort Run

Maintains possible maximum orbital motion without stopping, even if overload occurs.



Hot Top

If the surface temperature of the hot plate exceeds 50 degrees, the warning lamp turns on.



Computer Interface

PC communication via RS-232 / RS-485 / USB port.



Patented

Registered patent based on Jeio Tech's proprietary technology.



Temperature Calibration

More precise temperature control through temperature calibration.



Touch Display

Intuitive operation with color touch display.



External Sensor

Possible to check and control temperature of sample by connecting to device with an external sensor.



2 year Warranty

2 year warranty free A/S.

General Application

Hotplate: Acid and base digestions, trace metal analysis, sample drying, general reagent heating, evaporation.

Hotplate Stirrer: Accurate, stable, and highly-reproducible in various heating stirring experiments.

Magnetic Stirrer: Crystallization, solvent evaporation, chemical reaction, titration, distillation, media dissolution, food processing.

Overhead Stirrer: Petrochemical sector, beverage/hot water product/dairy product production, food processing.

Vortex Mixer: Powerful mixing of solutions.

	Description	Max. Speed (rpm)	Max. Capacity (L / cu ft)	Max. Temp. Range (°C / °F)	Model	Page	
Hotplate & Magnetic Stirrer	Digital	• Stable and precise control based on sample temperature	2000	20 / 0.71	350 / 662	TS	142
	Analog	• Easy to use and verified excellent durability	2000	20 / 0.71	350 / 662	TM	143
	Multi Type	• Differentiated multi-stirrer with temperature control	2000	0.25 x 15ea / 0.009 x 15ea	120 / 248	MS-MH	144
Hotplate	Digital	• Precise digital temperature control and timer function	N/A	N/A	350 / 662	T	145
Magnetic Stirrer	Digital	• Accurate and smooth control of sample viscosity changes	2000	5 / 0.18	N/A	MS-G	148
	Analog	• Fast stirring speed control with real time response	2500	5 / 0.18	N/A	MS-B	149
	Multi Type	• Optimized for diverse uses with differentiated stirring controls	2000	0.25 x 15ea / 0.009 x 15ea	N/A	MS-M	150
Overhead Stirrer	Prestige Touch	• Powerful torque motor and color touch display operation	2000	100 / 3.53	N/A	MSH	158
	High Performance	• Computer connection operation and stable control function	2000	20 / 0.71	N/A	MSD	160
	Easy Control	• Optimized for easy use overheating /overload safety function	2000	20 / 0.71	N/A	MSA	161
Vortex Mixer	Touch / Continuous	• True vortex mixing in innovative touch mode	3000	N/A	N/A	VM	168

※ The contents of the above and the contents of this catalog may differ depending on the specific model and conditions of use. For the information about the features and specifications that applying to each models, please check the information on the corresponding page of each models.



Hotplate & Magnetic Stirrer

Safe Hotplate & Magnetic Stirrer Proven durability and useful safety features

> **Variety of products optimized for purpose of use**

11 models available depending on the size and shape of the top plate, control method, and whether stirring is present.

> **Corrosion-resistant hotplates with white ceramic coating**

The white upper hotplate, which makes it easy to observe samples, is ceramic coated and has excellent corrosion resistance. (except for MS-MH)

> **Top plate with excellent heat transfer and durability**

Excellent durability and heat transfer ability as the heater is embedded in the top plate made of aluminum. (except for MS-MH)

> **Safety with spill-proof design**

Improved safety with structural design that prevents solution from flowing into the equipment.

> **Multiple overheating protection devices**

Equipment breakdown is prevented in advance with overheating prevention device for top plate, and main parts of the inside of the body, etc.

> **BLDC motor with excellent durability**

Excellent durability even for long-term use with BLDC motor. (except for T-series)

> **Maintains stirring ability even at high temperatures**

Using special permanent magnets that maintain strong magnetic force even at high temperatures. (except for T-Series)

> **Safety function of displaying top plate overheating**

A warning is displayed when the top plate is over 50°C even when the power is switched off. (except for MS-MH)



Digital feedback control. (left)
Analog type scale control. (right)



Top plate with excellent heat transfer.



Top plate overheating indicator light.



Stable & Powerful Magnetic Stirrer Stable Stirring with Differentiated Controllability

> Variety of products optimized for purpose of use

11 different models provided depending on the size and shape of the top plate, control method, and color, etc.

> BLDC motor with excellent durability

Excellent durability even during long-term repeated use with BLDC motor.

> Special magnet for maintaining stirring ability

Using a special permanent magnet to maintain powerful magnetic coupling.

> Safety with spill-proof design

Improved safety with structural design that prevents solution from flowing into the equipment.

> Slip prevention for glassware

Silicone top plate cover prevents accidents caused by slipping of glassware.

> Top plate cover for easy sample observation

White and black cover provided as standard makes it very convenient to observe changes according to the sample.

> Fast stirring speed control with real time response

Quick control of stirring speed (MS-B/T model) according to adjust with the of control knobs.

> Smooth and precise speed control

Smooth stirring control up to the settings values. Accurate feedback control. (MS-G/M model)

※ Some of the above contents are limited to specific models.



Silicone cover prevents slippage of glassware.



Easy to choose from various model configurations.



Selection of white or black cover according to sample.

Hotplate & Magnetic Stirrer

Digital type



Precise control based on sample temperature

Structural Functional Features


- Includes temperature probe (B class) as standard.
- Displays external temperature sensor errors.
- Top plate coated with white ceramic for excellent corrosion resistance.
- Excellent heat transfer and durability due to heater-integrated structure.
- Using special permanent magnets for maintaining stirring ability.

Use Convenience Features

- Precise temperature/stirring control with PID feedback control.
- Temperature control mode selection function. (Optimal/Fast/Slow/User/Point)
- Includes temperature auto-tuning and calibration functions.
- Count-down timer provided. (up to 99 hours and 59 minutes)
- Clamp Rod (option) can be added to the main body fix other experiment equipment.
- Convenient experiment observation with removable transparent shield. (option)
- Silicone Cover protects from the dust and overflow of samples. (TS/TM/T)

Outstanding Safety

- Spill-proof design minimizes inflow solution into device.
- Safety ensured by BLDC motor that does not generate sparks.
- Preventing equipment failure with multiple over temperature prevention devices.
- Warning displayed when top plate is 50°C or more.
- Prevents sample run-off thanks to smooth stirring start.
- The heating bath (option) is combined with the top plate to ensure safety. (except TS-QG)



TS-18QG

with 3 Prong Clamp, Clamp Rod, Magnetic Bar (option)



TS-17SG

Specification

Model		TS-14SG	TS-17SG	TS-18QG
Heating	Temperature range (°C / °F)	Top plate, Max, 350 / 662	Top plate, Max, 350 / 662	Top plate, Max, 350 / 662
	Control mode	PID Feedback (optimal, fast, slow)	PID Feedback (optimal, fast, slow)	PID Feedback (optimal, fast, slow)
	Display resolution (°C / °F)	0.1 / 32.18	0.1 / 32.18	0.1 / 32.18
Stirring	Speed range (rpm)	30 to 2000	30 to 2000	30 to 2000
	Capacity (L / cu ft, H ₂ O)	20 / 0.71	20 / 0.71	20 / 0.71
	Magnetic bar, Max. (L x Ø, mm / inch)	40 x 8 / 0.31 x 1.57	40 x 8 / 0.31 x 1.57	40 x 8 / 0.31 x 1.57
Material	Top plate	White ceramic coated aluminum alloy	White ceramic coated aluminum alloy	White ceramic coated aluminum alloy
	Body	Powder coated aluminum alloy	Powder coated aluminum alloy	Powder coated aluminum alloy
	Motor type	BLDC (Brushless Direct Current)	BLDC (Brushless Direct Current)	BLDC (Brushless Direct Current)
	Temperature probe	PT 100 (B class, Max. 250°C, 482°F)	PT 100 (B class, Max. 250°C, 482°F)	PT 100 (B class, Max. 250°C, 482°F)
Dimensions	Top plate (Ø or W x D, mm / inch)	140 / 5.51	170 / 6.69	180 x 180 / 7.08 x 7.08
	Exterior (W x D x H, mm / inch)	161 x 290 x 100 / 6.34 x 11.42 x 3.94	191 x 330 x 101 / 7.52 x 12.99 x 3.98	209 x 326 x 102 / 8.22 x 12.83 x 4.0
	Net weight (kg / lbs)	2.8 / 6.17	3.5 / 7.71	3.8 / 8.36
Timer (Heating and/or Stirring)		Max. 99 hrs 59 min.	Max. 99 hrs 59 min.	Max. 99 hrs 59 min.
Electrical requirements (230V, 50/60Hz, A)		3.0	4.0	4.0
Cat. No.		AAH34445K	AAH34475K	AAH34485K
Electrical requirements (120V/60Hz, A)		5.0	6.7	6.7
Cat. No.		AAH34443U	AAH34473U	AAH34483U

※ Except TS-18QG CSA certification

Accessories [Page 152](#) Temperature Probe, Heating Bath, Clamp Rod, Clamp, Transparent Shield, Silicone Cover

Hotplate & Magnetic Stirrer

Analog type



Easy to use with excellent durability

Structural Functional Features


- Three models provided according to the size and shape of hot plate.
- Top plate coated with white ceramic for excellent corrosion resistance.
- Excellent heat transfer and durability due to heater-integrated structure.
- Using special permanent magnets for maintaining stirring ability.

Use Convenience Features

- Easy adjustment of heating rate and stirring speed with individual control knobs.
- Maintains a constant stirring speed even when sample viscosity changes.
- Rapid stirring and stop function according to the user's needs.
- Clamp Rod (option) can be added to the main body to fix other experiment equipment.
- Convenient experiment observation with removable transparent shield. (option)
- Silicone Cover protects from the dust and overflow of samples. (TS/TM/T)

Outstanding Safety

- Spill-proof design minimizes inflow solution into device.
- Safety ensured by BLDC motor that does not generate sparks.
- Preventing equipment failure with multiple over temperature prevention devices.
- Warning displayed when top plate is 50°C or more.
- Prevents sample run-off thanks to smooth stirring start.
- The heating bath (option) is combined with the top plate to ensure safety. (except TM-QG)



TM-14SG

with Heating Bath, Clamp Holder, 3 Prong Clamp, Clamp Rod (option)



TM-18QG



Specification

Model		TM-14SG	TM-17SG	TM-18QG
Heating	Temperature range (°C / °F)	Top plate, Max, 350 / 662	Top plate, Max, 350 / 662	Top plate, Max, 350 / 662
	Control mode	Scale	Scale	Scale
	Heating rate	0 ~ 100% by 1%	0 ~ 100% by 1%	0 ~ 100% by 1%
Stirring	Speed range (rpm)	Max. 2000	Max. 2000	Max. 2000
	Capacity (L / cu ft, H ₂ O)	20 / 0.71	20 / 0.71	20 / 0.71
	Magnetic bar , Max. (L x Ø, mm / inch)	40 x 8 / 0.31 x 1.57	40 x 8 / 0.31 x 1.57	40 x 8 / 0.31 x 1.57
Material	Top plate	White ceramic coated aluminum alloy	White ceramic coated aluminum alloy	White ceramic coated aluminum alloy
	Body	Powder coated aluminum alloy	Powder coated aluminum alloy	Powder coated aluminum alloy
	Motor type	BLDC (Brushless Direct Current)	BLDC (Brushless Direct Current)	BLDC (Brushless Direct Current)
Dimensions	Top plate (Ø or W x D, mm / inch)	140 / 5.51	170 / 6.69	180 x 180 / 7.08 x 7.08
	Exterior (W x D x H, mm / inch)	161 x 290 x 100 / 6.34 x 11.42 x 3.94	191 x 330 x 101 / 7.52 x 12.99 x 3.98	209 x 326 x 102 / 8.22 x 12.83 x 4.0
	Net weight (kg / lbs)	2.8 / 6.17	3.5 / 7.72	3.8 / 8.37
Body Shape		Straight body		
Electrical requirements (230V, 50/60Hz, A)		3.0	4.0	4.0
Cat. No.		AAH34245K	AAH34275K	AAH34285K
Electrical requirements (120V, 60Hz, A)		5.0	6.7	6.7
Cat. No.		AAH34243U	AAH34273U	AAH34283U

※ Except TM-18QG CSA certification

Accessories [Page 152](#) Heating Bath, Clamp Rod, Clamp, Transparent Shield, Silicone Cover

Hotplate & Magnetic Stirrer Multi type



Differentiated multi-stirrer with temperature control

Structural Functional Features

- Integrated design with top plate and heater for fast heat transfer and excellent durability.
- Using special permanent magnets for maintaining stirring ability.
- Control by All/Column/Row/Point with four stirring control methods and timer function.

Use Convenience Features


- Separate heater knobs and operation LED display.
- Maintains accurate speed with feedback control even when sample viscosity and amount change.
- Count-down timer provided. (up to 99 hours and 59 minutes)
- Rapid stirring and stop function according to the user's needs.
- VFD display with excellent visual perception.



MS-53MH

Outstanding Safety

- The upper part of the top plate is composed of silicon pads, allowing for prevention of slippage of glassware and display of stirring point.
- Spill-proof design minimizes inflow solution into device.
- Safety ensured by BLDC motor that does not generate sparks.
- Safe by automatically shutting down the heater power when overheated.
- Prevents sample run-off thanks to smooth stirring start.



Specification

Model		MS-33MH	MS-53MH
Heating	Temperature range (°C / °F) ¹⁾	Top plate, Max, 120 / 248	Top plate, Max, 120 / 248
	Control mode	Scale	Scale
	Heating Power (W)	600	1000
Stirring	Speed range (rpm)	30 to 2000	30 to 2000
	Capacity per point at 2,000rpm (H ₂ O, mL / cu ft)	250 / 0.008	250 / 0.008
	Operating mode	4 (All, Column, Row, Point)	4 (All, Column, Row, Point)
	Position (Row x Column)	9 (3 x 3)	15 (5 x 3)
	Point distance (W x D, mm / inch)	117 x 90 / 4.61 x 3.54	117 x 90 / 4.61 x 3.54
	Magnetic bar , Max. (L x Ø, mm / inch)	30 x 8 / 1.18 x 0.31	30 x 8 / 1.18 x 0.31
	Load, Max (kg / lbs)	30 / 66.14	30 / 66.14
Material	Top plate	Aluminum	Aluminum
	Body	Powder coated steel	Powder coated steel
	Motor type	BLDC (Brushless Direct Current)	BLDC (Brushless Direct Current)
Dimensions	Top plate (W x D, mm / inch)	385 x 328 / 15.16 x 12.91	385 x 493 / 15.16 x 19.41
	Exterior (W x D x H, mm / inch)	395 x 450 x 104 / 15.55 x 17.72 x 4.09	395 x 614 x 104 / 15.55 x 24.17 x 4.09
	Net weight (kg / lbs)	10 / 22.05	15 / 33.06
Timer (Stirring)		Max. 99 hrs 59 min.	Max. 99 hrs 59 min.
Electrical requirements (230V, 50/60Hz, A)		4.1	5.8
Cat. No.		AAHK34015K	AAHK34025K
Electrical requirements (120V, 60Hz, A)		8.0	11.2
Cat. No.		AAHK34013U	AAHK34023U

1) This refers to the maximum temperature of the top plate. The temperature of the actual specimen may be reduced depending on the conditions of the experiment.

Accessories [Page 152](#) Magnetic Bar, Magnetic Retriever

Hotplate Digital type



Precise hotplate using digital method

Structural Functional Features

- Three models provided according to the size and shape of hot plate.
- Top plate coated with white ceramic for excellent corrosion resistance.
- Excellent heat transfer and durability due to heater-integrated structure.
- Possible to set the upper and lower limit of temperature control range.

Use Convenience Features

- Temperature control by adjustment of heating rate.
- Temperature control mode selection function. (Optimal/Fast/Slow/User/Point)
- Includes temperature auto-tuning and calibration functions.
- Count-down timer provided. (up to 99 hours and 59 minutes)
- Clamp Rod (option) can be added to the main body to fix other experiment equipment.
- Convenient experiment observation with removable transparent shield. (option)
- Silicone Cover protects from the dust and overflow of samples. (TS/TM/T)

Outstanding Safety

- Spill-proof design minimizes inflow solution into device.
- Preventing equipment failure with multiple over temperature prevention devices.
- Separate buttons for safe operating.
- Warning displayed when top plate is 50°C or more.
- Malfunction prevented by controller lock function.
- The heating bath (option) is combined with the top plate to ensure safety. (except T-QG)



T-14SG

with Clamp Holder, 3 Prong Clamp, Clamp Rod (option)



T-18QG



Specification

Model		T-14SG	T-17SG	T-18QG
Heating	Temperature range (°C / °F)	Top plate, Max, 350 / 662	Top plate, Max, 350 / 662	Top plate, Max, 350 / 662
	Control mode	PID Feedback (optional, fast, slow)	PID Feedback (optional, fast, slow)	PID Feedback (optional, fast, slow)
	Display resolution (°C / °F)	0.1 / 32.18	0.1 / 32.18	0.1 / 32.18
Material	Top plate	White ceramic coated aluminum alloy	White ceramic coated aluminum alloy	White ceramic coated aluminum alloy
	Body	Powder coated aluminium alloy	Powder coated aluminium alloy	Powder coated aluminium alloy
Dimensions	Top plate (Ø or W x D, mm / inch)	140 / 5.51	170 / 6.69	180 x 180 / 7.09 x 7.09
	Exterior (W x D x H, mm / inch)	161 x 290 x 100 / 6.34 x 11.41 x 3.94	191 x 330 x 101 / 7.52 x 12.99 x 3.98	209 x 326 x 102 / 8.23 x 12.83 x 4.01
	Net weight (kg / lbs)	2.2 / 4.85	2.9 / 6.39	3.4 / 7.49
Timer		Max. 99 hrs 59 min.	Max. 99 hrs 59 min.	Max. 99 hrs 59 min.
Electrical requirements (230V, 50/60Hz, A)		3.0	4.0	4.0
Cat. No.		AAH35045K	AAH35075K	AAH35085K
Electrical requirements (120V, 60Hz, A)		5.0	6.7	6.7
Cat. No.		AAH35043U	AAH35073U	AAH35083U

※Except T-18QG CSA certification

Accessories [Page 152](#) Heating Bath, Clamp Rod, Clamp, Transparent Shield, Silicone Cover

Heating Block for Round Flask and Vial

A safe solution to replace oil bath and heating mantle when heating the contents over 100°C

Using with Hotplate

- The heating block is mounted on a hotplate or hotplate magnetic stirrer to heat safely the contents of the 25ml - 2000ml round flask and vial.
- Excellent compatibility with various brands of hotplates. (flask dia. 135/140, Vial-135/180/150 x 150/180 x 180 (mm))
- Combination of inserts (option) allows the use of flasks of various capacities.

Excellent Safety Standards

- No risk of burns or noxious fumes due to heated oil that may occur when using an oil bath.
- Hot spots do not occur even when heated up to 300°C, resulting in excellent temperature accuracy and prevention of damage to glassware.

Use Convenience Features

- Economical and eco-friendly by eliminating the need to purchase and dispose of oil separately.
- Since oil is not used, it is easy to clean glassware and the product.



Convenient maintenance by an excellent durability.

Hard anodized surface is treated on aluminum block body to improve durability such as wear and corrosion resistance. Also cleaning is convenient.



Excellent thermal conductivity makes rapid heating.

Design to contact closely flask to the surface, no crack on a glassware and improve the thermal delivery.



Precise temperature control

Holes for a temperature probe make control the temperature precisely. (4 holes on a plate for a round flask, and 3 holes on an insert for the vial.)



Various sizes

4 kinds of round flask plate for 100ml to 2L flask and its insert, and a vial plate available to place its insert up to 4.



Good compatibility with hotplate

Possible to use not only a round hotplate that has 135 or 140mm diameter but also a square or round hotplate that has 170mm. (when removes support.)



Stable structure fixing the bottom

Fixing support designed fitting to a plate size of round and square hotplate to use in a stable. (Round flask plate is able to fix only for Φ135 / Φ140 hotplate.)



Selection guide

For Round flask

Plate					
100 ml (RBP-100)		500 ml (RBP-500)		1000 ml (RBP-1000)	2000 ml (RBP-2000)
Insert					
25 ml (RBI-25)	50 ml (RBI-50)	100 ml (RBI-100)	250 ml (RBI-250)	500 ml (RBI-500)	1000 ml (RBI-1000)

For Vial

Plate					
190 x 190 x 20 mm (VBP)					
Insert (Max. 4 inserts on the plate.)					
4 ml (1 dram) Φ15 x 17 mm, 16 holes (VBI-1D)	8 ml (2 dram) Φ17 x 23 mm, 16 holes (VBI-2D)	15 ml (4 dram) Φ21 x 30 mm, 9 holes (VBI-4D)	20 ml Φ28 x 22 mm, 4 holes (VBI-20S)	30 ml Φ28 x 32 mm, 4 holes (VBI-30S)	40 ml Φ28 x 42 mm, 4 holes (VBI-40S)

Magnetic Stirrer

Digital type



Accurate and smooth feedback speed control

Structural Functional Features

- Maintains accurate speed with feedback control even when sample viscosity and amount change.
- Polypropylene (PP) material with excellent chemical resistance.
- Minimized vibration during high-speed stirring with rubber feet.
- Using special permanent magnets for maintaining stirring ability.

Use Convenience Features

- Simultaneous display of set speed and actual speed.
- Timer (up to 99 hours 59 minutes) and program operation function. (2-step repeat pattern)
- Clamp Rod (option) can be added to the main body to fix other equipment such as burette, sensor, etc.
- 4 models provided according to the quantity and size of stirring.
- The two types of black and white covers are very useful for observing changes according to samples.

MS-22G



MS-17G

Outstanding Safety

- Safety ensured by BLDC motor that does not generate sparks.
- Spill-proof design minimizes inflow solution into device.
- Prevents samples splash with smooth operation at the start and stop of stirring.
- Provides silicon-based top plate cover to prevent slipping accidents of glassware.



Specification

Model		MS-17G	MS-22G
Stirring	Speed range (rpm)	30 to 2000	30 to 2000
	Capacity (L / cu ft, H ₂ O)	4 / 0.14	5 / 0.18
	Speed stability (± %)	1	1
	Magnetic bar, Max. (L x Ø, mm / inch)	40 x 8 / 0.31 x 1.57	50 x 10 / 1.97 x 0.39
	Load, Max. (kg / lbs)	20 / 44.09	25 / 55.12
Material	Top plate	Polypropylene with silicone cover	Polypropylene with silicone cover
	Body	Polypropylene	Polypropylene
	Motor type	BLDC (Brushless Direct Current)	BLDC (Brushless Direct Current)
Dimensions	Top plate	170 (Ø, mm) / 6.69 (inch)	220 (Ø, mm) / 8.66 (inch)
	Exterior (W x D x H, mm / inch)	210 x 280 x 72 / 8.27 x 11.02 x 2.83	260 x 355 x 72 / 10.23 x 13.98 x 2.83
	Net weight (kg / lbs)	1.4 / 3.09	1.8 / 3.97
Electrical requirements		AC 100 ~ 240V, 50/60Hz, 0.4A Max.	
Power consumption (W)		6	6
Cat. No.	KR plug	AAH331215BK	AAH331315BK
	US plug	AAH331215BU	AAH331315BU

Accessories [Page 153](#) Plate Cover, Clamp Rod, Clamp Holder, Magnetic Bar, Magnetic Retriever

Magnetic Stirrer

Analog type



Quick speed control of stirring with real time response

Structural Functional Features

- 3 models according to size, shape.
- Polypropylene (PP) material with excellent chemical resistance.
- Minimized vibration during high-speed stirring with rubber feet.
- Using special permanent magnets for maintaining stirring ability.

Use Convenience Features

- Quick control of stirring speed in real time according to degree of control knobs.
- The two types of black and white covers are very useful for observing changes according to samples.
- Convenient LED lamp for checking operation status.
- Lightweight makes handling convenient.

Outstanding Safety

- Safety ensured by BLDC motor that does not generate sparks.
- Spill-proof design minimizes inflow solution into device.
- Provides Silicon-based top plate cover to prevent slipping accidents of glassware.



MS-17B

Stirrer & Mixer

Specification

Model		MS-12B	MS-17B	MS-22B
Stirring	Speed range (rpm)	150 to 2500	150 to 2500	150 to 2500
	Capacity (L / cu ft, H ₂ O)	2 / 0.07	4 / 0.14	5 / 0.18
	Magnetic bar, Max. (L x Ø, mm)	30 x 8 / 1.18 x 0.31	40 x 8 / 1.57 x 0.31	50 x 10 / 1.97 x 0.39
	Load, Max. (kg / lbs)	15 / 33.07	20 / 44.09	25 / 55.12
Material	Top plate	Polypropylene with silicone cover	Polypropylene with silicone cover	Polypropylene with silicone cover
	Body	Polypropylene	Polypropylene	Polypropylene
	Motor type	BLDC (Brushless Direct Current)	BLDC (Brushless Direct Current)	BLDC (Brushless Direct Current)
Dimensions	Top plate (Ø or W x D, mm / inch)	120 / 4.72	170 / 6.69	220 / 8.66
	Exterior (W x D x H, mm / inch)	130 x 133 x 65 / 5.11 x 5.24 x 2.56	180 x 186 x 65 / 7.09 x 7.32 x 2.56	230 x 235 x 65 / 9.05 x 9.25 x 2.56
	Net weight (kg / lbs)	0.6 / 1.32	0.9 / 1.98	1.1 / 2.42
Electrical requirements		AC 100 ~ 240V, 50/60Hz with power adapter, 0.2A Max.		
Power consumption (W)		5	6	6
Cat. No.	KR plug	AAH330115BK	AAH330215BK	AAH330315BK
	US plug	AAH330115BU	AAH330215BU	AAH330315BU

Accessories Page 153 Plate Cover, Power Adapter, Magnetic Bar, Magnetic Retriever

Magnetic Stirrer

Multi Position type

Differentiated multi-stirring control for multi-purpose use

Structural Functional Features

- Maintains accurate speed with feedback control even when sample viscosity and amount change.
- Display of stirring position on the upper part of the stainless steel.
- Using special permanent magnets for maintaining stirring ability.

Use Convenience Features

- Simultaneous display of set speed and actual speed.
- Count-down timer provided. (up to 99 hours and 59 minutes)
- Six models provided according to the quantity and size of stirring.



MS-53M

Outstanding Safety

- Safety ensured by BLDC motor that does not generate sparks.
- Spill-proof design minimizes inflow solution into device.
- Prevents samples splash with smooth operation at the start and stop of stirring.
- Silicon pad (option) to prevent slip-related accidents of glassware.

Differentiated Multi Control

- Four kinds of stirring control methods and timer function allow for performance of various stirring experiments simultaneously. (MS-51M has two control methods.)
- Display LED point window for easy setting and confirmation.

All Mode

Simultaneous stirring with same speed.

Column Mode


Stirring with different speeds for each column.

Row Mode

Stirring with different speeds for each row.

Point Mode

Stirring with different speeds for each point





Stirrer & Mixer

Specification

Model		MS-51M	MS-23M	MS-32M	MS-33M	MS-52M	MS-53M
Stirring	Speed range (rpm)	30 to 2000	30 to 2000	30 to 2000	30 to 2000	30 to 2000	30 to 2000
	Capacity per point at 2000rpm (H ₂ O, mL)	250	250	250	250	250	250
	Operating mode	2 (All, Point)	4 (All, Column, Row, Point)	4 (All, Column, Row, Point)	4 (All, Column, Row, Point)	4 (All, Column, Row, Point)	4 (All, Column, Row, Point)
	Postion (Row x Column)	5 (5 x 1)	6 (2 x 3)	6 (3 x 2)	9 (3 x 3)	10 (5 x 2)	15 (5 x 3)
	Point distance (W x D, mm / inch)	90 / 3.54	117 x 170 / 4.61 x 6.69	117 x 90 / 4.61 x 3.54	117 x 90 / 4.61 x 3.54	117 x 90 / 4.61 x 3.54	117 x 90 / 4.61 x 3.54
	Magnetic bar , Max. (L x Ø, mm / inch)	30 x 8 / 1.18 x 0.31	30 x 8 / 1.18 x 0.31	30 x 8 / 1.18 x 0.31	30 x 8 / 1.18 x 0.31	30 x 8 / 1.18 x 0.31	30 x 8 / 1.18 x 0.31
	Load, Max. (kg / lbs)	30 / 66.14	30 / 66.14	30 / 66.14	30 / 66.14	30 / 66.14	30 / 66.14
Material	Top plate	Stainless steel	Stainless steel	Stainless steel	Stainless steel	Stainless steel	Stainless steel
	Body	Powder coated steel, PET	Powder coated steel, PET	Powder coated steel, PET	Powder coated steel, PET	Powder coated steel, PET	Powder coated steel, PET
	Motor type	BLDC (Brushless Direct Current)	BLDC (Brushless Direct Current)	BLDC (Brushless Direct Current)	BLDC (Brushless Direct Current)	BLDC (Brushless Direct Current)	BLDC (Brushless Direct Current)
Dimensions	Top plate (W x D, mm / inch)	124 x 485 / 4.88 x 19.09	362 x 318 / 14.25 x 12.52	242 x 318 / 9.53 x 12.52	362 x 318 / 14.25 x 12.52	246 x 485 / 9.68 x 19.09	367 x 485 / 14.45 x 19.09
	Exterior (W x D x H, mm / inch)	131 x 603 x 82 / 5.56 x 23.74 x 3.23	377 x 436 x 82 / 14.84 x 17.16 x 3.23	254 x 436 x 82 / 10 x 17.16 x 3.23	344 x 436 x 81 / 13.54 x 17.16 x 3.19	254 x 602 x 82 / 10 x 23.7 x 3.23	377 x 599 x 82 / 14.84 x 23.58 x 3.23
	Net weight (kg / lbs)	4 / 8.81	6 / 13.23	5 / 11.02	7 / 15.43	8 / 17.64	11 / 24.25
Timer	Max. 99 hrs 59 min.	Max. 99 hrs 59 min.	Max. 99 hrs 59 min.	Max. 99 hrs 59 min.	Max. 99 hrs 59 min.	Max. 99 hrs 59 min.	
Electrical requirements		AC 100 ~ 240V, 50/60Hz, 11.0A Max.					
Power consumption (W)		45	45	45	54	85	135
Cat. No.	KR plug	AAH332415K	AAH332115K	AAH332215K	AAH332315K	AAH332515K	AAH332615K
	US plug	AAH332415U	AAH332115U	AAH332215U	AAH332315U	AAH332515U	AAH332615U

Accessories Page 153 Silicone Pad, Magnetic Bar, Magnetic Retriever

Accessories for Hotplate & Magnetic Stirrer



Temperature Probe

With external temperature sensor, it is possible to check and control the temperature of the sample when it is connected to the equipment. (TS model)



Heating Bath

Top plate combines concave and convex structure, resulting in no slippage.



Clamp Rod

Two pieces can be fixed on the left and right of the back of the main body.



3 Prong Clamp

Various types of instruments can be fixed.



Clamp Holder

Clamp rod and clamp can be fixed.



Transparent Shield

Safely observe experiments.



Silicone Cover

- Protect from the dust and overflow of samples.
- Simple and easy installation with lightweight silicone materials.



Magnetic Bar [Page 154](#)

- Applied to a magnetic stirrer when stirring within the temperature range up to +280°C.
- Turbo type can be applied at high temperature of 400°C.



Magnetic Retriever [Page 155](#)

- Used when taking out magnetic bar.
- Uses ALNICO V magnets and samarium cobalt magnets.




Jeio Tech



Protective Silicone Cover for Hotplate and hotplate stirrer

Cat. No.	Description
AAA34501	Temperature Probe, B Class (Max. 250°C)
AAA34502	Temperature Probe, A Class (Max. 400°C)
00HPS0000012	Heating Bath (TS-14SG, TM-14SG, T-14SG)
00HPS0000015	Heating Bath (TS-17SG, TM-17SG, T-17SG)
00MTT0000132	Clamp Rod (Ø12, 400 mm)
BEA1000011	3 Prong Clamp (80 mm grip)
BEA1000012	3 Prong Clamp (60 mm grip)
BEA1000013	3 Prong Clamp (20 mm grip)
AAA37511	Clamp Holder (max. Ø27, C-20)
AAA37512	Clamp Holder (max. Ø17, C-10)
AAA37513	Clamp Holder (max. Ø12, C-7)
AAA37514	Clamp Holder (max. Ø12, C-5)
00HPS0000059	Transparent Shield (TS-14SG, TM-14SG, T-14SG)
00HPS0000039	Transparent Shield (TS-17SG, TM-17SG, T-17SG)
00HPS0000190	Silicone cover (TS-14SG, TM-14SG, T-14SG)
00HPS0000191	Silicone cover (TS-17SG, TM-17SG, T-17SG)
00HPS0000192	Silicone cover (TS-18SG, TM-18SG, T-18SG)



Accessories for Magnetic Stirrer



Plate Cover (Black, White)



Silicone Pad



Clamp Rod

Plate Cover

Cat. No.	Description	Model
00MTT0000002	White color, Silicone (Ø120)	MS-12B
00MTT0000003	White color, Silicone (Ø170)	MS-17B/G
00MTT0000004	White color, Silicone (Ø220)	MS-22B/G
00MTT0000005	Black color, Silicone (Ø120)	MS-12B
00MTT0000006	Black color, Silicone (Ø170)	MS-17B/G
00MTT0000007	Black color, Silicone (Ø220)	MS-22B/G

Silicone Pad

Cat. No.	Description (W x D x H, mm / inch)	Model
00MTS0000015	355 x 475 x 1 / 13.98 x 18.7 x 0.04	MS-53M
00MTS0000016	235 x 475 x 1 / 9.25 x 18.7 x 0.04	MS-52M
00MTS0000017	115 x 475 x 1 / 4.53 x 18.7 x 0.04	MS-51M
00MTS0000041	235 x 310 x 1 / 9.25 x 12.2 x 0.04	MS-32M
00MTS0000042	355 x 310 x 1 / 13.98 x 12.2 x 0.04	MS-23M / 33M

Power Adapter & Clamp Rod

Cat. No.	Description	Model
00MTT0000001	Power Adapter (100~240V, 50/60Hz)	MS-12B/17B
00MTT0000132	Clamp Rod (Ø12 x 400 / Ø0.47 x 15.75", M10)	MS-17G/22G

Clamp Holder



- Up to Ø27 mm rods mountable.
- Manufactured using aluminum die-casting and plastic for lightweight but strong fixation.
- Knob structure that is comfortable to handle when locking / unlocking.
- The contact surface with the stand support is stable and does not slip.
- C-20 and C-10 models are suitable for the overhead stirrer.
- The C-10, C-7 and C-5 models are suitable for use with hotplate and magnetic stirrer.

Model	C-20	C-10	C-7	C-5
Max. Ø (mm / inch)	27 / 1.06	17 / 0.66	12 / 0.47	12 / 0.47
Open Ø (mm / inch)	20 / 0.79	15 / 0.59	-	-
Body material	Aluminum	Aluminum	Aluminum	Polypropylene
Knob material	Polypropylene	Polypropylene	Polypropylene	Polypropylene
Bolt material	Stainless Steel	Stainless Steel	Stainless Steel	Stainless Steel
Suitable for	Overhead Stirrer	Overhead Stirrer	-	-
	-	Hotplate & Magnetic Stirrer	Hotplate & Magnetic Stirrer	Hotplate & Magnetic Stirrer
Cat. No.	AAA37511	AAA37512	AAA37513	AAA37514

Accessories Magnetic Bar

Features of Magnetic Bar

- Chemically stable with PTFE inert material.
- Completely sealed structure with no cracks or pores.
- Smooth finish to avoid absorption or cross contamination.
- Working Temperature Range: up to +280°C.
- Can be sterilized using chemical or physical (thermal) methods. (except γ Radiation)
- Magnet-ALNICO V, PTFE Coated.


※ The size of the magnetic bar is the average value, and there may be an error of $\pm 5\%$ in length and $\pm 10\%$ in diameter.



Cylindrical Magnetic Bar

- Smooth and round shape.
- The most commonly used magnetic bar for a variety of applications.


Cat. No.	D (mm / inch)	L (mm / inch)
BCW0011138	8 / 0.31	13 / 0.51
BCW0011156	6 / 0.23	15 / 0.59
BCW0011208	8 / 0.31	20 / 0.78
BCW0011258	8 / 0.31	25 / 0.98
BCW0011308	8 / 0.31	30 / 1.18
BCW001140	8 / 0.31	40 / 1.57
BCW001150	8 / 0.31	50 / 1.96



Octahedral Magnetic Bar

- Octagonal design for more efficient sample mixing.
- Used in containers having slight slopes or uneven surfaces.
- The center of the pivot ring selects the best position for stirring.
- Bigger turbulence than pivot ring is generated at lower speeds.
- Pivot ring reduces friction and noise.


Cat. No.	D (mm / inch)	L (mm / inch)
BCW001513	8 / 0.31	13 / 0.51
BCW001515	8 / 0.31	15 / 0.59
BCW001522	8 / 0.31	22 / 0.86
BCW001525	8 / 0.31	25 / 0.98
BCW001528	8 / 0.31	28 / 0.10
BCW001538	8 / 0.31	38 / 1.49
BCW001541	8 / 0.31	41 / 1.61
BCW001551	8 / 0.31	51 / 2.00
BCW00155110	10 / 0.39	51 / 2.00



Oval Magnetic Bar

- Designed to fit round-bottom flasks, the bar itself has a central axis.

Cat. No.	D (mm / inch)	L (mm / inch)
BCW001610	5 / 0.19	10 / 0.39
BCW001615	6 / 0.23	15 / 0.59
BCW001620	10 / 0.39	20 / 0.78
BCW00162510	10 / 0.39	25 / 0.98
BCW001630	16 / 0.63	30 / 1.18
BCW001635	16 / 0.63	35 / 1.38
BCW001640	20 / 0.78	40 / 1.57
BCW001650	20 / 0.78	50 / 1.96



Accessories Magnetic Retriever

Polypropylene

- Product with built-in ALNICO V magnet and coated with polypropylene.
- Convenient structure having ring part at the end for hanging after use.

PTFE

- Product with built-in ALNICO V magnet and coated with PTFE. (teflon)



Polypropylene



PTFE



Stirrer & Mixer

Cat. No.	L (mm)	Ø (mm)	Description
BCW0041300	300	10	Polypropylene, ALNICO V
BCW0041350	350	10	
BCW0041450	450	10	
BCW004150	150	10	PTFE, ALNICO V
BCW004250	250	10	
BCW004350	350	10	
BCW004450	450	10	




Overhead Stirrer

Possesses practical functions required in actual experiments, and has a beneficial and convenient configuration system

- > **Safe, spark-free motor**
Use of spark-free BLDC motor for safer experimental environment.
- > **Reliable control technology**
Real-time PID feedback control method maintains accurate rpm even when sample viscosity changes.
- > **Excellent structure for long-term continuous use**
Excellent durability of BLDC motor and structure allows for powerful and excellent heat dissipation.
- > **Patented best effort performance of superior level**
Automatic rpm control in case overload operation. Attempt to stir with set rpm according to viscosity change.
- > **Over-temperature and overload protection**
Built-in over-temperature and overload protection device with safety functions for continuous operation of the equipment.
- > **Free selection of rotation direction**
Rotation direction is selectable for various impeller uses. (except for MSA)
- > **Computer interface control**
Convenient remote control and data management through provided dedicated software. (except MSA)
Available to control remotely through an external controller and Bluetooth. (MSH)
- > **Comfortable use with smooth stirring start**
Progressively and smoothly reaches the set rpm, preventing accidents such as sample spill-out.



Convenient remote control and data management through dedicated software. (except MSA)



Dedicated software provided as standard.



Free selection of rotation direction. (except MSA)



Optimized model configuration and performance Jeio Tech's Overhead Stirrer

> Model configuration for the appropriate choice

Four models provided for optimum selection according to torque, rpm, function, and display.

> Prestige touch model (MSH)

Color touch screen and control
Knob control via external controller. (option)

> High performance model (MSD)

Clear FND & control knob.
Control via PC using USB connection.

> Easy control model (MSA)

Clear FND and control knob.
Over temperature and overload protection functions included.

> Slim & compact design

Design that minimizes interference with other devices even in the case of a complex set of experiments.

> Pass-through shaft design

Impeller shaft structure can be passed through, making it easy to control the height of the impeller at the time of experiment configuration.

> Easy to use with intuitive operation

Simple and intuitive operation environment.
Useful function of displaying device information.

> Controller lock to prevent erroneous operation

Controller lock to prevent errors caused by unexpected pressing of buttons. (except for MSA)

※ Some of the above contents are limited to specific models.



Model configuration according to torque, rpm, function, and display part.



Structure allows for convenient adjustment of impeller height.



Separate external controller for remote control. (MSH)

Overhead Stirrer

Prestige touch type

Prestige touch model with highest performance and stability

Structural Functional Features

- Possible for user to select clockwise (CW) / counter clock wise. (CCW)
- Slim & compact design.
- Convenient adjustment of height with pass-through impeller shaft.
- Chuck with excellent corrosion resistance. (max. Ø10 mm)
- 2 models provided according to torque and rpm.

Use Convenience Features

- Best effort function to maintain maximum set rpm by monitoring sample viscosity changes in real time.
- Accurate stirring with PID feedback control.
- Count-down timer provided. (up to 99 hours and 59 minutes)
- Operating elapsed time display function.
- Machine operation and data management with USB connection.
- Complete with accessories such as impeller, stand, clamp, etc.

Outstanding Safety


- Safety ensured by BLDC motor that does not generate sparks.
- Excellent heat dissipation structure makes use for a long time without any trouble.
- Continual use in safe by over temperature and overload protection.
- Malfunction prevented by controller lock function.
- Prevents sample run-off by smooth stirring start and stop function.
- Motor protection by displaying motor temperature in three colors.



MSH-0512

Easy and convenient color touch LCD controller

- Intuitive operation with touch LCD controller and control knobs.
- Displays set value, actual value, time, torque, and motor temperature.



PC control and data management with USB connection

- Useful for data management through PC control.
- Dedicated software and USB cable provided as standard.





Specialized External Controller

- When harmful fumes are generated through stirring, the reaction proceeds inside the fume hood and the device is safely controlled by external controller at outside.
- Checking control status in real time and controls start/stop or stirring speed changes.
- Color LCD display for intuitive and easy control.



Model	RMS
Exterior (W x D x H, mm / inch)	60 x 16 x 120 / 2.4 x 0.6 x 4.7
Wire length (mm / inch)	3000 / 118
Cat. No.	AAA37591

Specification

Model		MSH-0512	MSH-0520
Speed range (rpm)		50 to 1200	50 to 2000
Viscosity, max. (cP)¹⁾		30000 (~1200rpm) 50000 (~700rpm)	5000 (~2000rpm) 10000 (~1600rpm) 30000 (~600rpm) 50000 (~300rpm)
Stirring capacity (H₂O), max. (L / cu ft)		100 / 3.53 (~1200rpm)	60 / 2.12 (~2000rpm)
Rated torque (N-cm / inch)		41.6 / 16.38	20.8 / 8.19
Motor input / output (W)		71 / 48	71 / 48
Speed control		PID feedback control	PID feedback control
Material	Body	Powder coated aluminum	Powder coated aluminum
	Cover	Polypropylene	Polypropylene
	Motor type	BLDC (Brushless Direct Current)	BLDC (Brushless Direct Current)
Dimension	Chuck range (Ø, mm / inch)	3~10 / 0.12~0.39	3~10 / 0.12~0.39
	Exterior (W x D x H, mm / inch)	80 x 185 x 235 / 3.15 x 7.28 x 9.25	80 x 185 x 235 / 3.15 x 7.28 x 9.25
	Net weight (kg / lbs)	3.1 / 6.83	3.1 / 6.83
Electrical requirements		AC 100 to 240V, 50/60Hz with power adapter, 11.0A Max.	
Cat. No.	KR plug	AAH371415K	AAH371515K
	US plug	AAH371415U	AAH371515U

1) Tested with silicone oil at 26°C and 60% RH, and uses triple-bladed propeller impeller (Ø50).

Overhead Stirrer

High performance type



High performance model for faithful advanced features and safety

Structural Functional Features

- Possible for user to select clockwise (CW) / counter clock wise. (CCW)
- Slim & compact design.
- Convenient adjustment of height with pass-through impeller shaft.
- Chuck with excellent corrosion resistance. (max. Ø10 mm)

Use Convenience Features

- Best effort function to maintain maximum set rpm by monitoring sample viscosity changes in real time.
- Accurate stirring with PID feedback control.
- Operating elapsed time display function.
- Machine operation and data management with USB connection.
- Complete with accessories such as impeller, stand, clamp, etc.

Outstanding Safety

- Safety ensured by BLDC motor that does not generate sparks.
- Excellent heat dissipation structure makes use for a long time without any trouble.
- Continual use in safe by over temperature and overload protection.
- Malfunction prevented by controller lock function.
- Prevents sample run-off by smooth stirring start and stop function.



MSD-0420

Specification

Model		MSD-0420
Speed range (rpm)		80 to 2000
Viscosity, max. (cP) ¹⁾		5000 (~2000rpm) 10000 (~1200rpm) 30000 (~ 300rpm) 50000 (~ 200rpm)
Stirring capacity (H ₂ O), max. (L / cu ft)		20 / 0.71 (~2000rpm)
Rated torque (N·cm / inch)		19 / 7.48
Moter input / output (W)		70 / 40
Speed control		PID feedback control
Material	Body	Powder coated aluminum
	Cover	Polypropylene
	Motor type	BLDC (Brushless Direct Current)
Dimension	Chuck range (Ø, mm / inch)	3~10 / 0.12~0.39
	Exterior (W x D x H, mm / inch)	75 x 215 x 150 / 2.95 x 8.46 x 5.9
	Net weight (kg / lbs)	2.8 / 6.17
Electrical requirements (230V, 50/60Hz, A)		1.1 Max
Cat. No.		AAH373225K
Electrical requirements (120V/60Hz, A)		0.5
Cat. No.		AAH373225U

1) Tested with silicone oil at 26°C and 60% RH, and uses triple-bladed propeller impeller (Ø50).

Accessories [Page 162](#) Impeller, Stand & Support, Clamp, Magnetic Drive, Stirring Seals

Overhead Stirrer

Easy control type



Easy control model optimized for simple use

Structural Functional Features

- Excellent durability BLDC motor and dedicated controller.
- Slim & compact design.
- Convenient adjustment of height with pass-through impeller shaft.
- Chuck with excellent corrosion resistance. (max. Ø10 mm)

Use Convenience Features

- Accurate stirring with PID feedback control.
- Easy operating for convenient use.
- Displays the set speed and the current speed at the same time.
- Pause buttons make it easier to observe experiments.
- Motor output detection and display.
- Complete with accessories such as impeller, stand, clamp, etc.

Outstanding Safety

- Safety ensured by BLDC motor that does not generate sparks.
- Excellent heat dissipation structure makes use for a long time without any trouble.
- Continual use in safe by over temperature and overload protection.
- Over temperature/overload warning notification displayed.
- Prevents sample run-off by smooth stirring start and stop function.



MSA-0420

Stirrer & Mixer

Specification

Model		MSA-0420
Speed range (rpm)		80 to 2000
Viscosity, max. (cP) ¹⁾		5000 (~2000rpm) 10000 (~1200rpm) 30000 (~ 300rpm) 50000 (~ 200rpm)
Stirring capacity (H ₂ O), max. (L / cu ft)		20 / 0.71 (~2000rpm)
Rated torque (N-cm)		19 / 7.48
Motor input / output (W)		70 / 40
Speed control		PID feedback control
Material	Body	Powder coated aluminum
	Cover	Polypropylene
	Motor type	BLDC (Brushless Direct Current)
Dimension	Chuck range (Ø, mm)	3~10 / 0.12~0.39
	Exterior (W x D x H, mm / inch)	75 x 215 x 150 / 2.95 x 8.46 x 5.9
	Net weight (kg / lbs)	2.8 / 6.17
Electrical requirements (230V, 50/60Hz, A)		1.1 Max.
Cat. No.		AAH373325K
Electrical requirements (120V/60Hz, A)		0.5
Cat. No.		AAH373325U

1) Tested with silicone oil at 26°C and 60% RH, and uses triple-bladed propeller impeller (Ø50).

Accessories [Page 162](#) Impeller, Stand & Support, Clamp, Magnetic Drive, Stirring Seals

Accessories Impeller

Stainless steel impellers

Propeller Impeller

- Suitable for high-speed stirring of low/medium viscosity solutions.
- Excellent mixing properties for homogeneous and suspended solids.
- It can be mounted on an overhead stirrer that rotates clockwise.

3-Bladed

(Unit: mm / inch)

Cat. No.	Blade Ø	Shaft Ø	Length
AAA37521	50 / 1.97	8 / 0.31	400 / 19.69
AAA37522	70 / 2.76	8 / 0.31	400 / 19.69
AAA37523	100 / 3.94	8 / 0.31	400 / 19.69



4-Bladed

(Unit: mm / inch)

Cat. No.	Blade Ø	Shaft Ø	Length
AAA37531	50 / 1.97	8 / 0.31	400 / 15.75
AAA37532	70 / 2.76	8 / 0.31	400 / 15.75
AAA37533	100 / 3.94	8 / 0.31	400 / 15.75



Anchor Impeller

- Suitable for low-speed stirring or reaction of medium/high viscosity solutions.
- Used for polymer reactions or dispersion of large amounts of liquid minerals, etc.
- Select and use an impeller having similar size with the diameter of vessels.

(Unit: mm / inch)

Cat. No.	Blade, W x H	Shaft Ø	Length
AAA37571	45 x 45 / 1.77 x 1.77	8 / 0.31	300 / 11.81
AAA37572	60 x 60 / 2.36 x 2.36	8 / 0.31	500 / 19.69



Paddle Impeller

- Suitable for high-speed stirring of low/medium viscosity solutions.
- Provides smooth flow during heat exchange.

(Unit: mm / inch)

Cat. No.	Blade, W x H	Shaft Ø	Length
AAA37581	70 x 70 / 2.76 x 2.76	8 / 0.31	500 / 19.69

Dissolver Impeller

- Suitable for high/medium-speed stirring of / low/medium (<500 mPas) viscosity solutions.

(Unit: mm / inch)

Cat. No.	Blade Ø	Shaft Ø	Length
BEA0570031	50 / 1.97	8 / 0.31	300 / 11.81
BEA0570032	70 / 2.76	8 / 0.31	300 / 11.81




Turbine Impeller

- Suitable for high-speed stirring of medium/high viscosity solutions.

(Unit: mm / inch)

Cat. No.	Blade Ø	Shaft Ø	Length
AAA37561	57 / 2.24	8 / 0.31	500 / 19.69
AAA37562	90 / 3.54	8 / 0.31	500 / 19.69



Half-Moon Impeller

- Suitable for stirring medium viscosity solutions
- The blades are folded for use in narrow containers.
- The blades are spread by centrifugal force during stirring.

(Unit: mm / inch)

Cat. No.	Blade, W x H	Shaft Ø	Length
AAA37541	65 x 20 / 2.6 x 0.79	8 / 0.31	300 / 11.81
AAA37542	90 x 25 / 3.54 x 0.98	8 / 0.31	500 / 19.69

**Centrifugal Impeller**

- Suitable for high-speed stirring of low/medium viscosity solutions.
- The blades are folded for use in narrow containers.
- The blades are spread by centrifugal force during stirring.

(Unit: mm / inch)

Cat. No.	Blade, W x H	Shaft Ø	Length
AAA37551	90 x 15 / 3.54 x 0.59	8 / 0.31	500 / 19.69
AAA37552	150 x 15 / 5.9 x 0.59	8 / 0.31	500 / 19.69

**PTFE coated impellers****4-Bladed propeller**

(Unit: mm / inch)

Cat. No.	Blade Ø	Shaft Ø	Length
BEA0570081	50 / 1.97	8 / 0.31	500 / 19.69
BEA0570082	70 / 2.76	8 / 0.31	500 / 19.69



Propeller (4-bladed) / Turbine

Turbine

(Unit: mm / inch)

Cat. No.	Blade Ø	Shaft Ø	Length
BEA0570131	70 / 2.76	8 / 0.31	500 / 19.69



Half-moon / Centrifugal

Half-moon


(Unit: mm / inch)

Cat. No.	Blade, W x H	Shaft Ø	Length
BEA0570091	60 x 18 / 2.36 x 0.71	8 / 0.31	500 / 19.69

Centrifugal

(Unit: mm / inch)

Cat. No.	Blade, W x H	Shaft Ø	Length
BEA0570101	76 x 17 / 2.99 x 0.67	8 / 0.31	500 / 19.69




Anchor (Square / Round)

Anchor (Square¹⁾ / Round²)

(Unit: mm / inch)

Cat. No.	Blade, W x H	Shaft Ø	Length
BEA0570111 ¹⁾	80 x 50 / 3.15 x 1.97	8 / 0.31	500 / 19.69
BEA0570121 ²⁾	80 x 40 / 3.15 x 1.57	8 / 0.31	500 / 19.69



Paddle / Paddle A

Paddle (Paddle¹⁾ / Paddle A²)

(Unit: mm / inch)

Cat. No.	Blade Ø / W x H	Shaft Ø	Length
BEA0570141 ¹⁾	70 / 2.76	8 / 0.31	500 / 19.69
BEA0570151 ²⁾	78 x 80 / 3.07 x 3.15	8 / 0.31	500 / 19.69


Accessories for Overhead Stirrer

Stand & Support



Easy-to-use Dial Stand

- Fine adjustment of position up, down, back and forth.
- Possible to lock after finishing fine adjustment using dials.
- The structure allows for no additional clamp holder to be used.



Basic Stand 1



Basic Stand 2



Basic Stand 3

Stand

Cat. No.	Description (W x D x H, mm / inch)
AAA37501	Dial stand (400 x 400 x 700 / 15.75 x 15.75 x 27.56)
AAA37502	Basic stand 1 (400 x 400 x 700 / 15.75 x 15.75 x 27.56)
AAA37503	Basic stand 2 (300 x 400 x 700 / 11.81 x 15.75 x 27.56)
AAA37504	Basic stand 3 (250 x 350 x 700 / 9.84 x 13.78 x 27.56)

Support

Cat. No.	Description (Ø x L, mm / inch)
AAA37505	Support Rod, M14 (Ø15 x 500mm / Ø0.59 x 19.68)
AAA37506	Support Rod, M14 (Ø18 x 500mm / Ø0.71 x 19.68)
AAA37507	Support Rod, M14 (Ø20 x 500mm / Ø0.79 x 19.68)

Clamp & Seal



Utility Clamp 3 Prong



Utility Clamp Plat



3 Prong Clamp



2 Prong Clamp



Fixing Clamp



Flexible Coupling



PTFE Stirring Seal

Cat. No.	Description
BEA1000001	Utility clamp 3 prong
BEA1000002	Utility clamp plat
BEA1000011	3 Prong clamp (80 mm grip)
BEA1000012	3 Prong clamp (60 mm grip)
BEA1000013	3 Prong clamp (20 mm grip)
BEA1000014	2 Prong clamp (60 mm grip)
BEA0570181	Fixing clamp
BEA0570191	Flexible coupling
BEA0570161	PTFE stirring seals (24/40) with 8mm hole
BEA0570162	PTFE stirring seals (29/42) with 8mm hole
BEA0570163	PTFE stirring seals (34/45) with 8mm hole

Clamp Holder Page 153



C-20



C-10

- Made of aluminum die casting for light weight but strong fixing.
- Knob structure that is comfortable to handle when locking / unlocking.
- The contact surface with the stand support is stable and does not slip.
- C-20 and C-10 models are suitable for the overhead stirrer.

Magnetic Drive



MD-T6



- Mixing reaction samples in case pressurized or depressurized in container, processing of resin solvents, decomposition and synthesis using enzymes, and uniform mixing of liquid raw materials.
- Experiments that require large stirring force while maintaining full sealing.
- Connected with overhead stirrer.

- The outer magnet and inner magnet are rare-earth element, which form a strong magnetic force to achieve high torque. (rotational force)
- Low vibration and high speed rotation even in a sealed state.
- Includes cooling line base for easy connection with external cooling device. (for experiments over 70°C)
- The connection point to the ground joint of glassware is easy to separate from glassware in a sealed state.

Specification

Model		MD-T6-24	MD-T6-29	MD-T6-35
Material	Body	SS316 / Cr-Plate	SS316 / Cr-Plate	SS316 / Cr-Plate
	Bushing & seal	PTFE & viton	PTFE & viton	PTFE & viton
Max. torque (kgf.cm)		6	6	6
Driving	Pressure	1 x 10 ⁻⁴ mmHg to 5kg/cm ²	1 x 10 ⁻⁴ mmHg to 5kg/cm ²	1 x 10 ⁻⁴ mmHg to 5kg/cm ²
	Temp. (with cooling)	Max. 300°C / 572°F	Max. 300°C / 572°F	Max. 300°C / 572°F
Shaft size (Ø x H, mm / inch)		8 x 209 / 0.31 x 8.23	8 x 209 / 0.31 x 8.23	8 x 209 / 0.31 x 8.23
Ground joint		24 / 40	29 / 42	34 / 45
Application		Glassware	Glassware	Glassware
Cat. No.		BEA511001	BEA511002	BEA511003



Vortex Mixer

More convenient and useful functions Vortex Mixer Developed with Jeio Tech's Proprietary Technology

> **Touch vortex function maximizes convenience**

No matter where you touch the top of the platform, it's very easy to use with fast vortex formation.

> **Patented touch switch method**

Stable and sensitive Touch vortex using touch switch connected to top plate.
(Patented in Korea)

> **Low height to minimize burden on the wrist**

The structure minimizes the burden on the wrist during repeated use due to its low height.

> **Structure to maintain stability at high speeds**

The low-profile design and rubber feet make for slip-free operation of the equipment.

> **Multiple tubes can be combined simultaneously**

Multiple vortex tubes can be used simultaneously to reduce the required time for efficient testing.

> **Offers use environment of diverse tube shapes and capacities**

Most tubes such as centrifuge tubes, conical tubes, and test tubes can be used in Touch mode.

> **Powerful and fast vortex formation**

Runs up to 3,000 rpm for reliable and strong sample mixing.

> **Proven structural excellence**

BLDC motor with excellent durability.
Structure allows for low noise and low vibration even at 3,000rpm.



Rapid vortex formation even touching anywhere on top plate.



Available to handle the several tubes at the same time.



Available in various shapes.



Diverse vortex mixing in one device with optimized configuration can be used for various experiments

> Variety of convenient accessories

Tube holder accessories up to Ø25mm.
Convenient use with easily interchangeable structure.

> Includes microplate shaking function

Dedicated tray and stable speed adjustment function for optimum microplate mixing.

> Intuitive and convenient operation

Selects operation mode easily and intuitive speed control knobs.

> Safety features for protection of equipment

In case excessive pressure is applied to the top plate, operation is temporarily stopped to protect the equipment.

> Selectable operation mode according to application

Rapid vortex formation with touch mode when needed.
Run-L/H mode provides stable continuous mixing environment.

- ① **Touch mode (touch operation, Max. 3,000rpm)**
Rapid vortex formation with touch mode when needed.
Run-L/H mode provides stable continuous mixing environment.
- ② **Run-L mode (continuous operation, Max. 1,500rpm)**
Low speed mode to prevent contamination between samples.
Stable and continuous mixing. (except for VM-T model)
- ③ **Run-H mode (continuous operation, Max. 3,000rpm)**
Strong and continuous vortex formation up to 3,000rpm with excellent durability.

※ Some of the above contents are limited to specific models.



Various accessories provided according to the sample.



Easy attachment and detachment.



Low-profile design and non-slip feet for stability at high speeds.

Vortex Mixer



True vortex mixing with innovative touch



VM-96T (Persian Blue)



VM-96A




Specification

Model		VM-96A	VM-96T
Control	Speed range (rpm)	Max. 3000	Max. 3000
	Speed control	Scale	Scale
	Operating mode	Touch, Run-L, Run-H	Touch, Run
	Motion	Orbital	Orbital
	Orbit (dia., mm / inch)	4 / 0.16	4 / 0.16
Material	Platform	Silicone	Silicone
	Body	Polypropylene	Polypropylene
	Motor type	BLDC (Brushless Direct Current)	BLDC (Brushless Direct Current)
Dimensions	Platform (W x D, mm / inch)	76 x 76 / 3 x 3	76 x 76 / 3 x 3
	Exterior (W x D x H, mm / inch)	154 x 210 x 84 / 6.1 x 8.3 x 3.3	148 x 159 x 77 / 5.8 x 6.3 x 3
	Net weight (kg / lbs)	3.0 / 6.6	2.6 / 5.7
Protection class (DIN EN 60529)		IP 42	-
Electrical requirements		AC 100V to 240V, 50/60Hz, 3.3A Max.	
Cat. No.	KR plug	AAH360215K	AAH361215BK
	US plug	AAH360215U	AAH361215BU

Accessories for Vortex Mixer



Platform (pop-off)



Platform (standard)
Included as standard on equipment purchase



Microplate Tray



Cat. No.	Description
00VMS0000063	Platform (Pop-off)
00VMS0000058	Platform (standard)
00VMS000008JT	Microplate tray
AAA36520	Tube holder (up to $\Phi 50$ mm)
AAA36521	Tube holder ($\Phi 10$ mm x 16ea)
AAA36522	Tube holder ($\Phi 12$ mm x 9ea)
AAA36523	Tube holder ($\Phi 15$ mm x 9ea)
AAA36524	Tube holder ($\Phi 20$ mm x 5ea)
AAA36525	Tube holder ($\Phi 25$ mm x 4ea)
AAA36526	Tube holder ($\Phi 8$ mm x 18ea & $\Phi 10$ mm x 28ea)



Life Science Equipment



Satisfaction with the requirements of European regulations and guidelines.



Real-time equipment monitoring and control system using mobile app.



Compliance with Guidelines for Restriction of Hazardous Substances Use.



Orbital motion either clockwise or counterclockwise direction can be set.



Control of temperature / humidity / rpm / time, etc. through dedicated program.



More precise temperature control through temperature calibration.



Over temperature protection function.



PC communication via RS-232 / RS-485 / USB port.



Possible to set end time or start time of device operation.



2 year warranty free A/S.

General Application

Heating Block: Required in culture experiments such as enrichment, enzyme reaction, immunoassay, etc.

Electrophoresis System: Electrophoresis of nucleic acid and protein.

Power Supply for Electrophoresis System: Electrophoresis of nucleic acids and proteins.

	Description	Note	Model	Page	
Preparation Reaction Detection	Heating Block	<ul style="list-style-type: none"> Precise feedback temperature control up to 130°C Selection of aluminum blocks for various tubes 	Amb. +5 to 130°C / Amb +41 to 266°F	CHB	172
	Heating & Cooling Block	<ul style="list-style-type: none"> Convenient controlling for 10 programs & 10 steps Temperature control in 0.1°C and timer function included 	Amb. -20 to 95°C / Amb. 68 to 203°F	CCB	174
	Heating Shaker	<ul style="list-style-type: none"> High performance heating shaker at 100°C and 1500rpm Three specialized Interval Mix programs provided 	Max. 100°C Max. 1500rpm	CBS	175
	Horizontal Electrophoresis System	<ul style="list-style-type: none"> Integral structure solid, clean, and easy to maintain UV-transparent gel tray for convenient checking of results 	Max 198 samples per run	EP	178
	Vertical Electrophoresis System	<ul style="list-style-type: none"> Perform both electrophoresis and membrane transfer in the same tank Electrophoresis of up to 4 gels at once 	Max 4 gels per run	VEP	180
	Power Supply for Electrophoresis System	<ul style="list-style-type: none"> Safe power supply with CE certificate Four-set simultaneous output and stackable structure 	Constant Voltage Constant Current	EPS	182

※ The contents of the above and the contents of this catalog may differ depending on the specific model and conditions of use. For the information about the features and specifications that applying to each models, please check the information on the corresponding page of each models.

Vortex Mixer

Page 168



Biosafety Cabinet

Page 196



Incubated Shaker

Page 116



PCR Workstation

Page 203



Ice Rack

Page 276



Heating Block



Precise feedback temperature control up to 130°C

Structural Functional Features

- Optimal heat transfer is achieved by the tight coupling design of the main body.
- Transparent and durable polycarbonate lid protects the sample and provides improved temperature uniformity.
- Compact design makes it easy to move and use space.
- Polypropylene (PP) material with excellent chemical resistance.
- Two extractors for easy removal block are provided as standard.
- Aluminum blocks with excellent thermal conductivity. (option)
- Two models with different cover heights are provided depending on the block (tube size) used.
- Its built-in temperature limit setting feature (with max. 0.2°C overshooting) allows you to perform highly temperature-sensitive reactions such as isothermal amplifications.




CHB-350S

with Exchangeable Block
(option)

Use Convenience Features

- The temperature can be set accurately in units of 0.1°C.
- Feedback control enables precise temperature control.
- Calibration function minimizes temperature deviation.
- Easy to operate with clear VFD display, touch keys and control knobs.
- Provides operating elapsed time check.
- Count-down timer provided. (up to 99 hours and 59 minutes)
- User can choose between CHB-350S and CHB-350T depending on the type of tubes and blocks.



CHB-350T

with Exchangeable Block 2ea
(option)

Outstanding Safety

- Structure with minimal inflow of reagent or sample.
- Convenient and safe operation by notifying the status of heating condition.
- Overheating temperature prevention function. (Showing the error message if the temperature of the main body exceeds 140°C)
- Over-current prevention function. (Automatic power cutoff if the temperature of the main body exceeds 150°C)



with Extractor 2ea
(standard)



Selects operation mode easily and intuitive speed control knobs.



Polycarbonate lid ensures excellent durability and convenient sample checking.



Provide various kinds of heating block.

Specification

Model	CHB-350S	CHB-350T
Temperature		
Range (°C / °F)	Amb. +5 to 130 / Amb +41 to 266	Amb. +5 to 130 / Amb +41 to 266
Fluctuation at 80°C (±°C / ±°F)	0.1 / 0.18	0.1 / 0.18
Variation at 80°C (±°C / ±°F)	0.5 / 0.9	0.5 / 0.9
Display	Up to 130°C, 0.1°C resolution	Up to 130°C, 0.1°C resolution
Heat up time Ambient to 130°C (min)	15	15
Control	Feedback control	Feedback control
Timer	1min to 99hrs 59min	1min to 99hrs 59min
Safety device	Over temperature protection, Over current protection	
Material		
Body	PP, Powder coated steel	PP, Powder coated steel
Cover	PC	PC
Block	Black hard anodized Aluminum	Black hard anodized Aluminum
Dimensions		
Exterior (W x D x H, mm / inch)	249 x 330 x 125 / 9.8 x 12.99 x 4.92	249 x 330 x 250 / 9.8 x 12.99 x 9.84
Net weight (kg / lbs)	3.9 / 8.6	4.3 / 9.48
Electrical Requirements (230V, 50/60Hz, A)	3.5	3.5
Cat. No.	AAHJ1015K	AAHJ1115K
Electrical Requirements (120V, 60Hz, A)	6.7	6.7
Cat. No.	AAHJ1013U	AAHJ1113U

Accessories [Page 162](#) Heating Block

Heating Block for Round Flask and Vial

Detail [Page 146](#)

- ▶ The heating block is mounted on a hotplate or hotplate magnetic stirrer to heat safely (up to 300°C) the contents of the 25ml - 2000ml round flask or 1 dram ~ 40 ml vial.
- ▶ A safe solution to replace oil bath and heating mantle.
- ▶ Economical and eco-friendly by eliminating the need to purchase and dispose of oil separately.



Heating & Cooling Block



Convenient use thanks to easy temperature control with 10 programs and 10 steps

Structural Functional Features

- 10 Programs & 10 steps program control.
- Program can be set to 99 repetitions or unlimited repetitions.
- Optimal heat transfer is achieved by the tight coupling design of the main body.
- Transparent and durable polycarbonate lid protects the sample and provides improved temperature uniformity.
- Compact design makes it easy to move and use space.
- Polypropylene (PP) material with excellent chemical resistance.
- Two extractors for easy removal block are provided as standard.
- Aluminum blocks with excellent thermal conductivity. (option)

Use Convenience Features

- The temperature can be set accurately in units of 0.1°C.
- Feedback control enables precise temperature control.
- Calibration function minimizes temperature deviation.
- Easy to operate with clear VFD display, touch keys and control knobs.
- Provides operating elapsed time check.
- Count-down timer provided. (up to 99 hours and 59 minutes)

Outstanding Safety

- Structure with minimal inflow of reagent or sample.
- Convenient and safe operation by notifying the status of heating condition.
- Overheating temperature prevention function. (heating panel, control unit)
- Over-current prevention function.



CCB-350
with Exchangeable Block (option)



Specification

Model	CCB-350
Temperature	
Range (°C / °F)	Amb. -20 to 95 / Amb. -36 to 203
Fluctuation at 37°C (±°C / ±°F)	0.1 / 0.18
Variation at 37°C (±°C / ±°F)	0.5 / 0.9
Display	Up to 95°C, 0.1°C resolution
Control	Feedback control
Timer	1min to 99hrs 59min
Safety device	Over temperature protection, Over current protection
Material	
Body	PP, Powder coated steel
Cover	PC
Block	Black hard anodized Aluminum
Dimensions	
Exterior (W x D x H, mm / inch)	249 x 330 x 168 / 9.8 x 12.99 x 6.61
Net weight (kg / lbs)	5 / 11.02
Electrical Requirements (230V, 50/60Hz, A)	2
Cat. No.	AAHJ5015K
Electrical Requirements (120V, 60Hz, A)	4
Cat. No.	AAHJ5013U

Heating Shaker



Up to 100°C, 1500rpm, and program control

Structural Functional Features

- Three specialized interval mix programs are provided. (45/55/65°C & 1,000rpm, 1 minute operation/stop)
- 10-step temperature and shaking control set.
- Program can be set to 99 repetitions or infinite repetitions.
- Optimal heat transfer is achieved by the tight coupling design of the main body.
- Compact design makes it easy to move and use space.
- Polypropylene (PP) material with excellent chemical resistance.
- Improved durability with integrated shaking module. (registration KR 10-1429256)
- Aluminum block with excellent thermal conductivity and block cover (option) made of polypropylene.



CBS-350
with Block (option)

Use Convenience Features

- The temperature can be set accurately in units of 0.1°C.
- Feedback control enables precise temperature control.
- Calibration function minimizes temperature deviation.
- Protecting the sample by starting with soft shaking.
- Easy to operate with clear VFD display, touch keys and control knobs.
- Provides operating elapsed time check.
- Count-down timer provided. (up to 99 hours and 59 minutes)

Outstanding Safety

- Structure with minimal inflow of reagent or sample.
- Convenient and safe operation by notifying the status of heating condition.
- Overheating temperature prevention function. (heating panel, control unit)
- Over-current prevention function.

Specification

Model		CBS-350
Temperature		
Range (°C / °F)		Amb. +5 to 100 / Amb. +41 to 212
Fluctuation at 80°C (±°C / ±°F)		0.3 / 32.54
Variation at 80°C (±°C / ±°F)		0.7 / 33.26
Display		Up to 100°C, 0.1°C resolution
Control		Feedback control
Shaking System		
Motion type		Orbital
Speed range (rpm)	96-well tube block	150 to 1500
	0.5 ml tube block	150 to 1000
	1.5 ml/2.0 ml tube block	150 to 1000
	15 ml tube block	150 to 900
	Ø12 tube block	150 to 850
	Ø13 tube block	150 to 850
Amplitude, Orbital size (mm, dia.)		2
Timer		1min to 99hr 59min
Safety device		Over temperature protection, Over current protection
Material		
Body		PP, Powder coated steel
Cover		PP
Block		Black hard anodized Aluminum
Dimensions		
Exterior (W x D x H, mm / inch)		249 x 325 x 120 / 9.8 x 12.79 x 4.72
Net weight (kg / lbs)		8.25 / 18.19
Electrical Requirements (230V, 50/60Hz, A)		1.7
Cat. No.		AAHJ4015K
Electrical Requirements (120V, 60Hz, A)		3.4
Cat. No.		AAHJ4013U

Accessories **Page 163** Heating Block

Accessories

Accessories for Heating Block, Heating & Cooling Block

Image	Model	Description	Max. Mountable Tube	Dimensions (W x D x H, mm / inch)	Mountable Capacity of Blocks		
					CHB-350S	CHB-350T	CCB-350
	10CHB0000028	1.5 / 2.0 ml Tube block (micro-tube)	48	153 x 98 x 41 / 6.02 x 3.86 x 1.61	1	1	-
	10CHB0000045	1.5 / 2.0 ml Tube block (micro-tube)	30	98 x 76.5 x 41 / 3.86 x 3.01 x 1.61	2	2	1
	10CHB0000029	0.5 ml Tube block (micro-tube)	48	98 x 76.5 x 41 / 3.86 x 3.01 x 1.61	2	2	1
	10CHB0000030	0.2 ml Tube block (micro tube / strips / 96-well plate)	96	153 x 98 x 41 / 6.02 x 3.86 x 1.61	1	1	-
	10CHB0000031	15 ml Tube block (conical tube)	15	98 x 76.5 x 51 / 3.86 x 3.01 x 2.01	2*	2	1*
	10CHB0000032	50 ml Tube block (conical tube)	6	98 x 76.5 x 51 / 3.86 x 3.01 x 2.01	2*	2	1*
	10CHB0000033	50 ml Tube block (conical tube)	6	98 x 76.5 x 87 / 3.86 x 3.01 x 3.42	2*	2	1*
	10CHB0000034	Φ10 mm Tube block	35	98 x 76.5 x 51 / 3.86 x 3.01 x 2.01	2*	2	1*
	10CHB0000035	Φ12 mm Tube block	24	98 x 76.5 x 51 / 3.86 x 3.01 x 2.01	2*	2	1*
	10CHB0000036	Φ13 mm Tube block	24	98 x 76.5 x 51 / 3.86 x 3.01 x 2.01	2*	2	1*
	10CHB0000037	Φ15 mm Tube block	20	98 x 76.5 x 51 / 3.86 x 3.01 x 2.01	2*	2	1*
	10CHB0000038	Φ16 mm Tube block	16	98 x 76.5 x 51 / 3.86 x 3.01 x 2.01	2*	2	1*
	10CHB0000039	Φ18 mm Tube block	12	98 x 76.5 x 51 / 3.86 x 3.01 x 2.01	2*	2	1*
	100CHB0000040	Φ20 mm Tube block	12	98 x 76.5 x 51 / 3.86 x 3.01 x 2.01	2*	2	1*

※ It can be equipped when the lid is opened.

Accessories for Heating Shaker

Image	Model	Description	Max. Mountable Tube	Dimensions (W x D x H, mm / inch)	Mountable Capacity of Blocks	
					CBS-350	
	AAAJ4501	0.5 ml Tube block with block cover	48	98 x 76.5 x 30 / 3.86 x 3.01 x 1.18	1	
	AAAJ4502	1.5 / 2.0 ml Tube block with block cover	24	98 x 76.5 x 41 / 3.86 x 3.01 x 1.61	1	
	AAAJ4503	50 ml Tube block with block cover	6	98 x 76.5 x 51 / 3.86 x 3.01 x 2.01	1	
	AAAJ4504	0.2 ml Tube block with block cover	96	111 x 76.5 x 15 / 4.37 x 3.01 x 0.59	1	
	AAAJ4505	15 ml Tube block with block cover	15	98 x 76.5 x 51 / 3.86 x 3.01 x 2.01	1	
	AAAJ4506	Φ12 mm Tube block with block cover	24	98 x 76.5 x 51 / 3.86 x 3.01 x 2.01	1	
	AAAJ4507	Φ13 mm Tube block with block cover	24	98 x 76.5 x 87 / 3.86 x 3.01 x 3.42	1	


※ Both block cover for AAAJ4506 and AAAJ4507 is the same specification.

Horizontal Electrophoresis System

Differentiated ease of use electrophoresis system

Structural Functional Features

- Utilization of Nucleic Acid Experiments with Relatively Large Molecule Samples.
- Integral structure of transparent polycarbonate material ensures excellent durability and convenient sample checking.
- Cathode/anode and direction marking on the bottom of the tank makes it very convenient to use for checking the direction of electrophoresis.
- Easy to move with handles on the sides of the tank.
- There are rubber feet on the bottom of the tank and gel caster to prevent slip.
- Providing 31 well and 16 well Comb at the same interval as a multi-channel pipette, enabling fast loading of large quantities of samples. (EP-10, EP-18)




EP-10

with Power Cord (standard)

Use Convenience Features

- The gel tray is UV-transparent, and no need of separation for the Agarose Gel from the tray to check the results. (UV-A)
- Level marker is engraved on the gel tray to make it easy to make gel of constant thickness.
- Level meter included in gel caster for easy to uniformity of gel generating. (EP-10, EP-18)
- Power cord being connected/disconnected easily are divided into cathode/anode by black/red.
- Large quantities of samples can be loaded at once. (based on EP-18, total 204EA (including marker))
- Electrodes are easy to remove from tank and it's easy to clean the tank.




Comb spacing fitting to Multichannel pipettes (31, 16 well Comb).



Level Marker.



Gel caster level meter.



Anti-slip rubber feet.



Connection to Power Supply



EP-05



EP-10



EP-18

Model	EP-05	EP-10	EP-18
Electrophoresis Unit (Tank & Lid, Electrode, Power Cord)	1 set	1 set	1 set
Gel Caster	1 ea	1 ea	1 ea
Gel Tray	3 ea	2 ea (Up to 2 combs can be installed per tray)	1 ea (Up to 6 combs can be installed on tray)
Max. Sample Throughput	39 ea (Marker 1 ea)	132 ea (Marker 4 ea)	198 ea (Marker 6 ea)
Double-sided Comb Number of upper well / Number of lower well (thickness x quantity provided)	27/40 well* (1.0 mm x 1ea) (40 well: 15 + 15 + 10 27 well: 10 + 10 + 7 See the image on the bottom)	20/31 well* (1.0 mm) x 2ea 20/31 well* (1.5 mm) x 2ea 16/20 well* (1.0 mm) x 2ea 16/20 well* (1.5 mm) x 2ea 25/34 well* (1.0 mm) x 2ea 25/34 well* (1.5 mm) x 2ea	20/31 well* (1.0 mm) x 2ea 20/31 well* (1.5 mm) x 2ea 16/20 well* (1.0 mm) x 2ea 16/20 well* (1.5 mm) x 2ea 25/34 well* (1.0 mm) x 2ea 25/34 well* (1.5 mm) x 2ea

* Number of wells depending on the direction of use.

Comb Specification

Comb	27/40 well		20/31 well		16/20 well		25/34 well	
Well (ea)	27	40	20	31*	16*	20	25	34
Width (mm / inch)	4.8 / 0.19	3.0 / 0.11	4.8 / 0.19	3.0 / 0.11	6.0 / 0.24	4.8 / 0.19	4.5 / 0.18	3.0 / 0.11
Volume (µl)	30	20	48 / 72**	30 / 45**	60 / 90*	48 / 72**	45 / 68**	30 / 45**

* Equal spacing with multichannel pipettes makes it convenient for loading large quantities of samples.

** Volume per well is 1.0mm/1.5mm depending on the thickness of the comb

Specification

Model	EP-05	EP-10	EP-18	
Buffer capacity (L / cu ft)	0.5 / 0.017	1 / 0.035	1.8 / 0.063	
Material	Electrophoresis tank, Lid	PC (Polycarbonate)		
	Electrodes	Platinum wire		
Dimensions (W x D x H, mm/inch)	Electrophoresis tank	219 x 134 x 71 / 8.62 x 5.27 x 2.79	219 x 214 x 71 / 8.62 x 8.42 x 2.79	387 x 214 x 71 / 15.24 x 8.42 x 2.79
	Gel caster	254 x 87 x 18 / 10 x 3.42 x 0.79	167 x 264 x 30 / 6.57 x 10.39 x 1.18	167 x 264 x 30 / 6.57 x 10.39 x 1.18
	Gel tray	81 x 82 x 17 / 3.19 x 3.23 x 0.67 58 x 82 x 17 / 2.28 x 3.23 x 0.67	181 x 103 x 30 / 7.12 x 4.05 x 1.18	181 x 206 x 30 / 7.12 x 8.11 x 1.18
	Comb	225 x 45 / 8.86 x 1.77	158 x 46 / 6.22 x 1.81	158 x 46 / 6.22 x 1.81
	Electrophoresis unit weight (kg / lbs)	0.39 / 0.86	0.6 / 1.32	0.93 / 2.05
Cat. No.	AAAJ6011	AAAJ6021	AAAJ6031	

Accessories **Page 167** Electrophoresis Unit, Gel Caster Unit, Electrode Assembly, Comb

Vertical Electrophoresis System


Both Electrophoresis and Trans Blot in the Same Tank

Structural Functional Features

- Suitable for protein experiments which is relatively small molecule.
- Run up to 4 gels simultaneously.
- Compatible with electrophoresis module components from other brand.
- Stable with non-slip rubber feet.

User-Friendly Features

- Allows both electrophoresis and Trans blot in the same tank.
- Unique design for convenient, leak-free gel casting.
- Color-coded by comb thickness for easy identification.



VEP-TB

Product Configuration

Electrophoresis system (VEP-4G)



Tank & Lid with cables 1ea



Electrode assembly & running module 1ea



Gel Caster with comb (4ea), plate (4set), gasket (4ea)



Buffer dam 1ea

Electrophoresis system with trans blot module set (VEP-TB)


Standard Electrophoresis System (VEP-4G) with an added Trans blot module



Trans blot electrode assembly 1ea



Gel holder cassettes 2ea



Foam pad 4ea



Cooling unit 1ea



Gel press 1ea



Gel releaser 1ea



Selection Guide for Comb & Plate Set

Choose the appropriate comb thickness and number of wells based on the sample quantity and volume for analysis.

Comb thickness		1.5 mm	1.0 mm	0.75 mm
Well number	10-well	66 μ l	44 μ l	33 μ l
	11-well	60 μ l	40 μ l	30 μ l
	13-well	46 μ l	30 μ l	23 μ l
	15-well	40 μ l	26 μ l	20 μ l

Convenient, Leak-Free Gel Casting



Gel caster positioning
(31, 16 well comb)



Wrap the plate with a gasket



Attach the wrapped plate to the caster



Pour the gel and insert the comb

Specification

Model		VEP-TB	VEP-4G
Component		VERTICAL ELECTROPHORESIS System & Trans blot module set	VERTICAL ELECTROPHORESIS System
Material	Electrophoresis tank, Lid	PC	PC
	Plate	Glass	Glass
	Electrode assembly	PP	PP
	Gel gasket	Silicone	Silicone
	Trans blot electrode assembly	PP	-
Dimensions (W x D x H, mm)	Electrophoresis tank	172 x 130 x 127	
	Gel caster	115 x 75 x 98	
	Comb	83 x 21.5 x 2.1	
Max. number of gel running per tank		4ea	
Glass plate size (W x L, mm)		Short: 100 x 73 / Spacer: 100 x 83	
Total Buffer volume, 2gel / 4gel		2gel: 700 ml / 4gel: 1000ml	
Recommended power supply		EPS-300, 2~250V range	
Cat. No.		BEA30307001	BEA30307002

Power Supply for Electrophoresis System



Multipurpose Power Supply can be used Horizontal agarose gel, Vertical polyacrylamide gel running, various blotting experiments, etc.

Structural Functional Features

- Power supply optimized for electrophoresis system.
- Constant voltage and constant current, convenient use as two operation modes provided.
- It consists of four sets of output terminals, so four electrophoresis devices can be used at the same time.
- Space-efficient use with stackable design.
- Folding legs on the lower part make it convenient to operate and check the status.



EPS-300

Use Convenience Features

- Accurate current and voltage by microprocessor.
- Pause and restart functions.
- Count-down timer. (up to 99 hours and 59 minutes)
- Displaying the operating time.
- Possible to monitor experimental conditions such as voltage, current value, change the setting value, stopping and reactivating. (when purchased LC GreenBox)
- PC communication via RS-232 Port.



Safe use with CE certification for electrical safety.



Compatible with other companies' products and can be connected to up to 4 devices at the same time.

Outstanding Safety

- Rapid resistance change detection and detection of current leakage to prevent accidents by cutting power.
- Protection against over-current/over-voltage/overpower.
- Over-temperature protection.
- Made of nonflammable plastic being safe in a fire.



Efficient time management with a mobile monitoring system. (option)



Space efficient due to a structure that can be stacked in two or more layers without slipping.

Specification

Model	EPS-300
Output range	5~250V, in 1V increments 10mA~3000mA, in 10mA increments (Max. 300W)
Type of output	Constant voltage mode, constant current mode
Output terminals	4 pair recessed banana jack in parallel
Timer	1min to 99hr 59min
Pause/resume function	Yes
Display	7-segment 4 digit LED, 2 line
Operating conditions	0 to 40°C, 0 to 80% RH
Dimensions (W x D x H, mm / inch)	240 x 310 x 100 / 9.45 x 12.2 x 3.94
Weight (kg / lbs)	2.5 / 5.51
Electrical requirement (VAC, Hz)	100-240VAC, 50/60Hz
Current consumption (A)	1.5 ~ 3.6
Cat. No.	AAAJ6541

Accessories Page 167 LC GreenBox

Accessories



Electrophoresis Unit




Gel Caster Unit




Electrode Assembly




LC GreenBox (Mobile monitoring system)



Comb 27/40 well



Comb 20/31 well



Comb 16/20 well



Comb 25/34 well

Cat. No.	Description	Model
AAAJ6501	Electrophoresis unit (tank, lid, electrode assembly)	EP-05
AAAJ6502	Electrophoresis unit (tank, lid, electrode assembly)	EP-10
AAAJ6503	Electrophoresis unit (tank, lid, electrode assembly)	EP-18
AAAJ6511	Gel caster unit (gel caster, gel tray)	EP-05
AAAJ6512	Gel caster unit (gel caster, gel tray)	EP-10
AAAJ6513	Gel caster unit (gel caster, gel tray)	EP-18
AAAJ6521	Electrode assembly	EP-05
AAAJ6522	Electrode assembly	EP-10
AAAJ6523	Electrode assembly	EP-18
AAAJ6531	Comb 27/40 well, thickness 1 mm - (5ea)	EP-05
AAAJ6532	Comb 20/31 well, thickness 1 mm - (5ea)	EP-10/18
AAAJ6533	Comb 20/31 well, thickness 1.5 mm - (5ea)	EP-10/18
AAAJ6534	Comb 16/20 well, thickness 1 mm - (5ea)	EP-10/18
AAAJ6535	Comb 16/20 well, thickness 1.5 mm - (5ea)	EP-10/18
AAAJ6536	Comb 25/34 well, thickness 1 mm - (5ea)	EP-10/18
AAAJ6537	Comb 25/34 well, thickness 1.5 mm - (5ea)	EP-10/18
AAHQ1011K	LC GreenBox	EPS-300

Autoclave



Over temperature protection function.



2 year warranty free A/S.



Control of temperature / humidity / rpm / time, etc. through dedicated program.

General Application

Sterilization of experimental equipment/ materials, growth media, and glassware using saturated steam

Easy to use with optimized program function for each application.
Sterilization process monitoring ensures reliable sterilization.

> Optimized usage mode provided

Three (General, Standard, Program) modes are provided to allow a choice of sterilization processes optimized for the application.

General mode: Operate by setting sterilization conditions manually.

Standard mode: Select and operate four representative actions.

Program mode: The mainly used conditions are saved and used.

> Reliable sterilization performance monitoring


The device itself monitors the sterilization process to perform more reliable sterilization.

Low Heat: When the set temperature is not reached within the designated time.

Over/Low Temp.: When the set temperature is deviated from set value during the time set.

Sterilization Fail: When sterilization is incomplete due to power failure, and etc.

Operational mode	Temp. (°C)	Time (min)
General mode	110 to 123	1 ~ 999
Standard mode	Decontamination	121
	Solid sterilization	121
	Liquid sterilization	121
	Melting	100
Program mode	Sterilization 1	110 to 123
	Sterilization 2	110 to 123
	Liquid sterilization	110 to 123
	Melting	60 to 100



High-capacity, top loading method.



Displays information on the progress of the sterilization process and operation of the equipment.



Pressure gauge to check pressure inside of chamber.

Autoclave

Suitable for temperature characteristic tests such as heat treatment and drying

Proven Control Capability

- Verified according to KS P 6102, providing a more reliable sterilization environment.
- Verified stable temperature control performance by measuring the inside with 10 of temperature sensor.

Structural Functional Features


- Easy to open and close with circular handle.
- Including pressure gauge for internal pressure verification.
- Exhaust tank provided as standard.
- Including drain bucket for use cleanly.
- Upper part is made of stainless steel and is easy to maintain.
- Vertical, top loading system makes it easy to sterilize large volume samples.
- Casters for easy mobility during installation or relocation.

Use Convenience Features

- Easy operation environment allowing the whole sterilization process to be started with just one operation.
- Three operation modes provided according to usage.
- Precise temperature control with PID Feedback control and auto tuning.
- Automatic exhaust after the sterilization process.
- LED display lamp, clear VFD and sterilization step to check operation status and sterilization progress conveniently.
- Including sterilization elapsed time display.
- Additional temperature sensor and recorder for recording sterilization procedure. (option)
- Pressure for each step can be checked with pressure sensor. (option)
- Including manual valve for releasing internal pressure.

Outstanding Safety

- Preventing users from get burnt by using safety door cover made of plastic with low thermal conductivity.
- Automatic operation of safety valve when overpressure occurs.
- Over temperature protection. (over temperature limiter)
- Low level of solution alarm. (low level limiter)
- Overcurrent and short circuit protection.




Max. Loading Capacity

Media Bottle				
Bottle Capacity (ml)	ST-50G (ea x floor)	ST-65G (ea x floor)	ST-85G (ea x floor)	ST-105G (ea x floor)
250	20 x 3	20 x 3	20 x 5	20 x 6
500	14 x 2	14 x 3	14 x 4	14 x 4
1000	8 x 1	8 x 2	8 x 3	8 x 3
2000	5 x 1	5 x 1	5 x 2	5 x 3

Erlenmeyer Flask				
Bottle Capacity (ml)	ST-50G (ea x floor)	ST-65G (ea x floor)	ST-85G (ea x floor)	ST-105G (ea x floor)
250	14 x 3	14 x 4	14 x 5	14 x 6
500	8 x 2	8 x 3	8 x 4	8 x 4
1000	5 x 2	5 x 2	5 x 3	5 x 3
2000	3 x 1	3 x 2	3 x 2	3 x 3

※ The loading quantity is based on the theoretical maximum capacity, so the actual amount may vary depending on the flask used by the customer. Loading to the maximum capacity could pose a risk of sterilization failure.





Drain Bucket



Exhaust Tank



Basket

Recorder
(Thermal Line)

Pressure Sensor



Step Stool




Standard Accessories

Cat. No.	Description	Standard	Max. mounted No.	Model
00ATC0000002	Basket, Stainless steel (Ø380 x 180 mm)	2	2	ST-50G
		-	2	ST-65G
		-	3	ST-85G
		-	4	ST-105G
00ATC0000003	Basket, Stainless steel (Ø380 x 230 mm)	-	1	ST-50G
		2	2	ST-65G
		-	2	ST-85G
		-	3	ST-105G
00ATC0000004	Basket, Stainless steel (Ø380 x 330 mm)	-	1	ST-50G
		-	1	ST-65G
		2	2	ST-85G
		-	2	ST-105G
00ATC0000126	Basket, Stainless steel (Ø380 x 405 mm)	-	-	ST-50G
		-	1	ST-65G
		-	1	ST-85G
		2	2	ST-105G
AAAL1521	Exhaust Tank	1	1	ST-50G/65G/85G/105G
00ATC0000091	Drain Bucket (490 x 291 x 107, mm)	1	1	ST-50G/65G/85G/105G
AAAL1531	Step Stool, Powder Coated Steel (500 x 400 x 250, mm)	-	-	ST-105G

Optional Accessories

Cat. No.	Description
AAAL1501	Recorder set : Recorder (Dot, 6 Point) & Temp. sensor (Double)
AAAL1502	Recorder set : Recorder (Thermal Line) & Temp. sensor (Double)
00ATC0000108	Temperature sensor for Recorder
AAAL1511	Pressure sensor

※ The above option is only installed when initial order.



Specification


Model	ST-50G	ST-65G	ST-85G	ST-105G
Temperature				
Sterilization (°C / °F)	110 to 123 / 230 to 253.4	110 to 123 / 230 to 253.4	110 to 123 / 230 to 253.4	110 to 123 / 230 to 253.4
Melting (°C / °F)	60 to 100 / 140 to 212	60 to 100 / 140 to 212	60 to 100 / 140 to 212	60 to 100 / 140 to 212
Pressure				
Sterilization at 121°C (kg / cm ²)	1.1 to 1.2	1.1 to 1.2	1.1 to 1.2	1.1 to 1.2
Gauge range in analog (Mpa)	-0.1 to 0.6	-0.1 to 0.6	-0.1 to 0.6	-0.1 to 0.6
Permissible altitude*	Up to 800	Up to 800	Up to 800	Up to 800
Dimension				
Volume (L / cu ft)	50 / 1.76	65 / 2.29	85 / 3	105 / 3.7
Internal (Ø x H, mm / inch)	400 x 458 / 15.75 x 18.03	400 x 558 / 15.75 x 21.97	400 x 758 / 15.75 x 29.84	400 x 908 / 15.75 x 35.75
External (W x D x H, mm / inch)	624 x 672 x 870 / 24.57 x 26.46 x 34.25	624 x 672 x 971 / 24.57 x 26.46 x 38.23	624 x 672 x 1083 / 24.57 x 26.46 x 42.64	624 x 672 x 1233 / 24.57 x 26.46 x 48.54
External with exhaust tank (W x D x H, mm / inch)	722 x 672 x 870 / 28.42 x 26.46 x 34.25	722 x 672 x 971 / 28.42 x 26.46 x 38.23	722 x 672 x 1083 / 28.42 x 26.46 x 42.64	722 x 672 x 1233 / 28.42 x 26.46 x 48.54
Door max. opening (W x D x H, mm / inch)	1012 x 672 x 882 / 39.84 x 26.46 x 34.72	1012 x 672 x 983 / 39.84 x 26.46 x 38.7	1012 x 672 x 1095 / 39.84 x 26.46 x 43.11	1012 x 672 x 1245 / 39.84 x 26.46 x 49.01
Basket (Ø x H, mm / inch)	380 x 180 / 14.96 x 7.09	380 x 230 / 14.96 x 9.05	380 x 330 / 14.96 x 12.99	380 x 405 / 14.96 x 15.94
Weight (kg / lbs)	94.5 / 208.34	99.5 / 219.36	104 / 229.28	110 / 242.51
Electrical data & Ordering information				
230V, 50/60Hz, A	14.3	14.3	14.3	21.3
Cat. No.	AAHL1015K	AAHL1025K	AAHL1035K	AAHL1045K

* This product is not available for usage at an altitude higher than 800m.

Pumps

Select the optimum model between two different types of pump.

Both are suitable for rotary vacuum evaporation, decompression distillation, and vacuum filtering. Especially rotary vacuum pump provides more powerful pumping ideal for vacuum drying.



VE-11
with optional vacuum gauge / regulator



W2V-10

Electrical Aspirator Pump



- Portable dual-channel water-jet aspirator pumps
 - Creating a vacuum by means of venturi effect, aspirator pumps are suitable for rotary vacuum evaporation, decompression distillation, vacuum drying, and vacuum filtering.
- Built-in circulating pump and water tank make the unit portable.
- Circulating pump continuously pumps water quietly across a set of aspirators, and is therefore clean and efficient.
- Built-in check valve prevents backflow of water into the aspirator pump.
- Corrosion-resistant construction of all wetted parts with 304 stainless steel, polypropylene, silicone, or nickel-coated brass.
- Water tank features a drain port and a spigot for easy changing of water.

Accessories & Options

- Vacuum gauge with regulator: to monitor and control the vacuum pressure. (between 0.0267~0.0998MPa / 7.9~29.5" Hg)
- Cooling coil: to prevent vacuum from decreasing due to temperature increase.



Cooling coil



Aspirators

Cat. No.	Description
AAA71531	Vacuum gauge with regulator
AAA71501	Cooling coil
00EFA0009111	Aspirator (1ea)
AAH71122	Aspirator (10ea)

Model	VE-11		
Bath volume (L / cu ft)	9.5 / 0.3		
Max. vacuum (Mpa / mmHg)	0.0998 / 680		
Free-air capacity	18L / minute (0.6 cfm) x 2ea		
Material (bath / aspirator)	Polypropylene / Nickel-coated brass x 2ea		
Bath interior (W x L , D) (mm / inch)	298 x 210, 227 / 11.7 x 8.3 ,8.9		
Exterior (W x L x H) (mm / inch)	330 x 265 x 390 / 13 x 10.4 x 15.4		
Inhale nozzle	External Ø9.5 mm (0.4") x 2ea		
Drain cock, Over flow nozzle (Ø, mm / inch)	17 / 0.7		
Net weight (kg / lbs)	5.7 / 12.6		
Electrical requirements	230V, 50Hz / 0.7A	230V, 60Hz / 0.7A	120V, 60Hz / 1.4A
Cat. No.	AAH71015K	AAH71011K	AAH71016U

Rotary Vacuum Pump



- Rotary vacuum pump equipped with direct connect oil seal prevents air reflux, essentially used with vacuum ovens to achieve a higher level of vacuum pressure.
- Oil mist trap, an optional feature, prevents oil in-flow.

Accessories & Options

Cat. No.	Model	Suitable for
BEA700351	Oil Mist Trap, WOF-150 (PC)	W2V-10 / 20 / 40
BEA700355	Oil Mist Trap WOF-150 (Metal)	W2V-10 / 20 / 40
BEA700352	Oil Mist Trap (PA)	W2V-60
BEA700353	Vacuum Pump Oil (1L)	All
BEA700354	Vacuum Pump Oil (20L)	All

Model	W2V-10	W2V-20	W2V-40	W2V-60
Pumping speed (L/min, gal/min)	100 / 26.4	200 / 52.8	400 / 105.7	600 / 158.5
Ultimate pressure (Pa)	≤1.3 x 10 ⁻¹ (1 x 10 ⁻³ Torr)			
Motor power (W / Hp)	400 / 0.5	400 / 0.5	750 / 1.0	1500 / 2.0
Motor speed (rpm)	1700			
Oil capacity (L)	0.5	0.6	1.6	2.3
Inlet port (DN)	NW25			NW40 / 25
Outlet port (DN)	NW25			NW40 / 25
Weight (kg / lbs)	22.5 / 49.6	23.3 / 51.4	34.8 / 76.7	53.5 / 118.0
Dimensions (W x D x H, mm / inch)	150 x 392 x 250 / 5.9 x 15.4 x 9.8	150 x 423 x 250 / 5.9 x 16.6 x 9.8	170 x 484 x 287 / 6.7 x 18.9 x 11.3	206 x 609 x 313 / 8.1 x 24.0 x 12.3
Electrical requirements (230V)	50/60Hz			
Cat. No.	BEH73205K	BEH73215K	BEH73225K	BEH73235K

Ultrasonic Cleaners

Aggressive cleaning ability in a variety of output modes!

Industrial B.L.T type ultrasonic transducer (Bolt Clamped Langevin Type Transducer) generates high-frequency sound waves that provide maximum cleaning power.



UCP-10



UCS-20

UCP models (ABS, Membrane key pad, and SS tank)



- Durable / flame resistance ABS material used and cover with 304 stainless steel tanks.
- Over-temperature protection with alarm. (except for UCP-02)
- Digital LED display using splash-proof keypads: timer (1 to 99 min, Max 4 hour operating), temperature monitor. (except for UCP-02)
- High-frequency 40 kHz sound waves provide greater cleaning power and increased reliability.
- Multiple output mode (Normal/Pulse/Sweep) with power selection (Low/Middle/ High). (except for UCP-02)
- Before cleaning, dissolved gas removed by degas function. (except for UCP-02)
- Drain valve, basket, cover as standard offer. (except for UCP-02)

Model	UCP-02	UCP-10	UCP-20
Bath volume (L / cu ft)	2 / 0.07	10 / 0.35	20 / 0.70
Temperature range (°C / °F)	Amb. to 70 / Amb. to 158		
Ultrasonic power (W)	100	300	500
Dimensions			
Bath interior (W x L, D) (mm / inch)	142 x 132, 104 / 5.6 x 5.2, 4.1	295 x 240, 160 / 11.6 x 9.4, 6.3	500 x 290, 160 / 19.7 x 11.4, 6.3
Exterior (W x L x H) (mm / inch)	290 x 210 x 200 / 11.4 x 8.4 x 7.9	490 x 368 x 297 / 19.3 x 14.5 x 11.7	721 x 451 x 297 / 28.4 x 17.8 x 11.7
Net weight (kg / lbs)	2 / 4.4	11 / 24.3	17 / 37.5
Electrical requirements (230V, 50/60Hz, A)	0.6	1.3	2.2
Cat. No.	BDH532011K	BDH532012K	BDH532013K

※CE certified model: UCP-10, UCP-20.

UCS models (SS external and tank)



- External and tank parts are made of durable 304 stainless steel.
- Over-temperature protection with alarm.
- Digital LED display using splash-proof keypads: timer (1 to 99 min, Max 4 hour operating), temperature monitor.
- High-frequency 40 kHz sound waves provide greater cleaning power and increased reliability.
- Multiple output mode(Normal/Pulse/Sweep) with power selection (Low/Middle/ High).
- Drain valve, basket, cover as standard offer.

Model	UCS-05	UCS-10	UCS-20
Bath volume (L / cu ft)	5.7 / 0.20	10 / 0.35	20 / 0.70
Temperature range (°C / °F)	Amb. to 70 / Amb. to 158		
Ultrasonic power (W)	200	300	500
Dimensions			
Bath interior (W x L, D) (mm / inch)	295 x 160, 160 / 11.6 x 6.3, 6.3	295 x 240, 160 / 11.6 x 9.4, 6.3	500 x 290, 160 / 19.7 x 11.4, 6.3
Exterior (W x L x H) (mm / inch)	325 x 230 x 320 / 12.8 x 9.1 x 12.6	325 x 330 x 320 / 12.8 x 12.9 x 12.6	530 x 330 x 297 / 20.9 x 12.9 x 11.7
Net weight (kg / lbs)	8 / 17.6	11 / 24.3	15 / 33.1
Electrical requirements (230V, 50/60Hz, A)	0.9	1.3	2.2
Cat. No.	BDH532031K	BDH532032K	BDH532033K



Biosafety Cabinet & Clean Bench, Fume Hood



Biosafety Cabinet & Clean Bench, Fume Hood



Biosafety Cabinet & Clean Bench



General Application

Biosafety Cabinet: Handling pathogens and infectious samples, tissue culture, media dispensing, drug development, etc. BSL3 is mandatory, BSL2 is recommended.

Clean Bench: Non-toxic sap, DNA/RNA purification, plant tissue culture, solid medium division, electronic inspection, medical equipment assembly, pharmaceuticals, etc.

PCR workstation: Medical equipment assembly, PCR sample preparation, volatile staining, sterile medium division, tissue fixation/staining preparation.

UV-sterilization Cabinet: Sterilized to prevent DNA/RNA contamination.

	Description	Filter	Working Table	Model	Page
Protect user/ sample/ environment	Biosafety Cabinet <ul style="list-style-type: none"> • Class II A2 type: Protects the user, environment, and samples • Durable ECM motor: Ensures long-lasting performance. 	Exhaust ULPA (standard) Supply ULPA (standard)	Stainless Steel #304	JB-A	196
Protect sample	Clean Bench <ul style="list-style-type: none"> • Fully validated • Superior user convenience 	HEPA Pre filter (standard)	Stainless Steel #304	BC-H BC-B	200 202
	PCR Workstation <ul style="list-style-type: none"> • Convenient PCR working environment • Easy wind speed control 	HEPA Pre filter (standard)	Polycarbonate	PW	203
	UV-Sterilization Cabinet <ul style="list-style-type: none"> • Sterilized with UV light • Stainless steel work surface 	N/A	Stainless Steel #304	UVC	204

※ ※ The contents of the above and the contents of this catalog may differ depending on the specific model and conditions of use. For the information about the features and specifications that applying to each models, please check the information on the corresponding page of each models.



Satisfaction with the requirements of European regulations and guidelines.



Ensuring the reliability of electronic records and signatures to support GMP audit.



Bluetooth connection available for mobile control.



Registered patent based on Jeio Tech's proprietary technology.



If sash door is opened, automatically UV light is off and fluorescent light and fan are operated.



Differentiated diffusing muffler produces good laminar flow.



Equipped with UV Light.



Equipped with the same controller inside and outside for convenience.



Purifies and discharges harmful gas through internal filter.



2 year warranty free A/S.

Biosafety Cabinet

Safe Experiment Environment Certified by EN 12469


- Airflow design that protects users, samples, and the environment by separating indoor air from internal air of the equipment.
- By applying ULPA filter that can filter particles of 0.1~0.3 um by 99.999%, it provides more reliable and clean workspace than traditional HEPA filter.
- Maintains a stable airflow velocity with an inflow of 0.53 m/s and a downflow of 0.31 m/s when the sash height is 210 mm. (Certified by EN 12469)



JB-15A



Microbial Challenge Test Completed




ULPA filter that removes 99.999% of 0.1~0.3 um particles



Canopy (Option) that removes small quantities of volatile chemicals

Application of ECM Motor with High Durability and Excellent Energy Efficiency

- Embedded ECM (Electronically Commutated Motor) which generates low heat and maintains a constant airflow velocity under voltage fluctuations (unstable voltage) or heavy filter loads.
- Low noise compared to conventional AC motor-applied equipment can reduce the fatigue of the user during prolonged experiments.
- Realize office environment noise level. (60 dB)



Comparison of the motor control depending on voltage change (AC vs ECM)


Safety Function to Concentrate on Experiments

- Monitors velocity and temperature, and generates an alarm when it exceeds the ranges.
- Generates visual and audio alarms for abnormal power supply, abnormal sensor, and accessories replacement for immediate safety measures.
- Records the most recent 20 alarms, and separately manages replacement alarms of the main accessories for useful experiment feedback.
- A UV-protecting tempered glass is applied and resistant to shock. The glass breaks into small pieces in case of breakage to prevent injury to the user.
- Removes small quantities of volatile chemicals or radioactive isotopes using Canopy. (Option)



Smart Technologies

- Sash door, UV lamp, LED lamp, and blower are interconnected allowing for integrated control by opening/closing sash door. (Patented in Korea)
- Real-time push notification to mobile devices and live monitoring of temperature and airflow data with LC Connected. (Requires LC GreenBox, option)
- Offers both 'Auto Start' mode (fan and lamp start/stop based on sash position), and 'Manual' mode (manually control fan and lamp to customize user environment).
- Users can easily check status during experiments with their mobile devices and adjust settings without removing their hands. (LC Bluetooth, option)



Real-time mobile monitoring (Option)



Push notification to smartphone (Option)




Using LC Bluetooth helps to prevent cross-contamination (Option)

User Convenience


- Integrates with LC DataKeeper, compliant with 21 CFR Part 11, to support GMP audits (option).
- Safe experiment environment with Post-purge: fan runs for a set time to remove residual contaminants before shutdown.
- Two outlets for small equipment; outlet power can be enabled/disabled via touch controller.
- Events logged for easy tracking and separate records for filter, lamp, and UV lamp replacements for streamlined maintenance.
- MagPost™: internal steel plate in walls allows for convenient magnetic board use.



Data integrity software for GMP audit support (Option)



MagPost™ minimizes human error by attaching protocols with magnets



Consumable replacement alarm history

Ergonomic Structure

- Easy to clean working tray with handles, and convenient disposal of dirty water through drain valve.
- An ergonomic armrest designed for comfortable arm positioning during experiments.
- Seamless single-plate workstation prevents microbial contamination and allows easy cleaning.
- Built-in outlets and valves for suction pumps come as standard, with additional service valve provided for easy space utilization.



Easy-to-clean Working Tray




Arm Rest for Comfortable Experiments



Seamless structure to prevent contamination in fine gaps

Dimension

(Unit: mm)



- | | | | |
|--------------------------|-----------------------|--------------------------|-------------------|
| 1. Control panel | 6. Service cock valve | 11. Fixable caster | 16. Alarm contact |
| 2. Sash door & Handle | 7. Air slot | 12. Front cover | 17. RS-232 port |
| 3. Safe sash height | 8. Arm rest | 13. Exhaust filter cover | |
| 4. Internal power outlet | 9. Work surface | 14. Service nozzle inlet | |
| 5. Stand* | 10. Drain valve | 15. Canopy alarm port | |

* Stand size (W x D x H, mm / inch)

JB-12A : 1353 x 700 x 735 / 53.3 x 27.6 x 28.9 JB-15A : 1653 x 700 x 735 / 65.1 x 27.6 x 28.9 JB-18A : 1953 x 700 x 735 / 76.9 x 27.6 x 28.9

Specification

Models		JB-12A	JB-15A	JB-18A
Type		Class II A2		
Air	Airflow type	Vertical laminar flow		
	Filter	ULPA, typical efficiency of 99.999% at 0.1~0.3 μm		
	Downflow velocity (m/s)	0.31±0.025		
	Inflow velocity (m/s)	0.53±0.025		
	Exhaust volume (m ³ /h)	452	562	672
Material	Main body	Epoxy powder coated steel		
	Working surface	Stainless steel #304		
	Sash	UV absorbing tempered glass (5.0T)		
Dimensions	External (W x D x H, mm / inch)	1353 x 824 x 2257 / 53.3 x 32.4 x 88.9	1653 x 824 x 2257 / 65.1 x 32.4 x 88.9	1953 x 824 x 2257 / 76.9 x 32.4 x 88.9
	Working area (W x D x H, mm / inch)	1230 x 572 x 658 / 48.4 x 22.5 x 25.9	1530 x 572 x 658 / 60.2 x 22.5 x 25.9	1830 x 572 x 658 / 72.0 x 22.5 x 25.9
	Sash max opening (mm / inch)	515 / 20.3		
	Working opening (mm / inch)	210 / 8.3		
Noise level (dB)		< 57 db	< 58 db	< 60 db
LED lamp intensity (Lux)		More than 800		
Net weight with stand (kg / lbs)		305 / 672	347 / 765	390 / 860
Power Consumption (230V, 50/60Hz) (W)		308	374	418
Cat. No. (Body + Stand)		AAAB8001	AAAB8002	AAAB8003
Cat. No. (Body only)		AAAB8101K	AAAB8102K	AAAB8103K
Cat. No. (Stand only)		AAAB5101K	AAAB5102K	AAAB5103K

※ The unit comes standard with an ULPA filter, UV lamp, LED lamp, and 2 service valves.

Accessories

(Unit: mm)

Designation	Cat. No.	Model	Description
ULPA filter set (2 ULPA filters)	AAAB8511	JB-12A	670 x 457 x 120 and 1209 x 457 x 65
	AAAB8512	JB-15A	870 x 457 x 120 and 1509 x 457 x 65
	AAAB8513	JB-18A	1070 x 457 x 120 and 1809 x 457 x 65
Canopy	AAAB8551	JB-12A	708 x 441 x 285
	AAAB8552	JB-15A	905 x 441 x 285
	AAAB8553	JB-18A	1105 x 441 x 285
UV lamp	CHE0004423	JB-12A/JB-15A/JB-18A	Ø26*890
LED lamp	BSC0000071	JB-12A/JB-15A/JB-18A	Ø16*850
IV bar	AAAB8521	JB-12A	Ø15.8*1223
	AAAB8522	JB-15A	Ø15.8*1523
	AAAB8523	JB-18A	Ø15.8*1823
Service cock valve	AAAB8531	JB-12A/JB-15A/JB-18A	86 x 131 x 64
Suction pump	BEA1002901	JB-12A/JB-15A/JB-18A	160 x 120 x 150
LC-Bluetooth	AAAB8541	JB-12A/JB-15A/JB-18A	Supported from Android 9 and above, iOS 11 and above
Tablet PC for LC-Bluetooth	BEA0008542	JB-12A/JB-15A/JB-18A	Contact to salesperson
LC GreenBox	AAHQ1011K	JB-12A/JB-15A/JB-18A	Supported from Android 9 and above, iOS 11 and above
LC DataKeeper-Basic	SSA11001	JB-12A/JB-15A/JB-18A	FDA 21 CFR Part 11 compliant software, 1 year license
LC DataKeeper-Pro	SSA11002	JB-12A/JB-15A/JB-18A	FDA 21 CFR Part 11 compliant software, permanent license

※ The unit comes standard with an ULPA filter, UV lamp, LED lamp, and 2 service valves.
 ※ The Trial version of LC DataKeeper is available for 30 days and features same as Basic.

Biosafety Cabinet vs. Clean Bench vs. Fume Hood: Similar Appearance, Different Functions

Biosafety Cabinet

Protecting Users, Samples, and the Environment Simultaneously

Sample	User	Environment
✓	✓	✓

- Safe Zone
- Indoor Air
- Contaminated Air Passed Through Samples
- Clean Air Passed Through ULPA

Airflow of Biosafety Cabinet

Air entering through the sash door doesn't directly enter the chamber; instead, it moves to the top through a negative pressure plenum along with the air that has passed through the chamber. Some of this air is purified and exhausted through the exhaust ULPA filter, while some is purified by the supply ULPA filter to create vertical laminar airflow inside the chamber, ensuring the protection of users, samples, and the environment simultaneously.

Clean Bench

Preventing Cross-Contamination of Samples

Sample	User	Environment
✓	✗	✗

- Safe Zone
- Indoor Air
- Contaminated Air Passed Through Samples
- Clean Air Passed Through ULPA

Airflow of Clean Bench

Air entering from the top is filtered through a HEPA filter by a fan and provided creating vertical laminar airflow inside the chamber. The air passing through the chamber is discharged through the sash opening, preventing cross-contamination of samples.

Fume Hood

Protecting Users from Harmful Gases

Sample	User	Environment
✗	✓	✗

- Safe Zone
- Indoor Air
- Air Containing Harmful Gases such as Samples

Airflow of Fume Hood

Indoor air enters the chamber through the sash opening, and the air inside the chamber exits through the rear baffle into a central exhaust duct. It safely protects users from harmful gases that can affect respiratory health.

Clean Bench

Advanced type

Easy to use with dual controller and smart system

Structural Functional Features

- Dedicated controller and digital airflow sensor allow for operation at set wind speed even with fluctuations in sash height, minimizing power consumption and dramatically increasing HEPA filter life.
- Polyester fiber pre-filter extends the life of the HEPA filter.
- High-performance HEPA filter provides a clean environment. (99.99% removal of 0.3 μm particles, Class H14 according to EN 1822)
- Formation of reliable and uniform laminar flow through application of specialized diffusing muffler system.
- Easy to replace the filter on the front panel.
- Easy to move and fix the device.

Use Convenience Features

- Displays the temperature and humidity inside the workspace.
- Warning provided when the sash door is opened over the recommended height during operation.
- Easy and smooth sash door operation.
- Provides a quiet environment with minimal noise.
- Easy to operate with clear VFD and touch keys.
- UV lamp, fluorescent lamp, blower-dedicated button.
- Easy to check wind speed, temperature and humidity.

Outstanding Safety

- Smart door system protects the experimenter from UV light exposure.
- Toughened glass sash door blocks UV light.
- Protected from overcurrent.




BC-11H

with Stand, HEPA filter, gas cock (option)
Pre filter, UV Lamp, Fluorescent lamp Electrical socket (standard)

Completed internal cleanliness evaluation and microbiological testing compatibility test.

- Internal cleanliness rating (ISO 14644-1 standard)
Equivalent to ISO Class 2
(Experiment with sterile environment can be carried out from ISO Class 5, US Federal Class 100)
- Microbiological testing (ISO 14698 standard)
ISO Class 5 rating, satisfies the microbiological testing standards



Product Configuration For Easy Use

Easy operation during use with dual controller

- Equipped with the same controller inside and outside, it is convenient to operate the equipment even during work.

Built-in gas and air connection valves and power outlet

- Valves for gas and air use inside.
- Built-in electrical outlet with cover provided.




Smart Door System (Patent Granted)

- When the sash door is opened during UV sterilization, the UV lamp is automatically turned off to protect against UV exposure, and the fluorescent lamp and blower are automatically operated. (Korean Patent)



Smart Maintenance

- HEPA filter status checking**
Self-checking of HEPA filter status with a digital differential pressure sensor.
- UV lamp replacement notification**
UV irradiation intensity is monitored by a UV sensor, and an automatic replacement notification is provided when irradiation intensity drops below 80%.



Specification

Model		BC-01H	BC-11H	BC-21H
Air	Air flow type	Vertical Laminar Flow	Vertical Laminar Flow	Vertical Laminar Flow
	Airflow velocity (m/sec/fpm)	0.3 / 59	0.45 / 89	0.3 / 59
	Main filter	HEPA (99.99% at 0.3 μm)	HEPA (99.99% at 0.3 μm)	HEPA (99.99% at 0.3 μm)
	Pre-filter	Polyester Fibers	Polyester Fibers	Polyester Fibers
Illumination	U.V. lamp (254 nm, W)	25 x 1ea	30 x 1ea	25 x 2ea
	Fluorescent lamp (W)	30 x 2ea	30 x 2ea	32 x 2ea
	Fluorescent lamp intensity (Lux)	More than 650	More than 650	More than 650
Dimensions	Exterior with stand (W x D x H, mm / inch)	1135 x 620 x 1860 / 44.7 x 24.4 x 73.2	1435 x 620 x 1860 / 56.5 x 24.4 x 73.2	2035 x 620 x 1860 / 80.1 x 24.4 x 73.2
	Exterior without stand (W x D x H, mm / inch)	1135 x 620 x 1150 / 44.7 x 24.4 x 45.3	1435 x 620 x 1150 / 56.5 x 24.4 x 45.3	2035 x 620 x 1150 / 80.1 x 24.4 x 45.3
	Interior (W x D x H, mm / inch)	945 x 585 x 670 / 37.2 x 23.0 x 26.4	1245 x 585 x 670 / 49.0 x 23.0 x 26.4	1845 x 585 x 670 / 72.6 x 23.0 x 26.4
	Net weight (kg / lbs)	140 / 308.6	185 / 407.9	225 / 496
	Net weight (body + stand) (kg / lbs)	170 / 374.8	215 / 474	255 / 562.2
Noise level (dB)		Less than 65	Less than 65	Less than 65
Electrical requirements (230V, 50/60Hz, 1ph, A)		1.38	1.73	2.30
Cat. No.		AAHA5012K	AAHA5022K	AAHA5032K
Electrical requirements (120V, 60Hz, 1ph, A)		2.62	3.27	4.36
Cat. No.		AAHA5013U	AAHA5023U	AAHA5033U

Accessories **Page 205** HEPA Filter, Pre-filter, Stand, Gas Cock

Biosafety Cabinet & Clean Bench

Clean Bench

Basic type



Vertical laminar flow clean benches offer quick operation by simple adjustable analog control

Structural Functional Features

- High-performance HEPA filter provides a clean environment. (99.99% removal of 0.3 μm particles, Class H14 according to EN 1822)
- Polyester fiber pre-filter extends the life of the HEPA filter.
- Formation of reliable and uniform laminar flow through application of specialized diffusing muffler system.
- Easy to replace the filter on the front panel.
- Easy to move and fix the device.
- Internal cleanliness rating (ISO 14644-1 standard) Equivalent to ISO Class 2, Microbiological testing (ISO 14698 standard) ISO Class 5 rating, satisfies the microbiological testing standards (Experiment with sterile environment can be carried out from ISO Class 5)

Use Convenience Features

- Determines when to replace the filter with the built-in differential pressure gauge.
- Easy operation and wind speed control.
- Built-in gas and air connection valves and power outlet.
- Easy and smooth sash door operation.
- Provides a quiet environment with minimal noise.
- Height-adjustable caster for easy vertical installation.

Outstanding Safety

- When the sash door is opened in when UV is ON, UV is automatically turned off for user protection.
- Toughened glass sash door blocks UV light.
- Protects against UV exposure through three modes of illumination. (UV/Fluorescent/Off)



BC-11B

with stand, HEPA filter (option)
Pre filter, UV Lamp, Fluorescent lamp, Electrical socket
Differential pressure gauge (standard)



Controller



Differential pressure gauge

Specification

	Model	BC-01B	BC-11B	BC-21B
Air	Air flow type	Vertical Laminar Flow	Vertical Laminar Flow	Vertical Laminar Flow
	Airflow velocity (m/sec/fpm)	0.3 / 59	0.45 / 89	0.3 / 59
	Main filter	HEPA (99.99% at 0.3 μm)	HEPA (99.99% at 0.3 μm)	HEPA (99.99% at 0.3 μm)
	Pre-filter	Polyester Fibers	Polyester Fibers	Polyester Fibers
Illumination	U.V. lamp (254 nm, W)	25 x 1ea	30 x 1ea	25 x 2ea
	Fluorescent lamp (W)	30 x 2ea	30 x 2ea	32 x 2ea
	Fluorescent lamp intensity (Lux)	more than 650	more than 650	more than 650
Dimensions	Exterior with stand (W x D x H, mm / inch)	1135 x 620 x 1860 / 44.7 x 24.4 x 73.2	1435 x 620 x 1860 / 56.5 x 24.4 x 73.2	2035 x 620 x 1860 / 80.1 x 24.4 x 73.2
	Exterior without stand (W x D x H, mm / inch)	1135 x 620 x 1150 / 44.7 x 24.4 x 45.3	1435 x 620 x 1150 / 56.5 x 24.4 x 45.3	2035 x 620 x 1150 / 80.1 x 24.4 x 45.3
	Interior (W x D x H, mm / inch)	945 x 585 x 670 / 37.2 x 23.0 x 26.4	1245 x 585 x 670 / 49.0 x 23.0 x 26.4	1845 x 585 x 670 / 72.6 x 23.0 x 26.4
	Net weight (kg / lbs)	140 / 308.6	185 / 407.9	225 / 496
	Net weight (body + stand) (kg / lbs)	170 / 374.8	215 / 474	255 / 562.2
Noise Level (dB)	Less than 65	Less than 65	Less than 65	
Electrical requirements (230V, 50/60Hz, 1ph, A)	1.38	1.73	2.30	
Cat. No.	AAHA6012K	AAHA6022K	AAHA6032K	
Electrical requirements (120V, 60Hz, 1ph, A)	2.62	3.27	4.36	
Cat. No.	AAHA6013U	AAHA6023U	AAHA6033U	

Accessories [Page 205](#) HEPA Filter, Pre-filter, Stand, Gas Cock, Differential Pressure Gauge

Specially designed to minimize the sample contamination during PCR

Structural Functional Features

- High-performance HEPA filter provides a clean environment. (99.99% removal of 0.3 μm particles, Class H14 according to EN 1822)
- Convenient for use by mounting on an existing bench.
- Fluorescent lamp is included to make it convenient for conducting experiments.
- Fluorescent lamps are installed on the inside of the front to reduce eye fatigue.
- Effective sterilization with 254nm UV lamp. (replacement cycle: approx. 8,000 hours)
- Polyester fiber pre-filter extends the life of the HEPA filter.

Use Convenience Features

- Easy to replace the filter on the front panel.
- Front cover opening structure at the front for convenient inspection and maintenance.
- Convenient shelf with pipette holder to place pipette tip or vial rack
- Includes a hole on the right side to take out a power cable when using small device inside.
- Easy control of wind speed using the controller.
- Three model sizes depending on the bench size.
- UV timer function. (up to 30 minutes)

Outstanding Safety

- Polycarbonate work surface blocks UV by 99.9%.
- When the door is opened with UV light is on, the UV light is turned off automatically to protect the user.
- Protection from overcurrent.




PW-21

with Pre filter, HEPA filter, UV lamp, Fluorescent lamp (standard)

Completed internal cleanliness evaluation and microbiological testing compatibility test.

- Internal cleanliness rating (ISO 14644-1 standard) **Equivalent to ISO Class 4** (Experiment with sterile environment can be carried out from ISO Class 5, US Federal Class 100)
- Microbiological testing (ISO 14698 standard) **ISO Class 5 rating, satisfies the microbiological testing standards**



Specification

Model		PW-01	PW-11	PW-21
Air	Air flow type	Vertical	Vertical	Vertical
	Airflow velocity (m/sec/fpm)	0.4 / 78.7	0.4 / 78.7	0.4 / 78.7
	Main filter	HEPA (99.99% at 0.3 μm)	HEPA (99.99% at 0.3 μm)	HEPA (99.99% at 0.3 μm)
	Pre-filter	Polyester Fibers	Polyester Fibers	Polyester Fibers
Illumination	U.V. lamp (254 nm, W)	15	20	30
	Fluorescent lamp (W)	15	20	30
	Fluorescent lamp intensity (Lux)	More than 1,000	More than 1,000	More than 1,000
Material	Window (Front, Sides, Rear)	Polycarbonate (Clear)	Polycarbonate (Clear)	Polycarbonate (Clear)
	U.V. protection	99.90%	99.90%	99.90%
	Working table	Polycarbonate (Clear)	Polycarbonate (Clear)	Polycarbonate (Clear)
Dimensions	Exterior (W x D x H, mm / inch)	720 x 605 x 964 / 28 x 24 x 38	900 x 605 x 964 / 35.4 x 24 x 38	1200 x 605 x 964 / 47.2 x 24 x 38
	Interior (W x D x H, mm / inch)	700 x 585 x 602 / 27.5 x 23 x 23.7	880 x 585 x 602 / 34.6 x 23 x 23.7	1180 x 585 x 602 / 46.5 x 23 x 23.7
	Net weight (kg / lbs)	56.2 / 123.9	64.2 / 141.5	77.8 / 171.5
Noise level (dB)		Less than 60	Less than 60	Less than 60
Electrical requirements (230V, 50/60Hz, 1ph, A)		1.3	1.4	1.5
Cat. No.		AAHB3002K	AAHB3012K	AAHB3022K
Electrical requirements (120V, 60Hz, 1ph, A)		2.2	2.7	2.9
Cat. No.		AAHB3003U	AAHB3013U	AAHB3023U

Accessories Page 205 Pre-filter, HEPA Filter

UV Sterilization Cabinet



Clean environment due to sterilization and removal of sources of contamination with UV light

Structural Functional Features

- Built-in fluorescent lamp makes it easy to see inside.
- Compact design makes installation and movement convenient.
- VFD display allows for easy checking of operation status.

Use Convenience Features

- Convenient handling of equipment due to convenient front opening.
- Three sizes of models according to the cabinet size.
- Digital-type UV timer function.
(up to 59 minutes and 59 seconds)

Outstanding Safety

- Block UV 99.9% on the front/ left / right side with polycarbonate material.
- When the door is opened, the UV light is turned off automatically to protect the user and the fluorescent lamp is operated, making it convenient to use.
- Protection from overcurrent.



UVC-11

with UV lamp, Fluorescent lamp (standard)



Specification

Model		UVC-01	UVC-11	UVC-21
Illumination	U.V. Density ($\mu\text{W}/\text{cm}^2$)	300 \pm 10%	350 \pm 10%	300 \pm 20%
	U.V. lamp (254 nm, W)	8	15	20
	Fluorescent lamp (W)	8	15	20
	Fluorescent lamp intensity (Lux)	More than 800	More than 900	More than 1,000
Material	Window (Front, Sides)	Polycarbonate (Clear)	Polycarbonate (Clear)	Polycarbonate (Clear)
	Window (Rear)	Powder coated steel	Powder coated steel	Powder coated steel
	U.V. protection	99.9%	99.9%	99.9%
	Working table	Stainless steel	Stainless steel	Stainless steel
Dimensions	Exterior (W x D x H, mm / inch)	450 x 509 x 610 / 17.7 x 20 x 24	600 x 509 x 610 / 23.6 x 20 x 24	900 x 509 x 610 / 35.4 x 20 x 24
	Interior (W x D x H, mm / inch)	433 x 500 x 468 / 17.0 x 19.7 x 18.4	583 x 500 x 468 / 23 x 19.7 x 18.4	880 x 500 x 468 / 34.6 x 19.7 x 18.4
	Net weight (kg / lbs)	15 / 33.1	17.5 / 38.6	20 / 44.1
Electrical requirements (230V, 50/60Hz, 1ph, A)		0.1	0.2	0.3
Cat. No.		AAHB4002K	AAHB4012K	AAHB4022K
Electrical requirements (120V, 60Hz, 1ph, A)		0.2	0.4	0.6
Cat. No.		AAHB4003U	AAHB4013U	AAHB4023U

Accessories



HEPA Filter (option)

- 99.99% efficient HEPA filter, industry standard size for economical replacement.



Pre Filter (1EA included)

- Easy-to-replace pre-filter for trapping larger particles and extending the life of HEPA filters.



Stand with Casters (option)

- Stand with Casters : Mobile stand made of powder-coated steel ideal for easy relocation.



Differential Pressure Gauge (option-only for BC-B)

- Magnehelic differential pressure gauge measuring pressure drop across the HEPA filter to monitor its efficiency.

Filters

Description	BC-01H/B	BC-11H/B	BC-21H/B	PW-01	PW-11	PW-21
HEPA Filter	AAAB1601	AAAB1602	AAAB1603	00EDA0009219	00EDA0009220	00EDA0009221
Pre Filter	AAAB1611	AAAB1612	AAAB1613	AAAB3511	AAAB3512	AAAB3513

Others

Description	BC-01H/B	BC-11H/B	BC-21H/B
Stand with Casters	AAAB1621	AAAB1622	AAAB1623
Gas Cock	AAAB1561	AAAB1561	AAAB1561
Differential Pressure Gauge (Only BC-B)	AAAB1571	AAAB1571	AAAB1571
UV lamp	230V, 50/60Hz	AAAB1631	AAAB1633
	120V, 60Hz	AAAB1632	AAAB1634
			AAAB1635
			AAAB1636

Fume Hood



General Application

Hoods: chemical sampling, cosmetics, production, vapor-generating experiments, slide coating, welding.

- 1) **Ductless Fume Hood:** For spaces where duct installation is difficult but air purification is required.
- 2) **Fume Hood, Small & Mobile Fume Extractor:** For mobile laboratory setups or compact spaces where air purification is necessary.
- 3) **Arm Hood:** For areas with central exhaust systems and easy duct connection above the workbench.

	Description	Filter	Working Table	Model	Page
Safe removal of harmful gas	Ductless Fume Hood <ul style="list-style-type: none"> Structure without duct connection Various filters and workbenches 	(Option/Select up to 2ea) Complex Carbon Acid Formaldehyde Ammonia/Amine HEPA	(Option/Select 1ea) Stainless Steel (#304, #316), Phenolic Laminate, Polypropylene, Bakelite, PVC	DLH	208
	Fume Hood (small) <ul style="list-style-type: none"> Convenient small tabletop hood Polyethylene material 	N/A	Polyethylene	PMH POH	210 211
	Mobile Fume Extractor <ul style="list-style-type: none"> Local exhaust and air purifying type Wide range of choices provided High-efficiency filter box included 	(Option/Select up to 2ea) Complex Carbon Acid Formaldehyde Ammonia/Amine HEPA Pre filter (standard)	N/A	EFS PFS	212 213
	Arm Hood <ul style="list-style-type: none"> Mountable on the ceiling or wall Configurable with lab table 	Arm material: Polypropylene, Aluminum	N/A	AH	216

※ The contents of the above and the contents of this catalog may differ depending on the specific model and conditions of use. For the information about the features and specifications that applying to each models, please check the information on the corresponding page of each models.



Leakage minimized using an integrated mold.



Exhaust via duct connection.



Purifies and discharges harmful gas through internal filter.



2 year warranty free A/S.

Ductless Fume Hood

Built-in filter hood without duct connection

Structural Functional Features

- Easy to move and install due to there being no duct connection.
- Filtered air circulates through the room, dramatically reducing the air conditioning costs for room air.
- Structure that allows for selection of the appropriate filter and table material according to the chemical used.
- Polyester fiber pre-filter is included as standard to maximize the life of the main filter.
- Fluorescent lamp is included to make it convenient for conducting experiments.
- Convenient utility hole for connection with external equipment is included in the back of the work surface.
- Transparent acrylic material for easy internal observation.
- Easy to move and fix the device.

Use Convenience Features

- Convenient handling of equipment due to convenient front opening.
- Provides a quiet environment with minimal noise.
- Port for filter performance test provided on the front for convenient periodical validation.
- It can be used more safely by applying the fume hood monitoring system. (option)


Reliability Verification Completed

- Designed, manufactured and self-certified according to international standards to protect users and the environment from fumes generated during experiments.
- SPS-KACA002-132 (Korean Air Cleaning Association)
 - Applies verified composite carbon filter to remove harmful gases.
 - Completed deodorization efficiency test. (VOC, acid, base).
- ANSI / ASHRAE 110-1995, NF EN 14175-1 (U.S.A / EU)
 - Front suction wind speed measurement. (average 0.4m/s or more)
 - Maintains average surface wind speed $\pm 20\%$ at all measuring points.
 - Localization of air flow observation / observation of large-volume smoke flow. (no opening leakage)
 - Tracer-gas measurement, and user bodily inhalation pollution test. (less than 0.05 ppm)



DLH-01G

with Stainless Steel work surface & Stand with Casters (option)
Pre filter, Fluorescent lamp (standard)



Specification

Model		DLH-01G	DLH-11G
Air	Face velocity (m/s / fpm)	0.4 / 80 (Initial set point)	0.4 / 80 (Initial set point)
	Main filter	Max. 2ea available	Max. 2ea available
	Pre-Filter	Polyester fibers	Polyester fibers
Fluorescent lamp intensity (Lux)		More than 600	More than 600
Dimensions	Exterior with stand (mm / inch)	900 x 660 x 1985 / 35.4 x 26 x 78.1	1200 x 660 x 1985 / 47.2 x 26 x 78.1
	Exterior without stand (mm / inch)	900 x 660 x 1250 / 35.4 x 26 x 49.2	1200 x 660 x 1250 / 47.2 x 26 x 49.2
	Interior (mm / inch)	880 x 640 x 800 / 34.6 x 25.2 x 31.5	1180 x 640 x 800 / 46.5 x 25.2 x 31.5
	Net weight (body) (kg / lbs)	100 / 220.5	118 / 260.1
	Net weight (body + stand) (kg / lbs)	120 / 264.6	140 / 308.6
Noise level (dB)		Less than 55	Less than 55
Electrical requirements (230V, 50/60Hz, A)		0.7	0.7
Cat. No.		AAHB2002K	AAHB2012K
Electrical requirements (120V, 60Hz, A)		1.3	1.3
Cat. No.		AAHB2003U	AAHB2013U

※ Standard : Main Body, Pre-Filter, Fluorescent Lamp



Useful Dedicated Controller

- Easy use with Jeio Tech dedicated controller.
- Fan speed control function to adjust inlet wind speed.
- Available to check the filter usage time.
- Set filter replacement timing, and alarm function for the elapsed time.

The Front Design Makes Convenient Usage

- Equipped wind gauge on the front to check the airflow state intuitively.
- Easily access the top cover at the front to do easy maintenance such as the replacement of the filter or fluorescent lamp.



Dedicated Controller



Anemometer



Front Cover

Accessories

Filters

Description	Application	DLH-01G	DLH-11G
Complex carbon filter	Filters a wide range of chemicals, including acid, ammonia, and aldehyde	AAAB2543	AAAB2544
Acid filter¹⁾	Filters acid gases such as HCl, HF, H ₂ S, and CH ₃ COOH	AAAB2545	AAAB2546
Formaldehyde filter	Filters formaldehyde, acetic aldehyde, etc.	AAAB2547	AAAB2548
Ammonia / Amines filter	Filters basic gases such as NH ₃ , (CH ₃) ₃ N, etc.	AAAB2549	AAAB2550
HEPA filter	Filters 99.99% of toxic dust and other 0.3µm particles	00EDA0009193	00EDA0009194
Pre filter	Primarily filters particles to extend the life of other filters	00EDA0009196	00EDA0009198

1) Please do not use high percentage of reactive acid such as perchloric acid.

Work surfaces / Stands

Description	DLH-01G	DLH-11G
Work surface (SS #304)	AAAB2501	AAAB2502
Work surface (SS #316)	AAAB2503	AAAB2504
Work surface (phenolic laminate)	AAAB2513	AAAB2514
Work surface (polypropylene)	AAAB2507	AAAB2508
Work surface (bakelite)	AAAB2509	AAAB2510
Work surface (PVC)	AAAB2511	AAAB2512
Stand with casters	AAAB2521	AAAB2522

- Select an appropriate table for your purpose.
- Easy installation with the dedicated stand.



Gas detector / Gas detecting tubes

Description	Cat. No.
Gas detector (KITAGAWA AP-20)	AAAB2531
Gas detecting tube (benzene, GASTEC-121SP (0.5-10ppm))	00EAA0001550
Gas detecting tube (toluene, GASTEC-122L (1-100ppm))	00EAA0001551
Gas detecting tube (acetic acid, GASTEC-81L (0.125-25ppm))	00EAA0001552
Gas detecting tube (chloroform, GASTEC-137 (4-400ppm))	00EAA0001553
Gas detecting tube (formaldehyde, GASTEC-91L (0.05-1ppm))	00EAA0001554
Gas detecting tube (ammonia, GASTEC-3L (0.5-60ppm))	00EAA0001555

- Easily check filter saturation.




Fume Hood

Small type



Corrosion-resistant and compact design for multi-purposes



PMH-720
with Fume Extraction Module, Plinth (Option)



PMH-720
with Mobile Fume Extractor & Extraction Arm Hood (Option)

Structural Functional Features

- Lightweight, easy-to-move compact tabletop hood.
- Connected to the exhaust system of the facility or connected to a dedicated fume extractor.
- Easy to move and install if connected with exhaust system. (EM-33 or EFS-33)
- Sash door has 5 levels of height adjustment, transparent and made of solid polycarbonate material.
- Integral mold product that does not leak.
- Excellent space utilization through compact design.

Use Convenience Features

- No need for a separate exhaust system by using the dedicated Fume Extraction Module(EM-33) or Mobile Fume Extractor(EFS-33) where use a verified dedicated filter. (option)
- Excellent chemical resistance (polycarbonate material) and safe for handling strong acid and corrosive samples.
- Easy to clean as inner edges are curved.
- Bottom plate included. (standard)

Specification

Model		PMH-720
Required exhaust volume (m ³ /h) ¹⁾		100 / 170 ²⁾
Material	Body	PE (Polyethylene)
	Door / Guide	PC (Polycarbonate)
Dimension	Exterior (W x D x H, mm / inch)	720 x 472 x 535 / 28.3 x 18.6 x 21.1
	Front opening (W x H, mm / inch)	560 x 300 / 22.0 x 11.8
	Duct hole (O.D., mm / inch)	150 / 5.91
	Net weight (kg / lbs)	12 / 26.5
Cat. No.		AAAB5011

1) Air volume could be higher depending on the sash door height adjustment.
2) Exhaust volume measured based on airflow speed of 0.4 m/s for gas and 0.7 m/s for particle.

Accessories

Model	Cat.No.
Fume Extraction Module, EM-33	AAAF1536
Plinth for Fume Extraction Module	AAAF15363
Brackets for Extraction Arm Hood (Φ75)	AAAB6411
Extraction Arm Hood	AAAB6541
Mobile Fume Extractor with Extraction Arm Hood	AAHB6314



Molded Fume Hood

Small type (for observation)



Convenient internal observation and safe removal of generated fumes.



POH-720




Free of the harmful fume can be increased learning effects and focused observation

Transparent design for educational purpose

- The front and back side are made by wide transparent window (PC), useful for simulation and observing in school. (only POH models)
- Excellent space utilization through compact design.
- Free of the harmful fume can be increased learning effect by selecting the filter type.

Integral mold product with compact size

- Integral mold product that does not leak.
- Polyethylene material ensures excellent chemical resistance and safe handling of strong acids and corrosive samples.
- Lightweight, easy-to-move compact tabletop hood
- Connected to the exhaust system of existing facilities or connected to a dedicated harmful gas eliminator.



Secure both the front rail and the rear transparent window with plastic bolts



Select exhaust type depending on space (PMH-720 / POH-720)

Specification

Model		POH-720
Required exhaust volume (m ³ /h) ¹⁾		100 / 170 ²⁾
Material	Body	PE (Polyethylene)
	Door / Guide	PC (Polycarbonate)
Dimension	Exterior (W x D x H, mm / inch)	720 x 472 x 535 / 28.3 x 18.6 x 21.1
	Front opening (W x H, mm / inch)	560 x 300 / 22.0 x 11.8
	Duct hole (O.D., mm / inch)	150 / 5.91
	Net weight (kg / lbs)	12 / 26.5
Cat. No.		AAAB5021

1) Air volume could be higher depending on the sash door height adjustment.
2) Exhaust volume measured based on airflow speed of 0.4 m/s for gas and 0.7 m/s for particle.

Accessories



Model	Cat.No.
Fume Extraction Module, EM-33	AAAF1536
Plinth for Fume Extraction Module	AAAF15363

Filter information for EM-33 [Page 225](#)

Mobile Fume Extractor

Local exhaust type



Optimized configuration according to an usage environment



EFS-33-50D

with Extraction Arm & Hood (option)



EFS-33-75S

with Extraction Arm & Hood (option)



EFS-33-75S

with with Extraction Arm & PMH-720



EFS-33-75S

with Extraction Arm, Hood & Table Bracket Set (Option)

Structural Functional Features

- Excellent space utilization through compact design.
- The body has a low height (650 mm) allowing it to be placed on the bench and on the bottom part of the table.
- Easy to control wind and intuitive display of fan level. (maximum 240m³/h)
- Better safe use of chemicals with durability and spark-free BLDC motor.
- The appropriate filter is selected according to the experiment environment. (option)
- Single / Dual arm, arm diameter, and hood are selected according to the usage environment. (option)


Use Convenience Features

- Pre-filter built in to extend other filter life.
- Up to 2 filters can be installed depending on the material used. (option)
 - ※ Example of filter configuration
 - HEPA Filter for removing dust and carbon filter for removing/neutralizing harmful gas.
- Soft wheels for easy movement and fixing.
- Low fan speed notification when fan speed is low.

Examples Of Diverse Applications


Single Arm

General configuration for 100% suction performance.




Dual Arm

Space-efficient configuration with two hoods used.



Laboratory table integral structure

Maximal space utilization with placement of the lower part of the testing bench.



Specification

Model	EFS-33-50S	EFS-33-50D	EFS-33-75S	EFS-33-75D
Adaptable Arm	Φ50 mm / 2.0 inch 1ea	Φ50 mm / 2.0 inch 2ea	Φ75 mm / 3.0 inch 1ea	Φ75 / 3.0 inch 2ea
Airflow Rate	240 m ³ /h			
Noise Level	60 dB			
Dimension (W x D x H, mm / inch)	460 x 460 x 650 / 18.1 x 18.1 x 25.6			
Net Weight (kg / lbs)	65 / 143.3			
Electric Requirement (230V, 50/60Hz, A)	1.5			
Cat. No.	AAHB6312	AAHB6313	AAHB6314	AAHB6315
Electric Requirement (120V, 60Hz, A)	2.5			
Cat. No.	AAHB6312U	AAHB6313U	AAHB6314U	AAHB6315U

Accessories Page 214 Extraction Arm, Hood, Filter, Table Bracket, Gas Detector & Tubes

Mobile Fume Extractor

Air purifying type



Effective removal of harmful gases and dust in laboratories

Structural Functional Features

- Excellent space utilization through compact design.
- Made of 100% impregnated filter, it is more than twice as efficient as general filters.
- Easy to control wind and intuitive display of fan level. (maximum 240m³/h)
- Better safe use of chemicals with durability and spark-free BLDC motor.
- Verified a cleaning capacity of 2.6 m³/min in an area of 20.3m² according to SPS-KACA002-132 standard.

Use Convenience Features

- Pre-filter built in to extend other filter life.
- Up to 2 filters can be installed depending on the material used. (option)
- HEPA Filter for removing dust and carbon filter for removing/neutralizing harmful gas.
- Soft wheels for easy movement and fixing.
- Low fan speed notification when fan speed is low.



PFS-33



Fan level controller

Specification

Model	PFS-33
Airflow rate	240m ³ /h
Maximum coverage ¹⁾	20.3m ²
Noise level	≤ 60 dB
Fan	BLDC motor, 169W, 1.35A
Electric requirement	230V, 50/60Hz, 1.5A
Dimension (W x D x H, mm / inch)	460 x 460 x 650 / 18.1 x 18.1 x 25.6
Net weight (kg / lbs)	65 / 143.3
Material	Exterior: Epoxy powder coated steel / Interior: Stainless steel 304
Cat. No.	AAHB6411K

1) Measured at 21°C and 45%RH conditions.

Accessories [Page 215](#) Filter

Accessories

Extraction Arm

- Polypropylene with excellent chemical resistance.
- Diameter selected according to work type.
- Joint selected according to the working radius.

	Type	Description
Material	Polypropylene	Excellent chemical resistance
Diameter	50 mm / 2.0 inch	For fine work
	75 mm / 3.0 inch	For general work
Joint	3 Joint	More than 1.5 m working radius

※ Check out the selection guide at [Page 217](#) for details.

Hood

- $\Phi 50$ hood mounted on extraction arm as standard. (selection of white or clear)
- Shape selected according to the nature of the fumes and gases used.
- Material selected according to the work type.
- Strong and chemical resistant using PP and SBC materials.

	Type	Description
Shape	Square	Suction of substances heavier than air on the floor
	Dome	Suction of material lighter than air from the top
Material	Polypropylene, Opaque	Simple suction
	Styrene butadiene copolymer, Transparent	Suction as work progresses

※ Check out the selection guide at [Page 218](#) for details.

LC Magnet System

Convenient to replace

- Strong magnet attached on Arm end(Basic type) and corrosion resistant iron attached on the hood.
- No need the additional assembly.
- Sturdy attach with strong magnet securing the fume leakage.
- Provide each one of Dome type and Square type.

Various hood type

- In case the space is too small to place several arm hoods.
- Hood type can be selected depending on the properties of gas and the diffuse range of fume.



Convenient and quick replacement with strong magnet

Basic type (Arm end)

Table Bracket Set

- Bracket for fixing extraction arm to table and flexible duct hose for connecting to the main body.

Model	Dimension (W x D x H, mm / inch)	Cat. No.
Table Bracket Set ($\Phi 50$ mm / 2.0 inch)	165 x 254 x 110 / 6.5 x 10 x 4.3	AAAB6351
Table Bracket Set ($\Phi 75$ mm / 3.0 inch)	165 x 254 x 110 / 6.5 x 10 x 4.3	AAAB6352

Gas Detector & Tubes




- Determine the exchange period by checking the status of the currently installed filter.


	Item	Cat. No.
Gas Detector	KITAGAWA AP-20	AAAB2531
	Benzene, GASTEC-121SP (0.5-10 ppm)	00EAA0001550
Gas Detecting Tube	Toluene, GASTEC-122L (1-100 ppm)	00EAA0001551
	Acetic acid, GASTEC-81L (0.125-25 ppm)	00EAA0001552
	Chloroform, GASTEC-137 (4-400 ppm)	00EAA0001553
	Formaldehyde, GASTEC-91L (0.05-1 ppm)	00EAA0001554
	Ammonia, GASTEC-3L (0.5-60 ppm)	00EAA0001555

Accessories Filter

Jeio Tech Filter Type and Process

The exclusive filter series developed by Jeio Tech has been objectively verified by an accredited laboratory.


Type	Image	Magnification
Pre-Filter		
HEPA Filter		
Carbon Filter		



Pre-Filter
HEPA Filter
Carbon Filter
Clean air

- ※ The Jeio Tech's carbon filter is an excellent amorphous adsorbent with a high surface area, has fine pores, and is manufactured through a multi-level impregnation process.
- ※ Made of 100% impregnated filter, it is more than twice as efficient as general filters.
- ※ Has the efficiency of an excellent composite filter as it simultaneously performs physical adsorption and chemical adsorption.

Purifying Ability of Jeio Tech Filter for Harmful Gases




Authorized by KOREA CONFORMITY LABORATORIES

Filter	Description	Cat. No.
Complex Carbon	Filters a wide range of chemicals, including acid, ammonia, and aldehyde.	AAAB6354
Acid	Filters acid gases such as HCl, HF, H ₂ S and CH ₃ COOH	AAAB6356
Formaldehyde	Filters formaldehyde, acetate aldehyde, etc.	AAAB6358
Ammonia / Amines	Filters basic gases such as NH ₃ , (CH ₃) ₃ N, etc.	AAAB6359
HEPA	Filters 99.99% of toxic dust and other 0.3 μm particles	AAAB6353
Pre-Filter	Primarily filters particles to extend the life of other filters.	AAAB6360

- ※ Select a suitable filter according to the chemical used.
- ※ For the built-in pre-filter, the Particle is filtered first to extend the life of the other filter.

Extraction Arm Hood

Direct connection to external exhaust system to remove dust and harmful gases



Structural Functional Features

- Powerful exhaust possible with central exhaust system.
- The place to install, extraction arm, hood, column, and bracket are selected according to user's working environment and usage.
- Knobs used for freely adjusting the angle of each joint.
- Easy 360-degree rotation for convenient placement to the desired point.
- Equipped with a damper that can open and close the passage of the arm and adjust the wind speed.
- A dome hood is fitted as standard on the Φ50 arm.

Product Configuration Guide

- Place and surrounding environment desired for installation
 - Column and bracket selected to decide whether to fix arm to wall/ceiling/table in the workspace.
- Working radius
 - 3 Joint, up to a radius of 1.5 m
- Operation type
 - 75mm for general use and 50mm for fine work selected according to the volume to be subjected to suction.
 - In case of suction only, opaque hood / In case of inhalation during operation, transparent hood.
- Properties of used fumes and gases
 - Polypropylene with excellent chemical resistance
 - Dome hood and square hood selected considering the weight compared to air.



Polypropylene with excellent chemical resistance



Knob for angle adjustment



360° rotatable structure




Damper for controlling wind speed

Extraction Arm

	Type	Description
Material	Polypropylene	Excellent chemical resistance
Diameter	50 mm	For fine work
	75 mm ²⁾	For general work
Joint	3 Joint	Radius 1480 ~ 1780 mm / 58.3 ~ 70.1 inch

- 1) Weak against scratching due to the nature of the material.
2) 75mm radius used in the case of mist.



Specification

Model			Extraction Arm, Length (mm / inch)					Weight (kg / lbs)	Cat. No.
			A	B	C	D	Total		
75mm / 3.0 inch diameter, Polypropylene	3 Joints	AH-3160	250 / 9.8	640 / 25.2	530 / 20.9	260 / 10.2	1680 / 66.1	2.4 / 5.3	AAAB6541
	3 Joints	AH-3160C	250 / 9.8	640 / 25.2	530 / 20.9	260 / 10.2	1680 / 66.1	2.4 / 5.3	AAAB6542
50mm / 2.0 inch diameter, Polypropylene	3 Joints	AH5-3160	250 / 9.8	640 / 25.2	540 / 21.3	250 / 9.8	1680 / 66.1	1.4 / 3.1	AAAB6563
	3 Joints	AH5-3160C	250 / 9.8	640 / 25.2	540 / 21.3	250 / 9.8	1680 / 66.1	1.4 / 3.1	AAAB6566

Bracket & Column

- Column of suitable height is selected based on the ceiling height.
- Excellent durability due to powder coating.



Ceiling Bracket



Ceiling Column



Ceiling Column (for textured ceiling)



Wall Bracket

Specification

Model	Dimension (W x D x H, mm / inch)	Cat. No.
Ceiling Bracket	260 x 260 x 2 / 10.2 x 10.2 x 0.08	AAAB6508
	100 x 100 x 250 / 3.9 x 3.9 x 9.8	AAAB6512
Ceiling Column	100 x 100 x 500 / 3.9 x 3.9 x 19.7	AAAB6513
	100 x 100 x 750 / 3.9 x 3.9 x 29.5	AAAB6514
	100 x 100 x 1000 / 3.9 x 3.9 x 39.4	AAAB6515
	100 x 100 x 1250 / 3.9 x 3.9 x 49.2	AAAB6516
	100 x 100 x 1500 / 3.9 x 3.9 x 59.1	AAAB6517
	100 x 250 / 3.9 x 9.8 (∅ x H, mm / inch)	AAAB6518
	100 x 500 / 3.9 x 19.7 (∅ x H, mm / inch)	AAAB6519
Wall Bracket	158 x 259 x 82 / 6.2 x 10.2 x 3.2	AAAB6510

Hood

Image	Type	Material	Description
	Basic type	Polypropylene	<ul style="list-style-type: none"> · When the emission range of the source is narrow. · When the source can be collected at close range.
	Dome type (Clear)	Styrene butadiene copolymer	<ul style="list-style-type: none"> · When the emission range of the source is wide. · Used for capturing light gas.
	Dome type (White)	Polypropylene	
	Square type (Clear)	Styrene butadiene copolymer	<ul style="list-style-type: none"> · When heavy gas is collected, one side is supported on the table and used like a shield.
	Square type (White)	Polypropylene	
	Dome & Square hood of LC Magnet System (Clear)	Styrene butadiene copolymer	<ul style="list-style-type: none"> · Applied a magnet on the connection part, easy to change a hood depending on the gas properties.
	Hood type (White)	Polyethylene	<ul style="list-style-type: none"> · Compact total hood form, securely blocking the user's gas inhalation.



Dome type



Square type



Hood type

Specification

Diameter	Type	Dimension	Color	Material	Cat. No.
Φ75 mm / 3.0 inch	Basic	171 / 6.7 (Ø, mm / inch)	White	Polypropylene	Standard
		280 x 110 / 11 x 4.3 (Ø x H, mm / inch)	Clear	Styrene butadiene copolymer	AAAB6251
	Dome	400 x 140 / 15.7 x 5.5 (Ø x H, mm / inch)	White	Polypropylene	AAAB6221
			Clear	Styrene butadiene copolymer	AAAB6211
	Square	494 x 390 x 245 / 19.4 x 15.4 x 9.6 (W x D x H, mm / inch)	White	Polypropylene	AAAB6241
			Clear	Styrene butadiene copolymer	AAAB6231
LC Magnet System* (Dome & Square)	Dome: 400 x 140 / 15.7 x 5.5 (Ø x H, mm / inch) Square: 494 x 390 x 245 / 19.4 x 15.4 x 9.6 (W x D x H, mm / inch)	Clear	Styrene butadiene copolymer	AAAB6250	
Φ50 mm / 2.0 inch	Dome	294 x 180 / 11.6 x 7.1 (Ø x H, mm / inch)	White	Polypropylene	Clear Dome (Standard) (Optional White Dome)
			Clear	Styrene butadiene copolymer	

※ Ø75 Basic hood is provided as an integral part of the extraction arm and is used for suction of small volumes of material in confined spaces.

※ Ø50 Dome hood is provided as standard on the extraction arm and able to select material when ordering.

* LC Magnet System is a set of arm end, dome hood, and square hood with magnetic attachments so that users can easily replace the hood within seconds.

Recommended configuration of Mobile Fume Extractor for each application

Application	Details	Recommended Hood Type	Hood Position / Airflow Velocity(m/s)	Recommended Filter
Medical / Dentistry / Veterinary	Fumes from laser treatment	50Φ Basic	Side / 0.5	HEPA & Carbon
	Dentist dust and chemical uses	75Φ Basic or Chamber	Side / 0.5	HEPA & Carbon
	Aerosol from anesthetic	75Φ Dome or Square	Side / 0.5	Carbon or HEPA & Carbon
	Use of biopsy chemicals in clinical laboratory	Chamber	Surround / 0.4	Carbon or HEPA & Carbon
Pharmaceuticals	Dust and fume from odorous chemicals in the pharmaceutical laboratory	75Φ Basic	Side / 1.0 Above / 1.2	HEPA or HEPA & Carbon
		Chamber	Surround / 0.7	HEPA
Metal / Acrylic craft	Hazardous fumes and dust from molten metal and welding during the metal cutting process	Basic or Dome	Dust: Side / 1.2 Fume: Side / 0.5	HEPA or HEPA & Carbon
Electronic engineering	Fumes from brazing and soldering	75Φ Dome or 50Φ Basic	Side / 0.5	Carbon or HEPA & Carbon
Laboratory	Science experimental education	75Φ Dome	Dust: Side / 1.2 Fume: Side / 0.5	Carbon or HEPA & Carbon
		Chamber (For Formalin Users)	Surround / 0.4	Carbon or HEPA & Carbon
Chemicals/Biologicals / Agricultural	Hazardous chemicals from acid and base	75Φ Dome or Square	Dust: Side / 1.2 Fume: Side / 0.5	Carbon or HEPA & Carbon
		Chamber (For Formalin Users)	Surround / 0.4	


※ Applied based on the Korea Occupational Safety and Health Act



Safety Storage & Cabinet



Safety Storage & Cabinet



Safety Storage & Cabinet



	Description	Note	Model	Page
Safe-type	Flammable reagents <ul style="list-style-type: none"> • Safe storage of flammable/combustible materials in the event of a surrounding fire • Designed according to the strictest European standards EN 14470-1 	Type 30, 90	SC3-F	226
	Chemical waste <ul style="list-style-type: none"> • Safety funnel having a low possibility of spill during waste disposal • Various color lid to easy check what properties are in a container 	General / Anti-static	SF, WT	234
	Acid/Corrosive <ul style="list-style-type: none"> • Safely store acid/base reagents with corrosion-resistant interior materials • Completed a chemical test of SEFA standard 	Anti-corrosion interior finish	PCC SC-C	235 236
Filter-type	Filtering <ul style="list-style-type: none"> • The harmful gas is filtered by the filter included inside • Economical reagent cabinets that does not require external duct connection 	Closed internal circulation type	FSC	240
By Use	Gas cylinder cabinet <ul style="list-style-type: none"> • Safeguards gas cylinders from external impact • Stores gas cylinders so that they do not fall over 	Gas cylinder	DS-GC	243
General type	Steel window <ul style="list-style-type: none"> • Extremely convenient to check inside storage with glass window door • Convenient drawer-type structure with excellent chemical resistance due to PP material 	Stand type	DS-OW	244

※ The contents of the above and the contents of this catalog may differ depending on the specific model and conditions of use. For the information about the features and specifications that applying to each models, please check the information on the corresponding page of each models.



Satisfaction with the requirements of European regulations and guidelines.



30 minute fire resistance confirmed.



90 minute fire resistance confirmed.



Door locking system with key.



Excellence in acid/alkaline storage with a corrosion-resistant structure and exhaust system.



Leakage minimized using an integrated mold.



Exhaust via duct connection.



Purifies and discharges harmful gas through internal filter.



2 year warranty free A/S.

Rational storage according to purpose



> Storage of flammable / combustible materials

In case of an unexpected fire that may occur at any time, flammable and combustible materials must be stored in a suitable location. Jeio Tech's fire safety storage cabinets certified the strictest EN 14470-1 standard, is the best choice.




> Storage of corrosive substances such as strong acids and strong bases

The structure is optimized for long-term safe storage of corrosive substances such as strong acids and strong bases used in experiments and processes. It is composed of a solid steel structure and corrosion-resistant interior material, and provides differentiated safety performance due to its space specific exhaust structure.



> Reagent cabinet with built-in circulation filter

It is a reagent cabinet in which air flow is uniformly formed and the harmful gas is absorbed through the internal filter. It is safe from fine particles and gas as it is equipped with Jeio Tech's exclusive carbon filter and HEPA filter.



> Physically / chemically sturdy cabinet

It is a reagent cabinet consisting of a polypropylene drawer and a stable structure with a rigid steel frame coated with epoxy. It is equipped with a glass door for easy observation of stored items and a filter exhaust system to maintain a clean interior.



> Safe storage of gas cylinders

Cylindrical high-pressure containers must be prevented from falling and must be properly protected from unexpected external impact. Gas cylinders can be safely stored in Jeio Tech's exclusive solid steel storage cabinets.



blog

Jeio Tech



How to Choose the Right Cabinet for your Lab

> Extraction module (option)

This is harmful gas exhaust equipment (EM-33) equipped with Jeio Tech's exclusive filter. It can be installed and used on Jeio Tech's storage cabinets and reagent cabinet.




Cat. No.	Description (W x D x H, mm / inch)
AAAF1536	342 x 505 x 362.4 / 13.5 x 19.9 x 14.3, 230V, 60Hz

- It is installed at the upper part of the storage cabinet without a separate duct connection to filter and remove the harmful gas inside.
- Pre-filter is included as standard, and up to 2 filters can be added depending on the stored material. (for example: HEPA filter + carbon filter)
- Applicable to Jeio Tech's fire safety storage cabinets, acid/corrosive cabinets, general reagent cabinet series, small fume hood, etc.

> Jeio Tech's exclusive filter

The exclusive filter series developed by Jeio Tech has been objectively verified by an accredited laboratory.

- The Jeio Tech's carbon filter is an excellent amorphous adsorbent with a high surface area, has fine pores and is manufactured through a multi-level impregnation process.
- Made of 100% impregnated filter, it is more than twice as efficient as general filters.
- Has the efficiency of an excellent composite filter as it simultaneously performs physical adsorption and chemical adsorption.



- Electron microscopic photographs of Jeio Tech's carbon filter confirms development of numerous pores such as Macropores (diameter of more than 500Å), Mesopores (diameter of 20-500Å), and Micropores. (diameter of 20Å or less)

Filter	Description	Mounted on EM-33
Complex Carbon	Filters a wide range of chemicals, including acid, ammonia, and aldehyde.	AAAB6354
Acid	Filters acid gases such as HCl, HF, H ₂ S and CH ₃ COOH	AAAB6356
Formaldehyde	Filters formaldehyde, Acete aldehyde, etc.	AAAB6358
Ammonia / Amines	Filters basic gases such as NH ₃ , (CH ₃) ₃ N etc.	AAAB6359
HEPA	Filters 99.99% of toxic dust and other 0.3µ m particles	AAAB6353

- ※ The appropriate filter is selected according to the chemical used.
- ※ Built-in pre-filter primarily filters particles to extend the life of the other filters.

> Chemical test completed according to SEFA-8 standard

Jeio Tech's cabinets and reagent series uses materials that have passed the SEFA (The Scientific Equipment and Furniture Association) standard test, and the equipment has excellent acid resistance, alkali resistance, heat resistance and durability.



Alkali resistance

Exterior Epoxy Powder Coated Steel

The exterior has excellent durability and excellent chemical resistance due to epoxy powder coating.



Drawer Polypropylene (PP)

The drawer part made by solid mold injection is made of polypropylene with excellent chemical resistance.

Fire Safety Storage Cabinet

Type 30, 90

Safeguard flammable liquid from fire with proven fire resistance



Flammable Storage Cabinet
EN 14470-1 : 2004

SC3-90F-0806D1D-C
Standard : Bottom Sump 1ea

SC3-90F-0812D2D-C
Standard : Bottom Sump 2ea

SC3-90F-1906D1-C
Standard : Shelves 3ea & Bottom Sump 1ea
Option : Fume Extraction Module & Filters

SC3-90F-1912D2-C
Standard : Shelves 3ea & Bottom Sump 1ea
Option : Fume Extraction Module & Filters

Achievement Various Certifications

- > Obtained CE (EN 14470-1) certification.
- > Designated as a Excellent R&D Innovation Product in Korea. (K-Innovative Product)
- > Awarded the 2018 Safe Korean Technology Prize by the Ministry of the Interior and Safety.

* EN14470-1 verifies that the equipment meets the European safety standards for dangerous goods containers such as fire resistance, internal ventilation, shelf stability, and automatic door closure.



Awarded the 2018 Safe Korean Technology Prize by the Ministry of the Interior and Safety.

Fire resistance performance 90-minutes type

Type	Model
Tall	SC3-90F-1912D2-C
	SC3-90F-1906D1-C
Under bench	SC3-90F-0812D2D-C
	SC3-90F-0806D1D-C

Fire resistance performance 30-minutes type

Type	Model
Tall	SC3-30F-1912D2-C
	SC3-30F-1906D1-C
Under bench	SC3-30F-0812D2D-C
	SC3-30F-0806D1D-C



Safe Storage From Fire

> High fire resistance

Fire resistance performance is in accordance with the strictest European standard, with the internal temperature kept below 180K for up to 90 minutes in the event of fire.

> Fire response sealing system

When the temperature reaches 50°C, the door automatically closes. When the ambient temperature reaches 70±10°C, the exhaust holes and ventilation openings are automatically closed, totally blocking the flow of heat from the outside in the event of fire.

> Select the right model for your application

Various models are available to choose considering the natural ignition point of the handling hazardous material and the prediction time of fire suppression.

> Chemical resistant interior material

Excellent chemical resistance even for long-term use with the interior material coated with polyphenol resin.

> Under bench type height-adjustable

730 mm/820 mm height can be selected considering the work height when hood is installed at the bottom.

Solid and Secure Shelf Structure

> Solid shelf structure

It is a solid structure allowing for up to 80kg to be held per shelf and it can be fixed firmly to the left and right inside walls.

> Safe shelf design

With 30mm-deep shelf shape for safe storage, it is designed to be safe even if the stored contents are leaked.

> Right shelf selection

The built-in shelf is excellent in chemical resistance due to epoxy powder coating. Also, stainless steel shelves can be selected. (option)

Patented Safe Exhaust System

The optimized exhaust system that the outside flame does not inflow to the inside of the cabinet and removing the flammable vapor effectively usual.



Award Excellent Patent Award (in Korea)

Proven Fireproof Performance!

Items proven by testing


- > 30 minutes / 60 minutes / 90 minutes fire resistance performance test
- > For fire test according to standards of EN 1363-1, the time required for ΔT to be kept below 180°C is expressed as fire resistance. (EN 14470-1)



Real fire simulation test for fire resistance performance 90min type



Comparison of fire resistance performance based on certification



SC3-90F-1912D2-C fire test results

➤ **Convenient reagent In/Out**

Equipped with door closure to stop the door at the desired position for easy put flammable reagents in and out.



➤ **Automatically sealed ventilation**

When the ambient temperature reaches $70 \pm 10^\circ\text{C}$ in case of surrounding fire, the exhaust opening, and ventilation opening are designed to be automatically closed to completely prevent the inflow of external heat.




➤ **Anti-static charge through ground**

The ground to prevent static spark accidents caused by inside the cabinet where flammable vapors may be generated. For proper installation, each shelf should be connected with busbar inside of cabinet. Also the conduction bolt (located on upper part behind the cabinet) should be connected with distribution panelboard. (In case of ground can not be connected with distribution panelboard, please use grounded plug. (Option))



➤ **Safe ventilation of flammable vapor**

It is possible to discharge gas easily through the vent inlet and vent outlet to prevent accidents due to accumulation of flammable vapor, and to protect the health of the experimenter.



➤ **Configurable with lab safety workbench**

It is composed of a fume hood and a testing bench bottom, making it excellent for storing hazardous materials and liquid wastes.





Effluent reservoir

Safe handling in the event of stored material leakage and prevention of leakage out of the cabinets.



Key lock device

Key lock is provided as standard to block unauthorized access.



Safe stored materials handling

Ergonomic structure with easy access to shelves for user's safety and convenience.



Convenient moving and installation

Forklift space for convenient movement. Leveling feet for stable installation.





Extraction module (option)

Filters flammable gas inside the cabinet with a carbon filter without connecting a separate duct.







Select a suitable shelf

Providing the 3ea of shelf (up to 80kg storage) for tall type cabinet, and can be exchanged to drawer. (option)

	EN 14470-1	NFPA 30
Fire resistance performance	<p>O</p> <p>Fire test according to EN 1363-1, 15 min.- 90 min. (15 min. / 30 min. / 60 min. / 90 min.)</p>	<p>O</p> <p>NFPA 251 standard fire resistance test, 10 minute level</p>
Flammable vapor ventilation	<p>O</p> <p>Flammable vapor ventilation is recommended. Check negative pressure formation inside the cabinet during ventilation.</p> <div style="text-align: center;">  </div> <p>Keeps flammable vapors below the LEL concentration and removes smell(odors) through ventilation.</p>	<p>X</p> <p>Ventilation is not recommended. (Recommended to use level of container recommended in NFPA 30)</p> <div style="text-align: center;">  </div> <p>It seals the ventilation holes, and an smell(odors) occurs inside the storage cabinet. The risk of ignition is increased as the concentration of flammable vapor increases.</p>
Fire detection, automatic door and ventilation closure	<p>O</p> <p>In the case of a fire test, test with the ventilation hole opened. Automatic door and ventilation closure during fire test.</p>	<p>X</p> <p>In the case of a fire test, test with the ventilation hole blocked. Fire safety not considered when connecting to external ventilation facilities.</p>

Comparison of EN 14470-1 and NFPA 30 International Standards

Specification (Type 90)

Model	Tall type		Underbench type	
	SC3-90F-1912D2-C	SC3-90F-1906D1-C	SC3-90F-0812D2D-C	SC3-90F-0806D1D-C
Internal volume (L / cu ft)	763 / 26.9	357 / 12.6	244 / 8.6	110 / 3.9
Interior dimensions (W x D x H, mm / inch)	968 x 484 x 1629 / 38.1 x 19.1 x 64.1	453 x 484 x 1629 / 17.8 x 19.1 x 64.1	496 x 396 x 622 * 2ea / 19.5 x 15.6 x 24.5 * 2ea	448 x 396 x 622 / 17.6 x 15.6 x 24.5
Exterior dimensions (W x D x H, mm / inch)	1116 x 625 x 1952 / 43.9 x 24.6 x 76.9	601 x 625 x 1952 / 23.7 x 24.6 x 76.9	1200 x 560 x 730 (or 820) ¹⁾ / 47.2 x 22.0 x 28.7 (or 32.3)	600 x 560 x 730 (or 820) ¹⁾ / 23.6 x 22.0 x 28.7 (or 32.3)
Quantity of shelf (standard / max.) ²⁾	3 / 5	3 / 5	-	-
Bottom collecting sump ³⁾ (Volume : L / cu ft)	42 / 1.48	22 / 0.78	24 / 0.85 x 2ea	22 / 0.78
Max. load per shelf (kg / lbs)	80 / 176.3	80 / 176.3	40 * 2ea / 88.2 * 2ea ⁴⁾	40 / 88.2 ⁴⁾
Weight (kg / lbs)	385 / 848.8	270 / 595.2	220 / 485.0	134 / 295.4
CE certification	O	O	O	O
Cat. No.	AAAF390170	AAAF390140	AAAF390130	AAAF390120
Product Exterior				
Product interior				

- 1) The height of underbench type is 820mm(32.3 inch) if you select cabinet plinth set accessory. (easy to move with forklift and hand truck, option)
- 2) Epoxy powder coated shelves are provided(standard), able to change to stainless steel, drawer type or add quantity.(option)
- 3) At the bottom of the cabinet, a collecting sum set is given as standard (epoxy powder coated steel material), and able to change to stainless steel.
For one drawer underbench model (0806D1D), it is able to mount a container guard(option) to fit waste container.
- 4) Max. load per drawer.

Accessories (Type 90)

Accessory	Description (W x D x H, mm / inch)	SC3-90F-1912D2-C	SC3-90F-1906D1-C	SC3-90F-0812D2D-C	SC3-90F-0806D1D-C
Shelf	Epoxy powder coated steel	AAAF35014 (965 x 444 x 30 / 38.0 x 17.5 x 1.2)	AAAF35013 (451 x 444 x 30 / 17.8 x 17.5 x 1.2)	-	-
	Stainless steel	AAAF35004 (965 x 444 x 30 / 38.0 x 17.5 x 1.2)	AAAF35003 (451 x 444 x 30 / 17.8 x 17.5 x 1.2)	-	-
Drawer	Epoxy powder coated steel	AAAF35313 (798 x 452 x 100 / 31.4 x 17.8 x 3.9)	AAAF35312 (350 x 452 x 100 / 13.8 x 17.8 x 3.9)	-	-
	Stainless steel	AAAF35303 (798 x 452 x 100 / 31.4 x 17.8 x 3.9)	AAAF35302 (350 x 452 x 100 / 13.8 x 17.8 x 3.9)	-	-
Plinth (bottom) ¹⁾	Epoxy powder coated steel (Fireproof paint)	-	-	AAAF35412 (1200 x 510 x 90 / 47.2 x 20.1 x 3.5)	AAAF35411 (600 x 510 x 90 / 23.6 x 20.1 x 3.5)
Plinth (side) ¹⁾	340 x 480 x 350 / 13.4 x 18.9 x 13.8	-	-	AAAF15364	AAAF15364
Fume Extraction Module (EM-33)	342 x 505 x 362.4 / 13.5 x 19.9 x 14.3 (230V, 60Hz)	AAAF1536	AAAF1536	AAAF1536	AAAF1536
Container Guard	403 x 350 x 141 / 15.8 x 13.7 x 5.5	-	-	-	AAAF35454

※ Please refer to page 225 for details on EM-33 compatible filter.


- 1) "Plinth (bottom)" is an accessory that raises the height of the underbench type, "Plinth (side)" is an accessory for installing the exhaust module (EM-33) on the underbench type.

Dimension


(Unit: mm / inch)

	Model	A	B	C ¹⁾	D (Diameter)
Tall	SC3-90F-1912D2-C	170 / 6.7	160 / 6.3	1135 / 44.7	76 / 3.0
	SC3-90F-1906D1-C	170 / 6.7	160 / 6.3	1135 / 44.7	76 / 3.0
Underbench	SC3-90F-0812D2D-C	185 / 7.3	255 / 10.0	840 / 33.1	60 / 2.4
	SC3-90F-0806D1D-C	165 / 6.5	255 / 10.0	840 / 33.1	60 / 2.4


1) C is the maximum depth when the door/drawer is opened.




Tall:
Upper side



Underbench:
Upper side



Underbench:
Back side



Underbench equipped
with EM-33 and Plinth (side)

Configuration with Safety product of Jeio Tech

Ductless Fume Hood

- > Built-in filter hood that no need the duct connection. Easy to move and install, and designed by the international standard to satisfy performance.
- > Configuration with underbench type cabinet.



Detail [Page 208](#)




Safety Waste Disposal System

- > Safe waste disposal through anti-static and spill-proof design, also removal of internal gases.
- > Dual safety with the configuration of the Underbench type cabinet.



Detail [Page 234](#)

Specification (Type 30)

Model	Tall type		Underbench type	
	SC3-30F-1912D2-C	SC3-30F-1906D1-C	SC3-30F-0812D2D-C	SC3-30F-0806D1D-C
Internal volume (L / cu ft)	802 / 28.3	375 / 13.2	252 / 8.9	116 / 4.1
Interior dimensions (W x D x H, mm / inch)	968 x 507 x 1635 / 38.1 x 20.0 x 64.4	453 x 507 x 1635 / 17.8 x 20.0 x 64.4	506 x 415 x 602 * 2ea / 19.9 x 16.3 x 23.7 * 2ea	468 x 415 x 602 / 18.4 x 16.3 x 23.7
Exterior dimensions (W x D x H, mm / inch)	1101 x 630 x 1951 / 43.3 x 24.8 x 76.8	586 x 630 x 1951 / 23.1 x 24.8 x 76.8	1200 x 560 x 730 (or 820) ¹⁾ / 47.2 x 22.0 x 28.7 (or 32.3)	600 x 560 x 730 (or 820) ¹⁾ / 23.6 x 22.0 x 28.7 (or 32.3)
Quantity of shelf (standard / max.) ²⁾	3 / 5	3 / 5	-	-
Max. load per shelf (kg / lbs)	80 / 176.4	80 / 176.4	40*2ea / 88.2*2ea ³⁾	40 / 88.2 ³⁾
Bottom collecting sump ⁴⁾ (Volume : L / cu ft)	42 / 1.48	22 / 0.78	26 / 0.92 x 2ea	24 / 0.85
CE Certification	O	O	O	O
Weight (kg / lbs)	330 / 727.5	220 / 485.0	200 / 440.9	120 / 264.6
Cat. No.	AAAF330470	AAAF330440	AAAF330430	AAAF330420
Product Exterior				
Product interior				

1) The height of underbench type is 820mm(32.3 inch) if you select cabinet plinth set accessory. (easy to move with forklift and hand truck, option)

2) Epoxy powder coated shelves are provided(standard), able to change to stainless steel or drawer type, able to add quantity.(option)

3) Max. load per drawer.

4) At the bottom of the cabinet, a collecting sum set is given as standard (epoxy powder coated steel material), and able to change to stainless steel. For one drawer underbench model (0806D1D), it is able to mount a container guard(option) to fit waste container.

Accessories (Type 30)

Accessory	Description (W x D x H, mm / inch)	SC3-30F-1912D2-C	SC3-30F-1906D1-C	SC3-30F-0812D2D-C	SC3-30F-0806D1D-C
Shelf	Epoxy powder coated steel	AAAF35014 (965 x 444 x 30 / 38.0 x 17.5 x 1.2)	AAAF35013 (451 x 444 x 30 / 17.8 x 17.5 x 1.2)	-	-
	Stainless steel	AAAF35004 (965 x 444 x 30 / 38.0 x 17.5 x 1.2)	AAAF35003 (451 x 444 x 30 / 17.8 x 17.5 x 1.2)	-	-
Drawer	Epoxy powder coated steel	AAAF35313 (798 x 452 x 100 / 31.4 x 17.8 x 3.9)	AAAF35312 (350 x 452 x 100 / 13.8 x 17.8 x 3.9)	-	-
	Stainless steel	AAAF35303 (798 x 452 x 100 / 31.4 x 17.8 x 3.9)	AAAF35302 (350 x 452 x 100 / 13.8 x 17.8 x 3.9)	-	-
Plinth (bottom) ¹⁾	Epoxy powder coated steel (Fireproof paint)	-	-	AAAF35412 (1200 x 510 x 90 / 47.2 x 20.1 x 3.5)	AAAF35411 (600 x 510 x 90 / 23.6 x 20.1 x 3.5)
Plinth (side) ¹⁾	340 x 480 x 350 / 13.4 x 18.9 x 13.8	-	-	AAAF15364	AAAF15364
Fume Extraction Module (EM-33)	342 x 505 x 362.4 / 13.5 x 19.9 x 14.3 (230V, 60Hz)	AAAF1536	AAAF1536	AAAF1536	AAAF1536
Container Guard	423 x 370 x 141 / 16.6 x 14.5 x 5.5	-	-	-	AAAF35455

※ Please refer to page 225 for details on EM-33 compatible filter.


1) "Plinth (bottom)" is an accessory that raises the height of the underbench type, "Plinth (side)" is an accessory for installing the exhaust module (EM-33) on the underbench type.

Dimension


(Unit: mm / inch)

	Model	A	B	C ¹⁾	D (Diameter)
Tall	SC3-30F-1912D2-C	180 / 7.1	170 / 6.7	1095 / 43.1	76 / 3.0
	SC3-30F-1906D2-C	180 / 7.1	170 / 6.7	1120 / 44.1	76 / 3.0
Underbench	SC3-30F-0812D2D-C	160 / 6.3	235 / 9.3	840 / 33.1	60 / 2.4
	SC3-30F-0806D1D-C	160 / 6.3	235 / 9.3	840 / 33.1	60 / 2.4


1) C is the maximum depth when the door/drawer is opened.




Underbench: Upper side



Underbench: Back side



Tall: Upper side



Underbench equipped with EM-33 and Plinth (side)

Safety Waste Disposal System



Easy to collect and store various chemicals and liquid waste from the lab.

> Safe device to prevent leakage

- Safety funnel dedicated to waste fluid treatment, which is convenient for handling due to wide funnel inlet opening.
- 80% level indicator is equipped as standard to easily know when to replace the waste container, and can be changed to 70% level indicator if required by user. (container height: 380mm)

> Classified emission according to the type of chemical properties

- Different colors of the safety funnel lid depending on the chemical properties to clearly distinguish when treating liquid wastes. (SF-14P (Persian green) / -14B (Blue) / -14Y (Yellow) / 14R (Red))
- Safety labels can be attached on the top of the safety funnel to prevent the mixing of waste fluids with different properties.

> Pressure and gas treatment by reaction

- Hinge-type lid minimizes the discharge of vapors from the waste solution into the lab.
- Removes bad odors and harmful gas. (carbon filter can be replaced easily without remove existing attachment without the use of any tools) (option)
- Easy to check the replacement time for a filter with timestrip.

> Anti-static protection

- Anti-Static protection via special plastic material and ground cable. (SF model)
- Easy to organize ground wires when configured with a flammable material storage cabinet. (SC3-90F Drawer type)
- The ground to prevent static spark accidents caused by inside the cabinet where flammable vapors may be generated.




SF-14P & WT-20P

with Carbon Filter (option)

SF-14C & WT-20C


with Carbon Filter (option)



Comparison of safety waste disposal funnel.



Wide inlet of waste liquid funnel.



Waste liquid level notification function.



Classified color label and lid color according to chemical properties.



Hinge-type lid.



A dedicated flammable material storage cabinets to contain waste disposal system. (only for SC3-90F drawer type)



Grounded Plug.

Specification

Model	SF-14C	SF-14P	WT-20C	WT-20P
Description	Anti-static Safety funnel for collection of waste liquid	General use Safety funnel for collection of waste liquid	Anti-static liquid waste container	General-purpose liquid waste container
Dimension (W x D x H, mm / inch)	160 x 142 x 270 / 6.3 x 5.6 x 10.6	160 x 142 x 270 / 6.3 x 5.6 x 10.6	220 x 300 x 380 / 8.7 x 11.8 x 15.0	220 x 300 x 380 / 8.7 x 11.8 x 15.0
Inlet (ø, mm / inch)	135 / 5.3	135 / 5.3	45 / 1.8	45 / 1.8
Lid size (ø, mm / inch)	142 / 5.6	142 / 5.6	78 / 3.1	78 / 3.1
Volume (L / cu ft)	1 / 0.04	1 / 0.04	22 / 0.78	22 / 0.78
Material	PP	PP	HDPE	HDPE
Cat. No.	AAAM6001	AAAM6002	AAAM6011	AAAM6012

Accessories

Description	Dimension (ø x H, mm / inch)	Cat. No.
Carbon Filter	76 x 150 / 3.0 x 5.9	AAAM6051
Earth Bonding Plug	52 x 70 / 2.0 x 2.8	AAAF35453

※ In case of ground can not be connected with distribution panelboard, please use grounded plug. (option)



Acid Corrosive Cabinet Small type



Compact size enables efficient space utilization

- **Safe storage from corrosive substances**

 - Full injection molding with excellent chemical resistance for durability and corrosion resistance. (body: PE, door: PVC, hinge: PP)
 - Suitable for storage of small amounts of corrosive reagents. (Max. 18ea based on a 1L reagent basis)
 - Inside of cabinet is designed by with no steel structure for corrosion protection.
- **Convenient and safe ventilation**

 - Use duct connection of existing facility or dedicated hazardous gas removal Fume Extraction Module (EM-33 option) connection.
 - Dedicated the hazardous gas removal module(EM-33) can be selected the validated filter from Jeio Tech. (EM-33 and filter is option)



PCC-720

with Fume Extraction Module (EM-33)(option) + Plinth (option)

Specification

Model		PCC-720
Material	Body / Door	PE (Polyethylene)
	Hinge	PP (Polypropylene)
Dimension	Exterior (W x D x H, mm / inch)	720 x 472 x 535 / 28.3 x 18.6 x 21.1
	Front opening (W x D, mm / inch)	560 x 320 / 22.0 x 12.6
	Duct hole (O.D., mm / inch)	150 / 5.9
	Net weight (kg / lbs)	12.5 / 21.4
Cat. No.	AAAF61001	



Accessories

Model	Cat.No.
Plinth for Fume Extraction Module	AAAF15363
Brackets for Extraction Arm Hood (ø75)	AAAB6411
Fume Extraction Module, EM-33	AAAF1536
Bottom Plate	AAAF61501

For detail information of EM-33 and its filter.

[Page 225](#)

Acid/Corrosive Cabinet

Safety cabinets for long-term storage of acid/base substances



SC-C1912D4

with Shelves 4ea, Tray 8ea (standard)


SC-C1906D2

with Shelves 2ea, Tray 4ea (standard),
Fume Extraction Module & Filters (option)



Filter-type harmful gas exhaust (option)

Filters harmful gas inside the cabinet with a carbon filter without connecting a separate duct. (see page 225 for details)



Key lock device

Key lock is provided as standard to block unauthorized access.



Corrosion-resistant PP Tray

Spill tray to prevent leakage of stored materials.



Safe stored materials handling

Ergonomic structure with easy access to shelves for user's safety and convenience.



Convenient moving and installation

Forklift space for convenient movement.
Leveling feet for stable installation.



Identification pictogram stickers provided

Easy to classify stored materials with the provision of removable acid/base stickers.

Safe Storage from Corrosive Substances

- EN 14727 standards test completed**
 Reliable acid/base material container that has passed laboratory cabinet testing for performance and durability according to the strictest European standards.
- Completed SEFA standards chemical tests**
 Verified suitability for long-term storage of acid/base materials with testing of corrosion-resistant interior material, polypropylene shelf and exhaust outlet coated with polyphenol resin on the inside.
- Solid steel structure**
 Stable construction of steel frame prevents structural deformation, and it has excellent external chemical resistance due to epoxy coating.
- No steel structure inside storage cabinet**
 Stainless steel shelves, shelf racks, bolts, door catches, exhausts etc. are not used inside to prevent internal corrosion.
- Space-efficient internal compartments**
 The internal compartment structure allows materials stored in one reagent cabinet to be classified according to their characteristics, allowing for safe storage. (except SC-C-0706D1)

Strong and Secure Shelf Structure

- Strong shelf structure**
 It is a solid structure allowing for up to 50kg to be held per shelf and it can be fixed firmly to the left and right inside walls.
- Safe shelf design**
 The shelf tray with a minimum depth of 30mm is made of polypropylene and is designed to be safe even when storage spills out. The corresponding tray is provided as standard for each model.


Convenient and Safe Ventilation System

- Ventilation opening per inner compartment**
 It is possible to achieve maximum exhaust with there being a ventilation opening inside each compartment, preventing accumulation of gas.
- Connection of exhaust system (option)**
 The exhaust port can be conveniently connected to the exhaust system, making it easy to remove harmful gas.

Safe storage with separate gas discharge path according to storage characteristics

Outlets for contaminated air are designed for each internal space, so there is no heterogeneous gas mixing inside the storage compartment, and safe storage is possible as materials are stored according to their characteristics.

(Patented in Korea)







Patented

JEIO TECH



Specification (Tall Type)

Model	SC-C-1912D4	SC-C-1912D2	SC-C-1906D2	SC-C-1906D1
Internal volume (L / cu ft)	852 (213 x 4ea) / 30.1 (7.5 x 4ea)	894 (447 x 2ea) / 31.6 (15.8 x 2ea)	426 (213 x 2) / 15.0 (7.5 x 2ea)	445 / 15.7
Interior dimensions (W x D x H, mm / inch)	561 x 466 x 818 * 4ea / 22.1 x 18.3 x 32.2 * 4ea	561 x 466 x 1712 * 2ea / 22.1 x 18.3 x 67.4 * 2ea	561 x 466 x 818 * 2ea / 22.1 x 18.3 x 32.2 * 2ea	558 x 466 x 1712 / 22.0 x 18.3 x 67.4
Exterior dimensions (W x D x H, mm / inch)	1229 x 555 x 1900 / 48.4 x 21.9 x 74.8	1229 x 555 x 1900 / 48.4 x 21.9 x 74.8	655 x 555 x 1900 / 25.8 x 21.9 x 74.8	655 x 555 x 1900 / 25.8 x 21.9 x 74.8
Quantity of shelf (standard / max.)	4 / 8	6 / 12	2 / 4	3 / 6
Quantity of tray (standard / max.)	8 / 12	8 / 14	4 / 6	4 / 7
Distance of between shelves (mm / inch)	120 / 4.7	100 / 3.9	130 / 5.1	100 / 3.9
Max. load per shelf (kg / lbs)	50 / 110.2	50 / 110.2	50 / 110.2	50 / 110.2
Weight (kg / lbs)	284 / 626.1	257 / 566.6	160 / 352.7	146 / 321.9
Cat. No.	AAAF60043	AAAF60025	AAAF60023	AAAF60011
Product Exterior				
Product interior				

※ If combustible vapor is delayed without ventilation during storage of corrosive materials such as acid/base materials, this may affect the durability of the product, or the inside air may become contaminated. Therefore, it is recommended to connect to the ventilation system of the building or to a separate exhaust system.

※ There is an exhaust vent at the top of the cabinets, so it can be used by connecting to the forced exhaust or purifier.


Accessories




Tray (PP)



Extraction Module (option)



Underbench type with EM-33 and plinth (side)







Tall type with EM-33 and plinth (top)

Accessories	Description (W x D x H, mm / inch)	SC-C-1912D4	SC-C-1912D2	SC-C-1906D2	SC-C-1906D1
Tray	545 x 419 x 33 / 21.5 x 16.5 x 1.3 (Polypropylene)	00SCF0000037 (8 included as standard)	00SCF0000037 (8 included as standard)	00SCF0000037 (4 included as standard)	00SCF0000037 (4 included as standard)
Fume extraction module (EM-33)	342 x 505 x 362.4 / 13.5 x 19.9 x 14.3 (230V, 60Hz)	AAAF1536	AAAF1536	AAAF1536	AAAF1536
Plinth (Top)¹⁾	852 x 479 x 51 / 33.5 x 18.9 x 2.0	AAAF15362	AAAF15362	-	-

1) Plinth (Top) is required if two upper exhaust ports are exhausted to one Fume Extraction Module.

※ Please refer to page 225 for details on EM-33 compatible filter.

Specification (Underbench Type)

Model	SC-C-0812D2	SC-C-0809D2	SC-C-0806D1	SC-C-0706D1
Internal volume (L / cu ft)	312 (156 x 2ea) / 11.0 (5.5 x 2ea)	226 (113 x 2ea) / 8.0 (4.0 x 2ea)	146 / 5.2	115 / 4.1
Interior dimensions (W x D x H, mm / inch)	543 x 456 x 634 * 2ea / 21.4 x 18.0 x 25.0 * 2ea	393 x 456 x 634 * 2ea / 15.5 x 18.0 x 25.0	505 x 456 x 634 / 19.9 x 18.0 x 25.0	505 x 421 x 542 / 19.9 x 18.0 x 21.3
Exterior dimensions (W x D x H, mm / inch)	1200 x 548 x 820 / 47.2 x 21.6 x 32.3	900 x 548 x 820 / 35.4 x 21.6 x 32.3	600 x 548 x 820 / 23.6 x 21.6 x 32.3	600 x 510 x 650 / 23.6 x 20.1 x 25.6
Quantity of tray	2	2	1	1
Max. load per tray (kg / lbs)	50 / 110.2	50 / 110.2	50 / 110.2	50 / 110.2
Weight (kg / lbs)	110 / 242.5	92 / 202.8	65 / 143.3	52 / 114.6
Cat. No.	AAAF60122	AAAF60121	AAAF60113	AAAF60012
Product Exterior				
Product interior				

※ If combustible vapor is delayed without ventilation during storage of corrosive materials such as acid/base materials, this may affect the durability of the product, or the inside air may become contaminated. Therefore, it is recommended to connect to the ventilation system of the building or to a separate exhaust system.

Accessories


Accessory	Description (W x D x H, mm / inch)	SC-C-0812D2	SC-C-0809D2	SC-C-0806D1	SC-C-0706D1
Tray	PP (470 x 370 x 100 / 18.5 x 14.6 x 3.9)	00CBN0000010	-	00CBN0000010	00CBN0000010
	PVC (390 x 450 x 37 / 15.3 x 17.7 x 1.5)	-	00SCF0000100	-	-
Fume extraction module (EM-33)	342 x 505 x 362.4 / 13.5 x 19.9 x 14.3 (230V, 60Hz)	AAAF1536	AAAF1536	AAAF1536	AAAF1536
Plinth (side)	340 x 480 x 350 / 13.4 x 18.9 x 13.8	AAAF15364	AAAF15364	AAAF15364	AAAF15364

※ Please refer to the page 225 for details on EM-33 compatible filter.


Dimension

(Unit: mm / inch)


Model	A	A'	B	C	D (Dia.)	
Tall	SC-C-1912D4	315 / 12.4	315 / 12.4	75 / 3.0	1110 / 43.7	102 / 4.0
	SC-C-1912D2	315 / 12.4	315 / 12.4	75 / 3.0	1110 / 43.7	102 / 4.0
	SC-C-1906D2	325 / 12.8	-	75 / 3.0	1120 / 44.1	102 / 4.0
	SC-C-1906D1	325 / 12.8	-	75 / 3.0	1120 / 44.1	102 / 4.0
Underbench	SC-C-0812D2	130 / 5.1	520 / 20.5	605 / 23.8	1190 / 46.9	102 / 4.0
	SC-C-0809D2	130 / 5.1	370 / 14.5	605 / 23.8	940 / 37.0	102 / 4.0
	SC-C-0806D1	130 / 5.1	-	605 / 23.8	1060 / 41.7	102 / 4.0
	SC-C-0706D1	130 / 5.1	-	435 / 17.1	1060 / 41.7	102 / 4.0




Option: Plinth (side)+EM-33



Tall: Upper side



Underbench: Upper side



Underbench: Back side

Filtering Storage Cabinet

The internal carbon filter removes harmful gas in the reagent cabinet to protect the experimenter




FSC-140


with Drawer Shelves 5ea (standard)

Completed hazardous gas removal capability/capacity verification and door durability test.

- Ability to remove harmful gases / Capacity Test. (SPS-KACA 002-132)
Completed verification of high-efficiency hazardous gas rates.
Apply composite carbon filter. (VOC, acid, base)
- Door durability test. (KS G 2020 specification)
40,000 open/close test verification completed.



FSC-140



FSC-280



Easy to use as a drawer-type shelf
Drawer-type shelf for convenient handling of inner samples installed as standard.



Brightness control LED lamp
Easy to observe the inside with LED lamp, and possible to adjust the brightness to 5 levels.



Safe with effluent storage container
Safe handling in the event of stored material leakage and prevention of leakage out of the reagent cabinet.



Structure convenient for moving and installation
Structure allows it to be moved and fixed and it can be also be safely installed horizontally.



Key lock device
Key lock is provided as standard to block unauthorized access.



Possible to two filters installed at the same time
Two filters can be installed at the same time, except the pre-filter provided as standard.



▮ Suitable for Storing Sealed Reagents

- **No external exhaust system connection required**
A filter is installed inside the cabinet to remove harmful gas using internal air circulation.
- **Air-conditioning energy reduction**
It is a closed structure that does not allow the inside air to be discharged through the exhaust device, thus saving energy used for indoor cooling/heating.
- **Efficient fan speed control**
Convenient 0~5 step fan speed control function increases filter life and minimizes energy consumption.
- **Turbo mode for pleasant use**
In order to minimize exposure to harmful gas, the fan speed is temporarily maximized when opening the door.

▮ Outstanding Safety

- **Internal temperature rise warning**
High temperature warning alarm function in case of temperature rise inside the reagent.
- **Air circulation detection and warning**
Notification of internal air flow abnormality with detection of fan operation.
- **Sturdy epoxy powder coating**
Excellent chemical and corrosion resistance and durability even for long-term use.
- **External sensor connection possible**
Can be connected with external sensor (option) so that alarm is generated when the set limit is exceeded.

▮ Various Functions Convenient for Management

- **Efficient use with sleep mode**
Operates only at the set time (for example, working time period) to save energy and extend the filter life.
- **Displays various information**
Displays useful information such as VOC, temperature, humidity, filter check notification, etc. in real time.
- **Filter configuration customizable for application**
Various custom filter configurations supported depending on the chemical nature of the reagents to be stored. (see page 220 for details)
- **Structure allows for convenient filter management**
It is very convenient to check filter saturation and replace by opening the front part of the reagent cabinet.
- **Check filter saturation with Gas Detector**
Saturation of the complex carbon filter is easily checked using dedicated gas detector. (option)
- **Filter check notification function**
Filter maintenance notification function of 1-24 month cycle makes it easy to manage the filter according to the usage environment.
- **Filter history management function**
Supports efficient filter management by saving filter usage time and checking history.
- **Recording system for filter management**
Records the filter check time, material name and concentration value in the controller.



Gather information at a glance





Simple replacement



Checking filter saturation

Specification

Model		FSC-140	FSC-280
Internal volume (L / cu ft)		604 / 21.3	1317 / 46.5
Capacity (L or quantity of 1 L bottle)		140	280
Airflow rate (m ³ /min)		0.67	1.34
Noise level (dB)		< 50	< 55
Safety	Warning / Alarm	High Temperature, Door Open, Fan Fail, Sensor Error Check Filter, External Input Error	
	Protection	Over Temperature Protection, Over Current Protection	
Dimensions	External (W x D x H, mm / inch)	930 x 636 x 1890 / 36.6 x 24.8 x 74.4	1868 x 636 x 1890 / 73.5 x 24.8 x 74.4
	Internal (W x D x H, mm / inch)	795 x 578 x 1315 / 31.3 x 22.8 x 51.8	1733 x 578 x 1315 / 68.2 x 22.8 x 51.8
	Drawer (W x D x H, mm / inch)	613 x 478 x 109 / 24.1 x 18.8 x 4.3	613 x 478 x 109 / 24.1 x 18.8 x 4.3
	Spill tray (W x D x H, mm / inch)	670 x 526 x 34 (1ea) / 26.4 x 20.7 x 1.3 (1ea)	670 x 526 x 34 (2ea) / 26.4 x 20.7 x 1.3 (2ea)
Weight (kg / lbs)		300 / 661.4	495 / 1091.3
Quantity of drawer (standard / max.)		5 / 10	10 / 20
Max. load capacity per drawer (L or quantity of 1 L bottle)		28	28
Max. load per drawer (kg / lbs)		40 / 88.2	40 / 88.2
Electrical requirements (230V)		50/60Hz, 0.5A	50/60Hz, 1.0A
Cat. No		AAHF1122K	AAHF1132K
Electrical requirements (120V)		60Hz, 1.0A	60Hz, 2.0A
Cat. No		AAHF1123U	AAHF1133U
		 <p>Product appearance Product interior</p>	 <p>Product appearance Product interior</p>

Accessories

Jeio Tech's Filter Series

- The Jeio Tech's carbon filter is an excellent amorphous adsorbent with a high surface area, has fine pores and is manufactured through a multi-level impregnation process.
- Made of 100% impregnated filter, it is more than twice as efficient as general filters.
- Has the efficiency of an excellent composite filter as it simultaneously performs physical adsorption and chemical adsorption.



Filter	Description	Cat. No.
Complex Carbon	Filters a wide range of chemicals, including acid, ammonia, and aldehyde	AAAF1537
Acid	Filters acid gases such as HCl, HF, H ₂ S and CH ₃ COOH	AAAF1539
Formaldehyde	Filters formaldehyde, Acete aldehyde, etc.	AAAF1541
Ammonia / Amines	Filters basic gases such as NH ₃ , (CH ₃) ₃ N etc.	AAAF1542
HEPA	Filters 99.99% of toxic dust and other 0.3µm particles	00FSC0000083

- ※ Above filter should be ordered two when you apply the FSC-280.
- ※ The appropriate filter is selected according to the chemical used.
- ※ Built-in pre-filter primarily filters particles to extend the life of the other filters.

Gas Detector & Tube & Drawer Set

- Check the condition of the Complex Carbon filter for the currently used chemical and determine the replacement time of the filter.



	Item	Cat. No.
Gas Detector	KITAGAWA AP-20	AAAB2531
	Benzene, GASTEC-121SP (0.5~10 ppm)	00EAA0001550
Gas Detecting Tube	Toluene, GASTEC-122L (1~100 ppm)	00EAA0001551
	Acetic acid, GASTEC-81L (0.125~25 ppm)	00EAA0001552
	Chloroform, GASTEC-137 (4~400 ppm)	00EAA0001553
	Formaldehyde, GASTEC-91L (0.05~1 ppm)	00EAA0001554
	Ammonia, GASTEC-3L (0.5~60 ppm)	00EAA0001555
Drawer set	Drawer shelf (slide) set	AAAF1512

Gas Cylinder Cabinet



Safe storage of high-pressure cylinders

- > **Stable storage environment**
Solid steel cabinet prevents external impact and prevents cylinders from falling over.
- > **Consideration of gas piping connection**
The structure features a hole in the upper space for the connection of a gas regulator and piping.
- > **Door Lock**
Key lock is provided as standard to block unauthorized access.
- > **Convenient internal observation**
Inspection windows provided on the front to allow for checking of the gas valves and remaining amount.
- > **Plastic chain provided**
Plastic chain for fixing cylinders included as standard.
Plastic material prevents scratches to cylinders/product.
- > **Convenient cylinder replacement**
Includes folding ramp, which is convenient for replacing heavy cylinders.



DS-GC-18107N2

with Plastic Chain, Folding Ramp (standard)



DS-GC-18074N2

with Plastic Chain, Folding Ramp (standard)

Safety Storage & Cabinet

Specification

Model	DS-GC (Height: 1800 mm / 70.9 inch)		
	Depth : 550 mm / 21.7inch	Cat. No.	AAAF1070
Model		DS-GC-18074N2	DS-GC-18107N2
Width (mm / inch)		740 / 29.1	1070 / 42.1
Gas cylinder (Ø230 mm / 9.1 inch) storage quantity		2	3



Plastic chain for fixing cylinders



Folding ramp for convenient cylinder replacement

Multipurpose Cabinet

Window type

Door with glass window on the front for easy checking of stored materials



DS-OW-18108N2

with Drawer Set, Partition, Name Plate (standard)



Fume Extraction Module (option)

Filtering and remove the harmful gas inside without a separate duct connection. (Refer to page 225 for the details.)

Convenient and Durable Structure

- **Front glass window door**
Front glass window door for easy checking of stored materials inside.
Easy information management for stored materials due to provision of name plates for each drawer.
- **Solid steel structure**
Stable structure with steel frame prevents structural deformation, and it has excellent external chemical resistance due to epoxy coating.
- **Door sealing & fixing**
Safe and secure sealing and fixing of glass window door with door packing and by applying magnet fixtures on the upper and lower parts.
- **Exhaust system connection**
Fume extraction module (option) for connection to the exhaust system of the installation site.

Improved User Convenience

- **Solid shelf structure**
Each shelf can accommodate up to 40kg, and the polypropylene material provides excellent chemical resistance.
- **Convenient drawer-type shelf**
Drawer-type shelf structure with sliding rails for convenient handling of stored materials located inside.
- **Shelf interior partitions provided**
There are adjustable partitions inside each shelf to make it more convenient to organize stored materials. (4 partitions provided for each shelf as standard)



Key lock device

Key lock is provided as standard to block unauthorized access.



Floor anchor

Anchor (option) for fixing the cabinet to the floor to help prevent accidents.



Drawer-type shelf

Convenient storage possible even in the back of the shelf.


Specification

Model	DS-OW-18055N1	DS-OW-18108N2
Exterior width (mm / inch)	556 / 21.9	1083 / 42.6
Exterior depth (mm / inch)	550 / 21.7	550 / 21.7
Exterior height (mm / inch)	1800 / 70.9	1800 / 70.9
Quantity of drawer	6	12
Max. load per shelf (kg / lbs)	30 / 66.1	30 / 66.1
Quantity of name plate	6	12
Cat. No.	AAAF1051	AAAF1053

Accessories

Partition (4ea per drawer provided as standard.)


- Adjustable partition fitting to the size of stored substances.
- Polypropylene (PP) material with excellent chemical resistance.



Description	Cat. No.	Dimension (mm / inch)	
		W	H
Partition	00CBN0000019	450 / 17.7	110 / 4.3

Name Plate (1ea per tray provided as standard.)

- Easy to check the contents stored on each shelf of the cabinet.
- Fix it at holes in front of the tray.



Description	Cat. No.	Dimension (mm / inch)	
		W	H
Name plate	AAAF1523	320 / 12.6	115 / 4.5

Extraction Module

Cat. No.	Description	Model
AAAF1536	Extraction Module	EM-33

Anchor

Cat. No.	Description
00CBN0000036	Anchor

Dry Cabinet & Desiccator



General Application

- Storage of various moisture-sensitive reagents, samples and seeds.
- Storage of various electronic devices, cameras, lenses, photographs, and film literature.

	Description	Note	Model	Page
Dry cabinet	Auto desiccator <ul style="list-style-type: none"> • Storage of moisture-sensitive materials and samples • Convenient use with automatic electronic desiccator 	Auto-dry	DC2	248
Vacuum desiccator	Cubic <ul style="list-style-type: none"> • Convenient checking of vacuum degree with the vacuum gauge on the front door • Height-adjustable shelves and space-efficient stacking structure 	General type	VDC	250
		UV blocking type	VDC-U	250
	Round <ul style="list-style-type: none"> • Transparent and durable polycarbonate material • Vacuum up to 1 Torr for up to 72 hours 	General type	VDR	250
		UV blocking type	VDR-U	250
	Cylinder <ul style="list-style-type: none"> • Convenient structure for large-volume sample desiccator • 3-way valve included as standard for easy gas replacement 	General type	VDP	250
		UV blocking type	VDP-U	250
Desiccator	Cubic <ul style="list-style-type: none"> • Easy to check with hygrometer attached to the front door • Drawer-type shelves included as standard for convenient storage of small samples 	General type	SDC	252
		UV blocking type	SDC-U	252

※ The contents of the above and the contents of this catalog may differ depending on the specific model and conditions of use. For the information about the features and specifications that applying to each model, please check the information on the corresponding page of each model.



Compliance with Guidelines for Restriction of Hazardous Substances Use.



User safety by blocking potentially harmful UV light.



Leakage minimized using an integrated mold.



Door locking system with key.



2 year warranty free A/S.

Dry Cabinet

Auto desiccator

Storage of moisture-sensitive materials and samples



DC2-21L
with Shelves 4 ea (standard)



DC2-41
with Shelves 8 ea (standard)

> Automatic electronic desiccator unit

Convenient maintenance due to there being no need for a desiccant by automatically controlling the set humidity inside the storage cabinet.

> Convenient humidity setting and control

20~40%RH (Actual set-up humid.), 10~60%RH (Setting range, 1% scale)

Deviation notification and calibration for humidity & temperature allow more reliable control.

> Real-time temperature and humidity display

Displays the temperature and humidity inside the library in real time.



Current temperature value

Current humidity value

Humidity setting



Convenient internal observation

Convenient internal observation with large transparent window made of acrylic material.



Height-adjustable shelves & optional drawers

Shelves can be adjusted in height, making it convenient loading of various samples. If you install drawers (optional), it is easy to handle storage inside.



Easy to move and install

Casters which make it easy to move/fix are provided as standard for convenient moving and installation of the product. (DC2-41(A) model)







Door lock

Key lock is provided as standard to block unauthorized access.



Specification

Model		DC2-11(A) ³⁾	DC2-21L(A) ³⁾	DC2-41(A) ³⁾												
Control range ^{1), 2)}		20~60%RH (1% scale)	20~60%RH (1% scale)	20~60%RH (1% scale)												
Internal volume (L / cu ft)		120 / 4.24	240 / 8.48	480 / 16.95												
Internal dimensions (W x D x H, mm / inch)		470 x 530 x 485 / 18.5 x 20.9 x 19.1	470 x 530 x 1120 / 18.5 x 20.9 x 44.1	1015 x 530 x 1120 / 40.0 x 20.9 x 38.2												
Exterior dimensions (W x D x H, mm / inch)		540 x 680 x 650 / 21.6 x 26.8 x 25.5	540 x 680 x 1195 / 21.6 x 26.8 x 47	1085 x 680 x 1265 / 42.7 x 26.8 x 49.8												
Quantity of door		1	2	4												
Stackable type		Yes	No	No												
Quantity of shelves (standard/max.)		2 / 4	4 / 8	8 / 16												
Max. load per shelf (kg / lbs)		20 / 44	20 / 44	20 / 44												
Net Weight (kg / lbs)		42 / 93	65 / 143	118 / 260												
Electrical requirements (230V, 50/60Hz, 1P,A)		2	2	4												
Cat. No.	Basic	AAHF70125K	AAHF70225K	AAHF70425K												
	Anti-electricity	AAHF70135K	AAHF70235K	AAHF70435K												
Product appearance																
Product interior																
		<table border="1"> <thead> <tr> <th>Model</th> <th>A (mm / inch)</th> <th>B (mm / inch)</th> </tr> </thead> <tbody> <tr> <td>DC2-11(A)</td> <td>1094 / 43.1</td> <td>1202 / 47.3</td> </tr> <tr> <td>DC2-21L(A)</td> <td>1094 / 43.1</td> <td>1202 / 47.3</td> </tr> <tr> <td>DC2-41(A)</td> <td>2200 / 86.6</td> <td>1202 / 47.3</td> </tr> </tbody> </table> 			Model	A (mm / inch)	B (mm / inch)	DC2-11(A)	1094 / 43.1	1202 / 47.3	DC2-21L(A)	1094 / 43.1	1202 / 47.3	DC2-41(A)	2200 / 86.6	1202 / 47.3
Model	A (mm / inch)	B (mm / inch)														
DC2-11(A)	1094 / 43.1	1202 / 47.3														
DC2-21L(A)	1094 / 43.1	1202 / 47.3														
DC2-41(A)	2200 / 86.6	1202 / 47.3														

- 1) Internal humidity may temporarily increase when the dehumidifying element of the unit is dried.
- 2) Performance may vary depending on your installation environment.
- 3) DC2-A model is an anti-static product that prevents static electricity.

Accessories

Model		Dimension (W x D x H, mm)	Cat. No.
Shelves	Epoxy powder coated steel	460 x 525 x 60 / 18.1 x 20.7 x 2.4	AAAF7521
	Anti-static paint steel	460 x 525 x 60 / 18.1 x 20.7 x 2.4	AAAF7522
Drawer	Epoxy powder coated steel	420 x 500 x 100 / 16.5 x 19.7 x 3.9	AAAF7531
	Anti-static paint steel	420 x 500 x 100 / 16.5 x 19.7 x 3.9	AAAF7541



Shelf



Drawer

Vacuum Desiccator

Optimum choice for dehumidification in vacuum state

- > **1 Torr, maintained for 72 hours**
Up to 1 Torr vacuum maintained for 72 hours. Suitable for samples that need to be dehumidified in vacuum state.
- > **Solid integral structure**
Made of integral mold to minimize leakage and there is no warping even in repeated high vacuum conditions.
- > **Transparent observation environment**
The material of the body is polycarbonate.(PC)
- > **UV blocking model**
UV blocking model provided for vacuum storage of light-sensitive samples.
- > **Desiccant Tray**
Includes a dedicated tray made of PP as standard that the user can conveniently place a desiccant in.



VDC-31 (Cubic type)
with Vacuum Gauge (standard)



VDR-25G (Round type)
with Vacuum Gauge (option)



VDP-30G (Cylinder type)
with 2ea Perforated Shelf (standard),
Vacuum Gauge (standard)



Useful 3-way valve
Convenient and useful valve structure for vacuum formation, remove and gas substitution replacement.



Multi-level stacking (VDC series)
Stable multi-level loading to maximize space efficiency.



Height-adjustable shelves (VDC series)
Height-adjustable shelf structure for convenient loading of various samples.



Locker (VDR, VDP series)
Conveniently and safely fixed at normal pressure when used or transported.



RoHS eco-friendly products
Eco-friendly products that have been deemed acceptable in accordance with Restriction of the Use of Certain Hazardous Substances. (RoHS) (except for U model)



UV blocking model (U model)
Vacuum storage of light-sensitive samples.



Specification

Cubic Type

Model	Clear type	VDC-11	VDC-21	VDC-31	VDC-41
	UV protection	VDC-11U	VDC-21U	VDC-31U	VDC-41U
Internal volume (L / cu ft)		11 / 0.4	23 / 0.8	35 / 1.2	45 / 1.6
Vacuum	Gauge range (Analog) (Mpa)	0 to -0.1	0 to -0.1	0 to -0.1	0 to -0.1
	Max. permissible vacuum (Body) (Mpa)	1.33×10^{-4} (1 Torr)	1.33×10^{-4} (1 Torr)	1.33×10^{-4} (1 Torr)	1.33×10^{-4} (1 Torr)
Dimension	Internal (W x D x H, mm / inch)	248 x 254 x 238 / 9.8 x 10 x 9.4	346 x 365 x 246 / 13.6 x 14.4 x 9.7	355 x 375 x 345 / 14 x 14.8 x 13.6	355 x 374 x 445 / 14 x 14.7 x 17.5
	Overall (W x D x H, mm / inch)	322 x 285 x 271 / 12.7 x 11.2 x 10.7	420 x 397 x 281 / 16.5 x 15.6 x 11.1	420 x 397 x 381 / 16.5 x 15.6 x 15	420 x 397 x 491 / 16.5 x 15.6 x 19.3
	Perforated shelf (W x D x H, mm / inch)	202 x 204 x 6.5 / 7.95 x 8.0 x 0.3	289 x 304 x 11 / 11.4 x 12.0 x 0.4	289 x 304 x 11 / 11.4 x 12.0 x 0.4	289 x 304 x 11 / 11.4 x 12.0 x 0.4
	Desiccant tray (W x D x H, mm / inch)	158 x 194 x 24 / 6.2 x 7.6 x 0.9	218 x 268 x 26 / 8.6 x 10.6 x 1.0	218 x 268 x 26 / 8.6 x 10.6 x 1.0	218 x 268 x 26 / 8.6 x 10.6 x 1.0
	Nozzle dia. (mm / inch)	9.5 / 0.37	9.5 / 0.37	9.5 / 0.37	9.5 / 0.37
	Net weight (kg / lbs)	4.2 / 9.3	8.2 / 18.1	10.8 / 23.8	12.2 / 26.9
Shelves	Quantity of shelves (standard / max.)	2 / 4	2 / 4	3 / 5	3 / 6
	Max. load per shelf (kg / lbs)	3 / 6.6	5 / 11.0	5 / 11.0	5 / 11.0
Cat. No.	Clear type	AAAD4001	AAAD4021	AAAD4011	AAAD4031
	UV protection	AAAD4101	AAAD4121	AAAD4111	AAAD4131

Accessories [Page 253](#) Desiccant Tray, Perforated Shelf, Vacuum Seal, 3-way Valve, Vacuum Gauge

Round Type

Model	Clear type	VDR-20 (20G)	VDR-25 (25G)	VDR-30 (30G)
	UV protection	VDR-20U (20UG)	VDR-25U (25UG)	VDR-30U (30UG)
Internal volume (L / cu ft)		6 / 0.2	10 / 0.4	20 / 0.7
Vacuum	Gauge range (Analog) (Mpa)	0 to -0.1	0 to -0.1	0 to -0.1
	Max. permissible vacuum (Body) (Mpa)	1.33×10^{-4} (1 Torr)	1.33×10^{-4} (1 Torr)	1.33×10^{-4} (1 Torr)
Dimension	Overall (\emptyset x H, mm / inch) ^{1), 2)}	252 x 279 (252 x 354) / 9.5 x 11 (9.5 x 13.9)	316 x 325 (316 x 400) / 12.1 x 12.8 (12.1 x 15.7)	403 x 399 (403 x 475) / 15.2 x 15.7 (15.2 x 18.7)
	Internal (\emptyset x H, mm / inch)	183 x 173 / 7.2 x 6.8	242 x 206 / 9.5 x 8.1	300 x 266 / 11.8 x 10.5
	Desiccant tray (W x D, mm / inch)	159 x 31 / 6.3 x 1.2	220 x 29 / 8.7 x 1.1	276 x 34 / 10.9 x 1.3
	Nozzle dia. (mm / inch)	9.5 / 0.37	9.5 / 0.37	9.5 / 0.37
	Net weight (kg / lbs) ²⁾	1.5 (1.7) / 3.3 (3.7)	2.5 (2.6) / 5.5 (5.7)	4.2 (4.3) / 9.3 (9.5)
Shelves	Perforated shelf (\emptyset , mm / inch)	192 / 7.6	250 / 9.8	310 / 12.2
	Max. load per shelf (kg / lbs)	15 / 33.1	20 / 44.1	25 / 55.1
Cat. No.	Clear type, without vacuum gauge	AAAD2011	AAAD2111	AAAD2211
	Clear type, with vacuum gauge (G type)	AAAD2021	AAAD2121	AAAD2221
	UV protection, without vacuum gauge (U type)	AAAD2031	AAAD2131	AAAD2231
	UV protection, with vacuum gauge (UG type)	AAAD2041	AAAD2141	AAAD2241

1) The \emptyset values are without the 3-way valve. (When including the valve, the max. width is as follows: VDR-20: 286mm/11.3inch, VDR-25: 353mm/13.9inch, VDR-30: 427mm/16.8inch)

2) The figures within parentheses represent the value for the Vacuum Gauge-equipped models.

Accessories [Page 253](#) Desiccant Tray, Perforated Shelf, Vacuum Seal, 3-way Valve, Vacuum Gauge, Locker

Cylinder Type

Model	Clear type	VDP-25G	VDP-30G
	UV protection	VDP-25UG	VDP-30UG
Internal volume (L / cu ft)		16 / 0.6	24.9 / 0.9
Vacuum	Gauge range (Analog) (Mpa)	0 to -0.1	0 to -0.1
	Max. permissible vacuum (Body) (Mpa)	1.33×10^{-4} (1 Torr)	1.33×10^{-4} (1 Torr)
Dimension	Internal (\emptyset x H, mm / inch)	244 x 318 / 9.6 x 12.5	305 x 322 / 12.0 x 12.7
	Overall (\emptyset x H, mm / inch)	308 x 449 / 12.1 x 17.7	387 x 459 / 15.2 x 18.1
	Nozzle dia. (mm / inch)	9.5 / 0.37	9.5 / 0.37
	Net weight (kg / lbs)	2.8 / 6.2	4.1 / 9.0
Quantity of perforated shelves		2 / 5 (or 3) ¹⁾	2 / 5 (or 3) ¹⁾
Cat. No.	Clear type, with vacuum gauge (G type)	AAAD3111	AAAD3121
	UV protection, with vacuum gauge (UG type)	AAAD3311	AAAD3321

1) Perforated shelves can be placed in a maximum of 5 (spacing 40mm / 1.6 inch) or 3 (spacing 70mm / 12.8 inch) depending on the stacking position.

Accessories [Page 253](#) Perforated Shelf, Vacuum Seal, 3-way Valve, Vacuum Gauge, Locker

※ From the model name, the last letter U means 'UV blocking', and the letter G means 'Vacuum Gauge'.

Hygrometer and multi-level shelf structure for convenient dehumidified desiccator

> Hygrometer included

Very convenient as equipped with a hygrometer to check the internal humidity at the front.

> Height-adjustable drawer

Height of drawer can be adjusted according to the shape and size of stored material.

> Transparent observation environment

Highly transparent, shatterproof, non-corrosive, and easy to clean polycarbonate construction.

> UV blocking model

UV blocking model provided for vacuum storage of light-sensitive samples.

> Desiccant Tray

Includes a dedicated tray made of PP as standard that the user can conveniently place a desiccant in.



SDC-45



Easy humidity check with front mounted hygrometer.



Drawer provided as standard for the convenient storage of the back inside.



Stable multi-level loading allows for efficient utilization of space.



Minimizing damages or discoloration of light-sensitive samples. (-U model)

Specification

Model	Clear type	SDC-30 ¹⁾	SDC-45 ¹⁾
	UV protection type	SDC-30U ¹⁾	SDC-45U ¹⁾
Internal volume (L / cu ft)		35 / 1.2	50 / 1.8
Dimension	Internal (W x D x H, mm / inch)	322 x 367 x 302 / 12.7 x 14.4 x 11.9	322 x 367 x 422 / 12.7 x 14.4 x 16.6
	Overall (W x D x H, mm / inch)	360 x 410 x 326 / 14.1 x 16.1 x 12.8	360 x 410 x 447 / 14.1 x 16.1 x 17.6
	Drawer (W x D x H, mm / inch)	279 x 346 x 45 / 11.0 x 13.6 x 1.8	279 x 346 x 45 / 11.0 x 13.6 x 1.8
	Desiccant tray (W x D x H, mm / inch)	218 x 268 x 26 / 8.6 x 10.6 x 1.0	218 x 268 x 26 / 8.6 x 10.6 x 1.0
	Net weight (kg / lbs)	4.9 / 10.8	6.2 / 13.7
Drawers	Quantity (Standard / Max.)	2 / 4	3 / 6
	Max. load per drawer (kg / lbs)	5 / 11.0	5 / 11.0
Cat. No.	Clear type	AAAD5011	AAAD5021
	UV protection type	AAAD5111	AAAD5121

1) The deviation of the fitted hygrometers will be between 10 and 20%. The measuring of accurate humidity percentage should be needed a extra equipment.

Accessories for Desiccator

Product Configuration



Desiccant Tray



Perforated Shelf (VDC)



Perforated Shelf (VDP)



Vacuum Seal



Locker



3-way Valve (VDR)



3-way Valve (VDC)



Vacuum Gauge


Product	Cat. No.	VDC				VDR			VDP		SDC	
		VDC-11	VDC-21	VDC-31	VDC-41	VDR-20	VDR-25	VDR-30	VDP-25	VDP-30	SDC-30	SDC-45
Desiccant tray	00VDC0000037	●	-	-	-	-	-	-	-	-	-	-
	00VDC0000005	-	●	●	●	-	-	-	-	-	●	●
	AAAD1501	-	-	-	-	●	-	-	●	-	-	-
	AAAD1502	-	-	-	-	-	●	-	-	●	-	-
	AAAD1503	-	-	-	-	-	-	●	-	-	-	-
Perforated shelf	00VDC0000035	●	-	-	-	-	-	-	-	-	-	-
	00VDC0000003	-	●	●	●	-	-	-	-	-	-	-
	AAAD1511	-	-	-	-	●	-	-	-	-	-	-
	AAAD1512	-	-	-	-	-	●	-	-	-	-	-
	AAAD1513	-	-	-	-	-	-	●	-	-	-	-
	AAAD3511	-	-	-	-	-	-	-	●	-	-	-
	AAAD3512	-	-	-	-	-	-	-	-	●	-	-
Drawer	AAAD5521	-	-	-	-	-	-	-	-	-	●	●
Vacuum seal	00VDC0000032	●	-	-	-	-	-	-	-	-	-	-
	00VDC0000033	-	●	-	-	-	-	-	-	-	-	-
	00VDC0000015	-	-	●	-	-	-	-	-	-	-	-
	00VDC0000034	-	-	-	●	-	-	-	-	-	-	-
	00VDR0000001	-	-	-	-	●	-	-	-	-	-	-
	00VDR0000002	-	-	-	-	-	●	-	●	-	-	-
	00VDR0000055	-	-	-	-	-	-	●	-	●	-	-
Locker	00VDR0000015	-	-	-	-	●	-	-	-	-	-	-
	00VDR0000016	-	-	-	-	-	●	-	●	-	-	-
	00VDR0000017	-	-	-	-	-	-	●	-	●	-	-
3-way valve	AAAD4501	●	●	●	●	-	-	-	-	-	-	-
	AAAD1521	-	-	-	-	●	●	●	-	-	-	-
	AAAD3501	-	-	-	-	-	-	-	●	●	-	-
Vacuum Gauge	00EDA0004301	●	●	●	●	-	-	-	-	-	-	-
	00EDA0004302	-	-	-	-	●	●	●	●	●	-	-

※ This table provides product information for additional installation or replacement of components included in the product

Accessory

Cat. No.	Product Name	Description
BEA00078603	Silicagel (1kg)	Desiccant to remove moisture





Refrigerator & Freezer for Laboratory/Medical



Refrigerator & Freezer for Laboratory/ Medical



Refrigerator & Freezer for Laboratory/ Medical



General Application

Lab Refrigerators & Freezers: Chemical samples, various reagents, storage of culture media, temperature-sensitive materials, enzymes, preliminary cooling before freezing.

Medical refrigerators: Storage of medicines and vaccines in hospitals and pharmacies.

	Description	Temp. Range (°C/°F)	Volume (L / Cu ft)	Model	Page
Laboratory	Refrigerator <ul style="list-style-type: none"> • Verified by international standards • Temperature stabilization by defrost control 	2 to 10 / 35.6 to 50	73, 156, 303, 660, 1451 / 2.6, 5.5, 10.7, 23.3, 51.2	CLG CLG3	260
		Freezer <ul style="list-style-type: none"> • Model selected based on temperature range • Excellent insulation and sealing structure 	-25 to -15 / -13 to 5	154, 303, 647 / 5.4, 10.7, 22.8	FMG
	-35 to -25 / -31 to -13		154, 303, 647 / 5.4, 10.7, 22.8	FCG	
Medical	Pharmaceutical Refrigerator <ul style="list-style-type: none"> • Stable cold storage of medicines • Designed based on KS P 6108 (pharmaceuticals refrigerator conditions and tests) 	2 to 10 / 35.6 to 50	73 / 2.6	PSR3	264

※ The contents of the above and the contents of this catalog may differ depending on the specific model and conditions of use. For the information about the features and specifications that applying to each models, please check the information on the corresponding page of each models.

	Satisfaction with the requirements of European regulations and guidelines.		More convenient security system with electronic lock function.		Real-time temperature monitoring by thermal line recorder.
	Real-time equipment monitoring and control system using mobile app.		Registered patent based on Jeio Tech's proprietary technology.		Comply with the standards for manufacturing and quality management of medical supplies.
	Certificate of medical devices management system for the manufacturers of medical devices.		Eco-friendly product that has passed environmental impact assessment tests.		Automatic defrosting does not need the work of unnecessary frost elimination.
	In case of abnormal operation of the device, a voice message alarm is transmitted via telephone.		Warning alarm in case of deviation of temperature control.		It memorizes the state of power failure and operates automatically when power is restored.
	2 year warranty free A/S.		Door locking system with key.		



Refrigerator & Freezer for Laboratory/Medical

Improved reliability with strict international standard verification
Stable temperature control with excellent sealing and defrost control

> **Stable temperature control and maintenance**

More precise and stable temperature control due to Jeio Tech's exclusive controller.

> **Internal uniform temperature**

Optimized structural design for each model maintains uniform temperature inside.

> **Complete chamber structure**

Excellent thermal insulation and sealing to minimize heat loss for long-term use.

> **Clean use and maintenance**

Internal stainless steel material allows for clean maintenance and excellent corrosion resistance.

> **Satisfies the requirements for safe cold storage of pharmaceuticals**

Temperature fluctuation and variation meets the requirements of ICH Guideline Q1A(R2). (CLG model)

> **Licensed Pharmaceutical refrigerator as a Medical device (PSR3)**

Meet the standard for medicine storage. (KS P 6108) Uniform temperature even with stricter standards.

> **Application of defrost control technology of Jeio Tech**

Automatic defrosting or dual refrigeration without the need for defrosting minimizes temperature changes.

> **Dew-free exterior**

Provides a clean environment with a structure that does not form dew on the door and outer wall.



Solid steel drawer shelf. (option)



Drawer type of wire shelf with excellent permeability. (option)



Key lock door.



Provides optimized selection for storage purpose Convenient safety features for safe storage

- > **Mobile monitoring system (CLG3, PSR3)**
LC Connected system (option) to monitor the current status of the equipment through the mobile app and receipt of status messages.
- > **Door lock with key**
Door structure with key lock as standard provides a safe storage environment.
- > **Convenient mobility and installation**
Easy to move and install with casters. (except for 150 model)
- > **Convenient accessory selection supported**
Selection of various accessories such as temperature recorder as well as various shelves according to storage type.
- > **Temperature deviation notification function**
Warning when the actual temperature is out of a certain range of the set temperature. (High & Low Temperature)
- > **Lock function of controller**
Lock function for preventing change of controller settings prevents erroneous operation due to arbitrary changing.
- > **Self-protection system**
Door open alarm overcurrent protection and current leakage protection.
- > **Voice alarm system supported**
Alarm system (option) that transmits voice messages to designated users when device malfunction occurs. (Included as standard only for PSR3.)

※ Some of the above contents are limited to specific models.



Intuitive checking of device operation information with large VFD.



Storing samples in a stable through a real-time temperature recording. (Included as standard only for PSR3, for the other models are optional.)



Mobile remote monitoring system. (option) (CLG3/PSR3 model)

Laboratory Refrigerator

Excellent structure to meet stringent international standard



Structural Functional Features

- Stable temperature control with excellent insulation and sealing structure.
- Clean construction with no condensation on outer walls.
- Interior has excellent corrosion resistance as made of stainless steel.
- Magnetic packing door with smooth opening and closing.
- Casters for easy movement / installation. (except for 70, 150 model)
- Space-efficient configuration with double stackable. (only for 70, 150 model)
- Auto defrost function eliminates the need for separate defrosting.

Use Convenience Features

- Stable control with Jeio Tech's exclusive controller.
- Easy to use with clear VFD and touch key.
- Simultaneously displays set temperature and actual temperature.
- Calibration function minimizes temperature deviation.
- Height-adjustable and strong shelf provided as standard.
- Built-in fluorescent lamps. (except for 150 model)
- Glass door model for easy storage status checking.
- Convenient connecting external devices through Ø50 mm cable port. (option) (except for 70,150 model)
- No need for drainage management with a drain evaporation system. (only CLG3)



Intuitive checking of device operation information with large VFD




Door key lock



Outstanding Safety

- The 1400 model is equipped with two refrigeration systems to maintain temperature even in the event of a single refrigeration system failure.
- Temperature deviation notification. (High & Low Temperature)
- Malfunction prevented by controller lock function. (except for CLG3)
- Door lock with key as standard.
- Convenient history check through recording the recent event such as door open or power on/off. (only CLG3)
- Notification function when the door is open for a long time.
- Enhanced security with digital door lock. (only CLG3)



(left) Drawer shelf, (right) Digital lock drawer (option)

Satisfies the Requirements for Safe Cold Storage of Pharmaceuticals

- Temperature fluctuation and variation meets the requirements of ICH Guideline Q1A(R2).
- Stable temperature control for long-term operation.



Wire Shelf (Slide) (option)

Accessories Provided for Convenient Use

- Monitor via mobile app anytime, anywhere with LC Connected. (CLG3-70 only) (Requires LC GreenBox, option)
- Voice messages are transmitted to designated phone in case of an abnormality in operation. (option)
- Convenient use with drawer shelf or wire(slide) shelf. (option)
- Selection of various temperature recorders. (option) (Dot, Chart, Thermal Line)



Digital drawer lock (CLG3 model, option)

Specification

Model	Solid door	-	CLG-150S	CLG-300S	CLG-650S	CLG-1400S
	Glass door	CLG3-70	CLG-150G	CLG-300G	CLG-650G	CLG-1400G
Temperature range (°C / °F)		2 to 10 / 35.6 to 50	2 to 10 / 35.6 to 50	2 to 10 / 35.6 to 50	2 to 10 / 35.6 to 50	2 to 10 / 35.6 to 50
Chamber volume (L / cu ft)		73 / 2.6	156 / 5.5	303 / 10.7	660 / 23.3	1451 / 51.2
Dimension						
Interior (W x D x H, mm / inch)		418 x 440 x 400 / 16.5 x 17.3 x 15.7	550 x 580 x 490 / 21.7 x 22.8 x 19.3	550 x 550 x 1000 / 21.7 x 21.7 x 39.4	650 x 640 x 1450 / 25.6 x 25.2 x 57.1	1430 x 640 x 1450 / 56.3 x 25.2 x 57.1
Exterior (W x D x H, mm / inch)		522 x 600 x 725 / 20.9 x 23.6 x 28.5	683 x 791 x 870 / 26.9 x 31.1 x 34.3	670 x 770 x 1515 / 26.4 x 30.3 x 59.6	770 x 865 x 1965 / 30.3 x 34.1 x 77.4	1550 x 865 x 1965 / 61 x 34.1 x 77.4
Net weight (Kg / lbs)		67 / 147.7	100 / 220.5	130 / 286.6	210 / 463	350 / 771.6
Shelf						
Quantity of wire shelf (standard / max.)		2 / 3	2 / 6	3 / 11	3 / 17	6 / 34
Quantity of wire shelf (slide) (max.)		3	4	7	11	22
Quantity of drawer shelf (slide) (max.)		2	3	5	8	16
Max. Load per shelf (kg / lbs)		30 / 66.1	26 / 57.3	26 / 57.3	30 / 66.1	30 / 66.1
Fluorescent lamp		4.2W, 1ea	-	18W, 1ea	18W, 1ea	32W, 1ea
Electrical requirements (230V, 50Hz, A)		1.1	2	2.3	4	6
Cat. No.	Solid door	-	AAHE1002K	AAHE1012K	AAHE1022K	AAHE1042K
	Glass door	AAHE1302K	AAHE1102K	AAHE1112K	AAHE1122K	AAHE1142K
Electrical requirements (230V, 60Hz, A)		1.1	2	2.3	4	6
Cat. No.	Solid door	-	AAHE1001K	AAHE1011K	AAHE1021K	AAHE1041K
	Glass door	AAHE1301K	AAHE1101K	AAHE1111K	AAHE1121K	AAHE1141K
Electrical requirements (120V, 60Hz, A)		2.9	4	4.5	5	9.5
Cat. No.	Solid door	-	AAHE1003U	AAHE1013U	AAHE1023U	AAHE1043U
	Glass door	AAHE1303U	AAHE1103U	AAHE1113U	AAHE1123U	AAHE1143U

Accessories Page 266 Shelf, Drawer, Digital Drawer Lock, Recorder, Remote Alarm, Battery(Alarm), Stand with Caster, Electric Socket, LC GreenBox

Laboratory Freezer

Maintain stable temperature by direct cooling-type refrigeration system

FCG-150



FCG-150



FMG-650

with Recorder (option)

Structural Functional Features

- Stable temperature control with excellent insulation and sealing structure.
- Clean construction with no condensation on outer walls.
- Interior has excellent corrosion resistance as made of stainless steel.
- Magnetic packing door with smooth opening and closing.
- Casters for easy movement / installation. (except for 150 model)
- Space-efficient configuration with double stackable. (only for 150 model)
- Auto defrost function eliminates the need for separate defrosting.

Use Convenience Features

- Stable control with Jeio Tech's exclusive controller.
- Easy to use with clear VFD and touch key.
- Simultaneously displays set temperature and actual temperature.
- Calibration function minimizes temperature deviation.
- Height-adjustable and strong shelf provided as standard.


Outstanding Safety

- Temperature deviation notification. (High & Low Temperature)
- Malfunction prevented by controller lock function.
- Door lock with key as standard.
- Notification function when door is opened for a long time.



Direct Cooling-type Refrigeration System

- Cooling method with structure surrounding the internal sample space.
- The structure maintains stable temperature uniformity throughout the entire interior where the sample is placed.
- Efficient operation provided with minimal heat loss.



Accessories Provided for Convenient Use

- Voice alarm system transmits voice messages to designated phone if there is an abnormality in operation. (option)
- Convenient use with drawer type shelf. (option)
- Selection of various temperature recorders. (option)
(Dot, Chart, Thermal Line)

Specification

Model		FCG-150	FCG-300	FCG-650
		FMG-150	FMG-300	FMG-650
Temperature range (°C / °F)	FCG	-35 to -25 / -31 to -13	-35 to -25 / -31 to -13	-35 to -25 / -31 to -13
	FMG	-25 to -15 / -13 to 5	-25 to -15 / -13 to 5	-25 to -15 / -13 to 5
Chamber volume (L / cu ft)		154 / 5.4	303 / 10.7	647 / 22.8
Dimension				
Interior (W x D x H, mm / inch)		550 x 570 x 490 / 21.7 x 22.4 x 19.3	550 x 550 x 1000 / 21.7 x 21.7 x 39.4	700 x 675 x 1370 / 27.6 x 26.6 x 53.9
Exterior (W x D x H, mm / inch)		763 x 815 x 1055 / 30 x 32.1 x 41.5	765 x 760 x 1650 / 30.1 x 29.9 x 65	915 x 885 x 2020 / 36 x 34.8 x 79.5
Net weight (Kg / lbs)	FCG	123 / 271.2	155 / 341.7	238 / 524.7
	FMG	120 / 264.6	150 / 330.7	233 / 513.7
Shelf				
Quantity of shelf (standard / max.)		2 / 6	3 / 11	3 / 17
Max. Load per shelf (kg / lbs)		26 / 57.3	26 / 57.3	30 / 66.1
Electrical data & Ordering information				
230V, 50Hz, A	FCG	5.0	5.5	6.5
Cat. No.		AAHE2302K	AAHE2312K	AAHE2322K
Cat. No. (EU only)		AAHE2305K	AAHE2315K	AAHE2325K
230V, 60Hz, A	FMG	5.0	5.5	6.5
Cat. No.		AAHE2301K	AAHE2311K	AAHE2321K
230V, 50/60Hz, A		4.5	4.8	5.5
Cat. No.	AAHE2202K	AAHE2212K	AAHE2222K	
Cat. No. (EU only)	AAHE2205K	AAHE2215K	AAHE2225K	

Accessories [Page 266](#) Shelf, Drawer, Recorder, Remote Alarm, Stand with Caster

Pharmaceutical Refrigerator

Stable cold storage for pharmaceuticals provided

Comprehensive Temperature Management.

- Temperature recording in real-time by thermal line recorder.
- Verified the temperature uniformity under KS P 6108 and DIN 13221 standard.
- Temperature deviation alarm function. (Temp. upper & lower)



Includes Thermal line recorder as standard to monitor a real-time temperature.

Licensed as a Medical Device.

- Licensed manufacturer for a medical device. (No. 3333 registered by KFDA)
- Certified by GMP.



More Convenient Security System with Electronic Lock Function.

- The electronic lock system is provided as standard.
- An electronic lock drawer 1 set is provided as standard.
- Easy unlock on a touch display only for the authorized user.
- Convenient history check through recording the recent event such as door open or power on/off.



Digital door lock.

Safe Storage with the Temperature Deviation Alarm

- Built-in voice alarm system transmits voice messages to specified users (up to 7 phones) if there is an abnormality in the device.
- Monitor via mobile app anytime, anywhere with LC Connected. (Requires LC GreenBox, option)



LC Connected



PSR3-70


wire shelf 2ea, drawer 1ea, thermal line recorder, remote alarm system with battery (standard)



Intuitive and convenient use with 3.5inch color touch display.




Convenient storage with a drawer.



Wire shelf with good air permeability.

Dimensions



Description		Dimension (mm / inch)
A	External width	522 / 20.6
A'	Internal width	390 / 15.4
A''	Drawer width	349 / 13.7
B	Internal height	420 / 16.5
C	Between the bottom shelf and the top shelf	136 / 5.4
C'	Between the top shelf and the drawer	110 / 4.3
C''	Between the drawer and inner ceiling	117 / 4.6
D	External depth	600 / 23.6
E	Max. depth when door is opened	1075 / 42.3
F	Max. width when door is opened	909 / 35.8


Specification

Model		PSR3-70	
Internal volume (L / cu ft)		73 / 2.6	
Temperature	Range (°C / °F) (5°C set point)	2 to 10 / 35.6 to 50	
	Variation at 5°C (±°C)	1.0	
	Fluctuation at 5°C (±°C)	2.0	
Recorder		0	
Digital door lock		0	
Digital drawer lock		0	
Quantity of wire shelf		2	
Quantity of SS drawer		1	
Remote alarm system		0	
External dimension (W x D x H, mm / inch)		530 x 600 x 725 / 20.5 x 23.6 x 28.5	
Electrical data (230V, 1P, A)		50Hz, 2.0	60Hz, 2.0
Cat. No.		AAHE6302K	AAHE6301K
Electrical data (120V, 60Hz, 1P, A)		4.0	
Cat. No.		AAHE6303U	

Accessories Page 266 Shelf, Drawer, Recorder, Remote Alarm, Battery(Alarm), LC GreenBox

Accessories

Wire Shelf



- Provides uniform temperature conditions throughout the storage area thanks to space between wires.

Wire Shelf (Slide)



- Wire type with good airflow, and convenient to use as drawer type.

Drawer Shelf (Slide)



- The drawer shelf of a rigid steel plate makes easy to handle storage inside.

Electric Socket



- Power supply for AC 230V.
- Includes cover as standard.

Recorder

Recorder with Sensor (Thermal Line, Single Sensor)



- Graph / Digital mode
- Thermal-type temperature recorder.
- Single sensor type, and built-in record for Medical products as standard. (PSR3)

Recorder with Sensor (Thermal Line, Double Sensor)



- Graph/Digital mode.
- Thermal-type temperature recorder.
- Double sensor type, so available to record a temperature continuously even one sensor is out of work.

Recorder with Sensor (Chart)



- 24/7 circular chart recorder.
- Easy replacement of the recording paper.
- Independent sensor that are separate from the product.

Recorder with Sensor (Dot)



- Up to 6 channels supported.
- Set recording range and speed.
- Paper width: 1000 mm / 39.4 inch.
- Independent sensor that are separate from the product.

Cat. No.	Description
AAAE15631	Thermal Line Recorder Single sensor for CLG/CLG3/FMG/FCG
AAAE1503	Thermal Line Recorder Double sensor for CLG/CLG3/FMG/FCG
AAAE1501	Chart type Recorder(24hr/7day) for CLG/FMG/FCG
AAAE1562	Dot type Recorder for CLG/FMG/FCG

Communication

LC DataKeeper



- High-integrity software compliant with FDA 21 CFR Part 11
- Integrated with Jeio Tech equipment to undergo GMP audits

Cat. No.	Description
SSA11001	LC DataKeeper-Basic
SSA11002	LC DataKeeper-Pro

- ※ Applied model: CLG3 / PSR3
- ※ Trial version is available for 30 days, offering the same features as the Basic version, and can be requested through the website.

LC GreenBox



- Real-time monitoring the state of a device via mobile app under an internet connection.
- Available to connect 4 devices in max for 1 unit of LC GreenBox.

Cat. No.	Description	Dimension (W x D x H, mm / inch)
AAHQ1011K	LC GreenBox	156 x 94 x 34 / 6.14 x 3.7 x 1.34

- ※ Applied model: CLG3 / PSR3

Remote Alarm with Battery



- Activate/deactivate battery via settings.
- When the setting of battery is ON, a warning alarm occurs if the power is turned off during operation.
- The recorder also records power failure events.
- Included as standard only for PSR3.

Cat. No.	Description
AAAE15511	Remote alarm system
AAAE15301	Battery

- ※ Available for all models.

Model	Standard		Option							
	Wire Shelf	Drawer Shelf (Slide)	Wire Shelf (Slide)	Drawer Shelf (Slide)	Recorder (Chart)	Recorder (Dot)	Recorder (Thermal line)	Electric Socket	Stand with Caster	Remote Alarm
CLG3-70	AAAE15081	N/A	AAAE15091	AAAE15101 ¹⁾	N/A	N/A	AAAE15631	N/A	N/A	AAAE15511
CLG-150	AAAE1511	N/A	AAAE1514	AAAE1517	N/A	AAAE1562	AAAE1563	N/A	AAAE1542	AAAE1551
CLG-300 ³⁾	AAAE1511	N/A	AAAE1514	AAAE1517	AAAE1501	AAAE1562	AAAE1563	AAAE1532 ²⁾	N/A	AAAE1551
CLG-650 ³⁾	AAAE1512	N/A	AAAE1515	AAAE1518	AAAE1501	AAAE1562	AAAE1563	AAAE1532 ²⁾	N/A	AAAE1551
CLG-850 ³⁾	AAAE1513	N/A	AAAE1516	AAAE1519	AAAE1501	AAAE1562	AAAE1563	AAAE1532 ²⁾	N/A	AAAE1551
CLG-1400 ³⁾	AAAE1512	N/A	AAAE1515	AAAE1518	AAAE1501	AAAE1562	AAAE1563	AAAE1532 ²⁾	N/A	AAAE1551
FCG-150	AAAE1511	N/A	AAAE1514	AAAE1517	AAAE1501	AAAE1562	AAAE1563	N/A	AAAE1542	AAAE1551
FCG-300	AAAE1511	N/A	AAAE1514	AAAE1517	AAAE1501	AAAE1562	AAAE1563	N/A	N/A	AAAE1551
FCG-650	AAAE1512-1	N/A	AAAE1515-1	AAAE1518-1	AAAE1501	AAAE1562	AAAE1563	N/A	N/A	AAAE1551
FMG-150	AAAE1511	N/A	AAAE1514	AAAE1517	AAAE1501	AAAE1562	AAAE1563	N/A	AAAE1542	AAAE1551
FMG-300	AAAE1511	N/A	AAAE1514	AAAE1517	AAAE1501	AAAE1562	AAAE1563	N/A	N/A	AAAE1551
FMG-650	AAAE1512-1	N/A	AAAE1515-1	AAAE1518-1	AAAE1501	AAAE1562	AAAE1563	N/A	N/A	AAAE1551
PSR3-70 ⁴⁾	AAAE15081	AAAE15101	AAAE15091	AAAE15101 ¹⁾	N/A	N/A	AAAE15631	N/A	N/A	AAAE15511

1) CLG3-70, PSR3-70 offer 3 different drawer type : Digital type (AAAE15101), Key lock type (AAAE15102), Unlocked type (AAAE15103)

2) For 230V (Korean-style). To use a 120V (US-style) socket, select AAAE1531.

3) We provide the cable port (AAAE1521) only for CLG-300, 650, and 1400.

4) Remote Alarm System (AAAE15511) is included for PSR3-70 model.

Related Product

Cryo Box

- Excellent durability even at ultra-low temperature with stable PC material.
- Easily distinguish samples with numbers on the covers.
- Cover, tray and dividers are separated for easy cleaning.
- Pressurized sterilization possible. (autoclavable)



Model	CRB-25	CRB-81S	CRB-81T
Tube size (ml)	1.2 to 2.0	1.2 to 2.0	5
Max. Mountable tube	25 (5 x 5)	81 (9 x 9)	81 (9 x 9)
Dimension (W x D x H, mm / inch)	77 x 77 x 53 / 3 x 3 x 2.1	132 x 132 x 54 / 5.2 x 5.2 x 2.1	132 x 132 x 96 / 5.2 x 5.2 x 3.8
Cat. No.	AAAJ7001	AAAJ7002	AAAJ7003

※ 5 packs for one set. (each Cat. No.)

Storage Box

- Excellent chemical resistance due to polypropylene material.
- The circular grooves at the bottom ensure stable mounting of tubes.
- Easy to separate samples with numbers and letters on covers.
- Pressurized sterilization possible. (autoclavable)



Model	SRB-100
Tube size (ml)	1.2 to 2.0
Max. Mountable tube	100 (10 x 10)
Dimension (W x D x H, mm / inch)	140 x 140 x 53 / 5.5 x 5.5 x 2.1
Cat. No.	AAAJ7011

※ 5 packs for one set.



Lab Essentials



Lab Essentials



Liquid Handling

Disposable Beaker

Made by
JeioTech

✓ Economical beaker with minimal price burden

- Domestic products made by mold designed and manufactured by Jeio Tech's, minimizing the cost burden.

✓ Structure that can allow for accurate pouring without spillage

- The unique triangular spout structure allows even highly viscous solutions to be poured precisely without any dripping.

✓ Extended use area and excellent material

- Flat bottom shape for use on a magnetic agitator
- Convenient to use as it does not bend easily even when grasped in the hand.
- Polypropylene (PP) material with excellent chemical resistance.
- Autoclavable, possible to use up to 121°C.
- Four types of composition: 50 ml, 100 ml, 250 ml, and 500 ml.



Disposable Beaker
with Magnetic Stirrer (MS-T)

✓ Convenient design for actual users

- Stackable design that can store stacked beakers of the same capacity and save storage space.
- Easy to remove overlapping beakers.
- Semi-transparent material makes it easy to check the solution state from the side.
- Graduation are engraved and do not come off even after long periods of use.



Structure allows for pouring without spillage

The triangular spout allows for dripless pouring of even highly viscous solutions.



Can be used on magnetic stirrer

Flat-bottom structure allows for use on a magnetic stirrer by placing a magnetic bar.



Structure allowing for stacked storage

It is possible to store beakers of the same capacity by stacking and it is designed to allow the beakers to be removed easily.

Cat. No.	Model	Capacity (ml)	Graduation (ml)	Dimension (Ø x H, mm / inch)
00PER0000023	Disposable beaker, 50 ml	50	5	45 x 55 / 1.8 x 2.2
00PER0000024	Disposable beaker, 100 ml	100	5	56 x 69 / 2.2 x 2.7
00PER0000025	Disposable beaker, 250 ml	250	10	77 x 92 / 3.0 x 3.6
00PER0000026	Disposable beaker, 500 ml	500	20	96 x 116 / 3.8 x 4.6

Plastic Griffin Beaker (Transparent / Translucent)

Made by
JeioTech

▮ Selectable materials according to purpose of use

- Choice of two series: polymethylpentene (PMP) transparent material and polypropylene (PP) translucent material with excellent chemical resistance.
- PMP: 135°C / PP: 121°C autoclavable.

▮ Structure that can allow for accurate pouring without spillage

- Pointed spout allows for dripless pouring of even highly viscous solutions.

▮ Extended use area and excellent material

- Flat bottom shape for use on a magnetic stirrer.
- Structure that does not become easily deformed.
- Convenient to use as it does not bend easily even when grasped in the hand.
- Plastic griffin beaker conforming to ISO 7056 standards.
- Repetitive use possible with excellent chemical resistance.
- Five different capacity configurations from 100 ml to 2000 ml.

▮ Convenient design for actual users

- Stackable design that can store stacked beakers of the same capacity and save storage space.
- Easy to remove stacked beakers.
- Graduation are engraved and do not come off even after long periods of use.



Plastic Griffin Beaker (PMP)



PMP griffin beaker (transparent)



PP griffin beaker (translucent)



Structure allows for pouring without spillage

Pointed spout allows for dripless pouring of even highly viscous solutions.



Engraved graduation type

Graduation are engraved and do not come off even after long periods of use.



Structure allowing for stacked storage

It is possible to store beakers of the same capacity by stacking and it is designed to allow the beakers to be removed easily.

PMP Griffin Beaker (transparent)

Cat. No.	Capacity (ml)	Graduation (ml)	Dimension (Ø x H, mm / inch)
00PER0000032	100	5	61 x 72 / 2.4 x 2.8
00PER0000033	250	10	81 x 95 / 3.2 x 3.7
00PER0000034	500	20	101 x 119 / 3.9 x 4.7
00PER0000035	1000	50	125 x 148 / 4.9 x 5.8
00PER0000036	2000	100	157 x 186 / 6.2 x 7.3

PP Griffin Beaker (translucent)

Cat. No.	Capacity (ml)	Graduation (ml)	Dimension (Ø x H, mm / inch)
00PER0000027	100	5	61x72 / 2.4x2.8
00PER0000028	250	10	81x95 / 3.2x3.7
00PER0000029	500	20	101x119 / 3.9x4.7
00PER0000030	1000	50	125x148 / 4.9x5.8
00PER0000031	2000	100	157x186 / 6.2x7.3

Liquid Handling

Handle Beaker (Transparent / Translucent)

Made by JeioTech

Beaker with a convenient handle

- Handles are convenient for handling, and large capacities can be handled with just one hand.
- The inside of the handle is rounded to minimize pain of hand when holding.

Structure that can allow for accurate pouring without spillage

- Pointed spout allows for dripless pouring of even highly viscous solutions.

Extended use area and excellent material

- Flat bottom shape for use on a magnetic stirrer.
- Structure where beaker wall and floor thickness are both constant and do not become easily deformed.
- Plastic beaker conforming to ISO 7056 standards.
- Repetitive use possible with excellent chemical resistance.
- 500ml, 1000ml, and 2000ml capacity available.

Convenient design for actual users

- Stackable design that can store stacked beakers of the same capacity and save storage space.
- Easy to remove stacked beakers.
- Graduation are engraved and do not come off even after long periods of use.



Handle Beaker (PMP)



PMP Beaker (transparent)

Cat. No.	Capacity (ml)	Graduation (ml)	Dimension (W x Ø x H, mm / inch)
00PER0000042	500	10	148 x 95 x 134 / 5.8 x 3.7 x 5.3
00PER0000039	1000	20	183 x 121 x 166 / 7.2 x 4.8 x 6.5
00PER0000040	2000	50	221 x 149 x 212 / 8.7 x 5.9 x 8.3

PP Beaker (translucent)

Cat. No.	Capacity (ml)	Graduation (ml)	Dimension (W x Ø x H, mm / inch)
00PER0000041	500	10	148 x 95 x 134 / 5.8 x 3.7 x 5.3
00PER0000037	1000	20	183 x 121 x 166 / 7.2 x 4.8 x 6.5
00PER0000038	2000	50	221 x 149 x 212 / 8.7 x 5.9 x 8.3

Wash Bottle

Made by JeioTech

Wash Bottle Common Features

- Easily sprayed solution even when lightly pressed for convenience.
- Made of excellent chemical-resistant material. (Bottle: LDPE, Cap & Tube: PP)
- Easy to use with wide hole for supplying liquid.
- Double-helical structure increases tightness of both container and cap.
- Replaceable nozzle tip for effective maintenance.
- Engraved markings for convenient capacity checking
- Built-in safety labels for identification provided as standard.

Automatic Gas-vent Type

- Vapor pressure self-discharge structure inside the container eliminates the inconvenience of liquid flowing out.



Model	WB-N25	WB-N50	WB-G25	WB-G50
Description	Wash Bottle		Automatic Gas-Vent Wash Bottle	
Capacity (ml)	250	500	250	500
Cat. No.	AAAM2501	AAAM5001	AAAM2511	AAAM5011

Distilled water container (bottle with lower spigot), HDPE material


- Excellent chemical resistance, suitable for storage and dispensing of various solutions.
- Leak proof; No leakage even when using for a long time due to spigot.
- Adjustable dispensing capacity according to the rotation angle of the spigot handle.
- The spigot is rotatable, and uses a tubing connection with an adapter.
- Adapter can use tubing I.D. 6.35 - 7.94mm.
- The cap is made of silicone gasket, so it has excellent sealing.
- Easy to check the capacity with front-side scale.
- Space-efficient design for solution storage.
- DBR-20 model fits to Shaker Heavy Duty Holder. [Page 127](#)
- In case it is not able to connect the water direct system to Temperature & Humidity Chamber, Carboy can be used for the distilled water container. [Page 280](#)



Round Type



Rectangular Type



Model	DBR-10	DBR-20	DBC-05	DBC-10
Description	Round Type		Rectangular Type	
Volume (L / gal)	10 / 2.64	20 / 5.28	5 / 1.32	10 / 2.64
Dimension (W x D x H, mm / inch)	225 x 277 x 399 / 8.9 x 10.9 x 15.7	260 x 317 x 539 / 10.2 x 12.5 x 21.2	114 x 387 x 305 / 4.5 x 15.2 x 12	204 x 387 x 317 / 8 x 15.2 x 12.5
Cat. No.	AAAJ3021	AAAJ3031	AAAJ3111	AAAJ3121

Storage

Storage Box

- Excellent chemical resistance due to polypropylene material.
- Easy to separate samples with numbers and letters on covers.
- Pressurized sterilization possible. (autoclavable)
- Temperature range of use (°C / °F): -20 to 121 / -4 to 249.8



Model	SRB-100
Tube size (ml)	1.2 to 2.0
Max. Mountable tube	100 (10 x 10)
Dimension (W x D x H, mm / inch)	140 x 140 x 53 / 5.5 x 5.5 x 2.1
Cat. No.	AAAJ7011

※ 5 packs for one set.

Centrifuge Tube Rack, PC

- Excellent chemical resistance due to polycarbonate material
- 30ml centrifuge tube rack
- Pressurized sterilization possible. (autoclavable)



Model	30 ml Centrifuge Tube Rack
Tube size (ml)	30
Max. mountable tube	6
Dimension (W x D x H, mm / inch)	263 x 55 x 55 / 10.4 x 2.2 x 2.2
Cat. No.	AAAJ2111

Conical Tube Rack

- Possible to use independently or in microplate shaker. [Page 131](#) (except for 50ml conical tube rack, 6 holes)



Model	15 ml Conical Tube Rack	50 ml Conical Tube Rack	
Material	ABS (Acrylonitrile-butadienestyrene)	ABS (Acrylonitrile-butadienestyrene)	PC (Polycarbonate)
Dimension (W x D x H, mm / inch)	269 x 85 x 72 / 10.6 x 3.3 x 2.8	269 x 85 x 72 / 10.6 x 3.3 x 2.8	263 x 55 x 55 / 10.4 x 2.2 x 2.2
Max. mountable tube	30 (3 x 10)	12 (2 x 6)	6
Autoclavable	-	-	0
Cat. No.	00CPS000011	00CPS000012	AAAJ2211

Micro Tube Rack, PC

- Excellent chemical resistance due to polycarbonate material.
- Different size tubes available on both sides. (1.5ml/2.0ml micro tube on one side and 0.2ml/0.5ml micro tube on the other side)
- Available for independent use and for use in microplate shaker. [Page 131](#)
- Pressurized sterilization possible. (autoclavable)



Model	1.5 / 2.0 ml Micro Tube Rack
	0.2 / 0.5 ml Micro Tube Rack
Dimension (W x D x H, mm / inch)	269 x 90 x 50 / 10.6 x 3.5 x 2.0 (Including cover)
Max. mountable tube	80 (5 x 16), each side
Cat. No.	AAAJ2311

Cryo Box

Made by
JeioTech

- Excellent durability even at ultra-low temperature with stable PC material.
- Easily distinguish samples with numbers on the covers.
- Cover, tray and dividers are separated for easy cleaning. (Except for CRB-25)
- Pressurized sterilization possible. (autoclavable)
- Temperature range of use (°C / °F) :-196 to 121 / -320.8 to 249.8



Model	CRB-25	CRB-81S	CRB-81T
Tube size (ml)	1.2 to 2.0	1.2 to 2.0	5
Max. Mountable tube	25 (5 x 5)	81 (9 x 9)	81 (9 x 9)
Dimension (W x D x H, mm / inch)	77 x 77 x 53 / 3.0 x 3.0 x 2.1	132 x 132 x 54 / 5.2 x 5.2 x 2.1	132 x 132 x 96 / 5.2 x 5.2 x 3.8
Cat. No.	AAAJ7001	AAAJ7002	AAAJ7003

※ 5 packs for one set. (each Cat. No.)

Safety Waste Disposal System

More detail page 234

Made by
JeioTech

Easy to collect and store various chemicals and liquid wastes generated in the lab

- > Structure minimizes vapor emission**
 Hinge-type lid minimizes the discharge of vapors from the waste solution into the lab.
- > Exclusive carbon filter (option)**
 Removes bad odors and harmful gas. (carbon filter can be replaced easily without removing existing attachment without the use of any tools)
- > Safe anti-static function (SF-14C, WT-20C model)**
 Prevention of explosion due to static electricity by using special plastic material with anti-static performance, and grounding cable is provided as standard. (completed anti-static performance test by KEPCO)
- > Water level notification function (SF model)**
 Level bar displayed on the safety funnel when at 17.6L (80%) or 15.4L (70%) based on a 22L waste liquid container.



SF-14P & WT-20P
with Carbon Filter (option)



SF-14C & WT-20C
with Carbon Filter (option)



Accessories

Description	Dimension (ø x H, mm / inch)	Cat. No.
Carbon Filter	76 x 150 / 3.0 x 5.9	AAAM6051
Earth Bonding Plug	52 x 70 / 2.0 x 2.8	AAAF35453

Model	SF-14C	SF-14P	WT-20C	WT-20P
Description	Anti-static safety funnel for collection of waste liquid	General use safety funnel for collection of waste liquid	Anti-static liquid waste container	General-purpose liquid waste container
Dimension (W x D x H, mm / inch)	160 x 142 x 270 / 6.3 x 5.6 x 10.6	160 x 142 x 270 / 6.3 x 5.6 x 10.6	220 x 300 x 380 / 8.7 x 11.8 x 15.0	220 x 300 x 380 / 8.7 x 11.8 x 15.0
Inlet (ø, mm / inch)	135 / 5.3	135 / 5.3	45 / 1.8	45 / 1.8
Lid size (ø, mm / inch)	142 / 5.6	142 / 5.6	78 / 3.1	78 / 3.1
Volume (L / cu ft)	1 / 0.04	1 / 0.04	22 / 0.78	22 / 0.78
Cat. No.	AAAM6001	AAAM6002	AAAM6011	AAAM6012

Storage

On-Ice Rack

Made by JEIOtech

- Block made of aluminum with excellent thermal conductivity.
- Keeps the temperature of the sample constant when using on ice.
- Excellent thermal conductivity with aluminum block.



Model	On-Ice Rack	On-Ice Rack (96W)
Dimension (W x D x H, mm / inch)	157 x 152 x 21 / 6.18 x 5.98 x 0.83	157 x 152 x 20 / 6.18 x 5.98 x 0.79
Max. mountable tube	1.5 ml tube 50ea, 0.2 ml tube 64ea	1.5 ml tube 30ea, 96 well plate 1ea
Cat. No.	AAAJ2011	AAAJ2021


Flask Support

Made by JEIOtech

- Safe storage and handling of round flasks.
- Convenient handling of various shapes and simultaneous handling of several things.
- Excellent chemical resistance and heat resistance. (-30 to 120°C)
- Side handles for easy movement.



Various types of flasks can be used for various purposes simultaneous storage



Convenient movement with side handles

Easy to distinguish samples by name tags and numbers



Model	FT 010 13	FT 010 34	FT 025 12	FT 025 23	FT 050 12	FT 100 12	FT 200 12
Capacity (ml)	100	100	250	250	500	1000	2000
Row x column	1 x 3	3 x 4	1 x 2	2 x 3	1 x 2	1 x 2	1 x 2
Hole (dia. x H, mm/inch)	70 x 91 / 2.76 x 3.58	70 x 91 / 2.76 x 3.58	98 x 122 / 3.86 x 4.8	98 x 122 / 3.86 x 4.8	121 x 128 / 4.76 x 5.04	148 x 156 / 5.83 x 6.14	183 x 188 / 7.2 x 7.4
Dimension (W x D x H, mm/inch)	310 x 100 x 100 / 12.2 x 3.94 x 3.94	384 x 265 x 98 / 15.12 x 10.43 x 3.86	280 x 125 x 130 / 11.02 x 4.92 x 5.12	389 x 239 x 130 / 15.31 x 9.41 x 5.12	360 x 146 x 135 / 14.17 x 5.75 x 5.31	415 x 192 x 165 / 16.34 x 7.56 x 6.5	503 x 223 x 199 / 19.8 x 8.78 x 7.83
Cat. No.	AAA23706	AAA23701	AAA23707	AAA23702	AAA23703	AAA23704	AAA23705

Hood

More detail page 210

Made by
JeioTech



PMH-720



POH-720
with EM-33 (option)

Fume Hood, Small type

- Excellent chemical resistance due to polyethylene material.
- Top hole that can be connected to the exhaust system.
- The front and back side are made by wide transparent window (PC), useful for simulation and observing in school. (only POH models)

Model		POH-720	PMH-720
Material	Body	Polyethylene	Polyethylene
	Observation window	O	X
Dimension	Exterior (mm / inch)	720 x 472 x 535 / 28.3 x 18.6 x 21.1	720 x 472 x 535 / 28.3 x 18.6 x 21.1
	Duct hole (O.D., mm / inch)	150 / 5.91	150 / 5.91
	Net weight (kg / lbs)	12 / 26.5	12 / 26.5
Cat. No.		AAAB5021	AAAB5011

Desiccator

More detail page 250

Made by
JeioTech

Vacuum Desiccator, Cubic

- Max. 1 Torr vacuum maintained for 72 hours.
- Strong and transparent polycarbonate material.
- Includes vacuum gauge as standard.
- Height-adjustable shelves.
- Multi-level stacking structure for space efficiency.

Desiccator

- Hygrometer included.
- Transparent polycarbonate material makes for convenient internal observation.
- Two drawers provided.
- Includes desiccant tray as standard.
- UV blocking model provided for light-sensitive samples.



Cat. No.	Capacity (L/cu ft)
AAAD4001	11 / 0.39
AAAD4021	23 / 0.81
AAAD4011	35 / 1.24
AAAD4031	45 / 1.59

Cat. No.	Capacity (L / cu ft)
AAAD5011 (Clear)	35 / 1.24
AAAD5111 (UV protection)	35 / 1.24
AAAD5021 (Clear)	50 / 1.77
AAAD5121 (UV protection)	50 / 1.77

Connection

Tubing Connector

Made by JeioTech

- Excellent chemical resistance due to polypropylene material.
- 121°C, Autoclavable.
- Structure prevents tubes from dropping, adding convenience.
- T-shape, Y-shape, Cross Type, Straight Type.
- Various sizes from 3.2 mm to 12.7 mm.




O.D. (mm)	O.D. (inch)	Cat. No.			
		T-shape	Y-shape	Cross	Straight
3.2	1/8	00CJA0009042	00CJA0009057	00CJA0009048	N/A
4.0	5/32	00CJA0009041	00CJA0009058	00CJA0009047	N/A
4.7	3/16	00CJA0009044	00CJA0009046	00CJA0009050	00CJA0009051
6.0	1/4	00CJA0009043	00CJA0009045	00CJA0009049	00CJA0009054
8.0	5/16	00CJA0009055	00CJA0009059	00CJA0009061	N/A
9.5	3/8	00CJA0009036	00CJA0009037	00CJA0009062	00CJA0009052
12.7	1/2	00CJA0009056	00CJA0009060	00CJA0009063	00CJA0009053

Filter Adapter

Made by JeioTech

Filter Adapter

- Can be used between the filtering funnel and the filtering flask during vacuum filtration.
- Excellent sealing and adhesion with silicone material.
- Various uses for flasks as they come in various sizes.



Cat. No.	Description
AAA72200	Filter Adapter (5ea / pk)
AAA72210	Filter Adapter (10ea / pk)

Clamp

Made by
JeioTech

Clamp Holder

- Aluminum die-casting and plastic materials allow for lightweight and rigid fixation.
- C-20 and C-10 are suitable for overhead stirrers.
- C-10, C-7 and C-5 models are suitable for use with hotplates and magnetic stirrers.
- Applicable to 5/8" diameter support rod.



C-20



C-10



C-7



C-5

Model	C-20	C-10	C-7	C-5
Max. Ø (mm / inch)	27 / 1.06	17 / 0.67	12 / 0.47	12 / 0.47
Cat. No.	AAA37511	AAA37512	AAA37513	AAA37514

Utility Clamp



Utility Clamp 3 Prong



Utility Clamp Plat



3 Prong Clamp



2 Prong Clamp




Fixing Clamp

Cat. No.	Description	Cat. No.	Description
BEA1000001	Utility Clamp 3 Prong	BEA1000013	3 Prong Clamp (20 mm grip)
BEA1000002	Utility Clamp Plat (2P)	BEA1000014	2 Prong Clamp (60 mm grip)
BEA1000011	3 Prong Clamp (80 mm grip)	BEA0570181	Fixing Clamp
BEA1000012	3 Prong Clamp (60 mm grip)		



Environmental Test Equipment



Environmental Test Equipment



Temperature & Humidity Chamber



Satisfaction with the requirements of European regulations and guidelines.



Ensuring the reliability of electronic records and signatures to support GMP audit.



Real-time equipment monitoring and control system using mobile app.



Eco-friendly product that has passed environmental impact assessment tests.



Recognized as an excellent design product by the ministry of commerce, industry, and energy.



More precise temperature control through temperature calibration.



Electronic over temperature protection system.



Conventional over temperature protection, a kind of backup device against electronic over temperature protection failure.



Control of temperature / humidity / rpm / time, etc. through dedicated program.



PC communication via RS-232 / RS-485 / USB port.



User safety by blocking potentially harmful UV light.



Indicator to ascertain water level.



Alarm if the water is low level.



It memorizes the state of power failure and operates automatically when power is restored.



2 year warranty free A/S.

General Application

Temperature and humidity chamber: Test for reliability, durability, climatic, freezing resistance, quality assurance, thermal endurance.

Pharmaceuticals stability test chamber: Stability test in accordance with ICH Guideline Q1A (R2) and Q1B (Option 2)

	Description	Temperature Range (°C / °F)	Humidity Range (%RH)	Volume (L / cu ft)	Model	Page
Small type	Tabletop type <ul style="list-style-type: none"> • Compact design Temp. & Humid. chambers • Continuous innovative model 	-40 to 100 / -40 to 212	20 ~ 95	25, 65, 100 / 0.9, 2.3, 3.5	TH3-KE	286
		-20 to 100 / -4 to 212	20 ~ 95	25, 65, 100 / 0.9, 2.3, 3.5	TH3-ME	286
		0 to 100 / 32 to 212	30 ~ 95	25, 65, 100 / 0.9, 2.3, 3.5	TH3-PE	286
Medium-large size type	Ultra-low temperature type (-70°C) <ul style="list-style-type: none"> • Wide control range for temperature and humidity • Differentiated structural mechanisms 	-70 to 100 / -94 to 212	25 ~ 95	189, 302, 408 / 6.7, 10.7, 14.4	TH-KH	290
	Low temperature type (-40°C) <ul style="list-style-type: none"> • Proven reliability and reproducibility • Optimized control system 	-40 to 100 / -40 to 212	25 ~ 95	189, 302, 408, 800, 1000, 1450 / 6.7, 10.7, 14.4, 28.3, 35.3, 57.1	TH-G	290
	General type <ul style="list-style-type: none"> • Stable temperature and humidity control for even during long-term continuous use 	-5 to 100 / 23 to 212	10 ~ 95 25 ~ 95	180, 314, 419, 770, 997, 1540 / 6.4, 11.1, 14.8, 27.2, 35.2, 54.4	TH-TG	295
	Dual chamber <ul style="list-style-type: none"> • Space-efficient with two chambers. • Independent temperature and humidity control of each chamber 	-5 to 100 / 23 to 212	20 ~ 90	150, 288, 384 (x 2 chamber) / 5.3, 10.2, 13.6 (x 2 chamber)	TH3-DG	299
	Basic type <ul style="list-style-type: none"> • Dedicated controller with enhanced convenience • Improved water supply lines 	0 to 90 / 32 to 194	35 ~ 85	201, 398 / 7.1, 14.1	TH3-E	302
	Pharmaceuticals stability test <ul style="list-style-type: none"> • ICH Q1A (R2), Q1B (Option 2) • Dedicated lamp control system 	-5 to 85 / 23 to 185	20 ~ 85	300, 760 / 10.6, 26.8	TH-ICH	305
Room chamber type	Walk-in <ul style="list-style-type: none"> • Optimized temperature and humidity distribution with large room size walk-in chamber • Easy to use by optimizing the position of the control part according to the installation space and user requirements 					310
Custom type	Customized <ul style="list-style-type: none"> • Design and manufacture the customized Temp. & Humid. Chamber optimized by user requirements • Optimal customized equipment provided based on advanced technology and accumulated experience 					311

※ The contents of the above and the contents of this catalog may differ depending on the specific model and conditions of use. For the information about the features and specifications that applying to each models, please check the information on the corresponding page of each models.

Temperature & Humidity Chamber



Temperature & Humidity Chamber

Proven reliability and reproducibility Continuous innovative model

> Stable temperature and humidity distribution

Specialized internal circulation for fan/motor system structure.
Provides proven reproducibility and high reliability.

> Innovative external humidifier

Stable humidity control and maintains internal cleanliness with differentiated external humidification system.

> Complete chamber structure

Includes double-packing door.
Stable temperature control with excellent insulation and sealing.

> Minimal water consumption

Minimized refilling of water even for long-term high-humidity tests.
Efficient water supply and recirculation system.

Powerful & Useful Control Optimized precision control

> Intuitive operation with color touch display

Easy to understand operation of equipment and operation status with 5.6 inch color touch display controller.

> Powerful program operation

Provides diverse program operation environment with a maximum of 100 patterns and 2000 segments.

> Sophisticated custom control

PID zone subdivided into 6 zones for more precise control stability.

> Multi-functional dedicated software

Up to 32 devices can be controlled at same time.
Device operation and data processing function supported.



Ease of check the test status with LED display showing the operation condition.



Electronic over temperature shutdown system.



Controller for temperature and humidity control. (TH3-model)



Maximized usability with simple and useful functions

> Ease of check the test status

Convenient internal observation with LED lamp and looking glass window.
Figures out the operation state through displaying LED color even a long distance.

> Door lock with key lock

Door structure having well-sealed.
Included key lock features as standard.

> Provision of remote equipment monitoring service (TH3, Option)

Real-time monitoring of equipment operation status using smartphone.
(when purchased LC GreenBox)

> Convenient moving and installation

Easy to move and install with casters.

More convenient and safe structure applied the highest grade safety system

> Electronic over temperature shutdown system

Independent precision temperature sensor and controller.
Provides more secure and reliable over temperature shutdown.

> General over temperature shutdown system

Backup device for electronic system failure.
Implementation of over-temperature shutdown even for electronic errors.

> Self-protection system

Protection system as like water tank shortage, door open warning, etc.
Over-current and short circuit breaking device.

> Protection features for refrigeration system

Protection system for overload and over temperature on refrigerator.
Stops automatically in case refrigerant pressure (high/low pressure) is abnormal.

※ Some of the above contents are limited to specific models.



Tempered glass window with wire heater and LED lamp.



Door handle and locking device.



Ø50mm cable port and silicone plug included as standard.

Temperature & Humidity Chamber

Tabletop type

More convenience small type Temp. & Humid chamber with new advanced controller

Introduction video for TH chamber Tabletop type



TH3-ME

with Wire Shelves 2ea, Φ 50 Cable Port (standard), Stand (option)



Optimized Model Configuration

- 3 series (KE, ME, PE) according to temperature and humidity control range.
TH3-KE; -40 to 100°C (20~95%RH)
TH3-ME; -20 to 100°C (20~95%RH)
TH3-PE; 0 to 100°C (30~95%RH)
- Standard type of Temp. & Humid. Chamber to select the optimized model according to test conditions.
- Dedicated stand for convenient use and moving. (option)
- Recorder for maintaining test records. (option)

Structural Functional Features

- Proven reproducibility and reliability with structural design optimized for miniaturization.
- Compact design minimized foot space for excellent space saving.
- Accessible slim design to laboratories which have only 1-door entry.
- Stable humidity control ability with applying differentiated external humidification system and convenient maintenance for cleaning testing space.
- Excellent thermal insulation and sealing to ensure stable temperature and humidity reproducibility even for repeated test use.
- The door part has a double sealing structure and excellent airtightness.



Use Convenience Features

- Integrated with 21 CFR Part 11 compliant LC DataKeeper, suitable for GMP audits (option)
- Convenient operation with 5.6 inch color touch controller.
- Possible to check the operation status with graphs.
- Water level check function enables water replenishment on time.
- Water inlet is easy to use with easy cap.
- Efficient water supply and recirculation system minimizes water refilling even for long-term high-humidity tests.
- Operation and status monitoring from the front.
- Easy observation of samples with tempered glass door with built in heating wire and LED lamp.
- By displaying the operation status of the device with LED colors, convenient to identify chambers status from a far.
- The cable port and silicone plug included as standard make it convenient to connect external equipment.
- Easy to clean as it is easy to remove refrigerator condenser grill, making it convenient to maintain efficiency of refrigeration.

Multifunction Control System

- Storage and documentation of programs.
- Save program information and back up the set values.
- RS-232 port (default) and RS-485 port (option) supported.
- Control and data processing by connecting up to 32 devices to PC at same time via RS-485 port.
- Reliable control with temperature/humidity auto-tuning and temperature calibration.
- PID zone subdivided into 9 zones for more precise control stability.
- Monitor and control via mobile app anytime, anywhere with LC Connected. (Requires LC GreenBox, option)
- Convenient checking with useful graph display.
- Save data as excel and text files.

Description	No.
Pattern	100
Repeat time	999
Part repeat time	255
Max. segments / pattern	100
Available max. segments	2000
Programmable process time / segment	99 hr. 59 min.

* It is possible to set 100 segments per pattern, but the maximum number of segments is not 10,000 (100 pattern x 100 segment) but 2,000.



Color touch display controller (5.6 inches)



Indicator LED of operation condition






Tempered glass window with wire heater and LED lamp



Remote monitoring and control via LC Connected app

▮ Outstanding Safety

- Electronic over temperature shutdown system.
(advanced over temperature (A-OT) limiter)
Independent precision temperature sensor and controller for safer and more accurate over temperature shutdown.
- General over power shutdown system.
(B-OT, backup over temperature limiter)
A backup device for electronic system failure, mechanically preventing over temperature even when it occurs electronic errors.
- Ensures that the surface temperature is kept below 60°C.
- Low and empty water level alarm and automatic shut off.
- Door open warning and automatic shut off.
- Overcurrent and short circuit protection of device.
- Fuse installed for each heater for safety.
- Auto stop when operating current of compressor is overloaded.



Specification

Model	TH3-KE-025	TH3-KE-065	TH3-KE-100	TH3-ME-025	TH3-ME-065	TH3-ME-100	TH3-PE-025	TH3-PE-065	TH3-PE-100
Interior dimensions									
Chamber volume (L / cu ft)	25 / 0.9	65 / 2.3	100 / 3.5	25 / 0.9	65 / 2.3	100 / 3.5	25 / 0.9	65 / 2.3	100 / 3.5
Width (mm / inch)	350 / 13.8	400 / 15.7	500 / 19.7	350 / 13.8	400 / 15.7	500 / 19.7	350 / 13.8	400 / 15.7	500 / 19.7
Depth (mm / inch)	220 / 8.7	330 / 13.0	380 / 15.0	220 / 8.7	330 / 13.0	380 / 15.0	220 / 8.7	330 / 13.0	380 / 15.0
Height (mm / inch)	350 / 13.8	495 / 19.5	530 / 20.9	350 / 13.8	495 / 19.5	530 / 20.9	350 / 13.8	495 / 19.5	530 / 20.9
Quantity of shelves (standard/max.)	2 / 7	2 / 10	2 / 10	2 / 7	2 / 10	2 / 10	2 / 7	2 / 10	2 / 10
Distance of between shelves (mm / inch)	20 / 0.8	20 / 0.8	20 / 0.8	20 / 0.8	20 / 0.8	20 / 0.8	20 / 0.8	20 / 0.8	20 / 0.8
Max. load per shelf (kg / lbs)	25 / 55.1	25 / 55.1	25 / 55.1	25 / 55.1	25 / 55.1	25 / 55.1	25 / 55.1	25 / 55.1	25 / 55.1
Permitted total load (kg / lbs)	65 / 143.3	65 / 143.3	65 / 143.3	65 / 143.3	65 / 143.3	65 / 143.3	65 / 143.3	65 / 143.3	65 / 143.3
Cable port (Ø, mm / inch)	50 / 2.0	50 / 2.0	50 / 2.0	50 / 2.0	50 / 2.0	50 / 2.0	50 / 2.0	50 / 2.0	50 / 2.0
Water tank capacity (L / cu ft)	8 / 0.3	12 / 0.4	12 / 0.4	8 / 0.3	12 / 0.4	12 / 0.4	8 / 0.3	12 / 0.4	12 / 0.4
Weight (kg / lbs)	110 / 242.5	140 / 308.6	150 / 330.7	80 / 176.4	110 / 242.5	130 / 286.6	80 / 176.4	110 / 242.5	130 / 286.6
Temperature data									
Range without humidity (°C / °F)	-35 to 100 / -31 to 212 (50Hz)			-20 to 100 / -4 to 212 (50/60Hz)			0 to 100 / 32 to 212 (50/60Hz)		
	-40 to 100 / -40 to 212 (60Hz)								
Fluctuation (±°C / °F) ¹⁾	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54
Variation (±°C / °F) ¹⁾	0.5 / 0.9	0.5 / 0.9	0.5 / 0.9	0.5 / 0.9	0.5 / 0.9	0.5 / 0.9	0.5 / 0.9	0.5 / 0.9	0.5 / 0.9
Heating time (min) ²⁾	60	60	60	50	50	50	40	40	40
Cooling time (min) ³⁾	70	70	70	55	55	55	35	35	35
Humidity data (※ Ref.: Temp./ Humid. control range graph)									
Range (%RH)	20 to 95	20 to 95	20 to 95	20 to 95	20 to 95	20 to 95	30 to 95	30 to 95	30 to 95
Fluctuation (±%RH) ¹⁾	1	1	1	1	1	1	1	1	1
Variation (±%RH) ¹⁾	3	3	3	3	3	3	3	3	3
Electrical data & Ordering information									
Electrical requirements (AC230V, 1P, 50Hz, A)	12.0	13.9	15.8						
Cat. No.	AAH824132KAAH824232KAAH824332K			-	-	-	-	-	-
Cat. No. (EU only)	AAH824135KAAH824235KAAH824335K			-	-	-	-	-	-
Electrical requirements (AC230V, 1P, 60Hz, A)	12.0	13.9	15.8						
Cat. No.	AAH824131KAAH824231KAAH824331K			-	-	-	-	-	-
Electrical requirements (AC230V, 1P, 50/60Hz, A)				7.7	9.6	13.8	7.2	9.4	10.5
Cat. No.	-	-	-	AAH821132KAAH821232KAAH821332KAAH822132KAAH822232KAAH822332K					
Cat. No. (EU only)	-	-	-	AAH821135KAAH821235KAAH821335KAAH822135KAAH822235KAAH822335K					

※ The above specifications are tested according to IEC 60068 standard.


※ Permissible environmental conditions are 18 - 30°C (recommended 20°C), 85%RH or less, under 2000m above sea level.

※ The DI water condition is pH 6.2 - 7.2 and the electrical conductivity is less than 20µs/cm.

1) The set value is set under the conditions of 25°C, 60%RH.

2) Heating time condition value: TH3-KE (-40°C to 100°C), TH3-ME (-20°C to 100°C), TH3-PE (0°C to 100°C)

3) Cooling time condition value: TH3-KE (20°C to -40°C), TH3-ME (20°C to -20°C), TH3-PE (20°C to 0°C)



Dimension

Model	TH3-KE-025	TH3-KE-065	TH3-KE-100	TH3-ME-025	TH3-ME-065	TH3-ME-100	TH3-PE-025	TH3-PE-065	TH3-PE-100
W (mm / inch)	550 / 21.7	600 / 23.6	700 / 27.6	550 / 21.7	600 / 23.6	700 / 27.6	550 / 21.7	600 / 23.6	700 / 27.6
D (mm / inch)	880 / 34.6	1050 / 41.3	1110 / 43.7	880 / 34.6	1050 / 41.3	1110 / 43.7	880 / 34.6	1050 / 41.3	1110 / 43.7
H (mm / inch)	915 / 36.0	1060 / 41.7	1100 / 43.3	915 / 36.0	1060 / 41.7	1100 / 43.3	915 / 36.0	1060 / 41.7	1100 / 43.3
A (mm / inch)	211 / 8.3	276 / 10.9	276 / 10.9	211 / 8.3	276 / 10.9	276 / 10.9	211 / 8.3	276 / 10.9	276 / 10.9
B (mm / inch)	520 / 20.5	580 / 22.8	580 / 22.8	520 / 20.5	580 / 22.8	580 / 22.8	520 / 20.5	580 / 22.8	580 / 22.8
E (mm / inch)	170 / 6.7	170 / 6.7	170 / 6.7	170 / 6.7	170 / 6.7	170 / 6.7	170 / 6.7	170 / 6.7	170 / 6.7
F (mm / inch)	229 / 9.0	229 / 9.0	229 / 9.0	229 / 9.0	229 / 9.0	229 / 9.0	229 / 9.0	229 / 9.0	229 / 9.0

Accessories Page 308 Cable Port, Shelves, Signal Lamp, Recorder, Water Purifier, Direct Water System, Stand, LC GreenBox, LC DataKeeper

Temperature & Humidity Chamber

Low temperature type (-70 to 100°C / -40 to 100°C)

General Temp. & Humid chamber with proven reliability and reproducibility

- **TH-KH** Ultralow temp. type (-70°C)
- **TH-G** Low temp. type (-40°C)



TH-G-180

with Wire Shelves 2ea, Φ 50 Cable Port (standard), Recorder (option)

Optimized Model Configuration

- Various models optimized according to temperature and humidity control range and specialized functions.
- TH-KH Series (-70 to 100°C, 25~95%RH)
Provides wide temperature control range.
- TH-G Series (-40 to 100°C, 25~95%RH)
Representative model with proven reliability.

Structural Functional Features

- Stable humidity control ability with applying differentiated external humidification system and convenient maintenance for cleaning testing space.
- Excellent thermal insulation and sealing to ensure stable temperature and humidity reproducibility even for repeated test use.
- Proven reproducibility and reliability with structural design optimized for each capacity.
- The door part is a double-sealed structure that provides enhanced airtightness.
- Durability even after long-term repeated use has been verified by satisfying the door load test (SEFA standard).



Use Convenience Features

- Efficient water supply and recirculation system minimizes water refilling even for long-term high-humidity tests.
- Drawer-type water tank structure makes refill the water easily at the front part.
- Possible structure to refill the water even during operation.
- Direct water supply system. (option)
- Easy observation of samples with tempered glass door with built in heating wire and LED lamp.
- Open/close door with latches on either side.
- Door handle with built-in key lock.
- Cable port (Ø50 mm) is included as standard, making it convenient to connect external equipment.
- Easy to clean as it is easy to remove refrigerator condenser grill, making it convenient to maintain efficiency of refrigeration.



Drawer-type water tank



Door handle and locking device

Multifunction Control System

- 5.7-inch color touch display controller for convenient intuitive operation and checking of operation status.
- Save program information and back up the set values.
- Includes temperature and humidity auto-tuning.
- Highly reliable control through temperature calibration.
- PID zone subdivided into 6 zones for more precise control stability.
- RS-232 port (default) and RS-485 port (option) supported.
- Control and data processing by connecting up to 32 devices to PC at same time via RS-485 port.
- Convenient checking with useful graph display.



Cable port



Tempered glass window with wire heater and LED lamp

Description	No.
Pattern	120
Repeat time	999
Part repeat time	255
Max. segments / pattern	100
Available max. segments	1200
Programmable Process time / segment	999 hr. 59 min.

* It is possible to set 100 segments per pattern, but the maximum number of segments is not 12,000 (120 pattern x 100 segment) but 1,200.








Controller (5.7 inch)

Temp. & humid. graph

▮ Outstanding Safety

- Electronic over temperature shutdown system (advanced over temperature. (A-OT) limiter)
Independent precision temperature sensor and controller for safer and more accurate over temperature shutdown.
- General over power shutdown system. (B-OT, backup over temperature limiter)
A backup device for electronic system failure, mechanically preventing over temperature even when it occurs electronic errors.
- Ensures that the surface temperature is kept below 60°C.
- Low and empty water level alarm and automatic shut off.
- Door open warning and automatic shut off.
- Overcurrent and short circuit protection of device.
- Fuse installed for each heater for safety.
- Auto stop when operating current of compressor is overloaded.
- Auto stop when compressor is overheating.
- Auto shut off in case refrigerant pressure (high/low pressure) is abnormal.



Specification

Model	TH-KH-180	TH-KH-300	TH-KH-408
Interior dimensions			
Chamber volume (L / cu ft)	189 / 6.7	302 / 10.7	408 / 14.4
Width (mm / inch)	500 / 19.7	600 / 23.6	600 / 23.6
Depth (mm / inch)	630 / 24.8	630 / 24.8	850 / 33.5
Height (mm / inch)	600 / 23.6	800 / 31.5	800 / 31.5
Quantity of shelves (standard/max.)	2 / 9	2 / 10	2 / 10
Distance of between shelves (mm / inch)	60 / 2.4	60 / 2.4	60 / 2.4
Max. load per shelf (kg / lbs)	20 / 44.1	25 / 55.1	25 / 55.1
Permitted total load (kg / lbs)	50 / 110.2	65 / 143.3	65 / 143.3
Cable port (∅, mm / inch)	50 / 2.0	50 / 2.0	50 / 2.0
Water tank capacity (L / cu ft)	15 / 0.5	15 / 0.5	15 / 0.5
Weight (kg / lbs)	320 / 705.5	385 / 848.8	410 / 903.9
Temperature data			
Range without humidity (°C / °F)	-70 to 100 / -94 to 212	-70 to 100 / -94 to 212	-70 to 100 / -94 to 212
Fluctuation (±°C / °F) ¹⁾	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54
Variation (±°C / °F) ¹⁾	0.5 / 0.9	0.5 / 0.9	0.5 / 0.9
Heating time (-70°C to 100°C) (min)	70	70	70
Cooling time (20°C to -70°C) (min)	80	80	80
Humidity data (※ Ref.: Temp./ Humid. control range graph)			
Range (%RH)	25 to 95	25 to 95	25 to 95
Fluctuation (±%RH) ¹⁾	0.4	0.4	0.4
Variation (±%RH) ¹⁾	3	3	3
Electrical data & Ordering information			
Electrical requirements (380V, 3P, 50Hz, A)	10.1	10.2	10.2
Cat. No.	AAHC3014K	AAHC3114K	AAHC3214K
Cat. No. (EU only)	AAHC3010K	AAHC3110K	AAHC3210K

※ The above specifications are tested according to IEC 60068 standard.


※ Permissible environmental conditions are 18 - 30°C (recommended 20°C), 85%RH or less, under 2000m above sea level.

※ The DI water condition is pH 6.2 - 7.2 and the electrical conductivity is less than 20µs/cm.

1) The set value is set under the conditions of 25°C, 60%RH.

Dimension

Model	TH-KH-180	TH-KH-300	TH-KH-408
W (mm / inch)	1050 / 41.3	1150 / 45.3	1150 / 45.3
D (mm / inch)	1615 / 63.6	1615 / 63.6	1835 / 72.2
H (mm / inch)	1560 / 61.4	1760 / 69.3	1760 / 69.3
A (mm / inch)	550 / 21.7	550 / 21.7	670 / 26.4
B (mm / inch)	1225 / 48.2	1325 / 52.2	1325 / 52.2
E (mm / inch)	200 / 7.9	200 / 7.9	200 / 7.9
F (mm / inch)	330 / 13.0	330 / 13.0	330 / 13.0



Accessories Page 308 Cable Port, Shelves, Signal Lamp, Recorder, Water Purifier, Direct Water System, Inner Glass Door

Specification


Model	TH-G-180	TH-G-300	TH-G-408	TH-G-800	TH-G-1000	TH-G-1500
Interior dimensions						
Chamber volume (L / cu ft)	189 / 6.7	302 / 10.7	408 / 14.4	800 / 28.3	1000 / 35.3	1450 / 57.1
Width (mm / inch)	500 / 19.7	600 / 23.6	600 / 23.6	1000 / 39.4	1000 / 39.4	1100 / 43.3
Depth (mm / inch)	630 / 24.8	630 / 24.8	850 / 33.5	800 / 31.5	1000 / 39.4	1000 / 39.4
Height (mm / inch)	600 / 23.6	800 / 31.5	800 / 31.5	1000 / 39.4	1000 / 39.4	1450 / 57.1
Quantity of shelves (standard/max.)	2 / 9	2 / 10	2 / 10	2 / 14	2 / 14	3 / 18
Distance of between shelves (mm / inch)	40 / 1.6	40 / 1.6	40 / 1.6	40 / 1.6	40 / 1.6	40 / 1.6
Max. load per shelf (kg / lbs)	20 / 44.1	25 / 55.1	25 / 55.1	40 / 88.2	40 / 88.2	40 / 88.2
Permitted total load (kg / lbs)	50 / 110.2	65 / 143.3	65 / 143.3	100 / 220.5	100 / 220.5	100 / 220.5
Cable port (Ø, mm / inch)	50 / 2.0	50 / 2.0	50 / 2.0	50 / 2.0	50 / 2.0	50 / 2.0
Water tank capacity (L / cu ft)	15 / 0.5	15 / 0.5	15 / 0.5	15 / 0.5	15 / 0.5	15 / 0.5
Weight (kg / lbs)	320 / 705.5	385 / 848.8	410 / 903.9	620 / 1366.9	680 / 1499.1	860 / 1896.0
Temperature data						
Range without humidity (°C / °F)	-40 to 100 / -40 to 212 (60Hz) -35 to 100 / -31 to 212 (50Hz)					
Fluctuation (±°C / °F) ¹⁾	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54
Variation (±°C / °F) ¹⁾	0.5 / 0.9	0.5 / 0.9	0.5 / 0.9	0.5 / 0.9	0.5 / 0.9	0.5 / 0.9
Heating time (-40°C to 100°C) (min)	60	60	60	70	70	90
Cooling time (20°C to -40°C) (min)	60	60	60	75	75	90
Humidity data (※ Ref.: Temp./ Humid. control range graph)						
Range (%RH)	25 to 95	25 to 95	25 to 95	25 to 95	25 to 95	25 to 95
Fluctuation (±%RH) ¹⁾	0.4	0.4	0.4	1	1	2
Variation (±%RH) ¹⁾	3	3	3	3	3	3
Electrical data & Ordering information						
Electrical requirements (230V, 1P, 50Hz, A)	25.2	25.6	25.6			
Cat. No.	AAH80013K	AAH80023K	AAH80033K	-	-	-
Cat. No. (EU only)	AAH80015K	AAH80025K	AAH80035K	-	-	-
Electrical requirements (230V, 1P, 60Hz, A)	25.2	25.6	25.6			
Cat. No.	AAH80011K	AAH80021K	AAH80031K	-	-	-
Electrical requirements (380V, 3P, 50/60Hz, A)	9.1	9.2	9.2	12.2	12.2	15.5
Cat. No.	AAH80018K	AAH80028K	AAH80038K	AAH80048K	AAH80058K	AAH80078K
Cat. No. (EU only)	AAH80010K	AAH80020K	AAH80030K	AAH80040K	AAH80050K	AAH80070K

※ The above specifications are tested according to IEC 60068 standard.

※ Permissible environmental conditions are 18 - 30°C (recommended 20°C), 85%RH or less, under 2000m above sea level.

※ The DI water condition is pH 6.2 - 7.2 and the electrical conductivity is less than 20µs/cm.

1) The set value is set under the conditions of 40°C, 75%RH.



Dimension

Model	TH-G-180	TH-G-300	TH-G-408	TH-G-800	TH-G-1000	TH-G-1500
W (mm / inch)	1003 / 39.5	1103 / 43.4	1103 / 43.4	1503 / 59.2	1503 / 59.2	1503 / 59.2
D (mm / inch)	1139 / 44.8	1139 / 44.8	1359 / 53.5	1319 / 51.9	1519 / 59.8	2065 / 81.3
H (mm / inch)	1468 / 57.8	1668 / 65.7	1668 / 65.7	1968 / 77.5	1968 / 77.5	2007 / 79.0
A (mm / inch)	412 / 16.2	432 / 17.0	532 / 20.9	509 / 20.0	609 / 24.0	482 / 19.0
B (mm / inch)	1110 / 43.7	1246 / 49.1	1246 / 49.1	1465 / 57.7	1495 / 58.9	1500 / 59.1
E (mm / inch)	200 / 7.9	200 / 7.9	200 / 7.9	270 / 10.6	270 / 10.6	270 / 10.6
F (mm / inch)	330 / 13.0	330 / 13.0	330 / 13.0	330 / 13.0	330 / 13.0	330 / 13.0

Temperature & Humidity Chamber

General type (-5 to 100°C)

Stable control for long-term temperature and humidity tests



TH-TG-180

with Wire Shelves 2ea, Φ 50 Cable Port (standard), Recorder (option)



TH-TG-1500

with Wire Shelves 2ea, Φ 50 Cable Port (standard), Recorder, Signal Lamp (option)

Optimized For Long-term Testing

- Superior defrosting control technology ensures stable temperature and humidity control for long-term tests.
- Recommended in case no low temperature testing.
- 6 models up to 1500L depending on capacity.
- Satisfies ICH Guideline Q1A (R2) Stability Test.
- Recorder for maintaining test records. (option)
- Options available for convenient internal observation during long-term testing and access to samples under testing. (inner glass door and hand-in port)

Structural Functional Features

- Optimum temperature and humidity test can be carried out even with minimum refrigeration performance, it makes reduce the electric energy efficiently.
- Excellent thermal insulation and sealing to ensure stable temperature and humidity reproducibility even for repeated test use.
- Proven reproducibility and reliability with structural design optimized for each capacity.
- The door part is a double-sealed structure that provides enhanced airtightness.

Use Convenience Features

- Convenient operation with large color touch display controller.
- Efficient water supply and recirculation system minimizes water refilling even for long-term high-humidity tests.
- Drawer-type water tank structure makes refill the water easily at the front part.
- Possible structure to refill the water even during operation.
- Convenient checking the water level at the front of the chamber.
- Direct water supply system. (option)
- Easy observation of samples with tempered glass door with built in heating wire and LED lamp.
- Open/close door with latches on either side.
- Door handle with built-in key lock.
- Cable port (Ø50 mm) is included as standard, making it convenient to connect external equipment.
- Easy to clean as it is easy to remove refrigerator condenser grill, making it convenient to maintain efficiency of refrigeration.



Drawer-type water tank



Door handle and locking device

Multifunction Control System

- Save program information and back up the set values.
- Includes temperature and humidity auto-tuning.
- Highly reliable control through temperature calibration.
- PID zone subdivided into 9 zones for more precise control stability.
- RS-232 port (default) and RS-485 port (option) supported.
- Control and data processing by connecting up to 32 devices to PC at same time via RS-485 port.
- Convenient checking with useful graph display.



Cable port



Tempered glass window with wire heater and LED lamp

Description	No.
Pattern	100
Repeat time	999
Part repeat time	255
Max. segments / pattern	100
Available max. segments	2000
Programmable process time / segment	99 hr. 59 min.

* It is possible to set 100 segments per pattern, but the maximum number of segments is not 10,000 (100 pattern x 100 segment) but 2,000.






Controller



Outstanding Safety

- Electronic over temperature shutdown system (advanced over temperature (A-OT) limiter)
Independent precision temperature sensor and controller for safer and more accurate over temperature shutdown.
- General-type over power shutdown system. (B-OT, backup over temperature limiter)
A backup device for electronic system failure, mechanically preventing over temperature even when it occurs electronic errors.
- Ensures that the surface temperature is kept below 60°C.
- Low and empty water level alarm and automatic shut off.
- Door open warning and automatic shut off.
- Overcurrent and short circuit protection of device.
- Fuse installed for each heater for safety.
- Auto stop when operating current of compressor is overloaded.
- Auto stop when compressor is overheating.
- Auto shut off in case refrigerant pressure (high/low pressure) is abnormal.



Specification


Model	TH-TG-180	TH-TG-300	TH-TG-408	TH-TG-800	TH-TG-1000	TH-TG-1500
Interior dimensions						
Chamber volume (L / cu ft)	180 / 6.4	314 / 11.1	419 / 14.8	770 / 27.2	997 / 35.2	1540 / 54.4
Width (mm / inch)	500 / 19.7	748 / 29.4	748 / 29.4	748 / 29.4	748 / 29.4	1100 / 43.3
Depth (mm / inch)	600 / 23.6	600 / 23.6	800 / 31.5	920 / 36.2	1190 / 46.9	1000 / 39.4
Height (mm / inch)	600 / 23.6	700 / 27.6	700 / 27.6	1120 / 44.1	1120 / 44.1	1400 / 55.1
Quantity of shelves (standard/max.)	2 / 7	2 / 8	2 / 8	2 / 17	2 / 17	2 / 17
Distance of between shelves (mm / inch)	60 / 2.4	60 / 2.4	60 / 2.4	60 / 2.4	60 / 2.4	60 / 2.4
Max. load per shelf (kg / lbs)	20 / 44.1	25 / 55.1	25 / 55.1	40 / 88.2	40 / 88.2	40 / 88.2
Permitted total load (kg / lbs)	50 / 110.2	65 / 143.3	65 / 143.3	100 / 220.5	100 / 220.5	100 / 220.5
Cable port (Ø, mm / inch)	50 / 2.0	50 / 2.0	50 / 2.0	50 / 2.0	50 / 2.0	50 / 2.0
Water tank capacity (L / cu ft)	15 / 0.5	15 / 0.5	15 / 0.5	15 / 0.5	15 / 0.5	15 / 0.5
Weight (kg / lbs)	220 / 485.0	250 / 551.2	280 / 617.3	360 / 793.7	410 / 903.9	480 / 1058.2
Temperature data						
Range without humidity (°C / °F)	-5 to 100 /23 to 212					
Fluctuation (±°C / °F) ¹⁾	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54
Variation (±°C / °F) ¹⁾	0.5 / 0.9	0.5 / 0.9	0.5 / 0.9	0.5 / 0.9	0.5 / 0.9	0.5 / 0.9
Heating time (-5°C to 100°C) (min)	50	50	50	55	55	55
Cooling time (100°C to -5°C) (min)	120	120	120	145	145	145
Humidity data (※ Ref.: Temp./ Humid. control range graph)						
Range (%RH)	10 to 95	10 to 95	10 to 95	25 to 95	25 to 95	25 to 95
Fluctuation (±%RH) ¹⁾	2	2	2	2	2	2
Variation (±%RH) ¹⁾	3	3	3	3	3	3
Electrical data & Ordering information						
Electrical requirements (230V, 1P, 50/60Hz, A)	11.8	13.9	13.9	14.9	14.9	24.0
Cat. No.	AAHC2002K	AAHC2012K	AAHC2022K	AAHC2032K	AAHC2042K	AAHC2052K
Cat. No. (EU only)	AAHC2005K	AAHC2015K	AAHC2025K	AAHC2035K	AAHC2045K	AAHC2055K

※ The above specifications are tested according to IEC 60068 standard.

※ Permissible environmental conditions are 18 - 30°C (recommended 20°C), 85%RH or less, under 2000m above sea level.

※ The DI water condition is pH 6.2 - 7.2 and the electrical conductivity is less than 20µs/cm.

1) The set value is set under the conditions of 40°C, 60%RH.



Dimension

Model	TH-TG-180	TH-TG-300	TH-TG-408	TH-TG-800	TH-TG-1000	TH-TG-1500
W (mm / inch)	952 / 37.5	1201 / 47.3	1201 / 47.3	1201 / 47.3	1201 / 47.3	1554 / 61.2
D (mm / inch)	907 / 35.7	907 / 35.7	1107 / 43.6	1228 / 48.3	1498 / 59.0	1872 / 73.7
H (mm / inch)	1503 / 59.2	1525 / 60.0	1525 / 60.0	1993 / 78.5	1993 / 78.5	1935 / 76.2
A (mm / inch)	321 / 12.6	366 / 14.4	366 / 14.4	385 / 15.2	520 / 20.5	447 / 17.6
B (mm / inch)	1004 / 39.5	1016 / 40.0	1016 / 40.0	1335 / 52.6	1335 / 52.6	1384 / 54.5

Accessories Page 308 Cable Port, Shelves, Signal Lamp, Recorder, Water Purifier, Direct Water System, Inner Glass Door

Temperature & Humidity Chamber

Dual type

Space efficiency model that controls two upper and lower chambers individually



TH3-DG-150

with Wire Shelves 2ea Φ 50 Cable Port (standard), Recorder (option)

Optimized For Long-term Testing

- Independent control of the temperature and humidity with dual chambers divided into upper and lower.
- Suitable for use when various test are required.
- Superior defrosting control technology ensures stable temperature and humidity control for long-term tests.
- Recommended in case no low temperature testing.
- Provide 3 kinds of model depending on capacity.
- Satisfies ICH Guideline Q1A (R2) Stability Test.
- Recorder for maintaining test records. (option)

Structural Functional Features

- Optimum temperature and humidity test can be carried out even with minimum refrigeration performance, it makes reduce the electric energy efficiently.
- Excellent thermal insulation and sealing to ensure stable temperature and humidity reproducibility even for repeated test use.
- Proven reproducibility and reliability with structural design optimized for each capacity.
- The door part is a double-sealed structure that provides enhanced airtightness.

Use Convenience Features

- Integrated with 21 CFR Part 11 compliant LC DataKeeper, suitable for GMP audits (option)
- Convenient operation with 5.6 inch color touch display controller.
- Possible to check the operation status with graphs.
- Water level check function enables water replenishment at the right time.
- Efficient water supply and recirculation system minimizes water refilling even for long-term high-humidity tests.
- Easy observation of samples with tempered glass door with built in heating wire and LED lamp.
- By displaying the operation status of the device with LED colors, convenient to identify chambers status from a far.
- Easy to clean as it is easy to remove refrigerator condenser grill, making it convenient to maintain efficiency of refrigeration.

Multifunction Control System

- Storage and documentation of programs.
- Save program information and back up the set values.
- RS-232 port (default) and RS-485 port (option) supported.
- Control and data processing by connecting up to 32 devices to PC at same time via RS-485 port.
- Reliable control with temperature/humidity auto-tuning and temperature calibration.
- PID zone subdivided into 9 zones for more precise control stability.
- Monitor and control via mobile app anytime, anywhere with LC Connected. (Requires LC GreenBox, option)
- Convenient checking with useful graph display.
- Save data as excel and text files.

Description	No.
Pattern	100
Repeat time	999
Part repeat time	255
Max. segments / pattern	100
Available max. segments	2000
Programmable process time / segment	99 hr. 59 min.

* It is possible to set 100 segments per pattern, but the maximum number of segments is not 10,000 (100 pattern x 100 segment) but 2,000.

Outstanding Safety


- Electronic over temperature shutdown system. (advanced over temperature (A-OT) limiter)
Independent precision temperature sensor and controller for safer and more accurate over temperature shutdown.
- General over power shutdown system. (B-OT, backup over temperature limiter)
A backup device for electronic system failure, mechanically preventing over temperature even when it occurs electronic errors.
- Ensures that the surface temperature is kept below 60°C.
- Low and empty water level alarm and automatic shut off.
- Door open warning and automatic shut off.
- Overcurrent and short circuit protection of device.
- Fuse installed for each heater for safety.
- Auto stop when operating current of compressor is overloaded.
- Auto stop when compressor is overheating.
- Auto shut off in case refrigerant pressure (high/low pressure) is abnormal.





Door handle & locking device



Cable Port



Controller





Specification

Model	TH3-DG-150	TH3-DG-300	TH3-DG-400
Interior dimensions			
Chamber volume (L / cu ft)	150 x 2 / 5.3 x 2	288 x 2 / 10.2 x 2	384 x 2 / 13.6 x 2
Width (mm / inch) ²⁾	600 / 23.6	800 / 31.5	800 / 31.5
Depth (mm / inch) ²⁾	500 / 19.7	600 / 23.6	800 / 31.5
Height (mm / inch) ²⁾	500 / 19.7	600 / 23.6	600 / 23.6
Quantity of shelves (standard/max.) ²⁾	2 / 7	2 / 8	2 / 8
Distance of between shelves (mm / inch)	60 / 2.4	60 / 2.4	60 / 2.4
Max. load per shelf (kg / lbs)	25 / 55.1	25 / 55.1	25 / 55.1
Permitted total load (kg / lbs) ²⁾	75 / 165.3	75 / 165.3	75 / 165.3
Cable port (∅, mm / inch) ²⁾	50 / 2.0	50 / 2.0	50 / 2.0
Water tank capacity (L / cu ft)	20 / 0.7	20 / 0.7	20 / 0.7
Weight (kg / lbs)	420 / 926	516 / 1137.6	560 / 1234.6
Temperature data			
Range without humidity (°C / °F)	-5 to 100 / 23 to 212	-5 to 100 / 23 to 212	-5 to 100 / 23 to 212
Fluctuation (±°C / °F) ¹⁾	0.5 / 0.9	0.5 / 0.9	0.5 / 0.9
Gradient (±°C / °F) ¹⁾	< 2.0 / 3.6	< 2.0 / 3.6	< 2.0 / 3.6
Heating time (-5°C to 95°C) (min)	130	120	130
Cooling time (95°C to -5°C) (min)	140	130	140
Humidity data (※ Ref.: Temp./ Humid. control range graph)			
Range (%RH)	20 to 90	20 to 90	20 to 90
Fluctuation (±%RH) ¹⁾	1.5 (2.0 at the point of 24hr)	1.5 (2.0 at the point of 24hr)	1.5 (2.0 at the point of 24hr)
Variation (±%RH) ¹⁾	3	3	3
Electrical data & Ordering information			
Electrical requirements (230V, 1P, 50/60Hz, A)	13.0 x 2ea	16.8 x 2ea	16.8 x 2ea
Cat. No.	AAHC6302K	AAHC6312K	AAHC6322K

※ The above specifications are tested according to IEC 60068 standard.

※ Permissible environmental conditions are 18 - 25°C (recommended 20°C), 85%RH or less, under 2000m above sea level.


※ The DI water condition is pH 6.2 - 7.2 and the electrical conductivity is less than 20µs/cm.

1) The set value is set under the conditions of 40°C, 75%RH.

2) Each chamber.

Dimension

Model	TH3-DG-150	TH3-DG-300	TH3-DG-400
W (mm / inch)	1055 / 41.5	1255 / 49.4	1255 / 49.4
D (mm / inch)	1320 / 52.0	1470 / 57.9	1670 / 65.9
H (mm / inch)	1725 / 67.9	1925 / 75.8	1925 / 75.8
A (mm / inch)	710 / 28.0	760 / 29.9	760 / 29.9
B (mm / inch)	1450 / 57.1	1600 / 63.0	1600 / 63.0
C (mm / inch)	345 / 13.6	385 / 15.2	485 / 19.1
E (mm / inch)	330 / 13.0	330 / 13.0	330 / 13.0
F (mm / inch)	200 / 7.9	200 / 7.9	200 / 7.9



Humidity & Temperature Chamber

Accessories Page 308

Cable Port, Shelves, Signal Lamp, Recorder, Water Purifier, Direct Water System, Inner Glass Door, LC GreenBox, LC DataKeeper

Temperature & Humidity Chamber

Basic type

Applied the new developed controller specialized for Temp. & Humid chambers

▣ Differentiated Systems

- Convenient operation with 5.6 inch color touch display controller.
- Possible to check the operation status with graphs.
- Minimize the occurrence of faults by simplifying water direct line and electric circuit.
- Water level check function enables water replenishment at the right time.

▣ Structural Functional Features



- Excellent thermal insulation and sealing to ensure stable temperature and humidity reproducibility even for repeated test use.
- The door part is a double-sealed structure that provides enhanced airtightness.
- Maintains stable temperature and humidity by circulating air in the chamber at the optimum speed.
- Energy-saving system that can perform optimum temperature and humidity tests efficiently, reducing power consumption during operation.

▣ Outstanding Safety

- Electronic over temperature shutdown system.
(advanced over temperature (A-OT) limiter)
Independent precision temperature sensor and controller for safer and more accurate over temperature shutdown.
- General over power shutdown system.
(B-OT, backup over temperature limiter)
A backup device for electronic system failure, mechanically preventing over temperature even when it occurs electronic errors.
- Ensures that the surface temperature is kept below 60°C.
- Low and empty water level alarm and automatic shut off.
- Door open warning and automatic shut off.
- Overcurrent and short circuit protection of device.
- Fuse installed for each heater for safety.
- Auto stop when operating current of compressor is overloaded.
- Auto stop when compressor is overheating.
- Auto shut off in case refrigerant pressure (high/low pressure) is abnormal.



TH3-E-200
with Wire Shelves 2ea (standard)





Use Convenience Features

- Integrated with 21 CFR Part 11 compliant LC DataKeeper, suitable for GMP audits (option)
- Efficient water supply and recirculation system minimizes water refilling even for long-term high-humidity tests.
- Direct water supply system. (option)
- Easy observation of samples with tempered glass door with built in heating wire and LED lamp.
- A structure that allows the door to be opened and closed easily.
- Easy to clean as it is easy to remove refrigerator condenser grill, making it convenient to maintain efficiency of refrigeration.

Multifunction Control System

- Reliable control with temperature/humidity auto-tuning and temperature calibration.
- Control and data processing by connecting up to 32 devices to PC at same time via RS-485 port.
- PID zone subdivided into 9 zones for more precise control stability.
- Monitor via mobile app anytime, anywhere with LC Connected. (Requires LC GreenBox, option)
- Convenient checking with useful graph display.

Description	No.
Pattern	100
Repeat time	999
Part repeat time	255
Max. segments / pattern	100
Available max. segments	2000
Programmable process time / segment	99 hr. 59 min.

* It is possible to set 100 segments per pattern, but the maximum number of segments is not 10,000 (100 pattern x 100 segment) but 2,000.



Large LCD touch screen controller (5.6 inches)



Tempered glass window with wire heater and LED lamp



Remote monitoring and control via LC Connected app

Specification

Model	TH3-E-200	TH3-E-400
Interior dimensions		
Chamber volume (L / cu ft)	201 / 7.1	398 / 14.1
Width (mm / inch)	500 / 19.7	600 / 23.6
Depth (mm / inch)	670 / 26.4	830 / 32.7
Height (mm / inch)	600 / 23.6	800 / 31.5
Quantity of shelves (standard/max.)	2 / 9	2 / 10
Distance of between shelves (mm)	40 / 1.6	40 / 1.6
Max. load per shelf (kg / lbs)	20 / 44.1	25 / 55.1
Permitted total load (kg / lbs)	50 / 110.2	65 / 143.3
Cable port (∅, mm / inch)	50 / 2.0	50 / 2.0
Water tank capacity (L / cu ft)	20 / 0.7	20 / 0.7
Weight (kg / lbs)	200 / 440.9	270 / 595.2
Temperature data		
Range without humidity (°C / °F)	0 to 90 / 32 to 194	0 to 90 / 32 to 194
Fluctuation (±°C / °F) ¹⁾	0.5 / 0.9	0.5 / 0.9
Variation (±°C / °F) ¹⁾	0.8 / 1.44	0.8 / 1.44
Heating time (0°C to 90°C) (min)	80	80
Cooling time (90°C to 0°C) (min)	90	90
Humidity data (※ Ref.: Temp./ Humid. control range graph)		
Range (%RH)	35 to 85	35 to 85
Fluctuation (±%RH) ¹⁾	1	1
Variation (±%RH) ¹⁾	3	3
Electrical data & Ordering information		
Electrical requirements (230V, 1P, 50/60Hz, A)	15	16
Cat. No.	AAH841102K	AAH841112K

※ The above specifications are tested according to IEC 60068 standard.


※ Permissible environmental conditions are 18 - 30°C (recommended 20°C), 85%RH or less, under 2000m above sea level.

※ The DI water condition is pH 6.2 - 7.2 and the electrical conductivity is less than 20µs/cm.

1) The set value is set under the conditions of 40°C, 75%RH.

Dimension

Model	TH3-E-200	TH3-E-400
W (mm / inch)	655 / 25.8	755 / 29.7
H (mm / inch)	1495 / 58.9	1695 / 66.7
D (mm / inch)	1162 / 45.7	1322 / 52.0
A (mm / inch)	375 / 14.8	455 / 17.9
B (mm / inch)	985 / 38.8	1135 / 44.7



Accessories Page 308

Cable Port, Shelves, Signal Lamp, Recorder, Water Purifier, Direct Water System, Inner Glass Door, LC GreenBox, LC DataKeeper

Pharmaceutical Stability Test Chamber

Satisfies ICH Guideline Q1A (R2) & Q1B(Option 2) Test





TH-ICH-800

with Wire Shelves 2ea,
Φ50 Cable Port (standard),
Recorder (option)

TH-ICH-300

with Wire Shelves 2ea,
Φ50 Cable Port (standard),
Recorder (option)



▣ Stability Test (ICH Guideline Q1A (R2))

- Ease of data analysis with useful graphical shows.
- Superior defrosting control technology ensures stable temperature and humidity control for long-term tests.
- Drawer-type water tank structure makes refill the water easily at the front part.
- Possible structure to refill the water even during operation.
- Direct water supply system. (option)
- Recorder for maintaining test records. (option)
- Inner glass door allows for convenient internal observation during long-term testing and sample access during testing. (option)

▣ Photostability Test

- Built-in cool white fluorescent lamp and UVA lamp.
- Visible light sensor and UV sensor included as standard.
- Dedicated controller for light cumulative dose setting/control.
- Convenient management by displaying the usage time for each lamp.

▣ Structural Functional Features

- Optimum temperature and humidity test can be carried out even with minimum refrigeration performance, it makes reduce the electric energy efficiently.
- Excellent thermal insulation and sealing to ensure stable temperature and humidity reproducibility even for repeated test use.
- Proven reproducibility and reliability with structural design optimized for each capacity.
- The door part is a double-sealed structure that provides enhanced airtightness.

Use Convenience Features

- Convenient operation with 3.7 inch color touch display controller.
- A structure that allows the door to be opened and closed easily.
- Door handle with built-in key lock.
- Cable port (Ø50 mm) is included as standard, making it convenient to connect external equipment.
- Easy to clean as it is easy to remove refrigerator condenser grill, making it convenient to maintain efficiency of refrigeration.



Built-in cool white fluorescent lamp and UVA lamp

Structural Functional Features

- Fix control and program control are selectable.
- Process time displayable.
- Includes temperature and humidity auto-tuning.
- Highly reliable control through temperature calibration.
- PID zone subdivided into 9 zones for more precise control stability.
- RS-232 port (default) and RS-485 port (option) supported.
- Control and data processing by connecting up to 32 devices to PC at same time via RS-485 port.



Drawer type of water tank

Outstanding Safety

- Electronic over temperature shutdown system.
(advanced over temperature (A-OT) limiter)
Independent precision temperature sensor and controller for safer and more accurate over temperature shutdown.
- General over power shutdown system.
(B-OT, backup over temperature limiter)
A backup device for electronic system failure, mechanically preventing over temperature even when it occurs electronic errors.
- Ensures that the surface temperature is kept below 60°C.
- Low and empty water level alarm and automatic shut off.
- Door open warning and automatic shut off.
- Overcurrent and short circuit protection of device.
- Fuse installed for each heater for safety.
- Auto stop when operating current of compressor is overloaded.
- Auto stop when compressor is overheating.
- Auto shut off in case refrigerant pressure (high/low pressure) is abnormal.



UV / VIS sensor



Lamp controller

Description	No.
Pattern	100
Repeat time	999
Part repeat time	255
Max. segments / pattern	100
Available max. segments	2000
Programmable process time / segment	99 hr. 59 min.

* It is possible to set 100 segments per pattern, but the maximum number of segments is not 10,000 (100 pattern x 100 segment) but 2,000.



Controller



Specification

Model	TH-ICH-300	TH-ICH-800
Interior dimensions		
Chamber volume (L / cu ft)	300 / 10.6	760 / 26.8
Width (mm / inch)	750 / 29.5	750 / 29.5
Depth (mm / inch)	650 / 25.6	880 / 34.6
Height (mm / inch)	650 / 25.6	1160 / 45.7
Quantity of shelves (standard/max.)	2 / 8	2 / 17
Distance of between shelves (mm / inch)	60 / 2.4	60 / 2.4
Max. load per shelf (kg / lbs)	25 / 55.1	40 / 88.2
Permitted total load (kg / lbs)	75 / 165.3	100 / 220.5
Cable port (Ø, mm / inch)	50 / 2.0	50 / 2.0
Water tank capacity (L / cu ft)	15 / 0.5	15 / 0.5
Weight (kg / lbs)	280 / 617.3	400 / 881.8
Temperature data		
Range without humidity (°C / °F)	-5 to 85 / 23 to 185	-5 to 85 / 23 to 185
Fluctuation (±°C / °F) ¹⁾	0.3 / 0.54	0.3 / 0.54
Variation (±°C / °F) ¹⁾	0.5 / 0.9	0.7 / 1.26
Heating time (20°C to 85°C) (min)	25	30
Cooling time (20°C to -5°C) (min)	35	40
Humidity data (※ Ref.: Temp./ Humid. control range graph)		
Range (%RH)	20 to 85	20 to 85
Fluctuation (±%RH) ¹⁾	3	3
Variation (±%RH) ¹⁾	5	5
Light data (ICH Q1B option 2)		
Cool white (Lux)	6000	6500
Ultraviolet-A (W/m ²)	4.5	5
Uniformity of Visible light (±%)	10	15
Uniformity of Ultraviolet-A light (±%)	15	15
Typical time 1.2 million Lux-hr (hr)	200	184
Typical time 200 W-hr/m ² (hr)	50	40
Electrical data & Ordering information		
Electrical requirements (230V, 1P, 50/60Hz, A)	14.7	23.5
Cat. No.	AAHC1002K	AAHC1012K
Cat. No. (EU only)	AAHC1005K	AAHC1015K

※ The above specifications are tested according to IEC 60068 standard.


※ Permissible environmental conditions are 18 - 30°C (recommended 20°C), 85%RH or less, under 2000m above sea level.

※ The DI water condition is pH 6.2 - 7.2 and the electrical conductivity is less than 20µs/cm.

1) The set value is set under the conditions of 40°C, 60%RH.

Dimension

Model	TH-ICH-300	TH-ICH-800
W (mm / inch)	1130 / 44.5	1220 / 48.0
D (mm / inch)	985 / 38.8	1215 / 47.8
H (mm / inch)	1520 / 59.8	1955 / 77.0
A (mm / inch)	460 / 18.1	545 / 21.5
B (mm / inch)	905 / 35.6	1185 / 46.7



Accessories

Wire Shelf



- Excellent ventilation structure.
- Stainless steel.
- Includes 2ea as standard.

Perforated Shelf



- Suitable for high-load samples.
- Excellent maintenance as made of stainless steel.

Recorder (Paper)



- Up to 6 channels supported.
- Set recording range and speed.
- Paper width: 100 mm.

Recorder (Digital)



- Paperless. (5.7" LCD display)
- Up to 6 channels supported.
- Data stored in internal/external memory.

Cable Port



- Both cap and silicone plug are included as standard in Ø50 mm cable port.
- It is available to add Ø50 mm, Ø80 mm.

Signal Lamp



- Figuring out equipment status remotely.
- Run/Stand-by/Error displayed by color.

Water Purifier



- Water purifier composed of particle filter and carbon filter.
- Solid steel box construction.

Inner Glass Door



- Minimized temperature and humidity changes and observation of internal samples enabled.
- Hand-in port is provided as standard. (selectable when purchasing)

Stand (1 Stage)



- Dedicated stand for small appliances.
- Provides lower drawer and loading space.
- Casters provided as standard for moving/fixing.

Stand (2 Stage)



- Placing of two compact devices vertically.
- Minimized installation space for different temperature and humidity tests.

Direct Water System



- Direct connection of distilled water to direct water line of the equipment from a water supply pipe.
- Useful for long-term testing.

Model	Cable Port (ø50)	Cable Port (ø80)	Wire Shelf	Perforated Shelf	Signal Lamp	Recorder (Paper)	Recorder (Digital)	Water Purifier	Direct Water System	Inner Glass Door	Stand (1 stage)	Stand (2 stage)
TH3-KE-025	AAA8T612	N/A	00RTD0000679	AAA80602-6	AAA80550	AAA8T504	AAA8T506	AAA8T625	AAA80620	N/A	AAA80631	AAA80634
TH3-KE-065	AAA8T612	AAA8T611	00RTD0000680	AAA80602-7	AAA80550	AAA8T504	AAA8T506	AAA8T625	AAA80620	N/A	AAA80632	AAA80635
TH3-KE-100	AAA8T612	AAA8T611	00RTD0000681	AAA80602-8	AAA80550	AAA8T504	AAA8T506	AAA8T625	AAA80620	N/A	AAA80633	N/A
TH3-ME-025	AAA8T612	N/A	00RTD0000679	AAA80602-6	AAA80550	AAA8T504	AAA8T506	AAA8T625	AAA80620	N/A	AAA80631	AAA80634
TH3-ME-065	AAA8T612	AAA8T611	00RTD0000680	AAA80602-7	AAA80550	AAA8T504	AAA8T506	AAA8T625	AAA80620	N/A	AAA80632	AAA80635
TH3-ME-100	AAA8T612	AAA8T611	00RTD0000681	AAA80602-8	AAA80550	AAA8T504	AAA8T506	AAA8T625	AAA80620	N/A	AAA80633	N/A
TH3-PE-025	AAA8T612	N/A	00RTD0000679	AAA80602-6	AAA80550	AAA8T504	AAA8T506	AAA8T625	AAA80620	N/A	AAA80631	AAA80634
TH3-PE-065	AAA8T612	AAA8T611	00RTD0000680	AAA80602-7	AAA80550	AAA8T504	AAA8T506	AAA8T625	AAA80620	N/A	AAA80632	AAA80635
TH3-PE-100	AAA8T612	AAA8T611	00RTD0000681	AAA80602-8	AAA80550	AAA8T504	AAA8T506	AAA8T625	AAA80620	N/A	AAA80633	N/A
TH-KH-180	AAA8T610	AAA8T6111	00EDA0008212	AAA80602-1	AAA80550	AAA8T500	AAA8T505	AAA8T625	AAA80620	AAA81552	N/A	N/A
TH-KH-300	AAA8T610	AAA8T6111	00EDA0008213	AAA80602-2	AAA80550	AAA8T500	AAA8T505	AAA8T625	AAA80620	AAA81551	N/A	N/A
TH-KH-408	AAA8T610	AAA8T6111	00EDA0008214	AAA80602-3	AAA80550	AAA8T500	AAA8T505	AAA8T625	AAA80620	AAA81551	N/A	N/A
TH-G-180	AAA8T610	AAA8T611	00EDA0008212	AAA80602-1	AAA80550	AAA8T500	AAA8T505	AAA8T625	AAA80620	AAA81552	N/A	N/A
TH-G-300	AAA8T610	AAA8T611	00EDA0008213	AAA80602-2	AAA80550	AAA8T500	AAA8T505	AAA8T625	AAA80620	AAA81551	N/A	N/A
TH-G-408	AAA8T610	AAA8T611	00EDA0008214	AAA80602-3	AAA80550	AAA8T500	AAA8T505	AAA8T625	AAA80620	AAA81551	N/A	N/A
TH-G-800	AAA8T610	AAA8T611	00EDA0008215	AAA80602-4	AAA80550	AAA8T500	AAA8T505	AAA8T625	AAA80620	AAA81550	N/A	N/A
TH-G-1000	AAA8T610	AAA8T611	00EDA0008216	AAA80602-5	AAA80550	AAA8T500	AAA8T505	AAA8T625	AAA80620	AAA81550	N/A	N/A
TH-G-1500	AAA8T610	AAA8T611	00EDA0008216	AAA80602-5	AAA80550	AAA8T500	AAA8T505	AAA8T625	AAA80620	N/A	N/A	N/A
TH-TG-180	AAA8T610	AAA8T611	00RTD0001302	AAA80604-11	AAA80550	AAA8T500	AAA8T505	AAA8T625	AAA80620	AAAC2501	N/A	N/A
TH-TG-300	AAA8T610	AAA8T611	00RTD0001303	AAA80604-12	AAA80550	AAA8T500	AAA8T505	AAA8T625	AAA80620	AAAC2502	N/A	N/A
TH-TG-408	AAA8T610	AAA8T611	00RTD0001300	AAA80604-13	AAA80550	AAA8T500	AAA8T505	AAA8T625	AAA80620	AAAC2502	N/A	N/A
TH-TG-800	AAA8T610	AAA8T611	00RTD0001304	AAA80604-14	AAA80550	AAA8T500	AAA8T505	AAA8T625	AAA80620	AAAC2503	N/A	N/A
TH-TG-1000	AAA8T610	AAA8T611	00RTD0001301	AAA80604-15	AAA80550	AAA8T500	AAA8T505	AAA8T625	AAA80620	AAAC2503	N/A	N/A
TH-TG-1500	AAA8T610	AAA8T611	00RTD0001305	AAA80604-16	AAA80550	AAA8T500	AAA8T505	AAA8T625	AAA80620	N/A	N/A	N/A
TH3-DG-150	AAA8T610	AAA8T611	00EDA0008251	AAAC6501	AAA80550	AAA8T500	AAA8T505	AAA8T625	AAA80620	AAAC6511	N/A	N/A
TH3-DG-300	AAA8T610	AAA8T611	00RTD0001280	AAAC6502	AAA80550	AAA8T500	AAA8T505	AAA8T625	AAA80620	AAAC6512	N/A	N/A
TH3-DG-400	AAA8T610	AAA8T611	00EDA0008250	AAAC6503	AAA80550	AAA8T500	AAA8T505	AAA8T625	AAA80620	AAAC6513	N/A	N/A
TH3-E-200	AAA8T610	AAA8T611	00EDA0008212	AAA80602-1	AAA80550	AAA8T500	AAA8T505	AAA8T625	AAA80620	AAA81552	N/A	N/A
TH3-E-400	AAA8T610	AAA8T611	00EDA0008214	AAA80602-3	AAA80550	AAA8T500	AAA8T505	AAA8T625	AAA80620	AAA81551	N/A	N/A
TH-ICH-300	AAA8T610	AAA8T611	00RTD0001196	AAA80604-2	AAA80550	AAA8T500	AAA8T505	AAA8T625	AAA80620	AAA81554	N/A	N/A
TH-ICH-800	AAA8T610	AAA8T611	00RTD0001198	AAA80604-4	AAA80550	AAA8T500	AAA8T505	AAA8T625	AAA80620	AAA81555	N/A	N/A

LC GreenBox (only TH3 models)

- Monitor and control the operation status of the product in real time via mobile app with an internet connection.
- Up to 4 devices can be connected.

Description	Dimension (W x D x H, mm / inch)	Cat. No.
LC GreenBox	156 x 94 x 34 / 6.1 x 3.7 x 1.3	AAHQ1011K



LC DataKeeper (only TH3 models)

- High-integrity software compliant with FDA 21 CFR Part 11
- Integrated with Jeio Tech equipment for GMP audits

Description	Cat. No.
LC DataKeeper-Basic	SSA11001
LC DataKeeper-Pro	SSA11002



※ Trial version is available for 30 days, offering the same features as the Basic version, and can be requested through the website.

Temperature & Humidity Chamber

Walk-in type

Large room chamber with optimized structure and system configuration



▮ Structural Functional Features

- Excellent thermal insulation and sealing to ensure stable temperature and humidity reproducibility even for repeated test use.
- The structure panel has excellent durability, minimizing deformation of the structure.
- Clean form without leakage and water gaps through application of optimal drain system design.
- Reinforced door structure ensures excellent durability and door sealing for long-term use.

▮ Use Convenience Features

- Customized design and production are optimized by thoroughly inspecting the internal capacity, temperature/temperature and humidity range, freezing type, specimen jig, convenience devices, safety functions, control system, installation place, etc. together with the user.
- Designed as a modular control system structure and can be placed the position as the customer's request.
- Basic safety features such as over temperature, door opening and water shortage, also additional safety features can be added such as for smoke and flame.

Temperature & Humidity Chamber

Customized type

Customized production based on accumulated technology and experience



Large temperature & humidity chamber

Custom-sized temp. & humid. chamber according to customer needs



Clean room temperature & Humidity chamber

Clean temp. & humid. chamber with HEPA filter

> Expert technical consultation

We provide expert technical consultations based on our accumulated skills and experience to conduct appropriate environmental reliability tests.

> Optimized proposal

Proposal of optimal temperature and temperature & humidity chamber configuration based on the customer's use.

> Reasonable customization

Customized chamber comprehensively considering test conditions, usage patterns, place of use, etc.

Expert Consultation Contact

International Sales Team

RM1005, 219, Gasan digital 1-ro, Geumcheon-gu, Seoul, 08501, Republic of Korea
Tel. +82-2627-3816
E-mail. overseas@jeiotech.com

Heating & Cooling Chamber



- | | | | | | |
|--|--|--|---|--|--|
| | <p>Satisfaction with the requirements of European regulations and guidelines.</p> | | <p>Ensuring the reliability of electronic records and signatures to support GMP audit.</p> | | <p>Real-time equipment monitoring and control system using mobile app.</p> |
| | <p>Eco-friendly product that has passed environmental impact assessment tests.</p> | | <p>More precise temperature control through temperature calibration.</p> | | <p>Electronic over temperature protection system.</p> |
| | <p>Conventional over temperature protection, a kind of backup device against electronic over temperature protection failure.</p> | | <p>Over temperature protection function.</p> | | <p>Control of temperature / humidity / rpm / time, etc. through dedicated program.</p> |
| | <p>PC communication via RS-232 / RS-485 / USB port.</p> | | <p>It memorizes the state of power failure and operates automatically when power is restored.</p> | | <p>2 year warranty free A/S.</p> |

General Application

Reliability, Climate/freezing, Quality assurance, Thermal resistance testing.

- Electrical and electronic components, Sensor.
- The machinery / Military / Aircraft equipment.
- Construction materials / Plastics / Textile industries.
- Metal industries such as metallurgy.
- Semiconductor, PCB, LCD&LED.
- Chemical / Oil industry.
- Automotive industry.
- Pharmaceutical industry.

Description		Temp. Range (°C / °F)	Airflow	Volume (L / cu ft)	Model	Page
Small type	Tabletop type <ul style="list-style-type: none"> • Compact design • Excellent insulation and sealing structure 	-40 to 100 / -40 to 212	Vertical airflow	25, 65, 100 / 0.9, 2.3, 3.5	TC3-KE	316
		-20 to 100 / -4 to 212	Vertical airflow	25, 65, 100 / 0.9, 2.3, 3.5	TC3-ME	316
Medium large size type	Ultra low temp. type (-70°C) <ul style="list-style-type: none"> • Verification completed according to international standards • Convenient safety system configuration 	-70 to 180 / -94 to 356	Horizontal airflow	125, 253, 420, 720, 1000 / 4.4, 8.9, 14.8, 25.4, 35.3	KBD	320
	Low temp. and horizontal type (-35°C) <ul style="list-style-type: none"> • 27 points temperature verification completed • High-load samples are available 	-35 to 180 / -31 to 356	Horizontal airflow	125, 253, 420, 720, 1000 / 4.4, 8.9, 14.8, 25.4, 35.3	KMV	324
	Low temp. and vertical type (-35°C) <ul style="list-style-type: none"> • Vertical airflow model series • Provide good permeable shelves 	-35 to 180 / -31 to 356	Vertical airflow	253, 420, 720, 990 / 8.9, 14.8, 25.4, 35.0	JMV	328
	Low temp. type (-25°C) <ul style="list-style-type: none"> • Excellent durability and performance • Dual overheating temperature limiter 	-25 to 100 / -13 to 212	Horizontal airflow	125, 250, 400, 700, 1000 / 4.4, 8.8, 14.1, 24.7, 35.3	PBV	332
	Ambient temperature type <ul style="list-style-type: none"> • Stable temperature control • Offering a variety of custom options 	-5 to 100 / 23 to 212	Horizontal airflow	125, 250, 400, 700, 1000 / 4.4, 8.8, 14.1, 24.7, 35.3	PMV	336
General type	Compact design <ul style="list-style-type: none"> • Space-saving vertical structure • Includes Dual model (2-in-1) 	-20 to 100 / -4 to 212	Horizontal airflow	150, 255, 485, 150 x 2 chamber / 5.3, 9.0, 17.1, 5.3 x 2 chamber	LCH	340
		0 to 100 / 32 to 212	Horizontal airflow	150, 255, 485, 150 x 2 chamber / 5.3, 9.0, 17.1, 5.3 x 2 chamber	LCH-G	340

※ The contents of the above and the contents of this catalog may differ depending on the specific model and conditions of use. For the information about the features and specifications that applying to each models, please check the information on the corresponding page of each models.

Heating & Cooling Chamber



Heating & Cooling Chamber

Proven reliability and reproducibility Continuous innovative model

› **Completion of verification according to strict international standards**

Increase in reliability of accurate temperature implementation and test results (27 point measurement data) according to DIN 12880 and IEC 60068.

› **Complete chamber structure**

Includes double-packing door.
Stable temperature control with excellent insulation and sealing.

› **Steady and uniform airflow formation**

Air flow optimized for uniform heat transfer ensures stable overall temperature uniformity.

› **Height-adjustable and strong shelf**

Durable wire-shaped shelves provided as standard can be placed at different height.

Powerful & Useful control Optimized precision control

› **Dedicated controller to operate temperature.**

Temperature auto-tuning function.
Highly-reliable control through calibration.

› **Powerful program operation**

Provides diverse program operation environment with a maximum of 100 patterns and 6000 segments.

› **Sophisticated custom control**

Auto-tuning by subdivided PID zone for more precise control stability.

› **Multi-functional dedicated software**

Up to 32 devices can be controlled at same time.
Device operation and data processing function supported.



Easy to understand running status with LED display showing the operation condition.



Electronic over temperature protection system.



Controller for temperature. (TC3-model)




Maximized usability with simple and useful functions

- > **Ease of check the test status**
Convenient internal observation with LED lamp and looking glass window.
Figures out the operation state through displaying LED even a long distance.
- > **Door lock with key lock**
Door structure having well-sealed.
Included key lock features as standard.
- > **Provision of remote equipment monitoring service (TC3, option)**
Real-time monitoring of equipment operation status using smartphone.
(when purchased LC GreenBox)
- > **Convenient moving and installation**
Easy to move and install with casters.

More convenient and safe structure applied the top-grade safety system

- > **Electronic over temperature shutdown system**
Independent precision temperature sensor and controller.
Provides more secure and reliable over temperature shutdown.
- > **General over temperature shutdown system**
Backup device for electronic system failure.
Implementation of over-temperature shutdown even for electronic errors.
- > **Self-protection system**
Protection system as like water tank shortage, door open warning, etc.
Over-current and short circuit breaking device.
- > **Protection features for refrigeration system**
Protection system for overload and over temperature on refrigerator.
Stops automatically in case refrigerant pressure (high/low pressure) is abnormal.

※ Some of the above contents are limited to specific models.



27 Points temperature verification completed.



Dedicated stand with drawer and load space. (TC3-PE, KE)



Viewing Window as standard.

Heating & Cooling Chamber

Tabletop type

Compact temperature test chamber by continuing technological innovation

- TC3-KE (-40°C)
- TC3-ME (-20°C)




TC3-KE

with 2ea wire shelves (standard)
Stand (option)



TC3-ME



Optimized Model Configuration

- 2 series (KE, ME) according to temperature control range.
TC3-KE ; -40 to 100°C
TC3-ME ; -20 to 100°C
- Standard type heating & cooling chamber for selecting the optimized model according to test conditions.
- Dedicated stand for convenient use and moving. (option)
- Recorder for maintaining test records. (option)

Structural Functional Features

- Proven reproducibility and reliability with structural design optimized for miniaturization.
- Compact design minimized foot space for excellent space saving.
- Accessible slim design to laboratories which have only 1-door entry.
- Excellent thermal insulation and sealing to ensure stable temperature reproducibility even for repeated test use.
- The door part has a double sealing structure and excellent airtightness.



Use Convenience Features

- Integrated with 21 CFR Part 11 compliant LC DataKeeper, suitable for GMP audits (option)
- Convenient operation with 5.6 inch color touch controller.
- Easy observation of samples with tempered glass door with built in heating wire and LED lamp.
- By displaying the operation status of the device with LED colors, convenient to identify chambers status from a far.
- The cable port (ø 50 mm) included as standard make it convenient to connect external equipment.
- The door structure can be easily opened or closed.
- Door handle with built-in key lock.
- Identification of the control and operating conditions at the front panel.
- Design by stainless steel internal chamber provides clean maintenance and excellent corrosion resistance.
- Easy to clean as it is easy to remove refrigerator condenser grill, making it convenient to maintain efficiency of refrigeration.
- Equipped with easy-to-move/install caster.



Viewing window



Indicator LED of operation condition.

Multifunction Control System

- Highly-reliable control through calibration.
- Auto-tuning by subdivided PID zone for more precise control stability.
- RS-232 port (default) and RS-485 port (option) supported.
- Control and data processing by connecting up to 32 devices to PC at same time via RS-485 port.
- Convenient checking with useful graph display.
- Monitor and control via mobile app anytime, anywhere with LC Connected. (Requires LC GreenBox, option)



Cable port



Door handles & Locking device

Description	No.
Pattern	100
Repeat time	999
Part repeat time	255
Max. segments / pattern	100
Available max. segments	2000
Programmable process time / segment	99 hour 59 min.

* It is possible to set 100 segments per pattern, but the maximum number of segments is not 10,000 (100 pattern x 100 segment) but 2,000.

▮ Outstanding Safety

- Electronic over temperature protection system.
(A-OT, advanced over temperature limiter)
Independent precision temperature sensor and controller for safer and more accurate over temperature shutdown.
- General over-temperature protection.
(B-OT, backup over temperature limiter)
A backup device for electronic system failure, mechanically preventing over temperature even when it occurs electronic errors.
- Door open warning and automatic shut off.
- Over-current and short circuit protection of device.
- Each heater has fuse for more being safe use.
- Automatic stop when operating current of refrigerator is overloaded.
- Automatic stop in case of over temperature of the compressor.
- Auto shut off in case refrigerant pressure (high/low pressure) is abnormal.



Specification

Model	TC3-KE-025	TC3-KE-065	TC3-KE-100	TC3-ME-025	TC3-ME-065	TC3-ME-100
Temperature data						
Range without humidity (°C / °F)	-35 to 100 / -31 to 212 (50Hz)			-20 to 100 / -4 to 212 (50/60Hz)		
	-40 to 100 / -40 to 212 (60Hz)					
Fluctuation (±°C / °F) ¹⁾	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54
Variation (±°C / °F) ²⁾	0.5 / 0.9	0.5 / 0.9	0.5 / 0.9	0.5 / 0.9	0.5 / 0.9	0.5 / 0.9
Heating time (min) ³⁾	60	60	60	50	50	50
Cooling time (min) ⁴⁾	70	70	70	55	55	55
Interior dimensions						
Chamber volume (L / cu ft)	25 / 0.88	65 / 2.3	100 / 3.53	25 / 0.88	65 / 2.3	100 / 3.53
Width (mm / inch)	350 / 13.8	400 / 15.7	500 / 19.7	350 / 13.8	400 / 15.7	500 / 19.7
Depth (mm / inch)	220 / 8.7	330 / 13	380 / 15	220 / 8.7	330 / 13	380 / 15
Height (mm / inch)	350 / 13.8	495 / 19.5	530 / 20.9	350 / 13.8	495 / 19.5	530 / 20.9
Quantity of shelves (standard/max.)	2/7	2/10	2/10	2/7	2/10	2/10
Distance of between shelves (mm / inch)	20 / 0.8	20 / 0.8	20 / 0.8	20 / 0.8	20 / 0.8	20 / 0.8
Max. load per shelf (kg / lbs)	25 / 55.1	25 / 55.1	25 / 55.1	25 / 55.1	25 / 55.1	25 / 55.1
Permitted total load (kg / lbs)	65 / 143.3	65 / 143.3	65 / 143.3	65 / 143.3	65 / 143.3	65 / 143.3
Cable port (Ø, mm / inch)	50 / 2	50 / 2	50 / 2	50 / 2	50 / 2	50 / 2
Weight (kg / lbs)	110 / 242.5	140 / 308.6	150 / 330.7	80 / 176.4	110 / 242.5	130 / 286.6
Electrical data & Ordering information						
Electrical requirements (230V, 1P, 50Hz, A)	8.5	9.5	10.6			
Cat. No.	AAH861132K	AAH861232K	AAH861332K	-	-	-
Cat. No. (EU only)	AAH861135K	AAH861235K	AAH861335K	-	-	-
Electrical requirements (230V, 1P, 60Hz, A)	8.5	9.5	10.6			
Cat. No.	AAH861131K	AAH861231K	AAH861331K	-	-	-
Electrical requirements (230V, 1P, 50/60Hz, A)				5.5	6.5	8.6
Cat. No.	-	-	-	AAH860132K	AAH860232K	AAH860332K
Cat. No. (EU only)	-	-	-	AAH860135K	AAH860235K	AAH860335K

※ The above specifications are tested according to DIN 12880 and IEC 60068 standards.


※ Permissible environmental conditions are 18 - 30°C (recommended 20°C), 85%RH or less, under 2000m above sea level.

1) The set value is set under the conditions of 100°C.

2) The set value is set under the conditions of 25°C.

3) TC3-KE : -40°C to 100°C(-40°F to 212°F), TC3-ME : -20°C to 100°C(-4°F to 212°F)

4) TC3-KE : 20°C to -40°C(68°F to -40°F), TC3-ME : 20°C to -20°C(68°F to -4°F)



Dimension

Model	TC3-KE-025	TC3-KE-065	TC3-KE-100	TC3-ME-025	TC3-ME-065	TC3-ME-100
W (mm / inch)	550 / 21.7	600 / 23.6	700 / 27.6	550 / 21.7	600 / 23.6	700 / 27.6
D (mm / inch)	880 / 34.6	1050 / 41.3	1110 / 43.7	880 / 34.6	1050 / 41.3	1110 / 43.7
H (mm / inch)	915 / 36	1060 / 41.7	1100 / 43.3	915 / 36	1060 / 41.7	1100 / 43.3
A (mm / inch)	211 / 8.3	276 / 10.9	276 / 10.9	211 / 8.3	276 / 10.9	276 / 10.9
B (mm / inch)	520 / 20.5	580 / 22.8	580 / 22.8	520 / 20.5	580 / 22.8	580 / 22.8
E (mm / inch)	170 / 6.7	170 / 6.7	170 / 6.7	170 / 6.7	170 / 6.7	170 / 6.7
F (mm / inch)	229 / 9	229 / 9	229 / 9	229 / 9	229 / 9	229 / 9

Accessories Page 342 Cable Port, Shelves, Signal Lamp, Recorder, Stand, LC Green Box, LC DataKeeper

Heating & Cooling Chamber

Ultra low temperature type (-70°C), horizontal airflow

Wide range of temperature control from -70 to 180°C




KBD-040

with 2ea wire shelves (standard), Recorder, Signal Lamp (option)


Completion of Temperature Verification According to Strict International Standard

- Temperature verification in accordance with DIN 12880 and IEC 60068-3-5 provides excellent reliability and reproducibility.
- Provided specification the fluctuation and gradient in typical temperature point.
- Provided heating / cooling time data to help select the right model for testing purposes.
- Tested recovery time in accordance with specifications to provide data that is useful for real-use.



Structural Functional Features

- Chamber structure with excellent sealing and insulation provides a stable long-term test and reduce costs by minimizing energy loss.
- According to international standards, even when operated at upper temperature, the surface temperature does not exceed 51°C. (EN 563 standard)
- Door consists of a double-sealed structure to satisfy DIN 58371. (door airtight guidance)
- Durability even after long-term repeated use has been verified by satisfying the door load test (SEFA standard).
- Door structure is a two-point hang structure, building more stable seal with less force.





Optimized Dedicated Control System

- 5.6-inch color touch display controller.
- Intuitive screen configuration for easy and convenient control.
- Graph display makes it easier to check operation.
- PID zone subdivided into 4 zones for more precise control stability.
- RS-232 port (default) and RS-485 port (option) supported.
- Control and data processing by connecting up to 32 devices to PC at same time via RS-485 port.
- Includes software for PC control.
- Saving data, convenient for reporting.

Use Convenience Features

- Convenient opening and closing door with either side latches.
- Door handle with built-in key lock.
- Identification of the control and operating conditions at the front panel.
- Design by stainless steel internal chamber provides clean maintenance and excellent corrosion resistance.
- Easy to clean as it is easy to remove refrigerator condenser grill, making it convenient to maintain efficiency of refrigeration.
- Equipped with easy-to-move/install caster.
- Perforated shelf for heavy-load sample. (option)
- Convenient and various options such as signal lamps, cable port, air cycle rate, gas purge, etc.


Description	No.
Pattern	100
Repeat time	999
Part repeat time	255
Max. segments / pattern	100
Available max. segments	2000
Programmable process time / segment	99 hr. 59 min.

* It is possible to set 100 segments per pattern, but the maximum number of segments is not 10,000 (100 pattern x 100 segment) but 2,000.


Outstanding Safety

- Electronic over temperature protection system. (A-OT, advanced over temperature limiter)
Independent precision temperature sensor and controller for safer and more accurate over temperature shutdown.
- General over-temperature protection. (B-OT, backup over temperature limiter)
A backup device for electronic system failure, mechanically preventing over temperature even when it occurs electronic errors.
- Door open warning and automatic shut off.
- Emergency stop button on the front of the product.
- Over-current and short circuit protection of device.
- When main power connections, users can quickly respond by notifying of electric phase sequence errors.
- Electrical instrumentation access warning and system shutdown.
- Each heater by fuse for more safety using.
- Auto stop when operating current of refrigerator is overloaded.
- Automatic stop in case of over temperature of the compressor.
- Auto shut off in case refrigerant pressure (high/low pressure) is abnormal.

Temperature Fluctuation



Recovery Time after Door Opening



Specification


Model	KBD-012	KBD-025	KBD-040	KBD-070	KBD-100
Temperature data					
Range (°C / °F)	-70 to 180 / -94 to 356				
Fluctuation (±°C / °F)^{1),2)}					
at -70°C	0.2 / 0.36	0.2 / 0.36	0.2 / 0.36	0.2 / 0.36	0.2 / 0.36
at -50°C	0.2 / 0.36	0.2 / 0.36	0.2 / 0.36	0.2 / 0.36	0.2 / 0.36
at -25°C	0.2 / 0.36	0.2 / 0.36	0.2 / 0.36	0.2 / 0.36	0.2 / 0.36
at -10°C	0.2 / 0.36	0.2 / 0.36	0.3 / 0.54	0.2 / 0.36	0.2 / 0.36
at 25°C	0.2 / 0.36	0.2 / 0.36	0.2 / 0.36	0.2 / 0.36	0.2 / 0.36
at 70°C	0.2 / 0.36	0.2 / 0.36	0.2 / 0.36	0.2 / 0.36	0.2 / 0.36
at 100°C	0.2 / 0.36	0.2 / 0.36	0.3 / 0.54	0.2 / 0.36	0.2 / 0.36
at 150°C	0.2 / 0.36	0.3 / 0.54	0.2 / 0.36	0.2 / 0.36	0.2 / 0.36
at 180°C	0.3 / 0.54	0.3 / 0.54	0.2 / 0.36	0.2 / 0.36	0.2 / 0.36
Gradient (±°C / °F)¹⁾					
at -70°C	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54	0.2 / 0.36	0.2 / 0.36
at -50°C	0.4 / 0.72	0.5 / 0.9	0.4 / 0.72	0.3 / 0.54	0.3 / 0.54
at -25°C	0.4 / 0.72	0.4 / 0.72	0.4 / 0.72	0.3 / 0.54	0.3 / 0.54
at -10°C	0.5 / 0.9	0.4 / 0.72	0.5 / 0.9	0.4 / 0.72	0.4 / 0.72
at 25°C	0.5 / 0.9	0.4 / 0.72	0.4 / 0.72	0.4 / 0.72	0.4 / 0.72
at 70°C	0.6 / 1.08	0.7 / 1.26	0.6 / 1.08	0.7 / 1.26	0.7 / 1.26
at 100°C	1 / 1.8	1 / 1.8	1 / 1.8	0.9 / 1.62	0.9 / 1.62
at 150°C	1 / 1.8	1.1 / 1.98	1 / 1.8	0.8 / 1.44	0.8 / 1.44
at 180°C	1.2 / 2.16	1.2 / 2.16	1.1 / 1.98	1 / 1.8	1 / 1.8
Heating time (min)²⁾					
from -70°C to 25°C	26	28	30	30	31
from -50°C to 25°C	23	23	24	25	27
from -25°C to 25°C	15	15	16	17	18
from -10°C to 25°C	11	10	11	13	13
from 25°C to 70°C	13	13	14	14	16
from 25°C to 100°C	26	25	25	26	27
from 25°C to 150°C	41	40	38	42	44
from 25°C to 180°C	54	53	55	56	58
from -50°C to 180°C	69	68	68	72	75
from -50°C to 120°C	55	56	56	61	62
from -25°C to 120°C	47	48	48	53	55
from -50°C to 5°C	18	17	18	21	21
from -70°C to 5°C	25	25	25	28	28
from 5°C to 125°C	43	42	45	47	48
from 125°C to 180°C	18	17	18	22	23
Cooling time (min)²⁾					
from 25°C to -70°C	75	75	77	74	73
from 25°C to -40°C	33	33	32	30	32
from 25°C to -25°C	22	22	20	20	21
from 25°C to -10°C	17	17	14	13	15
from 70°C to 25°C	23	20	23	20	22
from 100°C to 25°C	34	35	35	32	31
from 120°C to -50°C	83	80	81	80	80
from 5°C to -40°C	23	23	23	23	23
from 5°C to -70°C	67	72	72	70	70
Recovery time after door was opened for 30sec.(min.)²⁾					
at -70°C	12	12	13	9	10
at -50°C	8	8	8	6	6
at -25°C	6	6	6	5	5
at -10°C	5	5	6	5	5
at 70°C	6	6	6	6	6
at 100°C	5	5	5	6	6
at 150°C	7	7	10	10	10
at 180°C	7	8	10	10	10
Electrical data & Ordering information					
Electrical requirements (380V, 3P, 50/60Hz, A)	10.0	12.0	18.0	29.0	29.0
Cat. No.	AAHKC1114K	AAHKC2114K	AAHKC3114K	AAHKC4114K	AAHKC5114K
Cat. No. (EU only)	AAHKC1110K	AAHKC2110K	AAHKC3110K	AAHKC4110K	AAHKC5110K

1) According to IEC 60068-3-5

2) According to DIN 12880

※ Above specification value is recorded by 380V/50Hz.

※ Above specifications can be changed without prior notice.



Dimension

Model	KBD-012	KBD-025	KBD-040	KBD-070	KBD-100
Interior dimensions					
Chamber volume (L / cu ft)	125 / 4.41	253 / 8.93	420 / 14.83	720 / 25.43	1000 / 35.31
Width (mm / inch)	500 / 19.7	600 / 23.6	750 / 29.5	900 / 35.4	1000 / 39.4
Depth (mm / inch)	500 / 19.7	650 / 25.6	700 / 27.6	800 / 31.5	910 / 35.8
Height (mm / inch)	500 / 19.7	650 / 25.6	800 / 31.5	1000 / 39.4	1100 / 43.3
Quantity of shelves (standard/max.)	2 / 6	2 / 9	2 / 11	2 / 15	2 / 16
Distance of between shelves (mm / inch)	30 / 1.2	30 / 1.2	30 / 1.2	30 / 1.2	30 / 1.2
Max. load per shelf (kg / lbs)	15 / 33.1	20 / 44.1	30 / 66.1	40 / 88.2	45 / 99.2
Permitted total load (kg / lbs)	50 / 110.2	70 / 154.3	90 / 198.4	120 / 264.6	150 / 330.7
Cable port (Φ, mm / inch)	50 / 2.0	50 / 2.0	80 / 3.1	80 / 3.1	80 / 3.1
Exterior dimensions					
Width (mm / inch), W	1114 / 43.9	1240 / 48.8	1424 / 56.1	1575 / 62	1675 / 65.9
Depth / with handle (mm - inch), D/D'	1327 / 1360 52.2 / 53.5	1492 / 1525 58.7 / 60	1567 / 1600 61.7 / 63	1777 / 1810 70 / 71.3	1887 / 1920 74.3 / 75.6
Height (mm / inch), H	1208 / 47.6	1308 / 51.5	1509 / 59.4	1819 / 71.6	1919 / 75.6
Weight (kg / lbs)	170 / 374.8	260 / 573.2	320 / 705.5	390 / 859.8	440 / 970

Accessories [Page 342](#)

Cable Port, Shelves, Signal Lamp, Recorder, Gas Purge System

Heating & Cooling Chamber

Low temperature type (-35°C), horizontal airflow

Horizontal airflow and -35 to 180°C temperature range




KMV-040

with 2ea wire shelves (standard), Recorder, Signal lamp (option)


Completion of Temperature Verification According to Strict International Standard

- Temperature verification in accordance with DIN 12880 and IEC 60068-3-5 provides excellent reliability and reproducibility.
- Provided specification the fluctuation and gradient in typical temperature point.
- Provided heating / cooling time data to help select the right model for testing purposes.
- Tested recovery time in accordance with specifications to provide data that is useful for real-use.



Structural Functional Features

- Chamber structure with excellent sealing and insulation provides a stable long-term test and reduce costs by minimizing energy loss.
- According to international standards, even when operated at upper temperature, the surface temperature does not exceed 51°C. (EN 563 standard)
- Door consists of a double-sealed structure to satisfy DIN 58371. (door airtight guidance)
- Durability even after long-term repeated use has been verified by satisfying the door load test. (SEFA standard)
- Door structure is a two-point hang structure, building more stable seal with less force.





Optimized Dedicated Control System

- 3.5-inch color touch display controller.
- Intuitive screen configuration for easy and convenient control.
- Graph display makes it easier to check operation.
- PID zone subdivided into 4 zones for more precise control stability.
- RS-232 port (default) and RS-485 port (option) supported.
- Control and data processing by connecting up to 32 devices to PC at same time via RS-485 port.
- Includes software for PC control.
- Saving data, convenient for reporting.

Use Convenience Features

- Convenient opening and closing door with either side latches.
- Door handle with built-in key lock.
- Identification of the control and operating conditions at the front panel.
- Design by stainless steel internal chamber provides clean maintenance and excellent corrosion resistance.
- Easy to clean as it is easy to remove refrigerator condenser grill, making it convenient to maintain efficiency of refrigeration.
- Equipped with easy-to-move/install caster.
- Perforated shelf for heavy-load sample. (option)
- Convenient and various options such as signal lamps, cable port, air cycle rate, gas purge, etc.


Description	No.
Pattern	100
Repeat time	999
Part repeat time	255
Max. segments / pattern	100
Available max. segments	2000
Programmable process time / segment	99 hr. 59 min.

* It is possible to set 100 segments per pattern, but the maximum number of segments is not 10,000 (100 pattern x 100 segment) but 2,000.


Outstanding Safety

- Electronic over temperature protection system. (A-OT, advanced over temperature limiter)
Independent precision temperature sensor and controller for safer and more accurate over temperature shutdown.
- General over-temperature protection. (B-OT, backup over temperature limiter)
A backup device for electronic system failure, mechanically preventing over temperature even when it occurs electronic errors.
- Door open warning and automatic shut off.
- Emergency stop button on the front of the product.
- Over-current and short circuit protection of device.
- When main power connections, users can quickly respond by notifying of electric phase sequence errors.
- Electrical instrumentation access warning and system shutdown.
- Each heater by fuse for more safety using.
- Auto stop when operating current of refrigerator is overloaded.
- Automatic stop in case of over temperature of the compressor
- Auto shut off in case refrigerant pressure (high/low pressure) is abnormal.

Temperature Fluctuation



Recovery Time after Door Opening



Specification


Model	KMV-012	KMV-025	KMV-040	KMV-070	KMV-100
Temperature data					
Range (°C / °F)	-35 to 180 / -31 to 356				
Fluctuation (±°C / °F)^{1), 2)}					
at -40°C	0.5 / 0.9	0.5 / 0.9	0.3 / 0.54	0.2 / 0.36	0.3 / 0.54
at -25°C	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54	0.2 / 0.36	0.2 / 0.36
at -10°C	0.3 / 0.54	0.3 / 0.54	0.2 / 0.36	0.3 / 0.54	0.3 / 0.54
at 25°C	0.2 / 0.36	0.2 / 0.36	0.5 / 0.9	0.4 / 0.72	0.2 / 0.36
at 70°C	0.3 / 0.54	0.3 / 0.54	0.4 / 0.72	0.3 / 0.54	0.2 / 0.36
at 100°C	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54
at 150°C	0.2 / 0.36	0.2 / 0.36	0.3 / 0.54	0.2 / 0.36	0.2 / 0.36
at 180°C	0.4 / 0.72	0.4 / 0.72	0.3 / 0.54	0.3 / 0.54	0.2 / 0.36
Gradient (±°C / °F)¹⁾					
at -40°C	0.5 / 0.9	0.5 / 0.9	1.2 / 2.16	0.6 / 1.08	0.3 / 0.54
at -25°C	0.6 / 1.08	0.6 / 1.08	1.2 / 2.16	0.7 / 1.26	0.4 / 0.72
at -10°C	0.7 / 1.26	0.7 / 1.26	1.2 / 2.16	0.8 / 1.44	0.4 / 0.72
at 25°C	0.7 / 1.26	0.7 / 1.26	1.1 / 1.98	0.8 / 1.44	0.5 / 0.9
at 70°C	0.8 / 1.44	0.8 / 1.44	1.7 / 3.06	0.9 / 1.62	0.7 / 1.26
at 100°C	1 / 1.8	1 / 1.8	1.9 / 3.42	1.1 / 1.98	0.8 / 1.44
at 150°C	1 / 1.8	1 / 1.8	2.3 / 4.14	0.5 / 0.9	0.7 / 1.26
at 180°C	1.4 / 2.52	1.4 / 2.52	2.3 / 4.14	0.8 / 1.44	0.8 / 1.44
Heating time (min)²⁾					
from -40°C to 25°C	23	23	23	20	23
from -25°C to 25°C	20	20	19	20	18
from -10°C to 25°C	14	14	17	15	14
from 25°C to 70°C	17	17	21	19	21
from 25°C to 100°C	31	31	29	44	44
from 25°C to 150°C	61	61	48	82	90
from 25°C to 180°C	65	65	55	100	76
from -40°C to 180°C	92	92	79	102	203
from -25°C to 120°C	68	68	66	67	99
from -40°C to 5°C	15	15	13	10	19
from 5°C to 125°C	53	53	54	80	171
from 125°C to 180°C	24	24	12	12	13
Cooling time (min)²⁾					
from 25°C to -40°C	41	41	37	38	44
from 25°C to -25°C	20	20	20	17	20
from 25°C to -10°C	13	13	15	10	15
from 70°C to 25°C	16	16	16	13	15
from 100°C to 25°C	25	25	25	18	27
from 150°C to 25°C	81	81	110	103	290
from 180°C to 25°C	152	152	173	177	310
from 180°C to -40°C	180	180	221	212	361
from 120°C to -25°C	43	43	46	48	48
from 180°C to 125°C	113	113	150	153	292
from 125°C to 5°C	33	33	36	27	31
from 5°C to -40°C	34	34	35	32	38
Recovery time after door was opened for 30sec (min.)²⁾					
at -40°C	8	8	6.5	7	6
at -25°C	3.5	3.5	7.5	6	7
at -10°C	2.5	2.5	6.5	6.5	6.5
at 25°C	4	4	5.5	4.5	5
at 70°C	4	4	5.5	5	8
at 100°C	4	4	6.5	6	6
at 150°C	4	4	4.5	4.3	3.5
at 180°C	11.5	11.5	6	5	4
Electrical data & Ordering information					
Electrical requirements (230V, 50Hz, 1P, A)	19.2	26.0			
Cat. No.	AAHK9012K	AAHK9022K	-	-	-
Cat. No. (EU only)	AAHK9015K	AAHK9025K	-	-	-
Electrical requirements (230V, 60Hz, 1P, A)	19.2	26.0			
Cat. No.	AAHK9011K	AAHK9021K	-	-	-
Electrical requirements (380V, 50/60Hz, 3P, A)			12.5	18.0	19.4
Cat. No.	-	-	AAHK9038K	AAHK9048K	AAHK9058K
Cat. No. (EU only)	-	-	AAHK9030K	AAHK9040K	AAHK9050K

1) According to IEC 60068-3-5

2) According to DIN 12880

※ Above specification value is recorded by 230V/60Hz, 380V/50Hz.

※ Above specifications can be changed without prior notice.



Dimension

Model	KMV-012	KMV-025	KMV-040	KMV-070	KMV-100
Interior dimensions					
Chamber volume (L / cu ft)	125 / 4.41	253 / 8.93	420 / 14.83	720 / 25.43	1000 / 35.31
Width (mm / inch)	500 / 19.7	600 / 23.6	750 / 29.5	900 / 35.4	1000 / 39.4
Depth (mm / inch)	500 / 19.7	650 / 25.6	700 / 27.6	800 / 31.5	910 / 35.8
Height (mm / inch)	500 / 19.7	650 / 25.6	800 / 31.5	1000 / 39.4	1100 / 43.3
Quantity of shelves (standard/max.)	2 / 6	2 / 9	2 / 11	2 / 15	2 / 16
Distance of between shelves (mm / inch)	30 / 1.2	30 / 1.2	30 / 1.2	30 / 1.2	30 / 1.2
Max. load per shelf (kg / lbs)	15 / 33.1	20 / 44.1	30 / 66.1	40 / 88.2	45 / 99.2
Permitted total load (kg / lbs)	50 / 110.2	70 / 154.3	90 / 198.4	120 / 264.6	150 / 330.7
Cable port (Φ, mm / inch)	50 / 2.0	50 / 2.0	80 / 3.1	80 / 3.1	80 / 3.1
Exterior dimensions					
Width (mm / inch), W	1070 / 42.1	1170 / 46.1	1370 / 53.9	1495 / 58.9	1595 / 62.8
Depth / with handle (mm - inch), D/D'	740 / 775 29.1 / 30.5	890 / 925 35 / 36.4	940 / 975 37 / 38.4	1595 / 1625 62.8 / 64	1695 / 1725 66.7 / 67.9
Height (mm / inch), H	1385 / 54.5	1585 / 62.4	1780 / 70.1	1650 / 65	1800 / 70.9
Weight (kg / lbs)	170 / 374.8	260 / 573.2	320 / 705.5	390 / 859.8	440 / 970

Accessories Page 342

Cable Port, Shelves, Viewing Window, Signal Lamp, Recorder, Fan Speed Adjuster, Gas Purge System, Hour Meter

Heating & Cooling Chamber

Low temperature type (-35°C), vertical airflow

Vertical airflow and -35 to 180°C temperature range



JMV-070

with 2ea wire shelves (standard), Recorder, Signal Lamp (option)




Completion of Temperature Verification According to Strict International Standard


- Temperature verification in accordance with DIN 12880 and IEC 60068-3-5 provides excellent reliability and reproducibility.
- Provided specification the fluctuation and gradient in typical temperature point.
- Provided heating / cooling time data to help select the right model for testing purposes.
- Tested recovery time in accordance with specifications to provide data that is useful for real-use.

Structural Functional Features

- Chamber structure with excellent sealing and insulation provides a stable long-term test and reduce costs by minimizing energy loss.
- According to international standards, even when operated at upper temperature, the surface temperature does not exceed 51°C. (EN 563 standard)
- Door consists of a double-sealed structure to satisfy DIN 58371. (door airtight guidance)
- Durability even after long-term repeated use has been verified by satisfying the door load test. (SEFA standard)
- Door structure is a two-point hang structure, building more stable seal with less force.



Side View





Optimized Dedicated Control System

- 3.5-inch color touch display controller.
- Intuitive screen configuration for easy and convenient control.
- Graph display makes it easier to check operation.
- PID zone subdivided into 4 zones for more precise control stability.
- RS-232 port (default) and RS-485 port (option) supported.
- Control and data processing by connecting up to 32 devices to PC at same time via RS-485 port.
- Includes software for PC control.
- Saving data, convenient for reporting.

Use Convenience Features

- Convenient opening and closing door with either side latches.
- Door handle with built-in key lock.
- Identification of the control and operating conditions at the front panel.
- Design by stainless steel internal chamber provides clean maintenance and excellent corrosion resistance.
- Easy to clean as it is easy to remove refrigerator condenser grill, making it convenient to maintain efficiency of refrigeration.
- Equipped with easy-to-move/install caster.
- Perforated shelf for heavy-load sample. (option)
- Convenient and various options such as signal lamps, cable port, air cycle rate, gas purge, etc.


Description	No.
Pattern	100
Repeat time	999
Part repeat time	255
Max. segments / pattern	100
Available max. segments	2000
Programmable process time / segment	99 hr. 59 min.

* It is possible to set 100 segments per pattern, but the maximum number of segments is not 10,000 (100 pattern x 100 segment) but 2,000.


Outstanding Safety

- Electronic over temperature protection system. (A-OT, advanced over temperature limiter)
Independent precision temperature sensor and controller for safer and more accurate over temperature shutdown.
- General over-temperature protection. (B-OT, backup over temperature limiter)
A backup device for electronic system failure, mechanically preventing over temperature even when it occurs electronic errors.
- Door open warning and automatic shut off.
- Emergency stop button on the front of the product.
- Over-current and short circuit protection of device.
- When main power connections, users can quickly respond by notifying of electric phase sequence errors.
- Electrical instrumentation access warning and system shutdown.
- Each heater by fuse for more safety using.
- Auto stop when operating current of refrigerator is overloaded.
- Automatic stop in case of over temperature of the compressor.
- Auto shut off in case refrigerant pressure (high/low pressure) is abnormal.

Temperature Fluctuation



Recovery Time after Door Opening



Specification


Model	JMV-025	JMV-040	JMV-070	JMV-100
Temperature data				
Range (°C / °F)	-35 to 180 / -31 to 356			
Fluctuation (±°C / °F)^{1), 2)}				
at -40°C	0.4 / 0.72	0.2 / 0.36	0.2 / 0.36	0.2 / 0.36
at -25°C	0.4 / 0.72	0.3 / 0.54	0.2 / 0.36	0.3 / 0.54
at -10°C	0.5 / 0.9	0.2 / 0.36	0.2 / 0.36	0.2 / 0.36
at 25°C	0.3 / 0.54	0.3 / 0.54	0.2 / 0.36	0.2 / 0.36
at 70°C	0.6 / 1.08	0.3 / 0.54	0.2 / 0.36	0.2 / 0.36
at 100°C	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54	0.2 / 0.36
at 150°C	0.4 / 0.72	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54
at 180°C	0.3 / 0.54	0.4 / 0.72	0.4 / 0.72	0.5 / 0.9
Gradient (±°C / °F)¹⁾				
at -40°C	0.7 / 1.26	0.3 / 0.54	0.9 / 1.62	0.7 / 1.26
at -25°C	0.6 / 1.08	0.3 / 0.54	0.9 / 1.62	0.7 / 1.26
at -10°C	0.6 / 1.08	0.3 / 0.54	1.1 / 1.98	0.8 / 1.44
at 25°C	0.8 / 1.44	0.4 / 0.72	0.7 / 1.26	0.6 / 1.08
at 70°C	0.8 / 1.44	0.6 / 1.08	1.1 / 1.98	0.9 / 1.62
at 100°C	0.9 / 1.62	0.7 / 1.26	1 / 1.8	0.8 / 1.44
at 150°C	1.1 / 1.98	1.2 / 2.16	0.8 / 1.44	0.9 / 1.62
at 180°C	1.4 / 2.52	1.5 / 2.7	1 / 1.8	1.3 / 2.34
Heating time (min)²⁾				
from -40°C to 25°C	30	31	20	21
from -25°C to 25°C	20	17	18	16
from -10°C to 25°C	17	13	13	15
from 25°C to 70°C	26	18	25	20
from 25°C to 100°C	28	21	26	25
from 25°C to 150°C	50	39	39	51
from 25°C to 180°C	56	63	46	56
from -40°C to 180°C	61	81	62	73
from -25°C to 120°C	52	39	47	47
from -40°C to 5°C	7	8	6	7
from 5°C to 125°C	44	59	40	49
from 125°C to 180°C	10	14	16	17
Cooling time (min)²⁾				
from 25°C to -40°C	37	29	28	46
from 25°C to -25°C	22	20	22	24
from 25°C to -10°C	18	17	19	17
from 70°C to 25°C	20	22	23	23
from 100°C to 25°C	30	26	34	31
from 150°C to 25°C	95	92	110	160
from 180°C to 25°C	144	158	186	289
from 180°C to -40°C	172	198	192	297
from 120°C to -25°C	56	57	54	46
from 180°C to 125°C	110	133	144	228
from 125°C to 5°C	37	43	30	35
from 5°C to -40°C	25	22	18	34
Recovery time after door was opened for 30sec.(min.)²⁾				
at -40°C	6	6	6.5	9
at -25°C	5	6.5	7	6.5
at -10°C	6.5	6.5	8.5	6.5
at 25°C	5.5	4.5	2	7
at 70°C	8.5	8	9.5	7.5
at 100°C	5.5	4.5	6	5.5
at 150°C	6.5	7	9	3.5
at 180°C	6.5	6.5	8	4
Electrical data & Ordering information				
Electrical requirements (230V, 50Hz, 1P, A)	26.0			
Cat. No.	AAHKB022K	-	-	-
Cat. No. (EU only)	AAHKB025K	-	-	-
Electrical requirements (230V, 60Hz, 1P, A)	26.0			
Cat. No.	AAHKB021K	-	-	-
Electrical requirements (380V, 50/60Hz, 3P, A)		12.5	19.6	19.6
Cat. No.	-	AAHKB038K	AAHKB048K	AAHKB058K
Cat. No. (EU only)	-	AAHKB030K	AAHKB040K	AAHKB050K

1) According to IEC 60068-3-5

2) According to DIN 12880

※ Above specification value is recorded by 230V/60Hz, 380V/50Hz.

※ Above specifications can be changed without prior notice.



Dimension

Model	JMV-025	JMV-040	JMV-070	JMV-100
Interior dimensions				
Chamber volume (L / cu ft)	253 / 8.93	420 / 14.83	720 / 25.43	990 / 34.96
Width (mm / inch)	600 / 23.6	750 / 29.5	900 / 35.4	1000 / 39.4
Depth (mm / inch)	650 / 25.6	700 / 27.6	800 / 31.5	900 / 35.4
Height (mm / inch)	650 / 25.6	800 / 31.5	1000 / 39.4	1100 / 43.3
Quantity of shelves (standard/max.)	2 / 9	2 / 11	2 / 15	2 / 16
Distance of between shelves (mm / inch)	30 / 1.2	30 / 1.2	30 / 1.2	30 / 1.2
Max. Load per shelf (kg / lbs)	20 / 44.1	25 / 55.1	40 / 88.2	40 / 88.2
Permitted total load (kg / lbs)	70 / 154.3	75 / 165.3	100 / 220.5	100 / 220.5
Cable port (Φ, mm/inch)	50 / 2.0	80 / 3.1	80 / 3.1	80 / 3.1
Exterior dimensions				
Width (mm / inch), W	990 / 39	1140 / 44.9	1290 / 50.8	1390 / 54.7
Depth (mm / inch), D	1180 / 46.5	1330 / 52.4	1430 / 56.3	1530 / 60.2
Height (mm / inch), H	1350 / 53.1	1590 / 62.6	1940 / 76.4	2040 / 80.3
Weight (kg / lbs)	260 / 573.2	320 / 705.5	390 / 859.8	440 / 970

Accessories Page 342

Cable Port, Shelves, Viewing Window, Signal Lamp, Recorder, Fan Speed Adjuster, Gas Purge System, Hour Meter

Heating & Cooling Chamber

Low temperature type (-25°C), horizontal airflow

-25 to 100°C temperature range to meet basic test requirements




PBV-040

with 2ea wire shelves (standard), Recorder, Signal lamp (option)




Completion of Temperature Verification According to Strict International Standard

- Temperature verification in accordance with DIN 12880 and IEC 60068-3-5 provides excellent reliability and reproducibility.
- Provided specification the fluctuation and gradient in typical temperature point.
- Provided heating / cooling time data to help select the right model for testing purposes.
- Tested recovery time in accordance with specifications to provide data that is useful for real-use.



Structural Functional Features

- Chamber structure with excellent sealing and insulation provides a stable long-term test and reduce costs by minimizing energy loss.
- According to international standards, even when operated at upper temperature, the surface temperature does not exceed 51°C. (EN 563 standard)
- Door consists of a double-sealed structure to satisfy DIN 58371. (door airtight guidance)
- Durability even after long-term repeated use has been verified by satisfying the door load test. (SEFA standard)
- Door structure is a two-point hang structure, building more stable seal with less force.





Optimized Dedicated Control System

- 3.5-inch color touch display controller.
- Intuitive screen configuration for easy and convenient control.
- Graph display makes it easier to check operation.
- PID zone subdivided into 4 zones for more precise control stability.
- RS-232 port (default) and RS-485 port (option) supported.
- Control and data processing by connecting up to 32 devices to PC at same time via RS-485 port.
- Includes software for PC control.
- Saving data, convenient for reporting.

Use Convenience Features

- Convenient opening and closing door with either side latches.
- Door handle with built-in key lock.
- Identification of the control and operating conditions at the front panel.
- Design by stainless steel internal chamber provides clean maintenance and excellent corrosion resistance.
- Easy to clean as it is easy to remove refrigerator condenser grill, making it convenient to maintain efficiency of refrigeration.
- Equipped with easy-to-move/install caster.
- Perforated shelf for heavy-load sample. (option)
- Convenient and various options such as signal lamps, cable port, air cycle rate, gas purge, etc.


Description	No.
Pattern	100
Repeat time	999
Part repeat time	255
Max. segments / pattern	100
Available max. segments	2000
Programmable process time / segment	99 hr. 59 min.

* It is possible to set 100 segments per pattern, but the maximum number of segments is not 10,000 (100 pattern x 100 segment) but 2,000.


Outstanding Safety

- General over-temperature protection. (over temperature limiter)
- Door open warning and automatic shut off.
- Emergency stop button on the front of the product.
- Over-current and short circuit protection of device.
- When main power connections, users can quickly respond by notifying of electric phase sequence errors.
- Electrical instrumentation access warning and system shutdown.
- Each heater by fuse for more safety using.
- Auto stop when operating current of refrigerator is overloaded.
- Automatic stop in case of over temperature of the compressor.
- Auto shut off in case refrigerant pressure (high/low pressure) is abnormal.

Temperature Fluctuation



Recovery Time after Door Opening



Specification


Model	PBV-012	PBV-025	PBV-040	PBV-070	PBV-100
Temperature data					
Range (°C / °F)	-25 to 100 /-13 to 212				
Fluctuation (±°C / °F)^{1), 2)}					
at -25°C	0.2 / 0.36	0.2 / 0.36	0.2 / 0.36	0.2 / 0.36	0.2 / 0.36
at -15°C	0.2 / 0.36	0.2 / 0.36	0.2 / 0.36	0.2 / 0.36	0.2 / 0.36
at -10°C	0.3 / 0.54	0.3 / 0.54	0.2 / 0.36	0.2 / 0.36	0.2 / 0.36
at 25°C	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54	0.4 / 0.72	0.2 / 0.36
at 70°C	0.2 / 0.36	0.2 / 0.36	0.2 / 0.36	0.3 / 0.54	0.2 / 0.36
at 100°C	0.2 / 0.36	0.2 / 0.36	0.2 / 0.36	0.3 / 0.54	0.2 / 0.36
Gradient (±°C / °F)¹⁾					
at -25°C	0.3 / 0.54	0.3 / 0.54	0.4 / 0.72	0.4 / 0.72	0.4 / 0.72
at -15°C	0.3 / 0.54	0.3 / 0.54	0.4 / 0.72	0.4 / 0.72	0.4 / 0.72
at -10°C	0.3 / 0.54	0.3 / 0.54	0.4 / 0.72	0.4 / 0.72	0.4 / 0.72
at 25°C	0.3 / 0.54	0.3 / 0.54	0.4 / 0.72	0.4 / 0.72	0.4 / 0.72
at 70°C	0.4 / 0.72	0.4 / 0.72	0.6 / 1.08	0.6 / 1.08	0.6 / 1.08
at 100°C	0.6 / 1.08	0.6 / 1.08	0.7 / 1.26	0.8 / 1.44	0.8 / 1.44
Heating time (min)²⁾					
from -25°C to 25°C	19	19	29	22	24
from -15°C to 25°C	17	17	22	20	22
from -10°C to 25°C	15	15	21	18	20
from 25°C to 70°C	20	20	33	22	27
from 25°C to 100°C	34	34	52	35	38
from -25°C to 100°C	47	47	76	46	61
Cooling time (min)²⁾					
from 25°C to -25°C	36	36	26	25	34
from 25°C to -15°C	21	21	21	20	24
from 25°C to -10°C	17	17	19	18	21
from 70°C to 25°C	17	17	20	20	25
from 100°C to 25°C	28	28	30	29	35
from 100°C to -25°C	65	65	50	53	56
Recovery time after door was opened for 30sec.(min.)²⁾					
at -25°C	5	5	5	4	7
at -15°C	4.5	4.5	3.5	5.5	5
at -10°C	6	6	4.5	4	5
at 25°C	5	5	4	5.5	3.5
at 70°C	5	5	9.5	5	7.5
at 100°C	5.5	5.5	9	6	7
Electrical data & Ordering information					
Electrical requirements (230V, 50Hz, 1P, A)	16.3	22.8			
Cat. No.	AAHK5012K	AAHK5022K	-	-	-
Cat. No. (EU only)	AAHK5015K	AAHK5025K	-	-	-
Electrical requirements (230V, 60Hz, 1P, A)	16.3	22.8			
Cat. No.	AAHK5011K	AAHK5021K	-	-	-
Electrical requirements (380V, 50/60Hz, 3P, A)			10.0	15.4	16.8
Cat. No.	-	-	AAHK5038K	AAHK5048K	AAHK5058K
Cat. No. (EU only)	-	-	AAHK5030K	AAHK5040K	AAHK5050K

1) According to IEC 60068-3-5

2) According to DIN 12880

※ Above specification value is recorded by 230V/60Hz, 380V/50Hz.

※ Above specifications can be changed without prior notice.



Dimension

Model	PBV-012	PBV-025	PBV-040	PBV-070	PBV-100
Interior dimensions					
Chamber volume (L / cu ft)	125 / 4.41	250 / 8.83	400 / 14.13	700 / 24.72	1000 / 35.31
Width (mm / inch)	500 / 19.7	600 / 23.6	750 / 29.5	900 / 35.4	1000 / 39.4
Depth (mm / inch)	500 / 19.7	650 / 25.6	700 / 27.6	800 / 31.5	910 / 35.8
Height (mm / inch)	500 / 19.7	650 / 25.6	800 / 31.5	1000 / 39.4	1100 / 43.3
Quantity of shelves (standard/max.)	2 / 6	2 / 9	2 / 11	2 / 15	2 / 16
Distance of between shelves (mm / inch)	30 / 1.2	30 / 1.2	30 / 1.2	30 / 1.2	30 / 1.2
Max. Load per shelf (kg / lbs)	15 / 33.1	20 / 44.1	30 / 66.1	40 / 88.2	45 / 99.2
Permitted total load (kg / lbs)	50 / 110.2	70 / 154.3	90 / 198.4	120 / 264.6	150 / 330.7
Cable port (Φ, mm/inch)	50 / 2.0	50 / 2.0	80 / 3.1	80 / 3.1	80 / 3.1
Exterior dimensions					
Width (mm / inch), W	1060 / 41.7	1154 / 45.4	1304 / 51.3	1454 / 57.2	1720 / 67.7
Depth / with handle (mm - inch), D/D'	672 / 702 26.5 / 27.6	822 / 852 32.4 / 33.5	872 / 902 34.3 / 35.5	1510 / 1540 59.4 / 60.6	1620 / 1650 63.8 / 65
Height (mm / inch), H	1340 / 52.8	1540 / 60.6	1724 / 67.9	1628 / 64.1	1534 / 60.4
Weight (kg / lbs)	360 / 793.7	420 / 925.9	520 / 1146.4	640 / 1411	700 / 1543.2

Accessories [Page 342](#)

Cable Port, Shelves, Viewing Window, Signal Lamp, Recorder, Fan Speed Adjuster, Gas Purge System, Hour Meter

Heating & Cooling Chamber

Ambient temperature type (-5°C), horizontal airflow

-5 to 100°C temperature range to meet basic test requirements




PMV-040

with 2ea wire shelves (standard), Recorder, Signal lamp (option)




Completion of Temperature Verification According to Strict International Standard

- Temperature verification in accordance with DIN 12880 and IEC 60068-3-5 provides excellent reliability and reproducibility.
- Provided specification the fluctuation and gradient in typical temperature point.
- Provided heating / cooling time data to help select the right model for testing purposes.
- Tested recovery time in accordance with specifications to provide data that is useful for real-use.



Structural Functional Features

- Chamber structure with excellent sealing and insulation provides a stable long-term test and reduce costs by minimizing energy loss.
- According to international standards, even when operated at upper temperature, the surface temperature does not exceed 51°C. (EN 563 standard)
- Door consists of a double-sealed structure to satisfy DIN 58371. (door airtight guidance)
- Durability even after long-term repeated use has been verified by satisfying the door load test. (SEFA standard)
- Door structure is a two-point hang structure, building more stable seal with less force.






Optimized Dedicated Control System

- 3.5-inch color touch display controller.
- Intuitive screen configuration for easy and convenient control.
- Graph display makes it easier to check operation.
- PID zone subdivided into 4 zones for more precise control stability.
- RS-232 port (default) and RS-485 port (option) supported.
- Control and data processing by connecting up to 32 devices to PC at same time via RS-485 port.
- Includes software for PC control.
- Saving data, convenient for reporting.

Use Convenience Features

- Convenient opening and closing door with either side latches.
- Door handle with built-in key lock.
- Identification of the control and operating conditions at the front panel.
- Design by stainless steel internal chamber provides clean maintenance and excellent corrosion resistance.
- Easy to clean as it is easy to remove refrigerator condenser grill, making it convenient to maintain efficiency of refrigeration.
- Equipped with easy-to-move/install caster.
- Perforated shelf for heavy-load sample. (option)
- Convenient and various options such as signal lamps, cable port, air cycle rate, gas purge, etc.




Description	No.
Pattern	100
Repeat time	999
Part repeat time	255
Max. segments / pattern	100
Available max. segments	2000
Programmable process time / segment	99 hr. 59 min.

* It is possible to set 100 segments per pattern, but the maximum number of segments is not 10,000 (100 pattern x 100 segment) but 2,000.

Outstanding Safety

- General over-temperature protection. (over temperature limiter)
- Door open warning and automatic shut off.
- Emergency stop button on the front of the product.
- Over-current and short circuit protection of device.
- When main power connections, users can quickly respond by notifying of electric phase sequence errors.
- Electrical instrumentation access warning and system shutdown.
- Each heater by fuse for more safety using.
- Auto stop when operating current of refrigerator is overloaded.
- Automatic stop in case of over temperature of the compressor.
- Auto shut off in case refrigerant pressure (high/low pressure) is abnormal.



Specification


Model	PMV-012	PMV-025	PMV-040	PMV-070	PMV-100
Temperature data					
Range (°C / °F)	-5 to 100 / 23 to 212				
Fluctuation (±°C / °F)^{1), 2)}					
at -5°C	0.3 / 0.54	0.3 / 0.54	0.2 / 0.36	0.2 / 0.36	0.3 / 0.54
at 0°C	0.3 / 0.54	0.3 / 0.54	0.2 / 0.36	0.2 / 0.36	0.3 / 0.54
at 25°C	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54
at 40°C	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54
at 60°C	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54	0.4 / 0.72	0.3 / 0.54
at 80°C	0.3 / 0.54	0.3 / 0.54	0.2 / 0.36	0.2 / 0.36	0.3 / 0.54
at 100°C	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54	0.4 / 0.72	0.3 / 0.54
Gradient (±°C / °F)¹⁾					
at -5°C	0.4 / 0.72	0.4 / 0.72	0.4 / 0.72	0.4 / 0.72	0.3 / 0.54
at 0°C	0.4 / 0.72	0.4 / 0.72	0.4 / 0.72	0.4 / 0.72	0.3 / 0.54
at 25°C	0.2 / 0.36	0.2 / 0.36	0.3 / 0.54	0.4 / 0.72	0.3 / 0.54
at 40°C	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54	0.4 / 0.72	0.3 / 0.54
at 60°C	0.5 / 0.9	0.5 / 0.9	0.4 / 0.72	0.4 / 0.72	0.4 / 0.72
at 80°C	0.7 / 1.26	0.7 / 1.26	0.6 / 1.08	0.5 / 0.9	0.5 / 0.9
at 100°C	1 / 1.8	1 / 1.8	0.8 / 1.44	0.7 / 1.26	0.7 / 1.26
Heating time (min)²⁾					
from -5°C to 25°C	24	24	21	17	19
from 0°C to 25°C	20	20	17	15	16
from 25°C to 40°C	16	16	15	11	11
from 25°C to 60°C	30	30	26	23	25
from 25°C to 80°C	32	32	33	34	38
from 25°C to 100°C	45	45	48	52	56
from -5°C to 100°C	60	60	62	64	75
Cooling time (min)²⁾					
from 25°C to -5°C	36	36	35	31	31
from 25°C to 0°C	30	30	27	25	25
from 40°C to 25°C	17	17	16	14	16
from 60°C to 25°C	37	37	35	31	31
from 80°C to 25°C	52	37	48	43	45
from 100°C to 25°C	65	65	62	59	57
from 100°C to -5°C	97	97	92	87	91
Recovery time after door was opened for 30sec.(min.)²⁾					
at -5°C	8.0	8.0	6.5	5.0	5.0
at 0°C	6.0	6.0	6.0	6.0	8.0
at 25°C	0.5	0.5	4.0	6.0	6.0
at 40°C	2.5	2.5	2.5	2.5	2.5
at 60°C	4.0	4.0	4.5	5.0	6.0
at 80°C	5.0	5.0	5.0	5.5	6.0
at 100°C	7.0	7.0	7.0	6.5	8.0
Electrical data & Ordering information					
Electrical requirements (230V, 50Hz, 1P, A)			15.5	21.7	
Cat. No.	-	-	AAHK4032K	AAHK4042K	-
Cat. No. (EU only)	-	-	AAHK4035K	AAHK4045K	-
Electrical requirements (230V, 60Hz, 1P, A)			15.5	21.7	
Cat. No.	-	-	AAHK4031K	AAHK4041K	-
Electrical requirements (380V, 50/60Hz, 1P, A)	9.7	13.2			
Cat. No.	AAHK4012K	AAHK4022K	-	-	-
Cat. No. (EU only)	AAHK4015K	AAHK4025K	-	-	-
Electrical requirements (380V, 50/60Hz, 3P, A)					6.0
Cat. No.	-	-	-	-	AAHK4058K
Cat. No. (EU only)	-	-	-	-	AAHK4050K

1) According to IEC 60068-3-5

2) According to DIN 12880

※ Above specification value is recorded by 230V/60Hz, 380V/50Hz.

※ Above specifications can be changed without prior notice.



Dimension

Model	PMV-012	PMV-025	PMV-040	PMV-070	PMV-100
Interior dimensions					
Chamber volume (L / cu ft)	125 / 4.41	250 / 8.83	400 / 14.13	700 / 24.72	1000 / 35.31
Width (mm / inch)	500 / 19.7	600 / 23.6	750 / 29.5	900 / 35.4	1000 / 39.4
Depth (mm / inch)	500 / 19.7	650 / 25.6	700 / 27.6	800 / 31.5	910 / 35.8
Height (mm / inch)	500 / 19.7	650 / 25.6	800 / 31.5	1000 / 39.4	1100 / 43.3
Quantity of shelves (standard/max.)	2 / 6	2 / 9	2 / 11	2 / 15	2 / 16
Distance of between shelves (mm / inch)	30 / 1.2	30 / 1.2	30 / 1.2	30 / 1.2	30 / 1.2
Max. Load per shelf (kg / lbs)	15 / 33.1	20 / 44.1	30 / 66.1	40 / 88.2	45 / 99.2
Permitted total load (kg / lbs)	50 / 110.2	70 / 154.3	90 / 198.4	120 / 264.6	150 / 330.7
Cable port (Φ, mm/inch)	50 / 2.0	50 / 2.0	80 / 3.1	80 / 3.1	80 / 3.1
Exterior dimensions					
Width (mm / inch), W	1060 / 41.7	1154 / 45.4	1304 / 51.3	1454 / 57.2	1720 / 67.7
Depth / with handle (mm - inch), D/D'	672 / 702 26.5 / 27.6	822 / 852 32.4 / 33.5	872 / 902 34.3 / 35.5	1510 / 1540 59.4 / 60.6	1620 / 1650 63.8 / 65
Height (mm / inch), H	1340 / 52.8	1540 / 60.6	1724 / 67.9	1628 / 64.1	1534 / 60.4
Weight (kg / lbs)	360 / 793.7	420 / 925.9	520 / 1146.4	640 / 1411	700 / 1543.2

Accessories Page 342

Cable Port, Shelves, Viewing Window, Signal Lamp, Recorder, Fan Speed Adjuster, Gas Purge System, Hour Meter

Heating & Cooling Chamber

General type (-20°C, 0), horizontal airflow

Includes space-saving vertical structure and dual chamber model

Optimized Model Configuration

- Two series (LCH, LCH-G) according to temperature control range.
LCH : -20 to 100°C
LCH-G : 0 to 100°C
- Space-saving vertical structural design.
- Independent control of the temperature with dual chambers divided into upper and lower.
- Standard temperature test chamber to select optimized model according to test conditions.

Structural Functional Features

- Chamber structure with excellent sealing and insulation provides a stable long-term test and reduce costs by minimizing energy loss.
- Ensures that the surface temperature is kept below 60°C.
- Cable port (Ø80 mm) is included as standard, making it convenient to connect external equipment.
- Easy to clean as it is easy to remove refrigerator condenser grill, making it convenient to maintain efficiency of refrigeration.
- Airflow optimized for uniform heat transfer.

Use Convenience Features

- Microprocessor PID method for precise temperature control.
- Optimized temperature control with temperature auto-tuning.
- Calibration function minimizes temperature deviation.
- 9 Steps Program Control. (200 times repeat function)
- Wait On/Off timer. (up to 99 hours 59 minutes)
- Save and use 3 frequently used temperatures.
- Identification of the control and operating conditions at the front panel.
- Internal chamber made of stainless steel provides clean maintenance and excellent corrosion resistance.
- Equipped with easy-to-move/install caster.

Outstanding Safety

- Over temperature limit function. (over temperature limiter)
- Malfunction prevented by controller lock function.
- Over-current and short circuit protection of device.
- Door open warning function.




LCH-21

with 2ea wire shelves(standard)
Recorder (option)

LCH-11G-2C


with 2ea wire shelves (standard)
per chamber



Specification

Model	LCH-11	LCH-21	LCH-31	LCH-11-2C	LCH-11G	LCH-21G	LCH-31G	LCH-11G-2C
Temperature data								
Range (°C / °F)	-20 to 100 / -4 to 212	0 to 100 / 32 to 212	0 to 100 / 32 to 212	0 to 100 / 32 to 212	0 to 100 / 32 to 212			
Fluctuation (±°C / °F) ¹⁾	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54
Variation at 80°C (±°C / °F) ¹⁾	0.8 / 1.44	0.8 / 1.44	0.8 / 1.44	0.8 / 1.44	0.8 / 1.44	0.8 / 1.44	0.8 / 1.44	0.8 / 1.44
Heating time (20°C to 100°C, min)	80	80	80	80	80	80	80	80
Cooling time (100°C to 40°C, min)	60	60	60	60	60	60	60	60
Interior dimensions								
Chamber volume (L / cu ft)	150 / 5.3	255 / 9	485 / 17.1	150x2/53x2	150 / 5.3	255 / 9	485 / 17.1	150x2/53x2
Width (mm / inch)	600 / 23.6	600 / 23.6	700 / 27.6	600 / 23.6 ²⁾	600 / 23.6	600 / 23.6	700 / 27.6	600 / 23.6 ²⁾
Depth (mm / inch)	500 / 19.7	500 / 19.7	680 / 26.8	500 / 19.7 ²⁾	500 / 19.7	500 / 19.7	680 / 26.8	500 / 19.7 ²⁾
Height (mm / inch)	500 / 19.7	850 / 33.5	1020 / 40.2	500 / 19.7 ²⁾	500 / 19.7	850 / 33.5	1020 / 40.2	500 / 19.7 ²⁾
Quantity of shelves (standard/max.)	2 / 6	2 / 12	2 / 14	2 / 6	2 / 6	2 / 12	2 / 14	2 ²⁾ / 6 ²⁾
Distance of between shelves (mm / inch)	30 / 1.2	30 / 1.2	30 / 1.2	30 / 1.2	30 / 1.2	30 / 1.2	30 / 1.2	30 / 1.2
Max. load per shelf (kg / lbs)	15 / 33.1	20 / 44.1	20 / 44.1	15 / 33.1	15 / 33.1	20 / 44.1	20 / 44.1	15 / 33.1
Permitted total load (kg / lbs)	50 / 110.2	70 / 154.3	70 / 154.3	50 / 110.2 ²⁾	50 / 110.2	70 / 154.3	70 / 154.3	50 / 110.2 ²⁾
Cable port (Ø, mm / inch)	80 / 3.1	80 / 3.1	80 / 3.1	80 / 3.1 ²⁾	80 / 3.1	80 / 3.1	80 / 3.1	80 / 3.1 ²⁾
Weight (kg / lbs)	125 / 275.6	185 / 407.9	225 / 496	195 / 429.9	120 / 264.6	180 / 396.8	220 / 485	190 / 418.9
Electrical data & Ordering information								
Electrical requirement (230V, 50Hz, 1P, A)	6.3	9.2	11.1	12.6	7.1	10.6	13.8	14.2
Cat. No.	AAHK6112K	AAHK6212K	AAHK6312K	AAHK6412K	AAHK6122K	AAHK6222K	AAHK6322K	AAHK6422K
Electrical requirement (230V, 60Hz, 1P, A)	7.0	10.5	13.0	14.0	6.2	9.1	11.5	12.4
Cat. No.	AAHK6111K	AAHK6211K	AAHK6311K	AAHK6411K	AAHK6121K	AAHK6221K	AAHK6321K	AAHK6421K
Electrical requirement (120V, 60Hz, 1P, A)	11.6	-	-	23.2	13	-	-	26
Cat. No.	AAHK6113U	-	-	AAHK6413U	AAHK6123U	-	-	AAHK6423U

- 1) According to IEC 60068-3-5
2) Based on inner chamber 1ea.



Dimension

Model	LCH-11G LCH-11	LCH-21G LCH-21	LCH-31G LCH-31
W (mm / inch)	740 / 29.1 794 / 31.3	740 / 29.1 794 / 31.3	840 / 33.1 894 / 35.2
D (mm / inch)	1032 / 40.6 1070 / 42.1	1032 / 40.6 1070 / 42.1	1212 / 47.7 1250 / 49.2
H (mm / inch)	1280 / 50.4	1630 / 64.2	1800 / 70.9
A (mm / inch)	355 / 14	355 / 14	445 / 17.5
B (mm / inch)	845 / 33.3	1050 / 41.3	1155 / 45.5

Dimension


Model	LCH-11G-2C
	LCH-11-2C
W (mm / inch)	740 / 29.1 794 / 31.3
D (mm / inch)	1032 / 40.6 1070 / 42.1
H (mm / inch)	1900 / 74.8
A (mm / inch)	355 / 14
B1 (mm / inch)	875 / 34.4
B2 (mm / inch)	1495 / 58.9

Accessories



Cable Port

- Both cap and silicone plug are included as standard in Ø50 mm cable port.
- It is available to add Ø50 mm, Ø80 mm.



Wire Shelf, Perforated Shelf

- Excellent corrosion resistance as made of stainless steel.
- 2 wire shelves are given as standard.
- Perforated shelf suitable for high load samples. (optional)



Viewing Window

- Constantly observing the samples under test in anytime.
- Tempered glass window with wire heater and LED lamp.



Signal Lamp

- Figuring out equipment status remotely.
- Run/Stand-by/Error displayed by color.



LC DataKeeper (only TC3 models)

- Cat. No.: SSA11001 (1yr), SSA11002 (permanent)
- High-integrity software compliant with FDA 21 CFR Part 11
- Integrated with Jeio Tech equipment for GMP audits



LC GreenBox

- Monitor and control the operation status of the device in real time via mobile app with simple internet connection.
- One LC GreenBox can connect up to 4 devices. (only TC3 model (Cat. No.: AAHQ1011K))

Model	Standard	Cable Port ø50	Cable Port ø80	Wire Shelf	Perforated Shelf	Viewing Window	Signal Lamp
TC3-KE-025	ø50 Cable Port	AAA8T612	N/A	00RTD0000679	AAA80602-6	N/A	AAA80550
TC3-KE-065	ø50 Cable Port	AAA8T612	AAA8T611	00RTD0000680	AAA80602-7	N/A	AAA80550
TC3-KE-100	ø50 Cable Port	AAA8T612	AAA8T611	00RTD0000681	AAA80602-8	N/A	AAA80550
TC3-ME-025	ø50 Cable Port	AAA8T612	N/A	00RTD0000679	AAA80602-6	N/A	AAA80550
TC3-ME-065	ø50 Cable Port	AAA8T612	AAA8T611	00RTD0000680	AAA80602-7	N/A	AAA80550
TC3-ME-100	ø50 Cable Port	AAA8T612	AAA8T611	00RTD0000681	AAA80602-8	N/A	AAA80550
KBD-012	ø50 Cable Port	AAAK8511	AAAK8512	00LTV0000070	AAAK1501	N/A	AAA80550
KBD-025	ø50 Cable Port	AAAK8511	AAAK8512	00LTV0000053	AAAK1502	N/A	AAA80550
KBD-040	ø80 Cable Port	AAAK8511	AAAK8512	00LTV0000029	AAAK1503	N/A	AAA80550
KBD-070	ø80 Cable Port	AAAK8511	AAAK8512	00LTV0000030	AAAK1504	N/A	AAA80550
KBD-100	ø80 Cable Port	AAAK8511	AAAK8512	00LTV0000057	AAAK1505	N/A	AAA80550
KMV-012	ø50 Cable Port	AAAK8511	AAAK8512	00LTV0000070	AAAK1501	AAAK8501	AAA80550
KMV-025	ø50 Cable Port	AAAK8511	AAAK8512	00LTV0000053	AAAK1502	AAAK8501	AAA80550
KMV-040	ø80 Cable Port	AAAK8511	AAAK8512	00LTV0000029	AAAK1503	AAAK8501	AAA80550
KMV-070	ø80 Cable Port	AAAK8511	AAAK8512	00LTV0000030	AAAK1504	AAAK8601	AAA80550
KMV-100	ø80 Cable Port	AAAK8511	AAAK8512	00LTV0000057	AAAK1505	AAAK8601	AAA80550
JMV-025	ø50 Cable Port	AAAK8511	AAAK8512	00LTV0000053	AAAK1502	AAAK8501	AAA80550
JMV-040	ø80 Cable Port	AAAK8511	AAAK8512	00LTV0000029	AAAK1503	AAAK8501	AAA80550
JMV-070	ø80 Cable Port	AAAK8511	AAAK8512	00LTV0000030	AAAK1504	AAAK8601	AAA80550
JMV-100	ø80 Cable Port	AAAK8511	AAAK8512	00LTV0000057	AAAK1505	AAAK8601	AAA80550
PBV-012	ø50 Cable Port	AAAK8511	AAAK8512	00LTV0000070	AAAK1501	AAAK8501	AAA80550
PBV-025	ø50 Cable Port	AAAK8511	AAAK8512	00LTV0000053	AAAK1502	AAAK8501	AAA80550
PBV-040	ø80 Cable Port	AAAK8511	AAAK8512	00LTV0000029	AAAK1503	AAAK8501	AAA80550
PBV-070	ø80 Cable Port	AAAK8511	AAAK8512	00LTV0000030	AAAK1504	AAAK8501	AAA80550
PBV-100	ø80 Cable Port	AAAK8511	AAAK8512	00LTV0000057	AAAK1505	AAAK8501	AAA80550
PMV-012	ø50 Cable Port	AAAK8511	AAAK8512	00LTV0000070	AAAK1501	AAAK8501	AAA80550
PMV-025	ø50 Cable Port	AAAK8511	AAAK8512	00LTV0000053	AAAK1502	AAAK8501	AAA80550
PMV-040	ø80 Cable Port	AAAK8511	AAAK8512	00LTV0000029	AAAK1503	AAAK8501	AAA80550
PMV-070	ø80 Cable Port	AAAK8511	AAAK8512	00LTV0000030	AAAK1504	AAAK8501	AAA80550
PMV-100	ø80 Cable Port	AAAK8511	AAAK8512	00LTV0000057	AAAK1505	AAAK8501	AAA80550
LCH-11	ø80 Cable Port	AAA8T610	AAA8T611	00EDA0008219	AAA22521	N/A	N/A
LCH-21	ø80 Cable Port	AAA8T610	AAA8T611	00EDA0008220	AAA22522	N/A	N/A
LCH-31	ø80 Cable Port	AAA8T610	AAA8T611	AAAK6531	AAAK6532	N/A	N/A
LCH-11-2C	ø80 Cable Port	AAA8T610	AAA8T611	00EDA0008219	AAA22521	N/A	N/A
LCH-11G	ø80 Cable Port	AAA8T610	AAA8T611	00EDA0008219	AAA22521	N/A	N/A
LCH-21G	ø80 Cable Port	AAA8T610	AAA8T611	00EDA0008220	AAA22522	N/A	N/A
LCH-31G	ø80 Cable Port	AAA8T610	AAA8T611	AAAK6531	AAAK6532	N/A	N/A
LCH-11G-2C	ø80 Cable Port	AAA8T610	AAA8T611	00EDA0008219	AAA22521	N/A	N/A



Recorder (Paper)

- Up to 6 channels supported.
- Set recording range and speed.
- Paper width: 100mm.



Recorder (Digital)

- Paperless. (5.7" LCD Display)
- Up to 6 channels supported.
- Data stored in internal/external memory.



Gas Purge System

- Gas system that can replace and purge the inside of chamber with nitrogen and CO₂ gas.



Hour Meter

- Cumulative usage time management. (No reset function)
- Minimum units of measurement: 0.1 hours. (6 minutes)
- Maximum display 99999.9 hours.



Stand (1 Stage)

- Dedicated stand for small appliances.
- Provides Lower drawer and loading space.
- Casters provided as standard for moving/fixing.





Stand (2 Stage)

- Placing of two compact devices vertically.
- Minimized installation space for different temperature and tests.

Model	Recorder (Paper)	Recorder (Digital)	Fan Speed Adjuster	Gas Purge System	Hour Meter	Stand (1 stage)	Stand (2 stage)
TC3-KE-025	AAA8T504	AAA8T506	N/A	N/A	N/A	AAA80631	AAA80634
TC3-KE-065	AAA8T504	AAA8T506	N/A	N/A	N/A	AAA80632	AAA80635
TC3-KE-100	AAA8T504	AAA8T506	N/A	N/A	N/A	AAA80633	N/A
TC3-ME-025	AAA8T504	AAA8T506	N/A	N/A	N/A	AAA80631	AAA80634
TC3-ME-065	AAA8T504	AAA8T506	N/A	N/A	N/A	AAA80632	AAA80635
TC3-ME-100	AAA8T504	AAA8T506	N/A	N/A	N/A	AAA80633	N/A
KBD-012	AAA8T500	AAA8T505	N/A	N/A	N/A	N/A	N/A
KBD-025	AAA8T500	AAA8T505	N/A	N/A	N/A	N/A	N/A
KBD-040	AAA8T500	AAA8T505	N/A	N/A	N/A	N/A	N/A
KBD-070	AAA8T500	AAA8T505	N/A	N/A	N/A	N/A	N/A
KBD-100	AAA8T500	AAA8T505	N/A	N/A	N/A	N/A	N/A
KMV-012	AAA8T500	AAA8T505	AAAK1531	AAA80691	AAAK1532	N/A	N/A
KMV-025	AAA8T500	AAA8T505	AAAK1531	AAA80691	AAAK1532	N/A	N/A
KMV-040	AAA8T500	AAA8T505	AAAK1531	AAA80691	AAAK1532	N/A	N/A
KMV-070	AAA8T500	AAA8T505	AAAK1531	AAA80691	AAAK1532	N/A	N/A
KMV-100	AAA8T500	AAA8T505	AAAK1531	AAA80691	AAAK1532	N/A	N/A
JMV-025	AAA8T500	AAA8T505	AAAK1531	AAA80691	AAAK1532	N/A	N/A
JMV-040	AAA8T500	AAA8T505	AAAK1531	AAA80691	AAAK1532	N/A	N/A
JMV-070	AAA8T500	AAA8T505	AAAK1531	AAA80691	AAAK1532	N/A	N/A
JMV-100	AAA8T500	AAA8T505	AAAK1531	AAA80691	AAAK1532	N/A	N/A
PBV-012	AAA8T500	AAA8T505	AAAK1531	AAA80691	AAAK1532	N/A	N/A
PBV-025	AAA8T500	AAA8T505	AAAK1531	AAA80691	AAAK1532	N/A	N/A
PBV-040	AAA8T500	AAA8T505	AAAK1531	AAA80691	AAAK1532	N/A	N/A
PBV-070	AAA8T500	AAA8T505	AAAK1531	AAA80691	AAAK1532	N/A	N/A
PBV-100	AAA8T500	AAA8T505	AAAK1531	AAA80691	AAAK1532	N/A	N/A
PMV-012	AAA8T500	AAA8T505	AAAK1531	AAA80691	AAAK1532	N/A	N/A
PMV-025	AAA8T500	AAA8T505	AAAK1531	AAA80691	AAAK1532	N/A	N/A
PMV-040	AAA8T500	AAA8T505	AAAK1531	AAA80691	AAAK1532	N/A	N/A
PMV-070	AAA8T500	AAA8T505	AAAK1531	AAA80691	AAAK1532	N/A	N/A
PMV-100	AAA8T500	AAA8T505	AAAK1531	AAA80691	AAAK1532	N/A	N/A
LCH-11	AAA8T500	AAA8T505	N/A	N/A	N/A	N/A	N/A
LCH-21	AAA8T500	AAA8T505	N/A	N/A	N/A	N/A	N/A
LCH-31	AAA8T500	AAA8T505	N/A	N/A	N/A	N/A	N/A
LCH-11-2C	AAA8T500	AAA8T505	N/A	N/A	N/A	N/A	N/A
LCH-11G	AAA8T500	AAA8T505	N/A	N/A	N/A	N/A	N/A
LCH-21G	AAA8T500	AAA8T505	N/A	N/A	N/A	N/A	N/A
LCH-31G	AAA8T500	AAA8T505	N/A	N/A	N/A	N/A	N/A
LCH-11G-2C	AAA8T500	AAA8T505	N/A	N/A	N/A	N/A	N/A

Heating Chamber



Real-time equipment monitoring and control system using mobile app.



Registered safety patent based on Jeio Tech's proprietary technology.



Recognized as an excellent design product by the Ministry of Commerce, Industry, and Energy.



More precise temperature control through 3-point temperature calibration.



More precise temperature control through temperature calibration.



Electronic over temperature protection system.



Conventional over temperature protection, a kind of backup device against electronic over temperature protection failure.



Warning alarm in case of deviation of temperature control.



Purifies and discharges harmful gas through internal filter.



Intuitive operation with color touch display.



Control of temperature / time, etc. through dedicated program.



Possible to set end time or start time of device operation.



PC communication via RS-232 / RS-485 / USB port.



It memorizes the state of power failure and operates automatically when power is restored.



2 year warranty free A/S.


General Application

Industrial oven: Heat up to 350°C depending on the industrial site, and size is selectable from 125L to 1000L.


Designed and tested based on 27 temperature measurement points in accordance with international standards

Horizontal airflow type with intake/damper included as standard.


- › **Completion of verification according to strict international standards**
High reliability by completing temperature verification at 27 points according to DIN 12880 standard.
- › **Includes air intake/damper as standard**
Adjustable air intake and damper are included as standard, making it easy to use.
- › **Rapid air exchange capacity**
Air intake and damper control as standard for quick air exchange and various tests.
- › **Enhanced safety thanks to over temperature double cut-off**
Independent precision-type electronic over temperature protection device.
Improved safety through double installation of mechanical devices.
- › **Steady and uniform horizontal airflow**
Horizontal airflow optimized for uniform heat transfer ensures stable overall temperature uniformity.
- › **Maintain precise temperature uniformity**
Air into Intake maintains precise control while supplying circulatory fresh air throughout the heater.
- › **Excellent chamber structure with sealing and insulation**
Excellent insulation and sealing to minimize energy loss and reduce test costs.
- › **Model provided based on temperature range and capacity**
10 standard models up to 350°C and 1000L.
Customized model can be provided according to user requirements.



Tested based on 27 temperature measurement points in accordance with DIN 12880





Uniform horizontal airflow



Structure allows for fresh air outside to be circulated through the heater

Heating Chamber

Suitable for test of temperature feature such as heat treatment and drying



Completion of Temperature Verification According to Strict International Standard

- Temperature verification in accordance with DIN 12880 and IEC 60068-3-5 provides excellent reliability and reproducibility.
- Provided specification with the fluctuation and gradient in typical temperature point.
- Provided heating / cooling time data to help select the right model for testing purposes.
- Tested recovery time in accordance with specifications to provide data that is useful for real-use.


Structural Functional Features

- Built-in air intake and damper that can be easily adjusted.
- There is a separate guard around the damper to protect the user from burns due to high temperature.
- Chamber structure with excellent sealing and insulation provides a stable long-term test and reduce costs by minimizing energy loss.
- According to international standards, even when operated at upper temperature, the surface temperature does not exceed 51°C. (EN-563 standard)
- Door consists of a double-sealed structure to satisfy DIN 58371. (door airtight guidance)
- Durability even after long-term repeated use has been verified by satisfying the door load test. (SEFA standard)
- Door locking device as standard offer.
- Equipped with easy-to-move/install caster.




LTV-012

with Wire Shelves 2ea (standard)
Recorder (option)




Uniform horizontal airflow

Temperature Fluctuation



LTV-070, 350°C

Recovery Time after Door Opening



LTV-070, 150°C



Optimized Dedicated Control System

- 3.5-inch color touch display controller applied.
- Intuitive screen configuration for easy and convenient control.
- Graph display makes it easier to check operation.
- PID zone subdivided into 4 zones for more precise control stability.
- RS-232 port (default) and RS-485 port (option) supported.
- Control and data processing by connecting up to 32 devices to PC at same time via RS-485 port.
- Includes software for PC control.
- Saving data, convenient for reporting.

Description	No.
Pattern	100
Repeat time	999
Part repeat time	255
Max. segments / pattern	100
Available max. segments	2000
Programmable process time / segment	99 hour 59 min.

* It is possible to set 100 segments per pattern, but the maximum number of segments is not 10,000 (100 pattern x 100 segment) but 2,000.




Use Convenience Features

- Convenient opening and closing door with either side latches.
- Door handle with built-in key lock.
- Identify all operations and operating conditions from the front panel.
- Design by stainless steel internal chamber provide clean maintenance and excellent corrosion resistance.
- Perforated shelf for heavy-load sample. (option)
- Convenient and various options such as signal lamps, cable port, air cycle rate, gas purge, etc.



Outstanding Safety

- Electronic over temperature protection system. (A-OT, advanced over temperature limiter)
Independent precision temperature sensor and controller for safer and more accurate over temperature shutdown.
- General over-temperature protection. (B-OT, backup over temperature limiter)
A backup device for electronic system failure, mechanically preventing over temperature even when it occurs electronic errors.
- Door open warning and automatic shut off.
- Emergency stop button on the front of the product.
- Over-current and short circuit protection of device.
- When main power connections, users can quickly respond by notifying of electric phase sequence errors.
- Electrical instrumentation access warning and system shutdown.
- Each heater by fuse for more safety using.



Structure allows for fresh air outside to be circulated through the heater.

Heating Chamber (Industry)

Specification

Model	Max. 250°C					Max. 350°C				
	LBV-012	LBV-025	LBV-040	LBV-070	LBV-100	LTV-012	LTV-025	LTV-040	LTV-070	LTV-100
Temperature data										
Range (intake/damper 100% close, °C / °F)	Amb.+45 to 250 / Amb.+81 to 482					Amb.+45 to 350 / Amb.+81 to 662				
Range (intake/damper 100% open, °C / °F)	Max. 80 / 176	Max. 80 / 176	Max. 120 / 248	Max. 80 / 176	Max. 140 / 284	Max. 115 / 239	Max. 115 / 239	Max. 180 / 356	Max. 150 / 302	Max. 210 / 410
Fluctuation (±C / °F)^{1), 2)}										
at 100°C / 212°F	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54	0.3 / 0.54	0.4 / 0.72	0.4 / 0.72	0.3 / 0.54	0.5 / 0.9	0.4 / 0.72
at 150°C / 302°F	0.3 / 0.54	0.3 / 0.54	0.4 / 0.72	0.4 / 0.72	0.4 / 0.72	0.4 / 0.72	0.4 / 0.72	0.3 / 0.54	0.4 / 0.72	0.6 / 1.08
at 200°C / 392°F	0.4 / 0.72	0.4 / 0.72	0.4 / 0.72	0.3 / 0.54	0.3 / 0.54	0.5 / 0.9	0.5 / 0.9	0.5 / 0.9	0.4 / 0.72	0.3 / 0.54
at 250°C / 482°F	0.5 / 0.9	0.5 / 0.9	0.5 / 0.9	0.6 / 1.08	0.4 / 0.72	0.6 / 1.08	0.6 / 1.08	0.6 / 1.08	0.4 / 0.72	0.4 / 0.72
at 300°C / 572°F	-	-	-	-	-	0.5 / 0.9	0.5 / 0.9	0.7 / 1.26	0.6 / 1.08	0.5 / 0.9
at 350°C / 662°F	-	-	-	-	-	0.7 / 1.26	0.7 / 1.26	0.6 / 1.08	0.6 / 1.08	0.8 / 1.44
Gradient (±C / °F)¹⁾										
at 100°C / 212°F	0.4 / 0.72	0.4 / 0.72	0.5 / 0.9	0.5 / 0.9	0.5 / 0.9	0.5 / 0.9	0.5 / 0.9	0.6 / 1.08	0.7 / 1.26	0.4 / 0.72
at 150°C / 302°F	0.7 / 1.26	0.7 / 1.26	1 / 1.8	0.7 / 1.26	0.8 / 1.44	0.7 / 1.26	0.7 / 1.26	1 / 1.8	1.1 / 1.98	0.6 / 1.08
at 200°C / 392°F	1.1 / 1.98	1.1 / 1.98	1.7 / 3.06	1.1 / 1.98	1.2 / 2.16	1.1 / 1.98	1.1 / 1.98	1.7 / 3.06	1.9 / 3.42	1.1 / 1.98
at 250°C / 482°F	1.7 / 3.06	1.7 / 3.06	2.6 / 4.68	1.5 / 34.7	1.7 / 3.06	1.6 / 2.88	1.6 / 2.88	2.5 / 4.5	2.7 / 4.86	1.6 / 2.88
at 300°C / 572°F	-	-	-	-	-	2.4 / 4.32	2.4 / 4.32	3.4 / 6.12	3.4 / 6.12	2.4 / 4.32
at 350°C / 662°F	-	-	-	-	-	3.4 / 6.12	3.4 / 6.12	3.4 / 6.12	4.1 / 7.38	3.2 / 5.76
Heating time (min, intake/damper 100% close)²⁾										
from Amb. to 100°C / 212°F	34	34	26	19	29	18	18	14	13	20
from Amb. to 150°C / 302°F	35	35	32	30	35	22	22	19	17	23
from Amb. to 200°C / 392°F	36	36	36	40	44	24	24	22	24	25
from Amb. to 250°C / 482°F	51	51	51	55	60	27	27	29	25	32
from Amb. to 300°C / 572°F	-	-	-	-	-	38	38	36	37	38
from Amb. to 350°C / 662°F	-	-	-	-	-	41	41	43	40	45
Cooling time (min, intake/damper 100% open)²⁾										
from 150°C / 302°F to 100°C / 212°F	31 / 87.8 ³⁾	31 / 87.8 ³⁾	21 / 69.8	18 / 64.4 ³⁾	31 / 87.8	19 / 66.2	19 / 66.2	9 / 48.2	11 / 51.8	15 / 59
from 200°C / 392°F to 100°C / 212°F	45 / 113 ³⁾	45 / 113 ³⁾	32 / 89.6	22 / 71.6 ³⁾	48 / 118.4	28 / 82.4	28 / 82.4	15 / 59	15 / 59	54 / 129.2
from 250°C / 482°F to 100°C / 212°F	54 / 129.2 ³⁾	54 / 129.2 ³⁾	34 / 93.2	34 / 93.2 ³⁾	62 / 143.6	41 / 105.8	41 / 105.8	21 / 69.8	16 / 60.8	59 / 138.2
from 300°C / 572°F to 100°C / 212°F	-	-	-	-	-	33 / 91.4	33 / 91.4	26 / 78.8	22 / 71.6	62 / 143.6
from 350°C / 662°F to 100°C / 212°F	-	-	-	-	-	43 / 109.4	43 / 109.4	30 / 86	27 / 80.6	62 / 143.6
Recovery time after door was opened (min.)²⁾										
at 100°C / 212°F	10.5	10.5	8	6.5	8	8.5	8.5	6	6	7.5
at 150°C / 302°F	8	8	9	8	9.5	8.5	8.5	9	7.5	9.5
at 200°C / 392°F	4.5	4.5	8.5	7	8.5	3.5	3.5	8.5	8.5	8.5
at 250°C / 482°F	5	5	9.5	5.5	9.5	5	5	11.5	8.5	9.5
at 300°C / 572°F	-	-	-	-	-	7.5	7.5	12.5	8.5	9
at 350°C / 662°F	-	-	-	-	-	8.5	8.5	12.5	10	10
Air circulation data										
Air change rate (approx, x/h)	213	266	173	233	193	213	266	173	233	193
Air circulation (approx, x/h)	374	590	294	348	220	374	590	294	348	220
Exhaust air volume flow (approx, L/min)	380	1240	1400	1910	1560	380	1240	1400	1910	1560
Air flow velocity (m/s)	0.2 to 0.6	0.2 to 0.6	0.2 to 0.6	0.2 to 0.6	0.2 to 0.6	0.2 to 0.6	0.2 to 0.6	0.2 to 0.6	0.2 to 0.6	0.2 to 0.6
Electrical data & Ordering information										
Electrical requirement (230V, 50/60Hz, 1P, A)	14	21.3				27.6				
Cat. No.	AAHK1012K	AAHK1022K	-	-	-	AAHK2012K	-	-	-	-
Electrical requirement (380V, 50/60Hz, 3P, A)			9.4	12	14.2		13.9	18.5	23.4	27.9
Cat. No.	-	-	AAHK1038K	AAHK1048K	AAHK1058K	-	AAHK2028K	AAHK2038K	AAHK2048K	AAHK2058K


1) According to IEC 60068-3-5

2) According to DIN 12880

3) Cooling time (min, Intake 50% close, Damper 100% open, According to DIN 12880)

※ Above specification value is recorded by 230V/60Hz, 380V/50Hz.

※ Above specifications can be changed without prior notice.



Dimension

Model	Max. 250°C					Max. 350°C				
	LBV-012	LBV-025	LBV-040	LBV-070	LBV-100	LTV-012	LTV-025	LTV-040	LTV-070	LTV-100
Interior dimension										
Chamber volume (L / cu ft)	125 / 4.4	253 / 8.9	420 / 14.8	720 / 25.4	1000 / 35.3	125 / 4.4	253 / 8.9	420 / 14.8	720 / 25.4	1000 / 35.3
Width (mm / inch)	500 / 19.7	600 / 23.6	750 / 29.5	900 / 35.4	1000 / 39.4	500 / 19.7	600 / 23.6	750 / 29.5	900 / 35.4	1000 / 39.4
Depth (mm / inch)	500 / 19.7	650 / 25.6	700 / 27.6	800 / 31.5	910 / 35.8	500 / 19.7	650 / 25.6	700 / 27.6	800 / 31.5	910 / 35.8
Height (mm / inch)	500 / 19.7	650 / 25.6	800 / 31.5	1000 / 39.4	1100 / 43.3	500 / 19.7	650 / 25.6	800 / 31.5	1000 / 39.4	1100 / 43.3
Quantity of shelves (standard/max.)	2 / 6	2 / 8	2 / 11	2 / 14	2 / 16	2 / 6	2 / 8	2 / 11	2 / 14	2 / 16
Distance of between shelves (mm / inch)	30 / 1.2	30 / 1.2	30 / 1.2	30 / 1.2	30 / 1.2	30 / 1.2	30 / 1.2	30 / 1.2	30 / 1.2	30 / 1.2
Max. load per shelf (kg / lbs)	15 / 33.1	20 / 44.1	30 / 66.1	40 / 88.2	45 / 99.2	15 / 33.1	20 / 44.1	30 / 66.1	40 / 88.2	45 / 99.2
Permitted total load (kg / lbs)	50 / 110.2	70 / 154.3	90 / 198.4	120 / 264.6	150 / 330.7	50 / 110.2	70 / 154.3	90 / 198.4	120 / 264.6	150 / 330.7
Damper (Ø, mm / inch)	96.4 / 3.8	146.4 / 5.8	146.4 / 5.8	146.4 / 5.8	146.4 / 5.8	96.4 / 3.8	146.4 / 5.8	146.4 / 5.8	146.4 / 5.8	146.4 / 5.8
Air Intake (Ø, mm / inch)	42.8 / 1.7	73.3 / 2.9	73.3 / 2.9	73.3 / 2.9	73.3 / 2.9	42.8 / 1.7	73.3 / 2.9	73.3 / 2.9	73.3 / 2.9	73.3 / 2.9
Exterior dimension										
Width (mm / inch), W	1100 / 43.3	1200 / 47.2	1350 / 53.1	1500 / 59.1	1600 / 63	1100 / 43.3	1200 / 47.2	1350 / 53.1	1500 / 59.1	1600 / 63
Depth / with handle (mm-inch), D/D'	714 / 745 28.1 / 29.3	965 / 995 38 / 39.2	1014 / 1045 39.9 / 41.1	1114 / 1145 43.9 / 45.1	1224 / 1255 48.2 / 49.4	714 / 745 28.1 / 29.3	965 / 995 38 / 39.2	1014 / 1045 39.9 / 41.1	1114 / 1145 43.9 / 45.1	1224 / 1255 48.2 / 49.4
Height / with damper (mm-inch), H/H'	1240 / 1362 48.8 / 53.6	1435 / 1559 56.5 / 61.4	1590 / 1712 62.6 / 67.4	1790 / 1912 70.5 / 75.3	1890 / 2012 74.4 / 79.2	1240 / 1362 48.8 / 53.6	1435 / 1559 56.5 / 61.4	1590 / 1712 62.6 / 67.4	1790 / 1912 70.5 / 75.3	1890 / 2012 74.4 / 79.2
Weight (kg / lbs)	170 / 374.8	260 / 573.2	320 / 705.5	390 / 859.8	440 / 970	170 / 374.8	260 / 573.2	320 / 705.5	390 / 859.8	440 / 970

Accessories

Model	LBV-012	LBV-025	LBV-040	LBV-070	LBV-100
	LTV-012	LTV-025	LTV-040	LTV-070	LTV-100
Wire shelf	00LTV0000070	00LTV0000053	00LTV0000029	00LTV0000030	00LTV0000057
Perforated shelf (Heavy load shelf)	AAAK1501	AAAK1502	AAAK1503	AAAK1504	AAAK1505
Viewing window (LBV Only)	AAAK1511		AAAK1512	AAAK1513	
Cable port (Φ80), SUS	AAAK1521				
Cable port (Φ50), SUS	AAAK1522				
Warning signal lamp	AAA80550				
Digital recorder (6 Channel)	AAA8T505				
Fan speed adjuster	AAAK1531				
Gas purge system	AAA80691				
Hour meter	AAAK1532				



Gas Purge System

Gas system that can replace and purge the inside of chamber with nitrogen and CO₂ gas.

APPENDIX

JeioTech 2025-2026 Catalog

A	
Acid/Corrosive Cabinet	236~239
Arm Hood	216
Autoclave (Vertical Type)	184
Aspirator	189
Aspirator [Water jet type]	189
Electric Aspirator	189

B	
Baths & Circulator	24
Chiller	40~47
Cold Trap Bath	57
General Heating Bath	54
Heating Bath Circulator	32
Heating Bath Circulator, Programmable	32
Refrigerated & Heating Bath Circulator, Programmable	28
Shaking Heating Bath	60
Thermal Bath	54~61
Thermal Bath Circulator	30
Ultrasonic Cleaner	190
Visible Bath	58
Water Bath	54~61
Water Bath Circulator	28, 32
Bench Top Fume Hood	210
Beaker	270~273
Biosafety Cabinet	196

C	
Cabinet	220
Acid/Corrosive Cabinet	236~240
Acid/Corrosive Cabinet-Small type	235
Biosafety Cabinet	196
Carboy	273
Clean Bench	200~202
Ductless Fume Hood	208
Electrophoresis System	178~183
Low temperature Incubator	100~103
Laboratory Refrigerator	260~261

Laminar Flow Cabinet	200~202
Mini Incubator	98
PCR Workstation	203
Pharmacy Refrigerator	264
Filtering Storage Cabinet	240
UV Sterilization Cabinet	204
Cell Culture	
Clean Bench	200~202
Heating Shaker	175
On-ice Rack	276
Incubator	90~103
Incubated Shaker	112~119
Low temperature Incubator	100~103
Laboratory Refrigerator	260~261
Laboratory Freezer	262
Mini Incubator	98
Open-air Shaker	126~133
Stirrer	148~151
Rotator	134
Chamber	
Environmental Test Chamber	286~311
Growth Chamber	104
Incubator	90~103
Laboratory Refrigerator	260~261
Laboratory Freezer	262
Mini Incubator	98
Oven	66~81
Pharmacy Refrigerator	264
Pharmaceutical Stability Chamber	305
Plant Growth Chamber	104
Temperature & Humidity Chamber	286~311
Temperature & Humidity Chamber, Dual	299
Chiller	36~49
Chiller	36~49
Circulator	
Heating Bath Circulator	32
Heating Bath Circulator, Programmable	32
Refrigerated & Heating Bath Circulator	30
Refrigerated & Heating Bath Circulator, Programmable	28

Clean Bench	200~205
Clean Bench	200~202
PCR Workstation	203
UV Sterilization Cabinet	204
Ultrasonic Cleaner	190
Cold Trap Bath	57

D	
Desiccator	246~253
Desiccator	252
Vacuum Desiccator (Cubic)	250
Vacuum Desiccator (Round)	250
Vacuum Desiccator (Cylindrical)	250
Drying	
Forced Convection Oven, Program	66
Forced Convection Oven, Advanced	68
Forced Convection Oven, Standard	72
Forced Convection Oven, Large	70
Forced Convection Oven, 4 Chambers	73
Natural Convection Oven, Advanced	74
Natural Convection Oven, Standard	76
Vacuum Oven	78
Ductless Fume Hood	208
Dry Cabinet	248

E	
Enclosure	
Biosafety Cabinet	196
Clean Bench	200~202
Electrophoresis System	178~183
Hood	208~219
PCR Workstation	203
UV Sterilization Cabinet	204
Fume Hood-Small type	210
Environmental Test (Chamber)	286~311
Temperature & Humidity Chamber Customized type	311

Index

Pharmaceutical Stability Chamber	305~307
Temperature & Humidity Chamber, Dual	299
Temperature & Humidity Chamber	286~311

F

Flask Support	276
Forced Convection	66~73
Forced Convection Oven, Program	66
Forced Convection Oven, Advanced	68
Forced Convection Oven, Standard	72
Forced Convection Oven, Large	70
Forced Convection Oven, 4 Chambers	73
Laboratory Refrigerator	260~261
Laboratory Freezer	262
Low Temperature Incubator (IL3)	100
Low Temp. Incubator, Multi Chambers	103
Low Temperature Incubator, Personal	102
Shaking Incubator	112~119
Pharmacy Refrigerator	305~307
Plant Growth Chamber	104
Freezer	262~263
Laboratory Freezer	262~263
Fume Hood	206~219
Ductless Fume Hood	208
Fume Hood-Small type	210
Mobile Fume Extractor, Local Exhaust	212
Mobile Fume Extractor, Air Purifying	213
Funnel	
Funnel Shaker	133
Furniture	220~245
Storage Cabinet	220~245

Filtering Storage Cabinet	240
Furnace	82
Muffle Furnace	82

G

Gravity Convection	
Natural Convection Oven, Advanced	74
Natural Convection Oven, Standard	76
Air-jacket Incubator, Advanced	90
Air-jacket Incubator, Standard	92
Low Temperature Incubator (IL3-A)	100
Growth Chamber	104
Plant Growth Chamber	104

H

Heating	
Bath	54~61
Bath Circulator	28~33
Heating Block	172~173
Heating & Cooling Block	174
Hotplate	145
Hotplate & Magnetic Stirrer	142~144
Incubator	90~103
Mini Incubator	98
Oven	66~81
Heating Shaker	175
Incubated Shaker	112~119
Temperature Chamber	312~349
HEPA Filtered	
Clean Bench	200~202
Temperature & Humidity Chamber, Customized type	311
Ductless Fume Hood	208
Arm Hood	216
PCR Workstation	203
Hood	
Arm Hoods	216~219
Clean Bench	200~202

Ductless Fume Hood	208
Fume Hood	208~219
Fume Hood-Small type	210
PCR Workstation	203
UV Sterilization Cabinet	204
Hotplate	142~145
Hotplate	145
Hotplate & Magnetic Stirrer	142~144
Humidity	
Plant Growth Chamber	104
Temperature & Humidity Chamber	286~311

I

ICH (guideline) Chamber	305
Incubator	90~103
Air-jacket Incubator	90~93, 100
Bench Top Shaking Incubator	112
General Purpose Incubator	90~99
Low Temperature Incubator	100~103
Low Temperature Incubator, Multi Chambers	103
Mini Incubator	98
Personal Low Temperature Incubator	102
Shaking Incubator	112~121
Shaking Incubator, Stackable type	112, 114
Shaking Incubator, Chamber type	112
Plant Growth Chamber	104

L

LC Connected	10
LC GreenBox	10
LC Datakeeper	8

M

Magnetic	
Hotplate & Magnetic Stirrer	142~144

Magnetic Stirrer	148~151
Multi-Position Magnetic Stirrer	144, 150
Medical Cooling	
Pharmacy Refrigerators	264
Mechanical Convection	
Forced Convection Oven, Program	66
Forced Convection Oven, Advanced	68
Forced Convection Oven, Standard	72
Forced Convection Oven, Large	70
Forced Convection Oven, 4 Chambers	73
Laboratory Refrigerator	260
Laboratory Freezer	262
Low Temperature Incubator (IL3)	100
Low Temp. Incubator, Multi Chambers	103
Low Temperature Incubator, Personal	102
Pharmacy Refrigerator	264
Shaking Incubator	112~121
Mixer	
Overhead Stirrer	156~165
Vortex Mixer	166
Muffle Furnace	82

N	
Natural Convection	
Natural(Gravity) Convection Oven, Advanced	74
Natural(Gravity) Convection Oven, Standard	76
Advanced Incubator	90
Standard Incubator	92
Low Temperature Incubator (IL3-A)	100
Network	
LC Connected	10
LC GreenBox	10
LC Datakeeper	8

O	
Orbital Shaker	
Open-air Shaker	126~133
Heating Shaker	175
Mini Shaker	130
Microplate Shaker	131
Shaking Incubator	112~121
Oven	
Forced Convection Oven	66~73
Forced Convection Oven, 4 Chambers	73
Gravity Convection Oven	74~77
High Temperature Oven	70, 344
Large Forced Convection Oven	70
Mechanical Convection Oven	66~73
Natural Convection Oven	74~77
Temperature Chamber	282~311
Vacuum Oven	78
Overhead Stirrer	156~165

P	
PCR	
PCR Workstation	203
UV Sterilization Cabinet	204
Pharmaceutical Stability Chamber	305
Pharmacy Refrigerator	264
Plant Growth Chamber	104
Plasticware	
Beaker	270~272
Flask Support	276
Cryo Box	275
Storage Box	274
Tubing connector	278
Pump	188
Aspirator	188
Rotary Vacuum Pump	188

R	
Reciprocating Shaker	
Open-air Shaker	126~133

Funnel Shaker	133
Reciprocating Shaking Bath	60
Shaking Heating Bath	60
Shaking Incubator	112~119
Recirculating Cooler	36~49
Refrigerated	
Laboratory Refrigerator	260
Low Temp. Incubator	100~103
Low Temp. Incubator, Multi Chambers	103
Low Temp. Incubator, Personal	102
Refrigerated & Heating Bath Circulator	30
Refrigerating & Heating Bath Circulator, Programmable	28
Shaking Incubator	112~119
Rocking Shaker	132
Rotary Vacuum Pump	189
Rotator	134

S	
Safety	
Flask Support	276
On-ice Rack	276
Shaker	122~137
Open-air Shaker	126~133
Funnel Shaker	133
Heating Shaker	175
Mini Shaker	130
Mixer	142~169
Microplate Shaker	131
Orbital/Reciprocating Shaker	126
Reciprocating Shaking Bath	60
Rocking Shaker	132
Shaking Incubator	112~119
Shaking Water Bath	60
Waving Shaker	132
Stirrer	142~165
Magnetic Stirrer	148~151
Multi-Position Magnetic Stirrer	144, 150
Hotplate & Magnetic Stirrer	142~144
Overhead Stirrer	158~165

Index

Stirring Hotplate	142~144
-------------------	---------

T

Table top	
Acid/Corrosive Cabinet-Small type	235
Bath	24~35, 50~61
Bath & Circulator	24~35
Circulator	24~35
Clean Bench (Vertical type Laminar Flow)	200~202
Desiccator	246~253
Electric Aspirator	189
Fume Hood-Small type	210~211
Heating Shaker	175
Hotplate Stirrer	143~144
Incubator	90~103
Magnetic Stirrer	148~155
Oven	66~81
Overhead Stirrer	158~165
Pump	188
Open-air Shaker	126~133
Shaking Incubator	112~119
Temperature Test Chamber	312~349
Temperature & Humidity Chamber	282~311
Temperature & Humidity Chamber	282~311
Temperature & Humidity Chamber, Dual	299
Floor Standing Model	290~309
Table top	286~289
Pharmaceutical Stability Chamber	305
Plant Growth Chamber	104
Thermal	
Thermal Bath	54~61
Thermal Bath Circulator	30
Tissue Culture	
Freezer	262
Heating Shaker	175
On-ice Rack	276
Incubated Shaker	112~119
Incubators	90~103
Open-air Shaker	126~133

Stirrer	142~165
Mixer	142~169
Multi-Position Magnetic Stirrer	144, 150
Orbital/Reciprocating Shaker	126
Overhead Stirrer	158~165
Rocking Shaker	132
Waving Shaker	132
Tubing connector	278

Ductless Fume Hood	208
Fume Hood-Small type	210
PCR Workstation	203
UV Sterilization Cabinet	204

U

Ultrasonic	
Ultrasonic Cleaner	190

V

Vacuum	
Cold Trap Bath	57
Electric Aspirator	188
Desiccator	246~253
Oven	66~81
Pump	188
Rotary Pump	188
Viscosity	
Visible Bath	58
Vortex Mixer	168


W

Washer	
Ultrasonic Cleaner	188
Water Bath	
Chiller	36~49
Circulator	24~35
General Bath	54
Shaking Water Bath	60
Waving Shaker	132
Workstation	
Bench Top Fume Hood	210~211
Biosafety Cabinet	196
Clean Bench	200~202

Build a New Lab!

Do you need help to set your new lab?
or Do you have a plan to move the lab?

Jeio Tech is your Lab Companion that provides what you need for your lab.
We always think about how we can help researchers to make those successful results.



One stop purchasing

You can find all what you need at once Jeio Tech.



Laboratory Equipment

- Oven & Incubator
- Furnace
- Shaker & Incubated Shaker
- Bath Circulator & Heating Bath
- Chiller (Recirculating Cooler)
- Stirrer & Mixer
- Autoclave



Biosafety Cabinet & Clean Bench, Fume Hood

- Biosafety Cabinet
- Clean Bench
- PCR Workstation
- UV Sterilization Cabinet
- Ductless Fume Hood
- Mobile Fume Extractor



Safety Storage & Cabinet

- Fire Safety Storage Cabinet
- Corrosive Cabinet
- Filtering Storage Cabinet
- Gas Cylinder Cabinet
- Multiple Drawer Cabinet
- Dry Cabinet & Desiccator



Refrigerator & Freezer

- Laboratory Refrigerator
- Laboratory Freezer
- Pharmaceutical Refrigerator



Lab Essentials

- Plastic Ware
- Safety Product



Environmental Test Equipment

- Temp. & Humid. Chamber
- Heating & Cooling Chamber
- Chamber
- Heating Chamber





Laboratory Equipment

Oven & Incubator
Furnace
Shaker & Incubated Shaker
Bath Circulator & Heating Bath
Chiller (Recirculating Cooler)
Stirrer & Mixer
Autoclave



Safety Storage & Cabinet

Fire Safety Storage Cabinet
Corrosive Cabinet
Filtering Storage Cabinet
Gas Cylinder Cabinet
Multiple Drawer Cabinet
Dry Cabinet & Desiccator



Biosafety Cabinet & Clean Bench, Fume Hood

Biosafety Cabinet
Clean Bench
PCR Workstation
UV Sterilization Cabinet
Ductless Fume Hood
Mobile Fume Extractor



Environmental Reliability Test Equipment

Temp. & Humid. Chamber
Heating & Cooling Chamber
Heating Chamber



Refrigerator & Freezer for Laboratory/Medical

Laboratory Refrigerator
Laboratory Freezer
Pharmaceutical Refrigerator



Lab Essentials

Plastic Ware
Safety Product

Jeio Tech

Head office & factory

Add. 153, Techno 2-ro, Yuseong-gu, Daejeon, 34025, South Korea

International sales office

Add. 10F-1005, 219, Gasandigital 1-ro, Geumcheon-gu, Seoul, 08501, South Korea

Tel. +82 2 2627 3816

Fax +82 2 3143 1824

E-mail overseas@jeiotech.com

Local offices in

Seoul, Suwon, Daejeon, Daegu, Busan

Jeio Tech, Inc. (US Branch)

19 Alexander Road, Ste. 7, Billerica, MA 01821-5094, U.S.A.

Tel +1 781 376 0700 Fax +1 781 376 0704

E-mail info@jeiotech.com

Jeio Tech Japan株式会社 (Japan Branch)

Luminas Minamidaira, Room 201, 2-70-60 Minamidaira, Hino City, Tokyo, JAPAN postal code 191-0041

Tel +81 42 506 751 Fax +81 42 506 7520