

MODEL U40-25

UPRIGHT ULTRA LOW FREEZER 0°C TO -40°C

Capacity [Cu. Ft.] 25 Cu. Ft.
Capacity [Liters] 710 Liters
Temperature Range [°C] 0°C to -40°C
Temperature Range [°F] 32°F to -40°F

Dimensions Width x Depth x Height
Exterior [Inches] 48.75" x 37" x 79.5"
Interior [Inches] 34" x 25" x 51.25"
Exterior [Centimeters] 123 x 93 x 201
Interior [Centimeters] 86 x 63 x 130

Refrigeration System Single Stage System
Compressors 1, Hermetic Compressor
Air Circulation Gravity Flow
Defrost Method Manual
Refrigerant HCFC & CFC Free
Exterior Powder Coated Steel
Interior Zinc Coated Galvanized Steel
Insulation Polyurethane (Ecomate) Foam
Inner Doors 5, Stainless Steel
Lockable Door Lockable with Padlock
Porthole Optional 1/2"
Casters 6, Swivel with Locks
Storage 4 Adjustable Shelves, with
1 Base Shelf
Shipping Weight 1125 Lbs.

Electrical

Standard Voltage 115 Volts, 60 Hertz, 1 Phase
Amps / Run Amps / Plug 20A | 7.3A | NEMA 5-20
(Optional Voltages) 208 Volts, 230 Volts
Certification(s) ETL Listed

Controller

Display Type Microprocessor, Digital
Alarms High/Low Temp, Power
Failure, Low Battery
Alarm Relay Yes, Dry Contacts
Automatic Restart Yes, Automatically Restarts
After Power Failure
Backup Battery Yes, Only Powers Alarm

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 1 Year Parts & Labor