




## Tube adapter

for 19400045 rotor



**Cat. No.** 19400041  
Swing-out rotor for DM0636 Centrifuge.  
(Order buckets/tubes separately)

Max. speed: 4000rpm  
Max. RCF: 2950×g  
Rotor material: stainless steel body and adapter



**Cat. No.** 19500025

15ml conical tube adapter, MOQ 8 Pcs (optional)



19400045



19400044



19400043



19400042

bucket	x 4										
tube	x 36		x 16			x 8	x 4	x 4	x 4	x 8	x 20
Nominal volume	3-7mL		10mL	12mL	15mL	50mL	100mL	85mL	50mL	15mL	3-7mL
φ×L mm	13×75	13×100	16×100	17×105	17×120	29.5×100	38×108	38×104	29.5×100	17×120 13×100 13×75	13×100 13×75
tube											
adapter	with plug					 A50V	 A100P80	 A100P50	 A50V	 A100P15	 A100P10
Package	36					8	4	4	4	8	4
Tubes per rotor	36	36	16	16	16	8	4	4	4	8	20
Radius(cm)	14	16.5	15.9	15.9	15.9	15.9	15.9	15.9	15.9	15.1	14.7
Max RCF (g-force)	2500	2950	2840	2840	2840	2840	2840	2840	2840	2700	2630