

## MODEL LRBFS12W1HC

### UNDERCOUNTER FMS REFRIGERATOR 1°C TO 15°C

<b>Capacity [ Cu. Ft. ]</b>	<b>11.2 Cubic Feet</b>
<b>Capacity [ Liters ]</b>	<b>318 Liters</b>
<b>Temperature Range [ °C ]</b>	<b>1°C to 15°C</b>
<b>Temperature Range [ °F ]</b>	<b>34°F to 59°F</b>
Dimensions	Width x Depth x Height
Exterior [ Inches ]	23.75" x 24.5" x 64.75"
Interior [ Inches ]	18.5" x 17.25" x 57"
Exterior [ Centimeters ]	59 x 61 x 164
Interior [ Centimeters ]	47 x 44 x 145
Compressor	Hermetic Compressor
Air Circulation	Forced Air
Defrost Method	Automatic
Primary Refrigerant(s)	R290
Exterior	Steel, Painted White
Interior	Molded Polystyrol
Insulation	Polyurethane Foam
Door	1 Door
Lockable Door	Keyed Lock
Door Material	Insulated Steel
Leveling Legs	4
Storage	5 Adjustable Glass Shelves
Shelving Weight Limit	40 Lbs.
Approximate Decibels	50 dB

#### Electrical:

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line   Running	15 Amp Dedicated   3.0A
NEMA Configuration	NEMA 5-15
Certification(s)	ETL Listed

#### Control:

Temperature Control	Analog
Temperature Display	Digital ( External )

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Interior is free from sources of ignition and meets the requirements for Class I, Division 2, for storage of flammable substances in closed containers. Product complies with the application requirements of NFPA 45 and NFPA 99.

**Warranty** 3 Years Parts & Labor