

MODEL MV4-6UCXR

HAZARDOUS LOCATION REFRIGERATOR 1°C TO 10°C

Capacity [Cu. Ft.]	5 Cu. Ft.
Capacity [Liters]	141 Liters
Temperature Range [°C]	1°C to 10°C
Temperature Range [°F]	34°F to 50°F
Dimensions	Width x Depth x Height
Exterior [Inches]	23.75" x 23.75" x 32.25"
Interior [Inches]	20.25" x 17.25" x 27"
Compressor	Hermetic Compressor
Air Circulation	Direct Cooling
Defrost Method	Cycle
Refrigerant	R600a
Exterior	Painted Steel
Insulation	EPA compliant H.D.U.F.
Lockable Door	Yes, Keyed
Porthole	None
Casters	No
Leveling Legs	2, Front
Storage	3 Adjustable Shelves
Shipping Weight	130 lbs.

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Electrical

Voltage:	115 Volts 60 Hertz 1 Phase Uses ATEX and non-sparking components
Amperage Line Running NEMA Configuration	15 Amp Dedicated 1.2A NONE Exposed Flying Leads, 3 Wires This unit comes without a plug and requires hardwiring in conduit, ensuring a secure connection in your high-risk environment
Certification(s)	ETL Listed Meets NFPA and OSHA for hazardous locations, Class 1, Division 2, Group C & D



Controller

Display Type	Microprocessor, Digital
Alarms	High/Low Temp
Min / Max Temperature	Yes, Display and Reset
Alarm Relay	Yes, Dry Contacts
Backup Battery	No

Warranty 2 Year Parts & Labor