

MODEL DH29-46SD

LAB AND PHARMACY FREEZER -15°C TO -23°C

Capacity [Cu. Ft.]	45 Cu. Ft.
Capacity [Liters]	1280 Liters
Temperature Range [°C]	-15°C to -23°C
Temperature Range [°F]	5°F to -13°F
Dimensions	Width x Depth x Height
Exterior [Inches]	52" x 34" x 82"
Interior [Inches]	48" x 25.75" x 62"
Exterior [Centimeters]	133 x 87 x 209
Interior [Centimeters]	121 x 65 x 157
Compressor	Hermetic Compressor
Air Circulation	Forced Air
Defrost Method	Adaptive defrost
Refrigerant	R290
Exterior	Stainless Steel
Interior	Aluminum
Insulation	Polyurethane Foam
Lockable Door	2 Doors, Keyed Lock
Door Material	Solid Stainless
Interior Light	LED
Porthole	1/2" Diameter Left Hand Side
Casters	4, Locking
Storage	6 Adjustable Shelves
Shipping Weight	525 lbs.

Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line Running	15 Amp Dedicated 10.0A
NEMA Configuration	NEMA 5-15
Certification(s)	UL Listed

Controller

Display Type	Microprocessor, Digital
Alarms	High/Low Temp
Alarm Relay	Dry Contacts

Warranty 3 Years Parts & Labor.

Additional 4 Years on Compressor Parts, No Labor.

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES

