

## MODEL MV20-2UCF

### LAB AND PHARMACY FREEZER -10°C TO -25°C

<b>Capacity [ Cu. Ft. ]</b>	<b>1.7 Cu. Ft.</b>
<b>Capacity [ Liters ]</b>	<b>48 Liters</b>
<b>Temperature Range [ °C ]</b>	<b>-15°C to -25°C</b>
<b>Temperature Range [ °F ]</b>	<b>5°F to -13°F</b>
Dimensions	Width x Depth x Height
Exterior [ Inches ]	23.75" x 24" x 20.5"
Interior [ Inches ]	17.75" x 13.75" x 16"
Exterior [ Centimeters ]	74 x 61 x 53
Interior [ Centimeters ]	45 x 34 x 40
Compressor	Hermetic Compressor
Air Circulation	Direct Cooling
Defrost Method	Manual
Refrigerant	R600a
Exterior	Painted Steel
Interior	Painted Aluminum
Insulation	US EPA and SNAP approved
Lockable Door	Yes, Keyed
Porthole	Yes, 3/8" Diameter
Casters	None
Leveling Legs	Yes, 4
Storage	2 Adjustable Shelves, Base Shelf
Shipping Weight	115 lbs.


#### Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line   Running	15 Amp Dedicated   1A
NEMA Configuration	NEMA 5-15
Certification(s)	UL Listed, C-UL, Energy Star

#### Controller

Display Type	Microprocessor, Digital
Alarms	High/Low Temp, Door Ajar Power Failure, Low Battery Sensor Failure, USB Failure
Data Logging	Yes, Adjustable Intervals
Data Download	N/A
Min / Max Temperature	Yes, Display and Reset
Alarm Relay	Yes, Dry Contacts
Backup Battery	N/A

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



**Warranty** 2 Year Parts & Labor