

ENGINEERS AND MANUFACTURERS OF ULTRA-LOW FREEZERS

## *SO-LOW EXPLOSION PROOF REFRIGERATOR MODEL DHH4-30SDRX MANUAL DEFROST*

*Meets NFPA and OSHA guidelines for Hazardous Locations, Class I, Division II, Group C & D.*

### **30 Cubic Ft. Upright +4°C / Solid Door**

#### **Specifications:**

* Capacity	30 cubic ft. / 849 liters
* Temperature	1°C to 10°C
* Dimensions	
* Exterior w x d x h	35 ¾" x 35" x 77 ½" 90.81 x 88.9 x 196.85 (cm)
* Interior w x d x h	32" x 22.5" x 65" 81.28 x 57.15 x 165.1 (cm)

#### **Features:**

*Thermostat	Exterior Adjustable Mechanical
*Doors	Single solid – key lock
*Legs	Four leveling
*Shelves	5 adjustable
*Interior Lighting	None
*Insulation:	Polyurethane foam
*Refrigerant:	R-134A
*Compressor	1/3 hp
*Air Circulation	Gravity
*Evaporator	Hot wall
*Defrost	Manual
*Voltage	120 VAC
*Running Amps	6
*Supply Plug	None / Unit must be hardwired
*Shipping weight	393 lbs.
*Approvals	UL Listed
*Warranty-(Domestic)	One year parts and labor.



Specifications subject to change without notice