

MYSTAIRE FILTRATION CHART (p.1 of 7)

MYSTAIRE®	Unit	Pre-filter Required	Main Filter Options	Safety Filter Options	Number of Filters Required
Isola Filtered Workstation					
	MY-ISL36 (-010) MY-ISL36P (-010) MY-ISL36P2 (-010)	MY-PRE-24	MY-4004	MY-5800	1
			MY-4007		1
			MY-4008		1
			MY-4009		1
			MY-4011		1
			MY-4014		1
			MY-4016		1
			MY-4017		1
			MY-4020		1
			MY-4016		1
	MY-ISL48 (-010) MY-ISL48P (-010) MY-ISL48P2 (-010)	MY-PRE-20	MY-4104	MY-6500	2
			MY-4107		2
			MY-4108		2
			MY-4109		2
			MY-4111		2
			MY-4114		2
			MY-4116		2
			MY-4117		2
			MY-4120		2
			MY-4116		2
	MY-ISL72 (-010) MY-ISL72P (-010) MY-ISL72P2 (-010)	MY-PRE-24	MY-4004	MY-5800	2
			MY-4007		2
			MY-4008		2
			MY-4009		2
			MY-4011		2
			MY-4014		2
			MY-4016		2
			MY-4017		2
			MY-4020		2
			MY-4016		2

MYSTAIRE FILTRATION CHART (p.2 of 7)

MYSTAIRE®	Unit	Pre-filter Required	Main Filter Options	Safety Filter Options	Number of Filters Required
Latitude Balance Enclosure					
	MY-LBE 48 MY-LFH 48	MY-PRE-24	MY-5904	N/A	2
			MY-5907		2
			MY-5908		2
			MY-5909		2
			MY-5911		2
			MY-5914		2
			MY-5916		2
			MY-5917		2
			MY-5920		2
	MY-LBE 72 MY-LFH 72	MY-PRE-36	MY-5904	N/A	3
			MY-5907		3
			MY-5908		3
			MY-5909		3
			MY-5911		3
			MY-5914		3
			MY-5916		3
			MY-5917		3
			MY-5920		3

MYSTAIRE FILTRATION CHART (p.3 of 7)

MYSTAIRE®	Unit	Pre-filter Required	Main Filter Options	Safety Filter Options	Number of Filters Required
Ductless Fume Hood - Aura Elite					
	MY-AU30 (-010)	MY-PRE24	MY-3804	MY-5804	1
			MY-3807	MY-5807	1
			MY-3808	MY-5808	1
			MY-3809	MY-5809	1
			MY-3811	MY-5811	1
			MY-3814	MY-5814	1
			MY-3816	MY-5816	1
			MY-3817	MY-5817	1
			MY-3820	MY-5820	1
	MY-AU30S (-010)	MY-PRE24	MY-3816	MY-5817	1
	MY-AU42 (-010)	MY-PRE24	MY-2804	MY-5904	2
			MY-2807	MY-5907	2
			MY-2808	MY-5908	2
			MY-2809	MY-5909	2
			MY-2811	MY-5911	2
			MY-2814	MY-5914	2
			MY-2816	MY-5916	2
			MY-2817	MY-5917	2
			MY-2820	MY-5920	2
	MY-AU42S (-010)	MY-PRE24	MY-2816	MY-5917	2
	MY-AU54 (-010)	MY-PRE24	MY-3804	MY-5804	2
			MY-3807	MY-5807	2
			MY-3808	MY-5808	2
			MY-3809	MY-5809	2
			MY-3811	MY-5811	2
			MY-3814	MY-5814	2
			MY-3816	MY-5816	2
			MY-3817	MY-5817	2
			MY-3820	MY-5820	2
	MY-AU54S (-010)	MY-PRE24	MY-3816	MY-5817	2

MYSTAIRE FILTRATION CHART (p.4 of 7)

MYSTAIRE®

Unit

Pre-filter Required

Main Filter
Required

Safety Filter
Required

Number of Filters
Required

Cyanoacrylate Fuming Chamber (initial filters included)



CA-3000 (-010)

FE-2013

MY-2818

N/A

1



CA-6000 (-010)

FE-2013

MY-2818

N/A

1



CA-9000 (-010)

FE-2013

MY-2818

N/A

2

MYSTAIRE FILTRATION CHART (p.5 of 7)

MYSTAIRE®

Unit

Pre-filter Required

Main Filter
Required

Safety Filter
Required

Number of Filters
Required

Downflow Workstation



FE-WS6 (-010)

FE-1013

FE-5004

N/A

1

FE-5007

N/A

1

FE-5014

N/A

1

FE-5020

N/A

1

PCR Prep Station (initial filters included)



MY-PCR24 (-010)

MY-PRE2

MY-HEPA18

N/A

1

MY-PCR32 (-010)

MY-PRELF

MY-HEPA

N/A

1

MY-PCR48 (-010)

MY-PRE7 (2)

MY-HEPA18 (2)

N/A

2

Ductless Forensic Workstation



FE-36 (-010)

FE-36PF

FE-36CF

N/A

1

FE-36HF

N/A

1

MYSTAIRE FILTRATION CHART (p.6 of 7)

MYSTAIRE®

Unit

Pre-filter Required

Main Filter
Required

Safety Filter
Required

Number of Filters
Required

Evidence Drying Cabinet - SecureDry™ (initial filters included)



MY-FDR32 (-010)

FE-1013 (1)
MY-PRE (1)

MY-1021S
MY-HEPA

N/A

1 HEPA
1 Carbon
1 FE1013
1MYPRE

MY-FDR36 (-010)

FE-1013 (1)
MY-PRE (1)

MY-1021S
MY-HEPA

N/A

1 HEPA
1 Carbon
1 FE1013
1MYPRE

MY-FDR48 (-010)

FE-1013 (1)
MY-PRE (1)

MY-1021S
MY-HEPA

N/A

1 HEPA
1 Carbon
1 FE1013
1MYPRE



MY-FDR64 (-010)

FE-1013 (2)
MY-PRE (2)

MY-1021S
MY-HEPA

N/A

2 HEPA
2 Carbon
2 FE1013
2MYPRE

MY-FDR64D (-010)

FE-1013 (2)
MY-PRE (2)

MY-1021S
MY-HEPA

N/A

2 HEPA
2 Carbon
2 FE1013
2MYPRE

MY-FDR72 (-010)

FE-1013 (2)
MY-PRE (2)

MY-1021S
MY-HEPA

N/A

2 HEPA
2 Carbon
2 FE1013
2MYPRE

MY-FDR72D (-010)

FE-1013 (2)
MY-PRE (2)

MY-1021S
MY-HEPA

N/A

2 HEPA
2 Carbon
2 FE1013
2MYPRE

MYSTAIRE FILTRATION CHART (p.7 of 7)

MYSTAIRE®	Unit	Pre-filter Required	Main Filter Required	Safety Filter Required	Number of Filters Required
Lab Partner Ductless Workstation					
	MY-2400 (-010)	MY-PRE4	MY-1404	N/A	1
			MY-1407	N/A	1
			MY-1408	N/A	1
			MY-1409	N/A	1
			MY-1414	N/A	1
			MY-1417	N/A	1
			MY-1420	N/A	1
	MY-3200 (-010)	MY-PRE1	MY-1104	N/A	1
			MY-1107	N/A	1
			MY-1108	N/A	1
			MY-1109	N/A	1
			MY-1114	N/A	1
			MY-1117	N/A	1
			MY-1120	N/A	1
	MY-4800 (-010)	MY-PRE4	MY-1404	N/A	2
			MY-1407	N/A	2
			MY-1408	N/A	2
			MY-1409	N/A	2
			MY-1414	N/A	2
			MY-1417	N/A	2
			MY-1420	N/A	2

MYSTAIRE®

Phone: +1 (919) 229-8511 • Fax: +1 (919) 528-1744 • www.mystaire.com