

MY-PCR Prep Station - 110V

MYSTAIRE®



Made in the USA



GENERAL SPECIFICATIONS

Construction Material	Work Surface	Polypropylene
	Shell	Polycarbonate
	Front Sash	Polycarbonate
	Filtration Housing	ABS
	Pipette Shelf	Polycarbonate
Air Quality		Class 100 / ISO 5
Controls		ON/OFF Switch
Blowers		Brushless-Sparkless
Lighting		Fluorescent
UV Light		254-nm Shortwave UV Light
UV Light Timer		15-Minute UV Light Timer
HEPA Filter Monitoring Package		Audible / Visual Alarm

DIMENSIONS

Product #	External Dimensions			Internal Dimensions			Sash Opening Dimensions		
	Width	Depth	Height	Width	Depth	Height	Width	Opening 1	Opening 2
MY-PCR24	24"	24"	30"	23"	23"	18"	21"	10"	20"
MY-PCR32	32"	24"	30"	31"	23"	18"	29"	10"	20"
MY-PCR48	48"	24"	30"	47"	23"	18"	45"	10"	20"

TECHNICAL SPECIFICATIONS

Product #	Electrical		Volume of Filtered Air	Prefilters Required	Main Filters Required
	Voltage	Power			
MY-PCR24	110V	130 W	1	1	55 lbs.
MY-PCR32	110V	130 W	1	1	70 lbs.
MY-PCR48	110V	230 W	2	2	95 lbs.

Unit weight is an approximation*

www.mystaire.com

Phone: +1 (919) 229-8511 • Toll Free: (877) 328-3912 • Email: contact@mystaire.com