

CleanPrep™ Dead Air Box - 110V

MYSTAIRE®



Made in the USA



GENERAL SPECIFICATIONS

Construction Material	Work Surface	Polypropylene
	Shell	Polypropylene
	Front Sash	Polycarbonate
	Pipette Shelf	Polycarbonate
Controls		ON/OFF Switch
Lighting		Fluorescent
UV Light		254-nm Shortwave UV Light
UV Light Timer		0-15 Minute UV Light Timer

DIMENSIONS

Product #	External Dimensions			Internal Dimensions			Sash Opening Dimensions		
	Width	Depth	Height	Width	Depth	Height	Width	Opening 1	Opening 2
MY-DB24	24.0"	24.0"	29.5"	23.5"	20.0"	29.0"	22.0"	10.5"	23"
MY-DB36	36.0"	24.0"	29.5"	35.5"	20.0"	29.0"	34.0"	10.5"	23"
MY-DB48	48.0"	24.0"	29.5"	47.5"	20.0"	29.0"	46.0"	10.5"	23"

TECHNICAL SPECIFICATIONS

Product #	Electrical		Unit Weight*
	Voltage	Power	
MY-DB24	110V	50 W	60 lbs.
MY-DB36	110V	50 W	75 lbs.
MY-DB48	110V	100 W	90 lbs.

Unit weight is an approximation*

www.mystaire.com

Phone: +1 (919) 229-8511 • Toll Free: (877) 328-3912 • Email: contact@mystaire.com