



SCIOLOGEX rotators give gentle but effective mixing of biological samples held in 1.5mL to 50mL tubes, and are used in a variety of applications including immune precipitations, prevention of blood coagulation, latex diagnostics, etc.

## SCI-RD-Pro LCD Digital Rotator

### Features

- Electronic adjustment of speed and time
- LCD display of speed and time
- Speed range from 10 to 70rpm
- Rolling action, gives gentle but effective mixing
- Adjustable mixing angle from 0 to 90°
- A wide variety of accessories are available

## SCI-RD-E Economical Rotator

### Features

- Provides gentle and effective mixing of biological samples
- Adjustable speed range from 0 to 80rpm
- A wide variety of accessories are available




LCD display



Adjustable angle



Wide range of accessories



### Specifications

	SCI-RD-Pro	SCI-RD-E
Motor type	DC motor	DC motor
Tilt angle	0-90°	-
Speed range	10-70rpm	0-80rpm
Speed display	LCD	Scale
Timer display	LCD	Scale
Timer	1-1199min	-
Operation mode	Continuous / timer	Continuous
Voltage	100-240V, 50/60Hz	100-240V, 50/60Hz
Power	40W	20W
Dimension [DxWxH ]	280x210x300mm	300x220x310mm
Permissible ambient temperature and humidity	5-40°C, 80%RH	5-40°C, 80%RH
Protection class	IP21	IP21

## SCI-RL-Pro LCD Digital Rotator

### Features

- Electronic adjustment of speed and time
- LCD display speed and time
- Speed range from 10 to 70rpm
- Provides gentle but effective rolling
- A wide variety of accessories are available



## SCI-RL-E Economical Rotator

### Features

- Provides gentle but effective rolling
- Adjustable speed range from 0 to 80 rpm
- A wide variety of accessories are available



LCD display



Wide range of accessories



### Specifications

	SCI-RL-Pro	SCI-RL-E
Motor type	DC motor	DC motor
Speed range	10-70rpm	0-80rpm
Speed display	LCD	Scale
Timer display	LCD	Scale
Timer	1-1199min	-
Operation mode	Continuous/timer	Continuous
Voltage	100-240V, 50/60Hz	100-240V, 50/60Hz
Power	40W	20W
Dimension[D×W×H]	220×510×260mm	150×530×190mm
Permissible ambient temperature and humidity	5-40°C, 80%RH	5-40°C, 80%RH
Protection class	IP21	IP21


# Tube Rotators | Accessories



## Disk Accessories Type-1

Used for SCI-RD-Pro  
for 1.5mLx60 centrifuge tubes holder

**Cat. No.**  
18900160



## Disk Accessories Type-2

Used for SCI-RD-Pro  
for 15mLx16 centrifuge tubes holder

**Cat. No.**  
18900161



## Disk Accessories Type-3

Used for SCI-RD-Pro  
for 50mLx8 centrifuge tubes holder

**Cat. No.**  
18900162



## Disk support rods

4pcs in a group  
use with disk accessories

**Cat. No.**  
18900140



## Rotisserie Accessory Type-1

Used for SCI-RL-Pro, for 1.5mLx48  
centrifuge tubes held horizontally

**Cat. No.**  
18900142



## Rotisserie Accessory Type-2

Used for SCI-RL-Pro, for 15mLx24  
centrifuge tubes held horizontally

**Cat. No.**  
18900143



## Rotisserie Accessory Type-3

Used for SCI-RL-Pro, for 50mLx24  
centrifuge tubes held horizontally

**Cat. No.**  
18900144



## Rotisserie Accessory Type-4

Used for SCI-RL-Pro, for 1.5mLx32  
centrifuge tubes held vertically

**Cat. No.**  
18900145



## Rotisserie Accessory Type-5

Used for SCI-RL-Pro, for 15mLx16  
centrifuge tubes held vertically

**Cat. No.**  
18900146



## Rotisserie Accessory Type-6

Used for SCI-RL-Pro, for 50mLx16  
centrifuge tubes held vertically

**Cat. No.**  
18900147



## Plastic Clips

Used with SCI-RL-E/sci-RD-E  
Disk/Rotisserie accessory, for  
50mL, 15mL and 1.5mL  
centrifuge tubes





## SCI-T6-Pro LCD Digital Tube Roller

### Features

- LCD displays the speed and time
- Speed range of 10-70rpm
- Six roller design for gentle rolling and easy to clean
- Continuous or timer operation

SCILOGEX tube rollers provide a gentle and highly efficient rolling motions, ideal for mixing blood samples, viscous substances and liquid-solid suspensions. Used in a variety of applications, such as the prevention of blood coagulation, immune precipitation, etc.



LCD display



Timer

## SCI-T6-S Tube Roller

### Features

- Speed range of 0-70rpm
- Six roller design for gentle rolling and easy to clean



### Specifications

	SCI-T6-Pro	SCI-T6-S
Amplitude	24mm	24mm
Motor type	DC motor	DC motor
Rocking motion	Rocking and rolling	Rocking and rolling
Max. load capacity	4kg	4kg
Number of rollers	6	6
Roller size[length]	280mm	280mm
Speed range	10 -70rpm	0-70rpm
Speed display	LCD	Scale
Timer	Yes	-
Timer display	LCD	-
Time setting range	1-1199min	-
Operation mode	Continuous/timer	Continuous
Voltage	100-240V,50/60Hz	100-240V,50/60Hz
Power	30W	25W
Dimension [D×W×H]	450×260×120mm	450×260×120mm
Weight	5.1kg	4.5kg
Permissible ambient temperature and humidity	5-40°C,80%RH	5-40°C,80%RH
Protection class	IP21	IP21